


# SILICONE HEMATURIA CATHETER

## Hematuria Catheter, Straight Whistle Tip:


REF. No.	Size	Balloon	Length	Description
1862-0318	18 FR	50-80 cc	420 mm	Straight whistle tip X-ray opaque line Attached plug Catheter spigot
1862-0320	20 FR			
1862-0322	22 FR	80-100 cc		
1862-0324	24 FR			



1862-03XX

## Hematuria Catheter, Coude Round Tip:


REF. No.	Size	Balloon	Length	Description
1863-0318	18 FR	50-80 cc	420 mm	Coude round tip X-ray opaque line Attached plug Catheter spigot
1863-0320	20 FR			
1863-0322	22 FR	80-100 cc		
1863-0324	24 FR			



1863-03XX

## Hematuria Catheter, Dufour:


REF. No.	Size	Balloon	Length	Description
1864-0318	18 FR	50-80 cc	420 mm	Dufour X-ray opaque line Attached plug Catheter spigot
1864-0320	20 FR			
1864-0322	22 FR	80-100 cc		
1864-0324	24 FR			



1864-03XX

## Hematuria Catheter, Couvelaire:

REF. No.	Size	Balloon	Length	Description
1865-0318	18 FR	50-80 cc	420 mm	Couvelaire X-ray opaque line Attached plug Catheter spigot
1865-0320	20 FR			
1865-0322	22 FR	80-100 cc		
1865-0324	24 FR			



1865-03XX



6 F., No. 29, Sec. 2, Jhongjheng E. Rd.,  
Danshuei Dist, New Taipei City. 251 Taiwan  
Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266  
www.fortunemed.com Mail:info@fortunemed.com

