

DAILY CLEANING SOLUTION KIT
REF 2118 - No Ammonium Biflouride
REF 7118 - with Ammonium Biflouride



INTENDED USE

The “Fluid Path” must be cleaned at the end of every working day to remove protein deposits. DAILY CLEANER? or DAILY RINSE? is the only daily maintenance required, and is essential to optimize trouble-free analyzer operation.

CONTENTS

REF 2118

Daily Cleaner Diluent

1 x 90 mL

HCl (0.1 N) and salts

Daily Cleaner Powder

6 bottles

pepsin (0.50 g)

REF 7118

Daily Cleaner Diluent

1 x 90 mL

HCl (0.2 N) and NH₅F₂ (0.05N) and salts

Daily Cleaner Powder

6 bottles

pepsin (0.50 g)