

SHOCKWAVE | M⁵

Intravascular Lithotripsy (IVL) for Safe and Effective Treatment of Vascular Calcium



MINIMIZE TRAUMA | **OPTIMIZE** OUTCOMES | **SIMPLIFY** PROCEDURES

Broad Applications in Calcified Peripheral Arteries

Common Femoral Disease

Enable Endovascular CFA Treatment



**Pre-Treatment
Angiogram**




**Final Angiogram
(Post IVL + DCB)**

Images courtesy of Bill Miller, MD

Femoropopliteal Disease

Safe and Effective Vessel Prep or Treatment



**Pre-Treatment
Angiogram**




**Final Angiogram
(Post IVL)**

Images courtesy of Marianne Brodmann, MD

Symptomatic Iliac Disease

Preservation of Future Treatment Options



**Pre-Treatment
Angiogram**



**Final Angiogram
(Post IVL)**


Images courtesy of Patricia Furey, MD

TAVI, TEVAR & EVAR Access

Facilitate Delivery of Large-Bore Devices



**Pre-Treatment
Angiogram**



**Final Angiogram
(Post IVL & 23Fr TEVAR)**

Images courtesy of David Caparrelli, MD



IVL GENERATOR AND CONNECTOR CABLE SPECS

Power	110-240 VAC; 50-60Hz; Single Phase, 15A service		IVL Generator CATALOG NUMBER: IVLGCC
Size	11" (28.0 cm) high x 6" (15.2 cm) wide x 11.5" (29.2 cm) deep		
Weight	15 pounds (6.8 kg)		
Output	Proprietary pulse delivery system. Output voltage 3000 volts peak, pulse frequency 1Hz		
Mobility	Product is designed to be mounted to an IV pole		
Length	5 ft (1.53m)		IVL Connector Cable CATALOG NUMBER: IVLCC
Compatibility	Proprietary male key distally designed to connect only to catheter.		
Operation	Lithotripsy pulsing is activated by pushing a button on the Connector Cable.		
Use	Re-usable		

M5 CATHETER SPECS

M5 Branding*	Catalog Number Original Branding*	Diameter (MM)	Length (MM)	Sheath Compatibility	Working Length	Pulses/ Cycle	Cycles	Pulses (Max)	Crossing Profile
M5IVL3560	M732LPBC3560DX1	3.5	60	6F	110	30	10	300	.054
M5IVL4060	M732LPBC4060DX1	4.0	60	6F	110	30	10	300	.057
M5IVL4560	M732LPBC4560DX1	4.5	60	6F	110	30	10	300	.058
M5IVL5060	M732LPBC5060DX1	5.0	60	6F	110	30	10	300	.062
M5IVL5560	M732LPBC5560DX1	5.5	60	6F	110	30	10	300	.064
M5IVL6060	M732LPBC6060DX1	6.0	60	6F	110	30	10	300	.066
M5IVL6560	M732LPBC6560DX1	6.5	60	7F	110	30	10	300	.068
M5IVL7060	M732LPBC7060DX1	7.0	60	7F	110	30	10	300	.073

*Products are equivalent; either Catalog Number may be used when ordering and either product may be received.

**Discover how you can treat calcium more effectively with the
Peripheral Intravascular Lithotripsy (IVL) System.**

Visit shockwavemedical.com or call +1877-77-LITHO (+1 877-775-4846) for more information.