

# ELECTROSURGERY

## NEONATAL

REF	DIMENSIONS	BAKING	ACTIVE SURFACE	REM/not REM	CABLE	CONNECTION
<b>F7805N</b>	<b>88x74 mm</b>	<b>FOAM</b>	<b>36 cm<sup>2</sup></b>	<b>Not REM</b>	<b>No</b>	<b>Tab Connection</b>
F7820N	88x74 mm	FOAM	33 cm <sup>2</sup>	REM	No	Tab Connection
F7820NW/V	88x74 mm	FOAM	33 cm <sup>2</sup>	REM	Yes 300 cm	Valleylab connector



## VERSATILE

REF	DIMENSIONS	BAKING	ACTIVE SURFACE	REM/not REM	CABLE	CONNECTION
F7305	164x117 mm	FOAM	114 cm <sup>2</sup>	Not REM	No	Tab Connection
F7305W/V	150x108 mm	FOAM	114 cm <sup>2</sup>	Not REM	Yes 300 cm	Valleylab connector
F7320	164x117 mm	FOAM	107 cm <sup>2</sup>	REM	No	Tab Connection
F7320W/V	150x108 mm	NOT WOVEN	107 cm <sup>2</sup>	REM	Yes 300 cm	Valleylab connector

