


Perfecting Performance



Occlutech Sizing Balloon




OSB

Occlutech Sizing Balloon

OSB

Changing the dimension

The Occlutech Sizing Balloon is intended for venous access. It is used without a guide catheter and is introduced directly via the guide wire.



Designations:

- (1) Usable length
- (2) X-ray markers
- (3) Soft tip
- (4) Inflation port

Product specifications

Article No.	Balloon Ø Nominal [mm]	Nominal Fill Vol. [ml]	Max. Fill Vol. [ml]	Balloon Ø Max. Fill Vol. [mm]	Balloon Length [mm]	Usable Length [mm]	Guide Wire [inch]
53SB254	25	25	35	27	45	700	0.035
53SB355	35	60	90	40	50	700	0.035



Occlutech International AB
Landskronavägen 2, SE-252 32 Helsingborg, Sweden
order@occlutech.com, www.occlutech.com

