





LOADING CARTS & ACCESSORIES FOR PLUSHER SERIES


PLUSHER WDS & PLUSHER WDC Series of instrument washer disinfectors and **PLUSHER WDB DS** series of bedpan washer disinfectors may be equipped by a large number of accessories. It is also possible to customize each basket depending on the demand in CSSD, laboratory, sluice room or ward. All carts and accessories are made of AISI 304 quality stainless steel in metal parts and necessary silicone or plastic parts resistant to chemicals used or heat if necessary to provide connection jacks.


	Ref. No	313001
	Description	Basic Trolley For Washing Large Objects For PLUSHER WDS 12
	Capacity	-
	Notes	Suitable for washing big objects Useful dimension for level: 550 x 490 x 560H mm. Baskets not included in the application


	Ref. No	313002
	Description	2 Level Washing Trolley W/ Support Grid For PLUSHER WDS 12
	Capacity	-
	Notes	Suitable for washing large objects Useful dimension for level: 490 x 500 x 290H mm. Baskets not included in the application


	Ref. No	313003
	Description	2 Level Washing Trolley W/ Support Grid For PLUSHER WDS 18
	Capacity	-
	Notes	Suitable for washing large objects Useful dimension 1st level: 590 x 720 x 290H mm. Useful dimension 2nd level: 590 x 720 x 170H mm. Baskets not included in the application


	Ref. No	313004
	Description	2 Level Washing Trolley For 4 DIN Baskets For PLUSHER WDC 8
	Capacity	4 x 1/1 DIN baskets
	Notes	Suitable for washing large quantities of standard surgical instruments Useful dimension for level: 500 x 480 x 200H mm. Basket not included in the application


	Ref. No	313005
	Description	2 Level Washing Trolley For 4 DIN Baskets For Plusher WDS 12
	Capacity	4 x 1/1 DIN baskets
	Notes	Suitable for washing large quantities of standard surgical instruments Useful dimension for level: 480 x 500 x 290H mm. Basket not included in the application


	Ref. No	313006
	Description	2 Level Washing Trolley For 6 DIN Baskets For PLUSHER WDS 18
	Capacity	6 x 1/1 DIN baskets
	Notes	Suitable for washing large quantities of standard surgical instruments Useful dimension for level: 590 x 720 x 300H mm. Basket not included in the application


	Ref. No	313007
	Description	3 Level Washing Trolley For 6 DIN Baskets For PLUSHER WDC 8
	Capacity	6 x 1/1 DIN baskets
	Notes	Suitable for washing large quantities of standard surgical instruments Useful dimension for level: 500 x 480 x 120H mm. Basket not included in the application


	Ref. No	313008
	Description	3 Level Washing Trolley For 6 DIN Baskets For PLUSHER WDS 12
	Capacity	6 x 1/1 DIN baskets
	Notes	Suitable for washing large quantities of standard surgical instruments Useful dimension for level: 480 x 500 x 180H mm. Basket not included in the application


	Ref. No	313009
	Description	3 Removable Level Washing Trolley For 6 DIN Baskets For PLUSHER WDS 12
	Capacity	6 x 1/1 DIN baskets
	Notes	Suitable for washing large quantities of standard surgical instruments Useful dimension for level: 480 x 500 x 170H-400H-630H mm. Basket not included in the application


	Ref. No	313010
	Description	3 Level Washing Trolley For 9 DIN Baskets For PLUSHER WDS 18
	Capacity	9 x 1/1 DIN baskets
	Notes	Suitable for washing large quantities of standard surgical instruments Useful dimension for level: 590 x 720 x 170H mm. Basket not included in the application


	Ref. No	313011
	Description	4 Level Washing Trolley For 8 DIN Baskets For PLUSHER WDC 8
	Capacity	8 x 1/1 DIN baskets
	Notes	Suitable for washing large quantities of standard surgical instruments Useful dimension for level: 500 x 480 x 80H mm. Basket not included in the application


	Ref. No	313012
	Description	4 Level Washing Trolley For 8 DIN Baskets Wo/ Side Support For PLUSHER WDS 12
	Capacity	8 x 1/1 DIN baskets
	Notes	Suitable for washing large quantities of standard surgical instruments Useful dimension for level: 480 x 550 x 110H mm. Basket not included in the application


	Ref. No	313013
	Description	4 Level Washing Trolley For 8 DIN Baskets For PLUSHER WDS 12
	Capacity	8 x 1/1 DIN baskets
	Notes	Suitable for washing large quantities of standard surgical instruments Useful dimension for level: 480 x 500 x 120H mm. Basket not included in the application


	Ref. No	313014
	Description	4 Removable Level Washing Trolley For 8 DIN Baskets For PLUSHER WDS 12
	Capacity	8 x 1/1 DIN baskets
	Notes	Suitable for washing large quantities of standard surgical instruments Useful dimension for level: 480 x 500 x 75H-210H-350H-490H-620H mm. Basket not included in the application


	Ref. No	313015
	Description	4 Level Washing Trolley For 12 DIN Baskets For PLUSHER WDS 18
	Capacity	12 x 1/1 DIN baskets
	Notes	Suitable for washing large quantities of standard surgical instruments Useful dimension for level: 590 x 720 x 120H mm. Basket not included in the application


	Ref. No	313016
	Description	4 Removable Level Washing Trolley For 12 DIN Baskets For PLUSHER WDS 18
	Capacity	12 x 1/1 DIN baskets
	Notes	Suitable for washing large quantities of standard surgical instruments Useful dimension for level: 480 x 500 x 110H-285H-460H-630H mm. Basket not included in the application


	Ref. No	313017
	Description	5 Level Washing Trolley For 10 DIN Baskets For PLUSHER WDS 12
	Capacity	10 x 1/1 DIN baskets
	Notes	Suitable for washing large quantities of standard surgical instruments Useful dimension for level: 480 x 500 x 90H mm. Basket not included in the application


	Ref. No	313018
	Description	5 Removable Level Washing Trolley For 10 DIN Baskets For PLUSHER WDS 12
	Capacity	10 x 1/1 DIN baskets
	Notes	Suitable for washing large quantities of standard surgical instruments Useful dimension for level: 480 x 500 x 75H-215H-355H-495H-630H mm. Basket not included in the application


	Ref. No	313019
	Description	5 Level Washing Trolley For 15 DIN Baskets For PLUSHER WDS 18
	Capacity	15 x 1/1 DIN baskets
	Notes	Suitable for washing large quantities of standard surgical instruments Useful dimension for level: 590 x 720 x 90H mm. Basket not included in the application


	Ref. No	313020
	Description	5 Removable Level Washing Trolley For 15 DIN Baskets For PLUSHER WDS 18
	Capacity	15 x 1/1 DIN baskets
	Notes	Suitable for washing large quantities of standard surgical instruments Useful dimension for level: 480 x 500 x 75H-210H-350H-480H-600H mm. Basket not included in the application


	Ref. No	313021
	Description	6 Level Washing Trolley For 12 DIN Baskets For PLUSHER WDS 12
	Capacity	12 x 1/1 DIN baskets
	Notes	Suitable for washing large quantities of standard surgical instruments Useful dimension for level: 480 x 500 x 60H mm. Basket not included in the application


	Ref. No	313022
	Description	6 Level Washing Trolley For 18 DIN Baskets For PLUSHER WDS 18
	Capacity	18 x 1/1 DIN baskets
	Notes	Suitable for washing large quantities of standard surgical instruments Useful dimension for level: 590 x 720 x 60H mm. Basket not included in the application


	Ref. No	313023
	Description	Container Washing Trolley For 2 Containers For PLUSHER WDC 8
	Capacity	2 containers max. dimensions 460 x 300 x 160H mm.
	Notes	Suitable for wash containers Basket not included in the application


	Ref. No	313024
	Description	Container Washing Trolley For 2 Containers For PLUSHER WDS 12
	Capacity	2 containers max. dimensions 300 x 600 x 300H mm.
	Notes	Suitable for wash containers Containers not included in the application


	Ref. No	313025
	Description	Container Washing Trolley For 6 Containers For PLUSHER WDS 18
	Capacity	6 containers max. dimensions 300 x 600 x 150H mm.
	Notes	Suitable for wash containers Containers not included in the application


	Ref. No	313097
	Description	Container Washing Trolley For 3 Containers With 3 Levels For PLUSHER WDS 18
	Capacity	3 containers max. dimensions 600x300x150H mm. & 2 washing levels
	Notes	Suitable for wash containers Useful dimension for level: 590 x 720 x 325H-95H-95H mm. Containers not included in the application Baskets not included in the application


	Ref. No	313026
	Description	Mini Invasive Surgery Washing Trolley For PLUSHER WDC 8
	Capacity	15 + 16 positions (8 luer lock + 8 silicone pipe + 15 guard connections for tubular hollow instruments Ø 3 - 7mm)
	Notes	Suitable for washing laparoscopy and ophthalmology instruments Useful dimension 1st level: 480 x 500 x 75H mm. Injectors included in the application Basket not included in the application


	Ref. No	313027
	Description	Mini Invasive Surgery Washing Trolley For PLUSHER WDS 12
	Capacity	18 + 18 position (9 luer lock + 9 silicon pipe + 18 guard connections for tubular hollow instruments Ø 3 - 7mm)
	Notes	Suitable for washing laparoscopy and ophthalmology instruments Useful dimension 2nd level: 480 x 500 x 80H mm. Useful dimension 1st level: 480 x 500 x 80H mm. Injectors included in the application Baskets not included in the application


	Ref. No	313028
	Description	Mini Invasive Surgery Washing Trolley For PLUSHER WDS 18
	Capacity	78 position (28 luer lock + 18 silicon pipe + 20 guard connections for tubular hollow instruments Ø 3÷7 mm + 12 nozzles for endotracheal pipes)
	Notes	Suitable for washing laparoscopy and ophthalmology instruments Useful dimension 2nd level: 610 x 750 x 80H mm. Useful dimension 1st level: 610 x 750 x 80H mm. Injectors included in the application Baskets not included in the application

	Ref. No	313029
	Description	Anesthetic Washing Trolley For PLUSHER WDC 8
	Capacity	6 connectors for flexible pipe, 6 balloons nozzles, 10 nozzles for endotracheal pipes, space for divers small items, hose clamps, etc.
	Notes	Suitable for washing anesthesia instruments Injectors included in the application Basket not included in the application


	Ref. No	313030
	Description	Anesthetic Washing Trolley For PLUSHER WDS 12
	Capacity	6 connectors for flexible pipe, 6 balloons nozzles, 10 nozzles for endotracheal pipes, space for divers small items, hose clamps, etc.
	Notes	Suitable for washing anesthesia instruments Injectors included in the application Basket not included in the application


	Ref. No	313031
	Description	Anesthetic Washing Trolley For PLUSHER WDS 18
	Capacity	16 connectors for flexible pipe, 12 balloons nozzles, 12 nozzles for endotracheal pipes, space for divers small items, hose clamps, etc.
	Notes	Suitable for washing anesthesia instruments Injectors included in the application Basket not included in the application


	Ref. No	313033
	Description	Shoe Washing Trolley (10 Pairs) For PLUSHER WDC 8
	Capacity	20 rubber shoes (10 pairs)
	Notes	Suitable for washing rubber OR shoes Max Temperature 55°C - Please check on the shoes' datasheet Useful dimension 1st level: 480 x 500 x 465H mm.


	Ref. No	313034
	Description	Shoe Washing Trolley (15 Pairs) For PLUSHER WDS 12
	Capacity	30 rubber shoes (15 pairs)
	Notes	Suitable for washing rubber OR shoes Max Temperature 55°C - Please check on the shoes' datasheet Useful dimension 1st level: 480 x 500 x 80H mm.


7


	Ref. No	313035
	Description	Shoe Washing Trolley (20 Pairs) For PLUSHER WDS 12
	Capacity	40 rubber shoes (20 pairs)
	Notes	Suitable for washing rubber OR shoes Max Temperature 55°C - Please check on the shoes' datasheet Insert 313103 included in the application.


	Ref. No	313036
	Description	Jar (Cylindrical Vessel) Washing Trolley For 4 Units For PLUSHER WDC 8
	Capacity	4 jars \varnothing_{\max} : 225 mm.
	Notes	Suitable for washing cylindrical vessels Useful dimension level: 450 x 500 x 440H mm. Injectors included in the application


	Ref. No	313037
	Description	Jar (Cylindrical Vessel) Washing Trolley For 6 Units For PLUSHER WDC 12
	Capacity	6 jars \varnothing_{\max} : 194 mm
	Notes	Suitable for washing cylindrical vessels Useful dimension 1st level: 490 x 500 x 150H mm. Injectors included in the application


	Ref. No	313038
	Description	Injection Washing Trolley For 70 Sleeve Injectors For Instruments For PLUSHER WDS 12
	Capacity	70 sleeve injectors for instruments \varnothing 3÷12 mm - H_{\max} : 580 mm
	Notes	Suitable for washing laparoscopy, ophthalmology and nail orthopedic surgery instruments Injectors included in the application


	Ref. No	313039
	Description	Inhalation Equipment Washing Trolley For PLUSHER WDC 8
	Capacity	15 injectors \varnothing 8 mm. h=150 mm., 15 injectors \varnothing 8 mm. h=95 mm., 10 injectors \varnothing 4 mm. h=50 mm., 6 cases \varnothing 40 mm. h=60 mm.
	Notes	Suitable for washing inhalation instruments Useful dimension 1st level: 480 x 500 x 180H mm. Injectors included in the application


	Ref. No	313098
	Description	Big Objects Washing Trolley With Removable Level FOR PLUSHER WDS 12
	Capacity	1 lower fix level free, 1 upper extractible level with spray arms.
	Notes	Suitable for washing of big objects, machines' covers, etc. Useful dimension 2nd level: 480 x 540 x 230H mm. Useful dimension 1st level: 480 x 540 x 295H mm. Injectors included in the application Baskets not included in the application


	Ref. No	313040
	Description	3 Level Washing Trolley For 6 Din Baskets W/ Bending For Dental & Ophthalmological Inserts For PLUSHER WDC 8
	Capacity	6 x 1/1 DIN baskets
	Notes	Suitable for washing large quantities of standard surgical instruments Useful dimension for level: 500 x 480 x 80H mm. Basket not included in the application


	Ref. No	313106
	Description	2 Level Dental Washing Trolley (28 Handpieces) For PLUSHER WDC 8
	Capacity	28 dental handpieces
	Notes	Suitable for washing of dental handpieces Useful dimension 1st level: 500 x 480 x 200H mm. Useful dimension 2nd level: 500 x 360 x 200H mm. Basket not included in the application



	Ref. No	313107
	Description	3 Removable Level Dental Washing Trolley For PLUSHER WDS 12
	Capacity	Coupling for turbines, handpieces, contra-angles & 2 x 1/1 DIN baskets
	Notes	Items positioned inside silicon adapters suitable for $\varnothing 14-17$ mm. There is a cartridge inside holder 2 levels with spray arms Basket not included in the application


	Ref. No	313041
	Description	2 Level Washing Trolley For Biberon (Total 4 Baskets Of 313054 & 313055) For PLUSHER WDC 8
	Capacity	2 baskets for biberons and 2 baskets teats or 4 baskets for biberon
	Notes	Suitable for washing of biberons and teats Useful dimension 2nd level: 500 x 480 x 200H mm. Useful dimension 1st level: 500 x 480 x 200H mm. Basket not included in the application


	Ref. No	313042
	Description	3 Level Washing Trolley For Biberon (Total 6 Baskets Of 313054 & 313055) For PLUSHER WDS 12 Sliding Door Models
	Capacity	4 baskets for biberon and 2 baskets teats
	Notes	Suitable for washing of biberons and teats 3 washing levels with spray arms ON 1 leve for teats and ferrules Injectors included in the application Baskets not included in the application


	Ref. No	313100
	Description	3 Level Washing Trolley For Biberon (Total 6 Baskets Of 313054 & 313055) For PLUSHER WDS 12 Manual Door Models
	Capacity	4 baskets for biberon and 2 baskets teats
	Notes	Suitable for washing of biberons and teats 3 washing levels with spray arms ON 1 leve for teats and ferrules Injectors included in the application Baskets not included in the application

	Ref. No	313043
	Description	3 Level Washing Trolley For Biberon (Total 9 Baskets Of 313054 & 313055) For PLUSHER WDS 18
	Capacity	6 baskets for biberon and 3 baskets teats
	Notes	Suitable for washing of biberons and teats Useful dimension 1st level: 590 x 700 x 240H mm. Injectors included in the application Baskets not included in the application


	Ref. No	313044
	Description	Lab. Glassware Washing Trolley W/ 2 Levels For 56 Positions (Wo/ Nozzles) For PLUSHER WDC 8
	Capacity	1 lower fix level with 56 injection positions, 1 upper extractible level with spray arms
	Notes	Suitable for washing glassware Useful dimension 2nd level: 470 x 490 x 160H mm. Useful dimension 1st level: 470 x 500 x 210H mm. To use 56 injectors type 313067 (standard configuration) Nozzles not included in the application
	Ref. No	313045
	Description	Lab. Glassware Washing Trolley W/ 2 Levels For 112 Positions (Wo/ Nozzles) For PLUSHER WDC 8
	Capacity	1 lower fix level with 56 injection positions, 1 upper extractible level with 56
	Notes	Suitable for washing glassware Useful dimension 2nd level: 470 x 490 x 170H mm. Useful dimension 1st level: 470 x 500 x 250H mm. To use 112 injectors type 313067 (standard configuration) Nozzles 1st-2nd level not included in the application


	Ref. No	313099
	Description	Lab. Glassware Washing Trolley W/ 2 Levels For 28 Positions (Wo/ Nozzles) & 1 DIN Basket For PLUSHER WDC 8
	Capacity	1 lower fix level with 28 injection positions and 1 DIN basket, 1 upper extractible level with spray arms
	Notes	Suitable for washing glassware and surgical instruments Useful dimension 2nd level: 470 x 490 x 165H mm. Useful dimension 1st level: 470 x 500 x 215H mm. To use 28 injectors type 313067 (standard configuration) Nozzles and DIN basket 1st level not included in the application


	Ref. No	313046
	Description	Lab. Glassware Washing Trolley W/ 3 Levels For 105 Positions (Wo/ Nozzles) For PLUSHER WDS 12
	Capacity	1 lower fix level with 35 injector positions, 2 upper extractible levels with 35 injector positions each
	Notes	Suitable for washing glassware Useful dimension 3th level: 490 x 510 x 110H mm. Useful dimension 2nd level: 490 x 510 x 130H mm. Useful dimension 1st level: 490 x 510 x 190H mm. To use injectors type 313068 x 35 + 313067 x 70 + 313070 x 105 (standard configuration) Nozzles not included in the application


	Ref. No	313047
	Description	Lab. Glassware Washing Trolley W/ 3 Levels For 35 Positions (Wo/ Nozzles) & 2 Levels For PLUSHER WDS 12
	Capacity	1 lower fix level with 35 injector positions, 2 upper extractible levels spray arms
	Notes	Suitable for washing glassware Useful dimension 3th level: 480 x 500 x 90H mm. Useful dimension 2nd level: 480 x 500 x 90H mm. Useful dimension 1st level: 480 x 500 x 160H mm. To use injectors type 313068 x 35 + 313070 x 35 (standard configuration) Nozzles not included in the application


11


	Ref. No	313048
	Description	Lab. Glassware Washing Trolley W/ 3 Levels For 70 Positions (Wo/ Nozzles) & 1 Level For PLUSHER WDS 12
	Capacity	1 lower fix level free, 2 upper extractible levels with 35 injector each
	Notes	Suitable for washing glassware Useful dimension 3th level: 480 x 500 x 130H mm. Useful dimension 2nd level: 480 x 500 x 150H mm. Useful dimension 1st level: 490 x 500 x 230H mm. To use injectors type 313068 x 70 + 313070 x 70 (standard configuration) Nozzles not included in the application


	Ref. No	313049
	Description	Lab. Large Glassware Washing Trolley W/ 3 Levels For 19 Positions (Wo/ Nozzles) For PLUSHER WDS 12
	Capacity	1 lower fix level with 19 injection positions
	Notes	Suitable for washing large glassware Useful dimension 1st level: 490 x 510 x 610H mm. To use injectors type 313069 x 9 + 313070 x 10 (standard configuration) Nozzles not included in the application


	Ref. No	313050
	Description	Lab. Glassware Washing Trolley W/ 3 Levels For PLUSHER WDS 12
	Capacity	1 lower fix level free, 2 upper extractible levels with spray arms
	Notes	Suitable for washing glassware Useful dimension 3th level: 480 x 540 x 90H mm. Useful dimension 2nd level: 480 x 540 x 110H mm. Useful dimension 1st level: 490 x 500 x 230H mm. Basket not included in the application


	Ref. No	313114
	Description	Da Vinci Instruments Si-Xi Series Washing Trolley For PLUSHER WDS 12
	Capacity	1 lower fix level free, 8 positions for Da Vinci instruments (Si series, Stapler for Si series, Xi series, Stapler for Xi series)
	Notes	Suitable for washing Si-Xi series of Da Vinci instruments Useful dimension 2nd level: 8 x Da Vinci instrument connectors (313116, 313117, 313118 & 313119) Useful dimension 1st level: 2 x 1/1 DIN baskets Baskets & connectors not included in the application


	Ref. No	313115
	Description	Da Vinci Instruments Si-Xi Series Washing Trolley For PLUSHER WDS 18
	Capacity	2 lower fix level free, 10 positions for Da Vinci instruments (Si series, Stapler for Si series, Xi series, Stapler for Xi series)
	Notes	Suitable for washing Si-Xi series of Da Vinci instruments Useful dimension 3rd level: 10 x Da Vinci instrument connectors (313116, 313117, 313118 & 313119) Useful dimension 2nd level: 3 x 1/1 DIN baskets Useful dimension 1st level: 3 x 1/1 DIN baskets Baskets & connectors not included in the application


	Ref. No	313051
	Description	Release Trolley For Manual Door Models Of PLUSHER WDS 12
	Capacity	-
	Notes	Suitable for the transport of washing trolley Loading trolley with block system to load washing material into the washer disinfectant Anti drop system to collect sinking water from load material Wheels able to be blocked Anti upside down system for security transport Ergonomic handle and design


	Ref. No	313052
	Description	Release Trolley For Sliding Door Models Of PLUSHER WDS 12
	Capacity	-
	Notes	Suitable for the transport of washing trolley Loading trolley with block system to load washing material into the washer disinfectant Anti drop system to collect sinking water from load material Wheels able to be blocked Anti upside down system for security transport Ergonomic handle and design


	Ref. No	313053
	Description	Release Trolley For PLUSHER WDS 18
	Capacity	-
	Notes	Suitable for the transport of washing trolley Loading trolley with block system to load washing material into the washer disinfectant Anti drop system to collect sinking water from load material Wheels able to be blocked Anti upside down system for security transport Ergonomic handle and design


	Ref. No	313054
	Description	Biberon Basket For 21 Bottles For PLUSHER Series
	Capacity	21 biberon bottles (\varnothing_{\max} 54 mm.)
	Notes	Suitable for washing biberon Basket dimension: 430 x 185 x 205H mm.


	Ref. No	313055
	Description	Teats Basket For 42 Teats Or Rings For PLUSHER Series
	Capacity	42 teats (\varnothing_{\max} 54 mm.) - 42 rings
	Notes	Suitable for washing teats and rings for biberon Basket dimension: 485 x 245 x 75H mm.


	Ref. No	313056
	Description	Specula Basket For 6 Speculus For PLUSHER Series
	Capacity	6 speculus
	Notes	Suitable for washing speculus Basket dimension: 280 x 175 x 155H mm.


	Ref. No	313057
	Description	Ophthalmic Basket W/ 2 Levels & Lid For PLUSHER Series
	Capacity	-
	Notes	Tray with double level inside and lid Basket dimension: 360 x 165 x 60H mm.


	Ref. No	313058
	Description	Ophthalmic Basket W/ 1 Level For PLUSHER Series
	Capacity	-
	Notes	Tray with only one level and without lid Basket dimension: 360 x 165 x 60H mm.


	Ref. No	313059
	Description	Tray Insert (30 Trays) For PLUSHER WDS 12
	Capacity	30 trays
	Notes	Suitable for washing trays Tray dimension max: 500 x 250 x 25H mm.


	Ref. No	313060
	Description	Tray Insert (4 Trays) For 313001
	Capacity	4 trays
	Notes	Suitable for washing trays Tray dimension max: 560 x 530 x 75H mm.


	Ref. No	313101
	Description	Tray Insert (49 Trays) For 313020 & 313097
	Capacity	28 trays 235x150x27H mm, 21 trays 250x200x51H mm.
	Notes	Suitable for washing trays on 313020 & 313097 Insert dimension: 515x700x185H mm.


	Ref. No	313061
	Description	Dental Insert 1/2 Din For 7 Handpieces For PLUSHER WDS 12 & PLUSHER WDS 18
	Capacity	7 handpieces
	Notes	Suitable for washing dental instruments into baskets


	Ref. No	313062
	Description	Dental Insert For 7 Handpieces (3 X W&H + 4 X Multiflex) For 313040, 313016 & 313020
	Capacity	7 handpieces (3 W&H + 4 Multiflex)
	Notes	Suitable for dental handpieces Dental insert for 313040, 313016 and 313020


	Ref. No	313102
	Description	Dental Insert For 7 Handpieces (3 X W&H + 4 X Nsk) For 313040, 313016 & 313020
	Capacity	7 handpieces (3 W&H + 4 NSK)
	Notes	Suitable for dental handpieces Dental insert for 313040, 313016 and 313020


	Ref. No	313108
	Description	Insert For Instruments and/or Dental Mirors For PLUSHER WDC 8 & PLUSHER WDS 12
	Capacity	32 Instruments / Mirors (8 x Large + 24 x Small)
	Notes	Insert Dimensions: 415x145x185H mm.


	Ref. No	313109
	Description	Tray Insert (9 Trays - 6 Steps of 55 Mm + 3 Steps of 25 Mm) For PLUSHER WDC 8 & PLUSHER WDS 12
	Capacity	9 Trays
	Notes	Suitable for washing trays Insert Dimensions: 455x245x127H mm.


	Ref. No	313063
	Description	Lab. Pipettes Insert For PLUSHER WDS 12
	Capacity	-
	Notes	Suitable for washing pipettes up to 600 mm. Cassette dimensions: 73 x 73 x 610H mm. Basket not included in the application


	Ref. No	313064
	Description	Lab. Insert For 30 Petri Dishes For PLUSHER WDC 8
	Capacity	30 Petri dishes
	Notes	Suitable for washing Petri dishes Basket dimension: 470 x 240 x 77H mm.


	Ref. No	313065
	Description	Lab. Test Tubes Basket (4 Baskets Per Level) For PLUSHER WDS 12
	Capacity	-
	Notes	Suitable for washing test tubes Up to 4 baskets for level Max test tubes height 150 mm.


	Ref. No	313066
	Description	Lab. Guard Connector For Hollow Instruments For PLUSHER Series
	Capacity	-
	Notes	Suitable for hollow instruments For hollow instruments with Ø 3 - 7 mm.


	Ref. No	313067
	Description	Lab. Injector (Nozzle) Size S Size For Lab. Trolleys For PLUSHER Series
	Capacity	-
	Notes	Glassware trolley insert Ø3 mm. L=90 mm.


	Ref. No	313068
	Description	Lab. Injector (Nozzle) Size M Size For Lab. Trolleys For PLUSHER Series
	Capacity	-
	Notes	Glassware trolley insert Possibility to adjust for different type of glassware Ø4 mm. L=160 mm.


	Ref. No	313069
	Description	Lab. Injector (Nozzle) Size L Size For Lab. Trolleys For PLUSHER Series
	Capacity	-
	Notes	Glassware trolley insert Ø6 mm. H=350 mm.


	Ref. No	313070
	Description	Lab. Cap For Lab. Trolleys For PLUSHER series
	Capacity	-
	Notes	-


	Ref. No	313071
	Description	Lab. Injector With Bottles Support For Lab. Trolleys For PLUSHER Series
	Capacity	-
	Notes	Glassware trolley insert with bottles support Ø6 mm. H=120 mm.


	Ref. No	313072
	Description	Lab. Bottles Support For Lab. Trolleys For PLUSHER Series
	Capacity	-
	Notes	Bottle support for glassware trolley insert H=125 mm.


	Ref. No	313073
	Description	Removable Spray Arm For 313016 & 313020
	Capacity	-
	Notes	Removable spray arm for 313016 - 313020


	Ref. No	313074
	Description	Removable DIN Basket Level For 313016 & 313020
	Capacity	-
	Notes	Removable level for DIN basket on 313016 - 313020


	Ref. No	313075
	Description	Removable Grid Level For 313016 & 313020
	Capacity	-
	Notes	Removable level with grid on 313016 - 313020


	Ref. No	313076
	Description	Removable Mini Invasive Insert For 313020 (Needs 313075 As Base)
	Capacity	-
	Notes	Mini invasive insert for 313020 (Needs of the level 313075 as base unit)


	Ref. No	313077
	Description	Shoe Insert (10 Pairs) For 313006, 313016 & 313020
	Capacity	20 rubber shoes (10 pairs)
	Notes	Shoes washing insert for 313006, 313016 and 313020


	Ref. No	313103
	Description	Shoe Insert (10 Pairs) For 313009, 313014 & 313018
	Capacity	20 rubber shoes (10 pairs)
	Notes	Shoes washing insert for 313009, 313014 and 313018


	Ref. No	313104	18
	Description	Containers Insert (3 Containers) For 313020 & 313097	
	Capacity	3 containers 600x300x150H mm. + 3 lids	
	Notes	Suitable for washing containers on 313020 and 313097 Insert dimension: 587x733x212H mm.	


	Ref. No	313105
	Description	Removable Mini Invasive Insert For 313009, 313014 & 313018
	Capacity	31 position (8 luer lock + 7 silicon pipe + 16 guard connections for tubular hollow instruments Ø 3-7 mm.)
	Notes	Suitable for washing laparoscopy instruments on 313009 -313014 - 313018


	Ref. No	313078
	Description	Removable Min Invasive Insert For 313016 (Needs 313075 As Base)
	Capacity	-
	Notes	Mini invasive insert for 313016 (Needs of the level 313075 as base unit)


	Ref. No	313079
	Description	Anesthetic Balloon Insert For 313020
	Capacity	-
	Notes	Anesthetic balloons insert for 313020


	Ref. No	313080
	Description	Anesthetic Balloon Insert For 313016
	Capacity	-
	Notes	Anesthetic balloons insert for 313016


	Ref. No	313081
	Description	Flexible Pipe Insert For 313016 & 313020
	Capacity	-
	Notes	Flexible pipe insert for 313016 and 313020


	Ref. No	313082
	Description	Removable Anesthetic Balloons Level For 313009 & 313014
	Capacity	-
	Notes	Anesthetic balloons level for 313009, 313014


	Ref. No	313083
	Description	Removable Flexible Anesthesia Pipe Insert Level For 313009 & 313014
	Capacity	-
	Notes	Flexible pipe level for 313009 and 313014


	Ref. No	313084
	Description	Removable Level For Flex Clean For 313009, 313034 & 313018
	Capacity	-
	Notes	Compatible with 313009 - 313034 - 313018


	Ref. No	313085
	Description	Breathing Mask Insert (8 Masks) For 313030 For PLUSHER WDS 12
	Capacity	8 breathing masks
	Notes	Suitable for washing breathing masks on 313030


	Ref. No	313086
	Description	Ophthalmological Basket Kit (Flex Clean Tray) W/ 1 Level & Wo/ Lid For PLUSHER Series
	Capacity	-
	Notes	Flex Clean tray with only one level and without lid Basket dimension: 360 x 165 x 60H mm.


	Ref. No	313087
	Description	Ophthalmological Basket Kit (Flex Clean Tray) W/ 2 Levels & Lid For PLUSHER Series
	Capacity	-
	Notes	Flex Clean tray with double level inside and lid Basket dimension: 360 x 165 x 60H mm.


	Ref. No	313116
	Description	Connection Kit For 1 Da Vinci Instrument Si Series For 313114 & 313115
	Capacity	1 Instrument
	Notes	-


	Ref. No	313117
	Description	Connection Kit For 1 Da Vinci Stapler Instrument Si Series For 313114 & 313115
	Capacity	1 Instrument
	Notes	-


	Ref. No	313118
	Description	Connection Kit For 1 Da Vinci Instrument Xi Series For 313114 & 313115
	Capacity	1 Instrument
	Notes	-


	Ref. No	313119
	Description	Connection Kit For 1 Da Vinci Stapler Instrument Xi Series For 313114 & 313115
	Capacity	1 Instrument
	Notes	-


	Ref. No	313088
	Description	Small Items Insert Basket For PLUSHER WDB DS
	Capacity	-
	Notes	Suitable for washing and disinfecting small objects and surgical instruments, with hinged cover. Size 380 x 250 x 130 H mm.


	Ref. No	313089
	Description	Slop Sink With Ringwash System For PLUSHER WDB DS (450X500X640/960 mm.)
	Capacity	-
	Notes	AISI 304 stainless steel sink, pressed in a conical en-bloc unit System is connected to bedpan washer disinfector's pump Dimensions of sink: 280 x 280 x 200H mm. Tilt-up stainless steel grate to sustain buckets Removable front panel with key Wall or floor drain trap connection DN 70 or DN90 Dimensions: 450 x 500 x 640/960H mm.


	Ref. No	313090
	Description	Slop Sink With Ringwash System For PLUSHER WDB DS (450X600X640/960 mm.)
	Capacity	-
	Notes	AISI 304 stainless steel sink, pressed in a conical en-bloc unit System is connected to bedpan washer disinfector's pump Dimensions of sink: 280 x 280 x 200H mm. Tilt-up stainless steel grate to sustain buckets Removable front panel with key Wall or floor drain trap connection DN 70 or DN90 Dimensions: 450 x 600 x 640/960H mm.


	Ref. No	313091
	Description	Slop Sink With Shelf & Ringwash System For PLUSHER WDB DS (450X500X640/960 mm.)
	Capacity	-
	Notes	AISI 304 stainless steel sink, pressed in a conical en-bloc unit System is connected to bedpan washer disinfector's pump Dimensions of sink: 280 x 280 x 200H mm. Tilt-up stainless steel grate to sustain buckets Removable front panel with key Wall or floor drain trap connection DN 70 or DN90 Dimensions: 450 x 500 x 640/960H mm.

	Ref. No	313092
	Description	Holder Designed For Japanese System For PLUSHER WDB DS
	Capacity	3 urinals or 1 urinal + 1 bedpan + cover or 1 commode + cover
	Notes	-

	Ref. No	313093
	Description	Holder Designed For Austrian System For PLUSHER WDB DS
	Capacity	3 urinals or 1 urinal + 1 bedpan + cover or 1 commode + cover
	Notes	-

	Ref. No	313094
	Description	Holder Designed For German System For PLUSHER WDB DS
	Capacity	3 urinals or 1 urinal + 1 bedpan + cover or 1 commode + cover
	Notes	-

	Ref. No	313095
	Description	Holder Designed For Italian System For PLUSHER WDB DS
	Capacity	3 urinals or 1 urinal + 1 bedpan + cover or 1 commode + cover
	Notes	-

	Ref. No	313096
	Description	Holder Designed For French System For PLUSHER WDB DS
	Capacity	3 urinals or 1 urinal + 1 bedpan + cover or 1 commode + cover
	Notes	-

Our company has various designs and several technical variants. There may be differences between the same product types and various versions. This document is intended for the use of Distributors/Dealers/Customers of TBT Inc. Reproduction in whole or in part by any other party is prohibited. Specifications are given for reference only and are subject to change without any previous notice. Accessory/Product images are for illustrative purposes only. rev.EN 005-2025 © 2015, TBT Inc. All rights reserved.