

Chromotrope 2R / Aniline Blue reagent

Blue counterstain reagent, component of BioGnost's Gomori Trichrome kit.

100 mL CABR-OT-100

Chromotrope 2R / Fast Green F.C.F. reagent

Alternative green counterstain reagent, usually used with Gomori Trichrome kit.

100 mL CFG-OT-100

Colloidal Iron, stock solution

Reagent for binding ferric ions to acid mucopolysaccharides and proteoglycans.

100 mL KOL-OT-100

Congo Red Puchtler reagent

For use with Congo Red Puchtler kit.

100 mL CRPR-OT-100

Congo Red reagent

For use with Congo Red Highman kit.

100 mL CR-OT-100



Cresyl Violet, solution

For use with Luxol Fast Blue kit.

100 mL CV-OT-100

Developer stock solution

For use with Bielschowsky kit.

50 mL STR-OT-50

Differentiation reagents based on acetic acid

For use with special staining kits in histopathology.

Reagent for differentiation in Gomori Trichrome kit	100 mL	RDG-OT-100
Acetic acid, 0.5% solution	100 mL	OK05-OT-100
Acetic acid, 1% solution	100 mL 500 mL	OK1-OT-100 OK1-OT-500
Acetic acid, Colloidal Iron (A) solution ●	250 mL	OKA-OT-250
Acetic acid, Colloidal Iron (B) solution	100 mL	OKB-OT-100

Differentiation reagents based on hydrochloric acid

For use with special staining kits in histopathology.

HCL reagent, A.F.O.G. ●	100 mL	HCLAF-OT-100
HCL reagent, Alcian Blue	100 mL	HCLAB-OT-100
HCL reagent, Colloidal Iron ●	100 mL	HCLK-OT-100
HCL reagent, Feulgen ●	100 mL	HCLF-OT-100
HCL reagent, HemoGnost Perls	100 mL	HCLH-OT-100
HCL reagent, HemoGnost Perls	250 mL	HCLH-OT-250
HCL reagent, P.A.S.	100 mL	HCLP-OT-100
HCL reagent, P.A.S.	500 mL	HCLP-OT-500
HCL reagent, Weigert-Van Gieson ●	100 mL	HCLW-OT-100

Ehrlich aldehyde reagent

Reagent for detecting urobilinogen in urine.

100 mL EAR-OT-100

Eosin-nigrosin vital

Fast detection (one-step detection) of sperm vitality and visualization of dead and living sperm cells with one reagent. Simple, easy and fast method for semen analysis.

30 mL ENV-30
100 mL ENV-100

Eosin Y, isotonic solution

Aqueous solution of Eosin Y used with Eosin and nigrosin vital kit.

100 mL EIO-OT-100

ErioGnost reagent

Acid fast reagent for nuclear staining, comparable to hematoxylin staining. Ecologically acceptable synthetic replacement for hematoxylin reagents.

100 mL ERS-OT-100
500 mL ERS-OT-500
1000 mL ERS-OT-1L

Fast Green F.C.F. reagent

Counterstain used with Masson-Goldner Trichrome and Paraldehyde Fuchsin staining kits.

100 mL FGR-OT-100
500 mL FGR-OT-500

Fast Green F.C.F. contrast reagent

Counterstain used in Grocott and Feulgen staining kits.

100 mL FGKR-OT-100

Field, solution A

Blue component of Field kit.

100 mL FDA-RTU
500 mL FDA-OT-500

Field, solution B

Red component of Field kit.

100 mL FDB-RTU
500 mL FDB-OT-500

Ferri reagent, Verhoeff B

For use with Verhoeff kit.

100 mL FRV-OT-100

Ferri reagent, Weigert B

Iron(III) nitrate solution in diluted hydrochloric acid, a component of Hematoxylin W kit.

100 mL FR-OT-100
500 mL FR-OT-500
1000 mL FR-OT-1L

Fuchsin Acid reagent

For use with Mallory Trichrome kit.

100 mL FA-OT-100