

# Sofia® Catheter

Access Devices



## Sofia® Catheter

Soft Torqueable Catheter Optimized For Intracranial Access

1 per box and includes shaping mandrel and introducer sheath

Product Code	Distal / Proximal ID		Distal OD			Proximal OD			Useable Lengths (cm)
	(in)	(mm)	(Fr)	(in)	(mm)	(Fr)	(in)	(mm)	
DA5115ST	0.055	1.40	5.0	0.0680	1.73	5.0	0.0680	1.73	115
DA5125ST	0.055	1.40	5.0	0.0680	1.73	5.0	0.0680	1.73	125

## Sofia® Plus Catheter

Soft Torqueable Catheter Optimized For Intracranial Access

1 per box and includes shaping mandrel and introducer sheath

Product Code	Distal / Proximal ID		Distal OD			Proximal OD			Useable Lengths (cm)
	(in)	(mm)	(Fr)	(in)	(mm)	(Fr)	(in)	(mm)	
DA6115ST	0.070	1.78	6.0	0.0815	2.07	6.0	0.0825	2.10	115
DA6125ST	0.070	1.78	6.0	0.0815	2.07	6.0	0.0825	2.10	125
DA6131ST	0.070	1.78	6.0	0.0815	2.07	6.0	0.0825	2.10	131

### INDICATIONS FOR USE:

The SOFIA® Plus and SOFIA® Distal Access Catheter is indicated for general intravascular use, including the neuro and peripheral vasculature. The SOFIA® Plus and SOFIA® Distal Access Catheter can be used to facilitate introduction of diagnostic or therapeutic devices. The SOFIA® Plus and SOFIA® Distal Access Catheter is not intended for use in coronary arteries.