ZONERAY AMALGAMATOR SERIES

Figure "8" rotation

The utility model patent No: ZL2005200149421.1



HL-AH(G10) Economical!



FEATURE:

- · With digital LCD display
- · Adjust time and speed freely (Speed: 2800~5000rpm, Time: 1~99s)
- · Time setting with memory function
- . When the lid is open, the motor will stop working
- Ergonomic design
- · Noiseless, very stable operation
- Automatic alarm to remind
- Easy maintenance

ZONERAY Amalgamator Comprehensive performance comparison

Model	G5	G6	G8 (High-end)	G10(Economical)
Picture				Ba
Speed	4200rpm	H:4200rpm L:2800rpm	2800~5000rpm	2800~5000rpm
Time range	4~60s	4~60s	1~99s	1~99s
G.W.	2.8KGS	2.8KGS	2.5KGS	2.8KGS
Usage	Well mix Dental Amalgam Capsule, Glass Ionomer Capsule, 3M Relyx Unicem Capsule, and so on.			

www.zoneray.com

www.zoneray.net

http://zoneray.en.alibaba.com





HANGZHOU ZHONGRUN MEDICAL INSTRUMENT CO., LTD

Declaration of Conformity

MANUFACTURER: HANGZHOU ZHONGRUN MEDICAL INSTRUMENT CO., LTD ADDRESS: ROOM 602-1QIUSHI BUILDING, NO. 6 JIAOGONG ROAD, HANGZHOU, 310012, CHINA TEL: 0086 571 88210205, FAX: 0086 571 88211825

We, the manufacturer, herewith declare that the products

CURING LIGHT UNIT AMALGAMATOR DENTAL SPOT WELDER ALGINATE MIXER

meet the provisions of Directive 93/42/EEC which apply to them.

The medical device has been assigned to class I according to Annex IX of the Directive 93/42/EEC. It bears the mark

 ϵ

Conformity assessment procedure: Annex VII of Directive 93/42/EEC

This Declaration of conformity is valid in connection with the release document for the respective batch of produced devices.

The above mentioned declaration of conformity is exclusively under the responsibility of

Company: HANGZHOU ZHONGRUN MEDICAL INSTRUMENT CO., LTD Address: 6B QIUSHI BUILDING, NO. 6 JIAOGONG BOAD, HANGZHOU, 310012, CHINA

HANGZHOU, August 28, 2018

Place, date

Stamp, Signature