

HOCATT™

Product of Signature Health LTD

Hyperthermic Ozone Carbonic Acid Transdermal Technology



**tablet and tablet holder not included.*

MODEL: HOCATT Plus

HOCATT Modalities:
<ul style="list-style-type: none"> • Whole Body Hyperthermia (Sauna): <ul style="list-style-type: none"> - Steam - Far Infrared • CO₂ / Carbonic Acid Sauna • Ozone (Transdermal + Auxiliary Functions) <ul style="list-style-type: none"> - Ozone Sauna (full-body transdermal) - Cupping + Breast Cupping (localized transdermal) - Vaginal Insufflation (systemic) - Ozone Water (for drinking) • Frequency Specific Microcurrents • Concentrated Oxygen Breathing <i>(for passive Exercise with Oxygen)</i> • Photon Light & Colour Technology • Essential Oil Infusions
External Components:
<ul style="list-style-type: none"> • 1 x 5 LPM Oxygen Concentrator (\pm 93 – 97% concentrated oxygen) • 1 x High-Volume Ozone Destructor <i>(for ozone destruction after session).</i>
Additional Features:
<ul style="list-style-type: none"> • Luxury Smith Chamber <i>(fiberglass chamber with special ABS coating)</i> • Auxiliary Ozone Functions can also be used separately <i>(outside of a sauna session)</i> • 1 x Adjustable Seat • 1 x Adjustable Foot Rest • 1 x No1 Ozone Generator <i>(for ozone sauna)</i> • 1 x No2 Ozone Generator <i>(for simultaneous auxiliary ozone functions)</i> • 1 x FSM Foot Plate Set and 1 x FSM Electrode Set • 1 x Electric CO₂ Regulator <i>(for a CO₂ cylinder - CO₂ cylinder is not included. This needs to be rented.)</i> • 1 x Water Resistant Heart-Rate Monitor • 1 x Full Colour Touch Screen Controller • 1 x Reusable Oxygen Cannula <i>(for concentrated oxygen breathing)</i> • 1 x Oxygen Humidifier <i>(for concentrated oxygen breathing)</i> • 4 x Flow Meters <i>(for setting the oxygen, ozone and CO₂ flow rates)</i>
Accessories and Consumables:
<ul style="list-style-type: none"> • 5 x Neck Capes • 3 x Glass Funnels (for cupping) • 3 x Breast Cupping Sets (S,M,L) • 1 x Bag of Disposable Vaginal Catheters (50 pcs) • 1 x HOCATT Silicone Tube Set • 3 x Water Dispensing Bottles

HOCATT™ EXTERNAL COMPONENTS

5 LPM Oxygen Concentrator:

Concentrates oxygen from the air into $\pm 93 - 97\%$ pure oxygen. Feeds concentrated oxygen to the HOCATT's two ozone generators and also passes oxygen through the humidifier into the oxygen breathing apparatus, so concentrated humidified oxygen may be inhaled during the HOCATT session.



High-Volume External Ozone Destructor:

This high-volume ozone destructor is unique to the HOCATT. It allows for ozone destruction after a HOCATT session by removing steam from the chamber and passing it through multiple activated carbon filters. This prevents a large cloud of the ozone-oxygen steam from being released into the room when the chamber doors open.

