

Electrotomes

for One-Stem Electrodes with Stabilizers

For use with HOPKINS® Forward-Oblique Telescopes 27005 BA/BIA/FA/FIA



Cutting by means of a spring

Movable thumb ring

In rest position the electrode is inside the sheath.



27050 EH

Electrotome

including:

27050 E

2x 27050

2x 27050

2x 277

280

Working Element (also for use with optical urethrotome)

Cutting Loop

Coagulation Electrode

High Frequency Cord


Protection Tube

27050 E

Cutting by means of a finger grip

Movable thumb ring

In rest position the electrode is outside the sheath.



27050 DH

Electrotome

including:

27050 D

2x 27050

2x 27050

2x 277

280

Working Element

Cutting Loop


Coagulation Electrode

High Frequency Cord

Protection Tube

27050 D

Cutting by means of rack and pinion



27050 CH

Electrotome

including:

27050 C

2x 27050

2x 27050

2x 277

280

Working Element

Cutting Loop

Coagulation Electrode

High Frequency Cord



Protection Tube

27050 C

High Frequency Cords see page 73

High Frequency Surgical Units see chapter 15, UNITS

21 Fr., working length 20 cm,
for use with HOPKINS® Straight Forward Telescope 27005AA




27068DK SACHSE **Urethrotome Sheath**, 21 Fr., with channel for FILIFORM bougies and 2 LUER-Lock cones

27068DO **Obturator**, for Urethrotome Sheath 21 Fr., 27068D

27068CD **Telescope Bridge**, 5 Fr., with channel for instruments up to 5 Fr.


27068F **Supplementary Sheath**, sides open, for introduction of a balloon catheter, to slip on Urethrotome Sheath 27068D

27068SK **Supplementary Sheath**, for continuous irrigation and suction, to slip on Urethrotome Sheath 27068D




27005AA

HOPKINS® Straight Forward Telescope 0°, enlarged view, diameter 4 mm, length 30 cm, **autoclavable**, fiber optic light transmission incorporated, color code: green




27005FA

HOPKINS® Telescope 12°, enlarged view, diameter 4 mm, length 30 cm, **autoclavable**, fiber optic light transmission incorporated, color code: black




27005BA

HOPKINS® Forward-Oblique Telescope 30°, enlarged view, diameter 4 mm, length 30 cm, **autoclavable**, fiber optic light transmission incorporated, color code: red



27005CA


HOPKINS® Lateral Telescope 70°, enlarged view, diameter 4 mm, length 30 cm, **autoclavable**, fiber optic light transmission incorporated, color code: yellow



27005EA

HOPKINS® Retrospective Telescope 120°, enlarged view, diameter 4 mm, length 30 cm, **autoclavable**, fiber optic light transmission incorporated, color code: white

21 Fr., working length 20 cm,
for use with HOPKINS® Straight Forward Telescope 27005AA



27050E Working Element



27069K	SACHSE Knife, straight
27069KS	Cold Knife, hook-shaped
27069L	Cold Knife, round
27069M	LUDVIK Knife, straight, serrated
27069P	Cold Knife, ring-shaped



280 Protection Tube, for sterilization and storage of electrodes, loops, curettes and knives

Please note: The knives displayed are not to be used with HF current.
Knives are delivered in packages of 6.