

Direct Ophthalmoscopes



Keeler Direct Ophthalmoscopes

A combination of optical perfection, superb ergonomics and versatile features make Keeler direct ophthalmoscopes the ideal choice. Whatever your requirement from basic examination to detailed diagnosis, there is a Keeler ophthalmoscope for you.

The next couple of pages will guide your choice but first read on to discover the benefits shared by all Keeler ophthalmoscopes as standard.

Key features

Keeler optical system

The more you see the better the diagnosis. Keeler's wide angle light beam and individual mirror alignment ensures that whatever the patient's pupil size, you will see the clearest view of the retina.

Comfort and fit

With curved surfaces and subtle patient side marking, the contoured, curved ophthalmoscope head guarantees an optimal viewing angle and the smooth line surfaces ensure a patient friendly non-threatening instrument.

Slimline handles

Keeler hand helds are supplied on a slimline, lightweight robust handle.

• Light, tough and built to last

Keeler ophthalmoscopes are built from resilient, lightweight materials to give you many years of trouble free service. High impact resistant covers combined with a strong, rigid chassis protect the carefully aligned optics.

A clearer view

Controllable light for bright, clear images

The bright white LED, Halogen or Xenon illumination is fully adjustable from zero. This ensures the comfort of your patients and the best possible diagnosis for you.

Reflex free images

Individually controlled optics reduce the effect of corneal reflex whilst minimising shadows and providing clearer retinal visualisation.





Apertures, graticules and filters

WIDE ANGLE - Illuminates the largest area of fundus for the best general diagnosis through a dilated pupil.

INTERMEDIATE - Easier access through undilated pupils and in peripheral examination. Ideal for paediatric examination.

MACULAR - For examination of the macular region of the fundus where a larger beam would create excessive pupillary reaction or patient discomfort.

SLIT - To determine retinal elevations and depressions, can also be used to assess anterior chamber depth.

GLAUCOMA - Projects a graticule onto the retina to assess the optic cup/disc ratio as an aid to glaucoma diagnosis.

FIXATION CROSS - Projects a graticule onto the retina for assessment of the degree and direction of eccentric fixation.

RED FREE AND FIXATION CROSS - For vessel examination with a fixation cross for assessment of eccentric fixation.

RED FREE - Used to examine blood vessels in fine detail. Blood vessels are silhouetted black against a dark green background.

COBALT BLUE - Use with fluorescein dye for the examination of corneal scars and abrasions.



LED Light Source

Keeler Direct Ophthalmoscopes are now available with a LED light source, giving bright, white illumination and 3 x longer lasting battery life.







3 x Longer lasting battery life

Battery Life	Battery will last 3 x longer with the LED light source than the existing filament bulb at maximum brightness before requiring a re-charge.	
Patch Quality	Homogenous, sharp, clear patch.	
Colour Temperature	3000K which equates to a neutral coloured patch.	
CRI	High-CRI therefore colours are rendered more accurately.	
Electronics	Uses the latest in LED technology.	
Compatibility	The LED light source is compatible with the Specialist, Professional, Practitioner and Standard Keeler Ophthalmoscope heads and the 3.6V lithium ion rechargeable handle and the 2.8V primary cell handle.	



Specialist

We recognise that illumination is of fundamental importance to good diagnosis. The re-engineered optics give 30% more illumination and because it is longer lasting, that means fewer bulb changes giving greater peace of mind.

- The World's finest Direct Ophthalmoscope
- Lightweight and well balanced

Almost 50% lighter, the Specialist combined with the slimline handle is perfectly balanced.

Expanded dioptre range

The most comprehensive lens range available in a direct ophthalmoscope, +44D to -45D lens range guarantees optimal diagnosis in any situation.



Key features

• 30% more illumination

Re-engineered optics give 30% more illumination to get through, for example, cloudy media and cataracts.

Slimmer optics

Slimmer and enclosed viewing optics enable you to get closer to the eye providing a wider field of view.

Premium hand made optics

The Morton Lens Track system, made up of individually ground glass lenses ensures an optimal view.

New Cobalt Blue filter

The new blue filter can be used in conjunction with fluorescein dye for the detection and examination of corneal scars and abrasions.

Pupillometer gauge

The pupillometer disc allows quick and easy estimation of pupil size.

Part Numbers

Specialist Ophthalmoscopes

1132-P-1157 3.6V Specialist Ophth (rechargeable) in case

1132-P-1149 2.8V Specialist Ophth in case

1132-P-1181 3.6V LED Specialist Ophth (rechargable) in case

1132-P-1173 2.8V LED Specialist Ophth in case

Accessories

1011-P-7034 Spec/Vista halogen bulb 3.6V Pkt2

1011-P-7042 Spec/Vista halogen bulb 2.8V Pkt2

1011-P-7229 Spec/Vista LED bulb 2.8V/3.6V Pk1



Professional

The Professional Ophthalmoscope has been designed to ensure ultimate control. Lenses, graticules and filters can all be introduced at the touch of a fingertip without the need to pull away from the patient.

Positive action maywheel and comprehensive lens range

The positive action, easy access, maywheel not only allows you to dial up the required lens quickly, but also gives you total control of where you are in the lens range making each exam quicker and easier.

The extensive range of +29D to -30D in single dioptre steps ensures easy examination of all ocular structures and covers a wide range of patient/user prescriptions. A quick step swing of +/- 20D helps speed up diagnosis.







Key features

Choice of illumination

The Professional now has a choice of Xenon or LED illumination for a brighter and whiter long lasting light. Giving you 3 x longer battery life.

Swing over filter

The red free filter for enhanced vessel examination can be used in conjunction with any of the 6 diagnostic beams.

Brow rest

The integral soft brow rest allows for easy positioning and protects your lenses if you are a spectacle wearer.



Part Numbers

Professional Ophthalmoscopes

1128-P-1014 3.6V Professional Ophth in case
1128-P-1015 3.6V LED Professional Ophth in case
1128-P-1013 2.8V Professional Ophth in case
1128-P-1008 2.8V LED Professional Ophth in case

Accessories

 1011-P-7114
 Std/Pract/Prof Xenon Bulb 3.6V Pkt2

 1011-P-7106
 Std/Pract/Prof Xenon Bulb 2.8V Pkt2

 1011-P-7210
 Std/Pract/Prof LED Bulb 3.6/2.8V Pkt1



Professional Retinoscope

Optimised for outstanding performance with either a Streak or Spot halogen beam with a simple bulb change. Plus the LED Professional Streak Retinoscope is now available!

• Neutralisation check (Patented)

The only Retinoscope allowing rapid confirmation of neutralisation. You can change the beam from divergent to convergent enabling you to cross check the accuracy of the neutral point, guaranteeing the assessment is correct every time.

Precision magnetic control (Patented)

Improved 360° access to a single control for both rotation and convergence adjustments. The precision movements are achieved by use of rare earth magnet technology.



Key features

Choice of illumination

Choose between 600 lux of brilliant halogen illumination or the LED Retinoscope. Whichever you choose, the Professional guarantees you enough light to refract even the most difficult eyes.

LED Professional Streak Retinoscope
 Gives you better contrast, image quality and
 longer battery life.

• Two position aperture

A simple finger movement allows a choice of 4mm or 1.7mm aperture to optimise brightness or resolution of the retinal reflex.

Part Numbers

Professional Retinoscopes

1302-P-1014 3.6V Professional Streak LED Ret and case1302-P-1011 3.6V Professional Streak Retinoscope in case1302-P-1009 3.6V Professional Spot Retinoscope in case

1302-P-1013 2.8V Professional Streak LED Ret and case

1302-P-1010 2.8V Professional Streak Retinoscope in case 1302-P-1008 2.8v Professional Spot Retinoscope in case

1302-P-1008 2.8v Professional Spot Retinoscope in ca

Accesso

1305-P-7000 Streak Retinoscope LED blister pack

Professional Combi Retinoscope Bulbs

1013-P-7009 Vista/Professional Streak bulb 3.6V Pkt2

1013-P-7011 Professional Spot bulb 3.6V Pkt2

1013-P-7010 Professional Spot bulb 2.8V Pkt2

Vista Retinoscope Bulbs

1013-P-7009 Vista/Professional Streak bulb 3.6V Pkt2

1013-P-7007 Vista Spot halogen bulb 3.6V Pkt2

1013-P-7006 Vista Spot halogen bulb 2.8V Pkt2

Ophthalmoscope/Retinoscope Sets

You can purchase your Ophthalmoscope and Retinoscope in a slimline double handle case and add in a Lithium Double Charger.



Specialist Ophthalmoscope/Retinoscope Set



Part Numbers

Specialist Ophthalmoscope/Retinoscope Sets

1427-P-104 3.6V Streak LED Ret and Specialist Ophth
1427-P-1063 Slimline 3.6V Streak Ret and Specialist Ophth
1427-P-1065 Slimline 3.6V Spot Ret and Specialist Ophth
1427-P-1103 2.8V Streak LED Ret and Specialist Ophth
1427-P-1062 Slimline 2.8V Streak Ret and Specialist Ophth
1427-P-1064 Slimline 2.8V Spot Ret and Specialist Ophth
1427-P-1116 3.6V Specialist & Combi LED only Student Set
1427-P-1119 3.6V Specialist Student Set LED only Volk

Professional Ophthalmoscope/Retinoscope sets

1427-P-102 3.6V Streak LED Ret and Professional Ophth
1427-P-1059 3.6V Streak/Professional Set with 2 handles
1427-P-1061 Slimline 3.6V Spot Ret and Professional Ophth
1427-P-1101 2.8V Streak LED Ret and Professional Ophth
1427-P-1058 Slimline 2.8V Streak Ret & Professional Ophth
1427-P-1060 Slimline 2.8V Spot Ret & Professional Ophth
1427-P-1071 3.6V Magnetic Streak Ret/Professional Ophth and handle
1427-P-1117 3.6V Professional LED only Student Set

1427-P-1118 3.6V Professional Student Set LED only Volk

Lithium Double Charger

1941-P-1368 Lithium Double Charger1999-P-7308 Double Charger Wall Mount Kit

Ophthalmoscope/Retinoscope Wall Units

Keeler's anti-theft designed Ophthalmoscopes and Retinoscopes have a lockable head and handle fitting so you'll always find them where you left them.



If you require a wall unit you will need to order the transformer unit which comes with the 2 x handles and each individual head.

Key features

• Tangle-free

The strong, coiled, tangle-free cable permits instrument use up to 3 metres from the transformer.

Durable

The robust construction is ideal for busy working conditions.

Automated

Automatic illumination when instruments are removed from the cradle. With LED confirmation.



Part Numbers

Wall Units

1964-P-2001 Double Cord Handle Power Supply
1132-P-1050 3.6V Specialist Ophth Head and Bulb
1128-P-1004 3.6V Professional Ophth in Carton
1132-P-1042 2.8V Specialist Ophth Head and Bulb
1128-P-1003 2.8V Professional Ophth in Carton
1302-P-1012 Professional Streak LED
Retinoscope Head
1302-P-1007 3.6V Streak Head in Carton
1302-P-1005 3.6V Spot Magnetic Retinoscope
1302-P-1004 2.8V Spot Magnetic Retinoscope
1302-P-1004 2.8V Spot Magnetic Retinoscope



Accessories

Colour Coded Grips

Add a bit of colour to your diagnostic area with colour coded grips. Available in four colours these new grips are not only fun but functional too!

Available individually or as an assorted pack of four, choose from pink, green, blue or black.





Near Point Cards for Retinoscopy

weak axis compare work	ing examples or
yesterday step turn arrow	s direction prin
accurate glasses testing	sbortcut summar
child error	effect evi
remind list	matured g
price system	flag flow
confirm lift	promote lo
justification television machi	ine optical obser
recognise examine envi	ronment calcul
1774-126H	#3Kee





• Keeler Slimline Handles

Spare handles for your instruments available in 3.6V Lithium rechargeable (red base) or 2.8V Battery (blue base) options.





Part Numbers

Colour Coded Grips

	1901-P-7028	Pink
	1901-P-7036	Green
	1901-P-7044	Blue
	EP29-05365	Black
1901-P-7052 Pack of four assorted		Pack of four assorted
		Colour Coded Grips

Near Point Cards for Retinoscopy

1321-P-7021 Near Point Cards for Retinoscopy

Keeler Slimline Handles

1911-P-1201 3.6V Slimline Handle **1901-P-1101** 2.8V Slimline Handle

Chargers and Upgrade Kits

A two hour fast charge and over 2 hours continuous use means your instrument is always ready for use. Choose between the Lithium-ion Double or Mini Chargers.

• Lithium-ion Double Charger

The compact charger takes two handles and fits neatly into your immediate testing environment, allowing easy access to your hand held instruments. It can be desk or wall mounted.

• Lithium-ion Mini Charger

Can be connected quickly and simply to the base of your Keeler handle, making it the perfect choice for domiciliary visits or practitioners on the move. Your instruments can still be used whilst the battery is charging.



Key features

• Over 2 hours of continuous use

Ideal for long examinations especially while first learning with your instrument.

• LED displays

LEDs indicate when battery life is low, charging or fully charged.

Upgrade Kits available

Keep your existing heads and upgrade to Lithium-ion.

The Double Lithium Upgrade Kit (pictured below) includes 2 x handles and a Lithium-ion Double Charger.

The Mini Charger Upgrade Kit includes 1 x handle and a Lithium-ion Mini Charger.



Part Numbers

Lithium-ion Chargers

1941-P-1368 Lithium Double Charger1941-P-1341 Lithium Mini Charger

Upgrade Kits

1911-P-1228 Lithium Duo Charger and 2 Slimline Handles

1911-P-1236 Slimline Handle and Mini Charger

Lithium Battery

EP39-18918 3.6V Lithium Battery

What's in the box?

Specialist

• Specialist Ophthalmoscope 1132-P-1157 with Halogen Bulb 3.6V SPECIALIST • Lithium Battery OPHTHALMOSCOPE • Spare Halogen Bulb (RECHARGEABLE) IN CASE • User Instructions • Specialist Ophthalmoscope 1132-P-1149 with Halogen Bulb 2.8V SPECIALIST • Spare Halogen Bulb OPHTHALMOSCOPE IN CASE • User Instructions 1132-P-1181 • Specialist Ophthalmoscope with LED Bulb 3.6V LED SPECIALIST OPHTHALMOSCOPE Lithium Battery (RECHARGABLE) Case IN CASE • User Instructions • Specialist Ophthalmoscope 1132-P-1173 2.8V LED SPECIALIST with LED Bulb **OPHTHALMOSCOPE** • Case IN CASE User Instructions

Professional

 Professional Ophthalmoscope with Halogen Bulb Lithium Battery Spare Halogen Bulb Case User Instructions
 Professional Ophthalmoscope with LED Bulb Lithium Battery Case User Instructions
 Professional Ophthalmoscope with Halogen Bulb Spare Halogen Bulb Case User Instructions
Professional Ophthalmoscope with LED Bulb Case User Instructions

Retinoscope

LED Streak RetinoscopeLithium BatteryEyebrow ExtensionCaseUser Instructions
 Halogen Streak Retinoscope Spare Halogen Bulb Lithium Battery Eyebrow Extension Case User Instructions
 Halogen Spot Retinoscope Spare Halogen Bulb Lithium Battery Eyebrow Extension Case User Instructions
LED Streak Retinoscope Eyebrow Extension Case User Instructions
 Halogen Streak Retinoscope Spare Halogen Bulb Eyebrow Extension Case User Instructions
 Halogen Spot Retinoscope Spare Halogen Bulb Eyebrow Extension Case User Instructions

Specialist Ophthalmoscope/ Retinoscope Sets

1427-P-1104 3.6V STREAK LED RETINOSCOPE AND SPECIALIST OPHTHALMOSCOPE	 LED Streak Retinoscope LED Specialist Ophthalmoscope Lithium Battery (x2) Eyebrow Extension Double Handheld Case User Instructions
1427-P-1063 SLIMLINE 3.6V STREAK RETINOSCOPE AND SPECIALIST OPHTHALMOSCOPE	Halogen Streak Retinoscope Halogen Specialist Ophthalmoscope Spare Halogen Streak Bulb Spare Halogen Specialist Bulb Lithium Battery (x2) Eyebrow Extension Double Handheld Case User Instructions
1427-P-1065 SLIMLINE 3.6V SPOT RETINOSCOPE AND SPECIALIST OPHTHALMOSCOPE	Halogen Spot Retinoscope Halogen Specialist Ophthalmoscope Spare Halogen Spot Bulb Spare Halogen Specialist Bulb Lithium Battery (x2) Eyebrow Extension Double Handheld Case User Instructions
1427-P-1103 2.8V STREAK LED RETINOSCOPE AND SPECIALIST OPHTHALMOSCOPE	 LED Streak Retinoscope LED Specialist Ophthalmoscope Eyebrow Extension Double Handheld Case User Instructions
1427-P-1062 SLIMLINE 2.8V STREAK RETINOSCOPE AND SPECIALIST OPHTHALMOSCOPE	Halogen Streak Retinoscope Halogen Specialist Ophthalmoscope Spare Halogen Streak Bulb Spare Halogen Specialist Bulb Eyebrow Extension Double Handheld Case User Instructions
1427-P-1064 SLIMLINE 2.8V SPOT RETINOSCOPE AND SPECIALIST OPHTHALMOSCOPE	Halogen Spot Retinoscope Halogen Specialist Ophthalmoscope Spare Halogen Spot Bulb Spare Halogen Specialist Bulb Eyebrow Extension Double Handheld Case User Instructions

Student LED Specialist Ophthalmoscope/ Retinoscope Sets

itemioscope sets	
1427-P-1116 3.6V SPECIALIST AND COMBI LED ONLY STUDENT SET	LED Streak Retinoscope LED Specialist Ophthalmoscope Spare Halogen Spot Bulb Lithium Battery (x2) Eyebrow Extension User Instructions Dynamic Targets Lithium Double Charger Lithium Mini Charger Lithium Mini Charger Keeler Lens Cloth Student How to Use Ophthalmoscope and Retinoscope Book Romanes Occulder Rule Fixation Bar Budgie Stick Keeler Reusable Pen Torch Double Handheld Soft Case
1427-P-1119 3.6V SPECIALIST STUDENT SET LED ONLY VOLK	LED Streak Retinoscope LED Specialist Ophthalmoscope Spare Halogen Spot Bulb Volk Digital Wild Field Lens Lithium Battery (x2) Eyebrow Extension User Instructions Dynamic Targets Lithium Double Charger Lithium Mini Charger Keeler Lens Cloth Student How to use Ophthalmoscope and Retinoscope Book Romanes Occulder Rule Fixation Bar Budgie Stick Keeler Reusable Pen Torch Double Handheld Soft Case Student City Bag

What's in the box?

Professional Ophthalmoscope/ Retinoscope Sets

• LED Streak Retinoscope 1427-P-1102 LED Professional Ophthalmoscope 3.6V STREAK LED • Lithium Battery (x2) **RETINOSCOPE AND** Eyebrow Extension **PROFESSIONAL** • Double Handheld Case OPHTHALMOSCOPE User Instructions • Halogen Streak Retinoscope • Halogen Professional Ophthalmoscope 1427-P-1059 • Spare Halogen Streak Bulb 3.6V STREAK/ • Spare Halogen Specialist Bulb PROFESSIONAL SET • Lithium Battery (x2) WITH 2 HANDLES • Eyebrow Extension • Double Handheld Case User Instructions • Halogen Spot Retinoscope • Halogen Professional 1427-P-1061 SLIMLINE Ophthalmoscope 3.6V SPOT • Spare Halogen Spot Bulb RETINOSCOPE AND • Spare Halogen Specialist Bulb **PROFESSIONAL** • Lithium Battery (x2) OPHTHALMOSCOPE Eyebrow Extension • Double Handheld Case User Instructions • LED Streak Retinoscope 1427-P-1101 LED Professional 2.8V STREAK LED Ophthalmoscope **RETINOSCOPE AND** Eyebrow Extension **PROFESSIONAL** • Double Handheld Case **OPHTHALMOSCOPE** User Instructions • Halogen Streak Retinoscope • Halogen Professional 1427-P-1058 Ophthalmoscope **SLIMLINE 2.8V** • Spare Halogen Streak Bulb STREAK RETINOSCOPE • Spare Halogen Specialist Bulb AND PROFESSIONAL • Eyebrow Extension OPHTHALMOSCOPE • Double Handheld Case • User Instructions • Halogen Spot Retinoscope • Halogen Professional 1427-P-1060 Ophthalmoscope **SLIMLINE 2.8V SPOT** • Spare Halogen Spot Bulb RETINOSCOPE AND • Spare Halogen Specialist Bulb **PROFESSIONAL** • Evebrow Extension **OPHTHALMOSCOPE**

• Double Handheld Case

User Instructions

• Halogen Professional Ophth

1427-P-1071
• 1 Handle shared between

Retinoscope and Ophth

STREAK RETINOSCOPE/ • Spare Halogen Streak Bulb

PROFESSIONAL
• Spare Halogen Specialist Bulb

OPHTHALMOSCOPE

AND HANDLE
• Eyebrow Extension
• Double Handheld Case
• User Instructions

• Halogen Streak Retinoscope

Student LED Professional Ophthalmoscope/ Retinoscope Sets

• LED Streak Retinoscope

• LED Professional Ophthalmoscope

* Spare Halogen Spot Bulb

Lithium Battery (x2)

Eyebrow Extension

User Instructions

Dynamic Targets

Lithium Double Charger

Lithium Mini Charger

Lithium Mini Charger

Keeler Lens Cloth

Student How to Use Opht
Retinoscope Book

Romanes Occulder Rule

Fixation Bar Budgie Stick

• Student How to Use Ophth and Retinoscope Book • Romanes Occulder Rule • Fixation Bar Budgie Stick • Keeler Reusable Pen Torch • Double Handheld Soft Case • Student City Bag • LED Streak Retinoscope • LED Professional Ophthalmoscope • Spare Halogen Spot Bulb • Volk Digital Wild Field Lens • Lithium Battery (x2) • Eyebrow Extension User Instructions 1427-P-1118 • Dynamic Targets 3.6V PROFESSIONAL • Lithium Double Charger STUDENT SET LED • Lithium Mini Charger **ONLY VOLK** • Keeler Lens Cloth • Student How to Use Ophth and Retinoscope Book • Romanes Occulder Rule • Fixation Bar Budgie Stick • Keeler Reusable Pen Torch • Double Handheld Soft Case • Student City Bag

Wall Units

• Handheld Wall Mount • 2 x Corded Handles 1964-P-2001 • 1.5mm Hexagon Wrench DOUBLE CORD • Ball Ended Hex Wrench 2mm HANDLE • Charger **POWER SUPPLY** • Wall Bracket, Wall Plugs & Screws • User Instructions 1132-P-1050 • Specialist Ophthalmoscope 3.6V SPECIALIST Head with Halogen Bulb OPHTHALMOSCOPE User Instructions **HEAD AND BULB** 1128-P-1004 3.6V PROFESSIONAL • Professional Ophthalmoscope **OPHTHALMOSCOPE** Head with Halogen Bulb IN CARTON 1302-P-1012 • Retinoscope Head with LED **PROFESSIONAL** Streak Bulb STREAK LED • Eyebrow Extension **RETINOSCOPE HEAD** • User Instructions • Retinoscope Head with 1302-P-1007 Streak Bulb 3.6V STREAK HEAD • Eyebrow Extension IN CARTON • User Instructions • Retinoscope Head with 1302-P-1005 3.6V SPOT Spot Bulb MAGNETIC • Eyebrow Extension RETINOSCOPE • User Instructions

Finoff Transilluminator

2903-P-7000
FINOFF • Finoff Transilluminator Head 2.8V
TRANSILLUMINATOR • User Instructions
2.8V

2903-P-7001
FINOFF • Finoff Transilluminator Head 3.6V
TRANSILLUMINATOR • User Instructions
3.6V

Specifications

Electrical Ratings

The ophthalmoscope/retinoscope/otoscope, the power supply (EP29-32777) with its charging dock (1941-P-5289 and 1941-P-5326) together constitute a Medical Electrical System as defined in EN 60601-1.

INPUT MAINS DATA	100-240V – 50/60Hz
POWER SUPPLY RATING	12V : 2.5amps
Operation	
RUN TIME	Maximum 15 minutes ON
KON TIIVIE	Minimum 10 minutes OFF

Instrument Heads and Handles

INPUT VOLTAGE (DC)

3V 2xAA Alkaline Batteries - BLUE 3.75V Lithium Ion Rechargeable Battery - RED (EP39-18918) 3.65V NiMH Rechargeable Battery - BLACK (1919-P-7149)

Dimensions and Weight

SPECIALIST	260 x 35 x 40mm (HxWxD) Weight 270g
PROFESSIONAL	212 x 43 x 28mm (HxWxD) Weight 260g
RETINOSCOPE	240 x 32 x 40mm (HxWxD) Weight 260g
FINOFF TRANSILLUMINATOR	203 x 44 x 29mm (HxWxD) Weight 217g
Classification	
EQUIPMENT	Class II
CLASSIFICATION	Type B protection against shock

Overview and intended use

The Keeler Ophthalmoscope is indicated for examining the posterior segment of the eye referred to as the fundus, in order to aid in the screening and diagnosis of retinal pathology including, but not limited to, diseases such as cataracts, papilloedema, glaucomatous disc cupping, diabetic retinopathy, hypertensive retinopathy and retinal detachments. When set to a high power and magnification, it can also be used to examine the anterior segment of the eye, which includes the eyelids, cornea, sclera, conjunctiva, iris, aqueous, crystalline lens and anterior vitreous.

The Keeler Retinoscope is indicated for the objective assessment of refractive state of the eye. By observing the retinal red reflex, it is also a means to providing information regarding the visual system, for example media and lens opacities, significant ocular aberrations and accommodative status.

The Keeler Finoff Transilluminator is indicated for the transillumination of tissue, and, in particular, scleral transillumination.



More Keeler Instruments

With over 100 years of industry experience, Keeler has become a household name among optometrists and ophthalmologists worldwide.

View our full range at www.keeler.co.uk | www.keelerusa.com

Got a question? - Call us

UK: +44 (0)1753 857177 India: +91 99 3031 1090

USA: +1 800 523 5620 Brazil: +55 11 4302 6053

China: +86-18512119109





