









Nellcor			1732A/B,24/5B; Nellcor N100, N200, N180, N6000, N190, N20, N3000		
	DEC-8	ADP-10017B	Datascope: Expert Datascope: Passport (Nellcor® SpO2);GE :Corometrics 188 Corometrics 129 Critikon Dinamap 9700 Critikon Dinamap 9710Critikon Dinamap 9720 Plus Critikon Dinamap Compact Critikon Dinamap MPS Critikon Dinamap Pro 300 Critikon Dinamap Pro 400 Critikon Dinamap Select;IVY 405T, 405P;M.D.E.:Escort II OPT30 Prism; Nellcor:N100 N200 N180 N-20 NPB-40 NPB-75; Nonin: 8500;8604; Welch Allyn Protocol: Propaq 206-EL, Protocol:Propaq 204-EL , Propaq 202-EL, Propaq 106EL - length 2,4m		DB9M-DB9F
	DEC-4	ADP-10017E	1.2 m		
Nellcor	DOC-10/2008773-001	ADP-10017C	DINAMAP™ ProCare with Nellcor OxiMax option DINAMAP™ ProSeries V2 with Nellcor OxiMax option DINAMAP™ Pro 1000 V3 with Nellcor OxiMax option; Fukuda Kontron 7100 seires; Philips SureSigns VS1; colin: S207 oximax module		3M14P-DB9F
	SCP-10; MC-10	ADP-10017D	DINAMAP™ ProSeries V1 and DINAMAP™ Pro 1000 V2 with Nellcor option		3M14P-DB9F
Nihon Kohden		ADP-10018A			Round 10P-DB9F
		ADP-10018B	Nihon Kohden AL-800PA, Module ,BSM -1101, 7100/7200 , OLV-1100A, liftscope,BSM-8200/1 Lifescope 12 , BSM-8301/2 Lifescope 9 , BSM-88000 Lifescope 14		Round 10P-DB9F
	JL-900P	ADP-10018C	BSM-2300 Lifescope,BSM-4100 Lifescope , BSM-5100 Lifescope ,BSM-9510 Lifescope , BSM-9800 Lifescope		Rectangle 14P-DB9F
	JL-302T	ADP-10018D	use with OPV-1500 Lifescope N; OPV-3100 patient monitor		AMP20 pins-DB9F
	JL-650P	ADP-10018E	Nellcor Oximax module; Use with Nellcor SpO2for modules AY-651/653P, cable length 3 meters,		AMP20 pins-DB9F