

Capsulorhexis

E2002 F

Utrata Capsulorhexis Forceps



SUB
2 mm



Delicate, fine angled tips ideal for grasping delicate capsule during 1.8 mm MICS procedures. And 12 mm shafts. Overall length 87 mm.



E2002 G

Gimbel Capsulorhexis Forceps



Similar in construction to the E2002 Utrata Capsulorhexis Forceps with jaw set to open only 3 mm. Light tension. Flat handle with dull finish. Overall length 80 mm.



E2002 GE

Utrata Capsulorhexis Forceps



Fine angled tips and 12 mm shafts. With iris stop. Overall length 87 mm.



E2002 M

Utrata MICS Capsulorhexis Forceps



SUB
2 mm



Fine angled tips ideal for grasping delicate capsule during 1.8 mm MICS procedures. Angled 12 mm shafts with iris stop allows ease of movement in the anterior chamber. Flat handle with dull finish. Overall length 80 mm.



E2002 R

O'Gawa-Utrata Capsulorhexis Forceps



Fine angled tips with angled 13 mm shafts and iris stop. Round handle for fingertip rotational control. Dull finish. Overall length 112 mm.



E2002 RC

Lehner-Utrata Capsulorhexis Forceps Round Handle



Utrata-like fine angled tips with gently curved shaft for performing capsulorhexis through tunnel incisions. Dull finish. Round Handle. Overall length 112 mm.



Fixation

E1793

Castroviejo Fixation Forceps



Straight shafts. 1 x 2 teeth 0.3 mm.
Overall length 108 mm.



E1795

Castroviejo Fixation Forceps



1 x 2 teeth 0.5 mm. Overall length 109 mm.



E1943 22

Pierse Colibri Forceps Micro



Colibri style 0.1 mm Pierse tips. Flat open handle with dull finish. Overall length 78 mm.



Intraocular Lens

ET2912 CE

IOL Loading Forceps



This IOL loading forceps is designed for use with BLIS injector. The forceps grasps and places the IOL into the IOL cartridge.



2-700

Shepard IOL Insertion Forceps



Curved 18 mm shafts with 1.91 mm wide grasping tip. Flat, serrated cross action handle with dull finish. Titanium. Overall length 108 mm.



65245

IOL Grasping Forceps



Delicate round jaws, angled 50°. Overall length 110 mm.



Forceps

Suturing

63120Y1

Suturing Forceps, Mod. Zuerich



35° angled. 1 x 2 teeth. Overall length 78 mm.



65150G

Barsewisch Anterior Chamber Forceps



45° angled. Tying platform 1.5 length 73 mm.



E1564

Osher Forceps Elschmig Teeth 0.12 mm



The Elschmig teeth of this forceps are ideal for grasping sclera during suturing. Delicate tying platform handles the finest suture material. 0.12 mm, 1 x 2 teeth, 3 mm platform. Overall length 104 mm.



E1686

Troutman Superior Rectus Forceps



Shafts angled 45°, 10 mm from tip to bend. 1 x 2 teeth. Serrated handle with polished finish. Overall length 106 mm.



E1796

Castroviejo Suturing Forceps 0.12 mm



Straight shafts with teeth and 5.5 mm tying platform. Wide serrated handle with polished finish. Teeth 0.12 mm. Overall length 108 mm.



E1796S

Sanders-Castroviejo Suturing Forceps



0.12 mm teeth. 5.5 mm tying platform. Overall length 109 mm.



E1797

Castroviejo Suturing Forceps




1 x 2 teeth. 0.3 mm. 6 mm tying platform. Overall length 109 mm.



Slit

LP7548ADB

Slit knife 2.2 mm angled double bevel




A = 2.2 mm, B = 2.3 mm.




SUB7548ADB

Slit knife 2.2 mm angled double bevel




B = 2.2 mm, C = 2.7 mm.



LP7549ADB

Slit knife 2.4 mm angled double bevel




A = 2.4 mm, B = 2.5 mm.




SUB7549ADB

Slit knife 2.4 mm angled double bevel




B = 2.4 mm, C = 2.9 mm.



LP7550A

Slit knife 2.5 mm angled bevel up




A = 2.5 mm, B = 3.3 mm.



LP7551A

Slit knife 2.65 mm angled bevel up




A = 2.65 mm, B = 3.7 mm.



Slit

LP7551ADB

Slit knife 2.65 mm angled double bevel




A = 2.65 mm, B = 2.8 mm.



LP7559A

Slit knife 2.75 mm angled bevel up




A = 2.75 mm, B = 3.4 mm.



SUB7559A

Slit knife 2.75 mm angled bevel up




B = 2.75 mm, C = 3.3 mm.



LP7552ADB

Slit knife 2.85 mm angled double bevel




A = 2.85 mm, B = 3.2 mm.



SUB7552ADB

Slit knife 2.85 mm angled double bevel




B = 2.8 mm, C = 3.4 mm.



LP7553A

Slit knife 2.85 mm angled bevel up




A = 2.85 mm, B = 3.5 mm.



Stab

LP7545

Stab knife 45 degree




A = 45°, B = 1.8 mm, C = 11.0 mm.



SUB7545

Stab knife 45 degree




B = 1.9 mm, C = 1.8 mm.



Trapezoidal

LP7599

Trapezoidal knife 1.5 - 1.7 mm angled




A = 1.5 mm, B = 1.98 mm, E = 2.20 mm.




SUB7599

Trapezoidal knife 1.5 - 1.7 mm angled

B = 1.7 mm, C = 1.9 mm, D = 1.5 mm,
E = 3.0 mm.

LP7596

Trapezoidal knife 1.5 - 2.0 mm angled



A = 1.5 mm, B = 2.0 mm, E = 2.0 mm.




Scissors

Capsule

24081

Vannas Capsulotomy Scissors, Straight




Overall length 82 mm.



24082

Vannas Capsulotomy Scissors, Curved




85 mm.



24083

Vannas Capsulotomy Scissors, Straight




Model Tuebingen. Overall length 86 mm.



24084

Vannas Capsulotomy Scissors, Curved



Model Tuebingen. Overall length 85 mm.



25030

Ong Capsule Scissors



Angled up. 11 mm from mid-screw to tip. 85 mm.



25040

Uthoff Capsule Scissors



Curved. 83 mm.



25030Y2

Ong Capsule Scissors




Slightly curved tip. Angled.



Capsule

E3386

Vannas Straight Capsulotomy Scissors




Sharp pointed tips. 8 mm from mid-screw to tip.
Overall length 84 mm.



E3387

Vannas Curved Capsulotomy Scissors



6mm blades with sharp pointed tips. 8 mm from mid-screw to tip. Serrated handle with polished finish.
Overall length 84 mm.



E3388

Vannas Angled Capsulotomy Scissors



6mm blades with sharp pointed tips. 9 mm from mid-screw to tip. Serrated handle with polished finish.
Overall length 81 mm.



E3388 GWL

Gills-Welsh Vannas Micro Scissors



Angled. Long model. 11 mm blades with sharp pointed tip. 13 mm from mid-screw to tip. Serrated handle with polished finish. Overall length 86 mm.



E3390 A

Angled Capsule Scissors



Delicate 10 mm blades angled 45°. 110 mm.



Conjunctival

E3315

Troutman Curved Conjunctival Scissors



Blunt tips. 17 mm from mid-screw to tip. With lock.
Overall length 112 mm.

