



Responsible, Efficient, Reliable
PanSW Med Corp
www.panswmed.com

Pansw EKG cable

- 1, Various products compatible with many marks and models;
- 2, effective customers support for compatibilities and after sales service;
- 3, Various connector styles available: banana, clip, DIN, snap, etc;

DIN 3.0 connector



banana 4.0 end



Snap end



Clips end



R

N

F

L

C1

C2

C3

C4

C5

C6

IEC color coding

RA

RL

LL

LA

V1

V2

V3

V4

V5

V6

AHA color coding

Note: (within this catalogue)
1, integrated means complete cable from machine to electrode: trunk cable means cable from machine to the yoke; leadwire means the wires part from yoke to the electrodes;
2, EKG cable: means 10 leads patient cable for ECG machine or defibrillator;

Schiller

Part No	OEM No	Description		Picture	Compatibility
---------	--------	-------------	--	---------	---------------

Part No	OEM No	Description		Connector	Patient group
EKG-10012A	21714 401	Suction Electrode use with DIN 3.0 connector; ball diameter: 26mm; metal electrode diameter: 21 mm;			adult
EKG-10012B		Suction Electrode use with banana 4.0 connector; ball diameter: 26mm; metal electrode diameter: 21 mm;			adult size
EKG-10012C		suction ECG electrode use with both DIN 3.0 and banana 4.0 connector; ball diameter: 26mm; metal electrode diameter: 21 mm;			adult size
EKG-10012D	20298 92-001; Philips 40418 A	Multi-Use Suction Electrode use for snap/banana/DIN/clip connectors; ball diameter: 26mm; metal electrode diameter: 21 mm;			adult size
EKG-10012E		Suction Electrode use with DIN 3.0 connector; ball diameter: 23mm; metal electrode diameter: 17mm;			pediatric size
EKG-10012F		Suction Electrode use with banana 4.0 connector; ball diameter: 23mm; metal electrode diameter: 17mm;			pediatric size
EKG-10012G		suction ECG electrode use with both DIN 3.0 and banana 4.0 connector; ball diameter: 23mm; metal electrode diameter: 17mm;			pediatric size
EKG-10012H		Multi-Use Suction Electrode use for snap/banana/DIN/clip connectors; ball diameter: 23mm; metal electrode diameter: 17mm;			pediatric size

Limb ECG Electrode

Part No	OEM No	Description		Connector	Compatibility
---------	--------	-------------	--	-----------	---------------

EKG-10013A		Limb Electrode for DIN 3.0 connector; clamps size: 30x130mm;			for adult use
EKG-10013B		Limb Electrode for banana 4.0 connector; clamps size: 30x130mm;			for adult use
EKG-10013C		limb Electrode for DIN 3.0 connector and banana 4.0 connector; clamps size: 30x130mm;			for adult use
EKG-10013D	20298 91-001; Philips 40494 E	Multi-Use Limb Electrode; for snap/DIN/banana/clip connectors; universal use; clamps size: 30x130mm;		 	for adult use
EKG-10013E		Multi-Use Limb Electrode; for snap/DIN/banana/clip connectors; universal use; clamps size: 20x85mm			for pediatric use

SIEMENS

EKG-10015A		10 leads; integrated; IEC, DB15M-Din3.0			
EKG-10015B		10 leads; integrated; AHA, DB15M-Din3.0			
EKG-10015C		10 leads; integrated; IEC, DB15M-Banana4.0			
EKG-10015D		10 leads; integrated; AHA, DB15M-Banana4.0			
EKG-10015E		10 leads; integrated; IEC, DB15M-Clip end			Cardiostat 31