



Navigation menu

Medical oxygen flow meter with klik adjustment

Flow meters Medical Gases

MEDICAL OXYGEN FLOW METER WITH CLIK ADJUSTMENT

Product Information:

SERIES RO.400

Medical oxygen flow meter with klik adjustment control for hospital use. Adjustments flow is fast and stable.



Information TECHNICAL SPECIFICATIONS

Performance

Gas	Medical Oxygen
Working pressure	3.5 bar
Flow rate	0-15 lit/min continuous flow
Intervals	1-1,5-2-3-4-6-9-12-15
Flow rate	0-6 lit/min (on request)
Intervals	0,5-1-1,5-2-2,5-3-4-5-6
Accuracy	+/-10% o +/- 0,5 lt.min.
Temperature range	-20° C +50°C

Materials

Body	CW617N-UNI EN 12165 (Brass Chrome plated)
------	--

Connections

Inlet	AFNOR/UNI
Outlet	9/16 UNF maschio (1/4" on request) con nut and hosetail

Size and weight

Dimension	170 x 100 x 110 mm.
Weight	0,780 kg.

DEVICE ACCESSORIES

Sterilisable humidifier

SAFETY DEVICES

Single-tested

TRACEABILITY

Serial number stamped on regulator

PRODUCT LIFETIME

10 years from date of manufacture

DOCUMENTATION

Declaration of conformity and instructions guide

CODE	DESCRIPTION	PACKAGING
RO.400.12039NF	Flow meter with snap adjustment flow regulation 1-15 litres/minute. AFNOR probe	1 Pz.
RO.400.12039UNI	Flow meter with snap adjustment flow regulation 1-15 litres/minute. UNI probe	1 Pz.

CONTACTS

ADDRESS:

Via Leonardo da Vinci, 168 - 20090
Trezzano sul Naviglio (MI)

TELEPHONE:

+39 02 4452075

FAX:

+39 02-4450634

EMAIL:
sanosub@sanosub.com

SAN-O-SUB MBB Srl - Via Leonardo da Vinci, 168 - 20090 Trezzano sul Naviglio (MI) | C.F. e P.I. 12321560968 | REA MI-2653762 | Capitale Sociale int. vers. € 100.000,00 Privacy | Cookies - Sito web realizzato da ExtraWeb