

NAILING SYSTEM
IMPLANTS & INSTRUMENTS





◦ 393 | PFN-2 Nail 130°

MATERIAL - Titanium / Stainless Steel
DIAMETER - 9mm, 10mm, 11mm, 12mm
LENGTH - 170mm, 200mm, 240mm



◦ 394 | PFN-2 Nail 135°

MATERIAL - Titanium / Stainless Steel
DIAMETER - 9mm, 10mm, 11mm, 12mm
LENGTH - 170mm, 200mm, 240mm



◦ 393 | Long PFN-2 Nail 130° (Left & Right)

MATERIAL - Titanium / Stainless Steel
DIAMETER - 9mm, 10mm, 11mm
LENGTH - 320mm, 340mm, 360mm, 380mm, 400mm, 420mm



◦ 394 | Long PFN-2 Nail 135° (Left & Right)

MATERIAL - Titanium / Stainless Steel
DIAMETER - 9mm, 10mm, 11mm
LENGTH - 320mm, 340mm, 360mm, 380mm, 400mm, 420mm



◦ 353 | 4.9mm Interlocking Bolt

MATERIAL - Titanium / Stainless Steel
LENGTH - 22mm to 60mm (2mm Increment)
65mm to 90mm (5mm Increment)



◦ 395 | PFN-2 Blade Screw

MATERIAL - Titanium / Stainless Steel
LENGTH - 75mm to 120mm (5mm Increment)



◦ 355.13 | End Cap For PFN-2 Nail

MATERIAL - Titanium / Stainless Steel



◦ 355.13A | End Cap For PFN-2 Nail Locking Blade

MATERIAL - Titanium / Stainless Steel



- **811.001** | Insertion Handle For PFNA



- **811.003** | Aiming Arm 130° For PFNA
- **811.004** | Aiming Arm 135° For PFNA



- **811.005A** | Aiming Arm For Dynamic Locking
- **811.005B** | Aiming Arm For Static Locking



- **811.006** | Connecting Bolt



- **811.007** | Buttress Compression Nut For PFNA Blade
- **811.008** | Protection Sleeve 16 x 11mm For PFNA Blade
- **811.009** | Drill Sleeve 11 x 3.2mm For PFNA Blade
- **811.010** | Trocar 3.2mm For PFNA Blade



- **811.011** | Protection Sleeve 11 x 8mm
- **811.012** | Drill Sleeve 8 x 4mm
- **811.013** | Trocar 4mm



- **811.014** | Protection Sleeve with Handle 20 x 17mm



- **811.015** | Sleeve 17 x 3.2mm For Guide Wire
- **811.039** | Pointer for Sleeve 17 x 3.2mm



• **811.016** | Inspector For PFNA Blade



• **811.017** | Extraction Screw For PFNA Blade



• **811.018** | Hammer Guide For PFNA Blade



• **811.019** | Pin Wrench L & Key Type




• **811.020** | Radiographic Ruler For Pfna130°,135°



• **811.021** | Direct Measuring Device For 3.2mm Guide Wire



• **811.022** | Key For PFNA Blade With Combination Wrench 11mm



• **811.023** | Measuring Device For Locking Bolt



- **804.010** | Hexagonal Screw Driver 3.5mm Tip Long for 4.9mm Bolt



- **811.025** | Starting Reamer 16.5mm For PFNA Nail



- **811.026** | Reamer 10.3mm For PFNA Blade Screw
- **811.028** | Fixation Sleeve Pfna for Reamer 10.3mm



- **811.027** | Drill bit 10.3mm For PFNA Blade Screw



- **811.029** | Guide Wire 2.8 x 460mm With Hook



- **811.030** | Guide Wire 3mm For PFNA Blade



- **811.031** | Guide Wire 3 x 450mm With Eye



- **811.032** | Hammer Guide



• **811.033** | Connector For Hummer Guide



• **811.034** | 8mm Hexagonal Shaft With Spherical Head



• **811.035** | Awl For PFNA



• **811.036** | Slotted Hammer with Handle



• **811.037** | T Handle Cannulated



• **811.038** | Connecting Bolt For PFNA Aiming Arm



• **724.001** | Flexible Reaming Shaft



• **724** | Bit for Flexible Reaming Shaft
LENGTH - 8mm, 8.5mm, 9mm, 9.5mm, 10mm,
10.5mm, 11mm, 11.5mm, 12mm



• 602.030 | Drill Bits 4.0 mm Length 300mm for Jacobs Chuck



• 811.040 | 8mm Hexagonal Shaft With Spherical Head



• 811.041 | PFN-2 Nail Extractor Attachment




• 804.012 | Stopper Guide Wire 2.5mm Length 950mm



• 804.013 | Guide Wire 2.5mm Length 950mm

AGC.48 PFN-2 Nail Instrument Set



Instruments

811.001	Insertion Handle For PFNA	1 Each
811.003	Aiming Arm 130° For PFNA	1 Each
811.004	Aiming Arm 135° For PFNA	1 Each
811.005A	Aiming Arm For Dynamic Locking	1 Each
811.005B	Aiming Arm For Static Locking	1 Each
811.006	Connecting Bolt	2 Each
811.007	Buttress Compression Nut For PFNA Blade	1 Each
811.008	Protection Sleeve 16 x 11mm For PFNA Blade	1 Each
811.009	Drill Sleeve 11 x 3.2mm For PFNA Blade	1 Each
811.010	Trocar 3.2mm For PFNA Blade	1 Each
811.011	Protection Sleeve 11 x 8mm	1 Each
811.012	Drill Sleeve 8 x 4mm	1 Each
811.013	Trocar 4mm	1 Each
811.014	Protection Sleeve With Handle 20 x 17mm	1 Each
811.015	Sleeve 17 x 3.2mm For Guide Wire	1 Each
811.016	Inspector For PFNA Blade	1 Each
811.017	Extraction Screw For PFNA Blade	1 Each
811.018	Hammer Guide For PFNA Blade	1 Each
811.019	Pin Wrench L & Key Type	1 Each
811.020	Radiographic Ruler For PFNA 130°, 135°	1 Each
811.021	Direct Measuring Device For 3.2mm Guide Wire	1 Each
811.022	Key For Pfna Blade With Combination Wrench 11mm	1 Each
811.028	Fixation Sleeve PFNA for Reamer 10.3mm	1 Each
811.023	Measuring Device For Locking Bolt	1 Each

Instruments

804.010	Hexagonal Screw Driver 3.5mm Tip Long for 4.9mm Bolt	1 Each
811.025	Starting Reamer 16.5mm For PFNA Nail	1 Each
811.026	Reamer 10.3mm For PFNA Blade Screw	1 Each
811.027	Drill bit 10.3mm For PFNA Blade Screw	1 Each
811.029	Guide Wire 2.8 x 460mm With Hook	1 Each
811.030	Guide Wire 3mm For PFNA Blade	1 Each
811.031	Guide Wire 3 x 450mm With Eye	1 Each
811.032	Hammer Guide	1 Each
811.033	Connector For Hammer Guide	1 Each
811.034	8mm Hexagonal Shaft With Spherical Head	1 Each
811.035	Cannulated Bone Awl	1 Each
811.036	Slotted Hammer with Handle	1 Each
811.037	T Handle Cannulated	1 Each
811.038	Connecting Bolt For PFNA Aiming Arm	1 Each
811.040	T Handle Hexagonal Wrench 8.0mm for universal	1 Each
811.041	PFN-2 Nail Extractor Attachment	1 Each
724.001	Flexible Reaming Shaft	1 Each
724	Bit for Flexible Reaming Shaft 8.0mm, 8.5mm, 9.0mm, 9.5mm, 10.0mm, 10.5mm, 11.0mm, 11.5mm, 12.0mm	1 Each
811.039	Pointer for Sleeve 17 x 3.2mm	1 Each
602.030	Drill Bits 4.0mm x 300mm for Jacobs Chuck	1 Each
804.012	Stopper Guide Wire 2.5mm Length 950mm	1 Each
804.013	Guide Wire 2.5mm Length 950mm	1 Each

AGC.49 PFN-2 Nail Implant Set



Implants

393	PFN-2 Nail 9.0mm, 10.0mm, 11.0mm, 12.0mm x 170mm Length 130°	1 Each
393	PFN-2 Nail 9.0mm, 10.0mm, 11.0mm, 12.0mm x 200mm Length 130°	1 Each
393	PFN-2 Nail 9.0mm, 10.0mm, 11.0mm, 12.0mm x 240mm Length 130°	1 Each
394	PFN-2 Nail 9.0mm, 10.0mm, 11.0mm, 12.0mm x 170mm Length 135°	1 Each

Implants

394	PFN-2 Nail 9.0mm, 10.0mm, 11.0mm, 12.0mm x 200mm Length 135°	1 Each
394	PFN-2 Nail 9.0mm, 10.0mm, 11.0mm, 12.0mm x 240mm Length 135°	1 Each
395	PFN-2 Blade Screw 75mm, 80mm, 85mm, 90mm, 95mm, 100mm, 105mm, 110mm, 115mm, 120mm	1 Each
353	4.9mm Interlocking Bolt 28mm, 30mm, 32mm, 34mm, 36mm, 38mm, 40mm, 42mm, 44mm, 46mm, 48mm, 50mm	1 Each

AGC.74 Long PFN-2 Nail 130° & 135° (Left) Set



Implants

393D	Long PFN-2 Nail 9.0mm x 320mm, 340mm, 360mm, 380mm, 400mm, 420mm 130°	1 Each
393E	Long PFN-2 Nail 10.0mm x 320mm, 340mm, 360mm, 380mm, 400mm, 420mm 130°	1 Each
393F	Long PFN-2 Nail 11.0mm x 320mm, 340mm, 360mm, 380mm, 400mm, 420mm 130°	1 Each
394D	Long PFN-2 Nail 9.0mm x 320mm, 340mm, 360mm, 380mm, 400mm, 420mm 135°	1 Each
394E	Long PFN-2 Nail 10.0mm x 320mm, 340mm, 360mm, 380mm, 400mm, 420mm 135°	1 Each
394F	Long PFN-2 Nail 11.0mm x 320mm, 340mm, 360mm, 380mm, 400mm, 420mm 135°	1 Each
395	PFN-2 Blade Screw 75mm, 80mm, 85mm, 90mm, 95mm, 100mm, 105mm, 110mm, 115mm, 120mm	1 Each
353	4.9mm Interlocking Bolt 22mm, 24mm, 26mm, 28mm, 30mm, 32mm, 34mm, 36mm, 38mm, 40mm, 42mm, 44mm, 46mm, 48mm, 50mm, 52mm, 54mm, 56mm, 58mm, 60mm, 65mm, 70mm, 75mm, 80mm	2 Each

AGC.75 Long PFN-2 Nail 130° & 135° (Right) Set



Implants

393D	Long PFN-2 Nail 9.0mm x 320mm, 340mm, 360mm, 380mm, 400mm, 420mm 130°	1 Each
393E	Long PFN-2 Nail 10.0mm x 320mm, 340mm, 360mm, 380mm, 400mm, 420mm 130°	1 Each
393F	Long PFN-2 Nail 11.0mm x 320mm, 340mm, 360mm, 380mm, 400mm, 420mm 130°	1 Each
394D	Long PFN-2 Nail 9.0mm x 320mm, 340mm, 360mm, 380mm, 400mm, 420mm 135°	1 Each
394E	Long PFN-2 Nail 10.0mm x 320mm, 340mm, 360mm, 380mm, 400mm, 420mm 135°	1 Each
394F	Long PFN-2 Nail 11.0mm x 320mm, 340mm, 360mm, 380mm, 400mm, 420mm 135°	1 Each
395	PFN-2 Blade Screw 75mm, 80mm, 85mm, 90mm, 95mm, 100mm, 105mm, 110mm, 115mm, 120mm	1 Each
353	4.9mm Interlocking Bolt 22mm, 24mm, 26mm, 28mm, 30mm, 32mm, 34mm, 36mm, 38mm, 40mm, 42mm, 44mm, 46mm, 48mm, 50mm, 52mm, 54mm, 56mm, 58mm, 60mm, 65mm, 70mm, 75mm, 80mm	2 Each



◦ 361 | Long Proximal Femoral Nail 130° (Left & Right)

MATERIAL - Titanium / Stainless Steel

DIAMETER - 9mm, 10mm, 11mm

LENGTH - 340mm, 360mm, 380mm, 400mm, 420mm, 440mm



◦ 391 | Long Proximal Femoral Nail 135° (Left & Right)

MATERIAL - Titanium / Stainless Steel

DIAMETER - 9mm, 10mm, 11mm

LENGTH - 340mm, 360mm, 380mm, 400mm, 420mm, 440mm



◦ 364 | Proximal Femoral Nail 130°

MATERIAL - Titanium / Stainless Steel

DIAMETER - 9.0mm, 10.0mm, 11.0mm, 12.0mm

LENGTH - 250mm



◦ 397 | Proximal Femoral Nail 135°

MATERIAL - Titanium / Stainless Steel

DIAMETER - 9.0mm, 10.0mm, 11.0mm, 12.0mm

LENGTH - 250mm



◦ 376 | Trochanteric Femoral Nail 130°

MATERIAL - Titanium / Stainless Steel

DIAMETER - 9.0mm, 10.0mm, 11.0mm, 12.0mm

LENGTH - 170mm



◦ 396 | Trochanteric Femoral Nail 135°

MATERIAL - Titanium / Stainless Steel

DIAMETER - 9.0mm, 10.0mm, 11.0mm, 12.0mm

LENGTH - 170mm



◦ 353 | 4.9mm Interlocking Bolt

MATERIAL - Titanium / Stainless Steel

LENGTH - 22mm to 60mm (2 Increment)
65mm to 80mm (5 Increment)



◦ 354 | 6.2mm Interlocking Bolt

MATERIAL - Titanium / Stainless Steel

LENGTH - 55mm to 120mm (5 Increment)



• **365 | 8.0mm Proximal Femoral Screw**

MATERIAL - Titanium / Stainless Steel

LENGTH - 55mm, 60mm, 65mm, 70mm, 75mm, 80mm,
85mm, 90mm, 95mm, 100mm, 105mm, 110mm,
115mm, 120mm



• **355.06 | End Cap For Proximal Femoral Nail**

MATERIAL - Titanium / Stainless Steel



- 807.016 | Insertion Handle For Pfn



- 807.018 | Aiming Arm 130° For PFN
- 807.019 | Aiming Arm 135° For PFN



- 807.020 | Aiming Arm For Dynamic Locking



- 807.021 | Threaded Conical Bolt For Pfn Nail



- 807.003 | Protection Sleeve 12.0/10.0mm Long
- 807.004 | Protection Sleeve 10.0/8.0mm Long



- 807.005A | Guide Wire Sleeve 10.0/2.3mm for 12mm Sleeve
- 807.027 | Guide wire pointer 8.0mm



- 807.006 | Guide Wire Sleeve 8.0/2.0mm for 10mm Sleeve
- 807.026 | Guide wire pointer 6.2mm



- 807.007 | Drill Sleeve 8.0/4.0mm Long



• **811.038** | Connecting Bolt For Pfna Aiming Arm




• **807.014** | Protection Sleeve with Handle



• **804.010** | Hexagonal Screw Driver 3.5mm Tip Long for 4.9mm Screw



• **807.012** | Canulated Screw Driver Long for P.F.N Screws



• **804C.008** | Depth Gauge for Proximal Femoral Nail



• **811.037** | T Handle Cannulated



• **807.025** | 15mm Starting Reamer for PFN Nail with Q/C



• **811.032** | Hammer Guide



- **807.024** | 6.2mm Step Drill Bits with Q/C for P.F.N Screw
- **807.023** | 8.0mm Step Drill Bits with Q/C for P.F.N Screw



- **807.022** | Radiographic Ruler For Pfn



- **811.034** | 8mm Hexagonal Shaft With Spherical Head



- **802.002** | Tissue Protector



- **811.036** | Slotted Hammer 400gm



- **804.011** | Medullary Tube



- **811.035** | Cannulated Bone Awl



- **811.033** | Connector For Hammer Guide



- **807.026** | Guide wire pointer 6.2mm
- **807.027** | Guide wire pointer 8.0mm



- **807.028** | Guide wire Sleeve for Protection Sleeve with Handle
- **807.029** | Pointer For Sleeve 15 x 2.3mm



- **811.022** | Key For Pfna Blade With Combination Wrench 11mm



- **724.001** | Flexible Reaming Shaft



- **724** | Bit for Flexible Reaming Shaft 8.0mm, 8.5mm, 9.0mm, 9.5mm, 10.0mm, 10.5mm, 11.0mm, 11.5mm, 12.0mm



- **807.015** | Guide Wire 2.0mm x 500mm With Threaded Tip



- **807.013** | Guide Wire Simple 2.0mm Length 405mm



- **602.030** | Drill Bits 4.0 mm Length 300mm for Jacobs Chuck



• **807A.013** | Guide Wire Simple 2.3mm Length 500mm
For 8.0mm screw



• **811.040** | T Handle Hexagonal Wrench 8.0mm for universal




• **804.013** | Guide Wire 2.5mm Length 950mm



• **804.012** | Stopper Guide Wire 2.5mm Length 950mm

AGC.47 Instrument for Proximal Femoral Nail 130° & 135° with short PFN and locking Bolt Set



Instruments

807.016	Insertion Handle For Pfn	1 Each
807.018	Aiming Arm 130° For Pfn	1 Each
807.019	Aiming Arm 135° For Pfn	1 Each
807.020	Aiming Arm For Dynamic Locking PFN	1 Each
807.021	Threaded Conical Bolt For Pfn Nail	2 Each
807.003	Protection Sleeve 12.0/10.0mm Long	1 Each
807.004	Protection Sleeve 10.0/8.0mm Long	1 Each
807.005A	Guide Wire Sleeve 10.0/2.3mm for 12mm Sleeve	1 Each
807.006	Guide Wire Sleeve 8.0/2.0mm for 10mm Sleeve	1 Each
807.007	Drill Sleeve 8.0/4.0mm Long	1 Each
811.038	Connecting Bolt For PFN Aiming Arm	1 Each
807.014	Protection Sleeve with Handle	1 Each
804.010	Hexagonal Screw Driver 3.5mm Tip Long for 4.9mm Bolt	1 Each
807.012	Canulated Screw Dirver Long for P.F.N Screws	1 Each
804C.008	Depth Gauge for Proximal Femoral Nail	1 Each
811.037	T Handle Cannulated	1 Each
811.032	Hammer Guide	1 Each
811.033	Connector For Hummer Guide	1 Each
811.036	Slotted Hammer with Handle	1 Each
807.022	Radiographic Ruler For PFN Nail	1 Each
807.023	8.0mm Step Drill Bits with Q/C for PFN Screw	1 Each
807.024	6.2mm Step Drill Bits with Q/C for PFN Screw	1 Each
807.025	15mm Starting Reamer for PFN Nail with Q/C	1 Each
807.026	Guide wire pointer 6.2mm	1 Each
807.027	Guide wire pointer 8.0mm	1 Each
807.028	Guide wire Sleeve for Protection Sleeve with Handle 2.5mm	1 Each
807.029	Pointer for sleeve 15 x 2.5mm	1 Each
811.034	8mm Hexagonal Shaft With Spherical Head	1 Each
811.040	T Handle Hexagonal Wrench 8.0mm for universal	1 Each
811.035	Cannulated Bone Awl	1 Each
802.002	Tissue Protector	1 Each
804.011	Medullary Tube	1 Each
807.015	Guide Wire 2.0mm x 405mm With Threaded Tip	1 Each
807.013	Guide Wire Simple 2.0mm Length 405mm	1 Each
807A.013	Guide Wire Simple 2.3mm Length 500mm For 8.0mm screw	1 Each
602.030	Drill Bits 4.0mm x 300mm for Jacobs Chuck	1 Each
811.022	Key For Pfna Blade With Combination Wrench 11mm	1 Each
804.012	Stopper Guide Wire 2.5mm Length 950mm	1 Each
804.013	Guide Wire 2.5mm Length 950mm	1 Each
724.001	Flexible Reaming Shaft	1 Each
724	Bit for Flexible Reaming Shaft 8.0mm, 8.5mm, 9.0mm, 9.5mm, 10.0mm, 10.5mm, 11.0mm, 11.5mm, 12.0mm	1 Each

Implant

364	Proximal Femoral Nail 9.0mm, 10.0mm, 11.0mm, 12.0 x 250mm 130°	1 Each
397	Proximal Femoral Nail 9.0mm, 10.0mm, 11.0mm, 12.0 x 250mm 135°	1 Each
376	Trochanteric Femoral Nail 9.0mm, 10.0mm, 11.0mm, 12.0mm x 170mm 130°	1 Each
396	Trochanteric Femoral Nail 9.0mm, 10.0mm, 11.0mm, 12.0mm x 170mm 135°	1 Each
353	4.9mm Interlocking Bolt 22mm, 24mm, 26mm, 28mm, 30mm, 32mm, 34mm, 36mm, 38mm, 40mm, 42mm, 44mm, 46mm, 48mm, 50mm, 52mm, 54mm, 56mm, 58mm, 60mm, 65mm, 70mm, 75mm, 80mm	2 Each
354	6.2mm Hip Locking Bolt 55mm, 60mm, 65mm, 70mm, 75mm, 80mm, 85mm, 90mm, 95mm, 100mm, 105mm, 110mm, 115mm, 120mm	1 Each
365	8.0mm Proximal Femoral Screw 55mm, 60mm, 65mm, 70mm, 75mm, 80mm, 85mm, 90mm, 95mm, 100mm, 105mm, 110mm, 115mm, 120mm	1 Each

AGC.10 Long Proximal Femoral Nail Set 130° & 135° (Left)



Implants

361D	Long Proximal Femoral Nail 9.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 440mm 130°	1 Each
361E	Long Proximal Femoral Nail 10.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 440mm 130°	1 Each
361F	Long Proximal Femoral Nail 11.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 440mm 130°	1 Each
391D	Long Proximal Femoral Nail 9.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 440mm 135°	1 Each
391E	Long Proximal Femoral Nail 10.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 440mm 135°	1 Each
391F	Long Proximal Femoral Nail 11.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 440mm 135°	1 Each

AGC.11 Long Proximal Femoral Nail Set 130° & 135° (Right)



Implants

361D	Long Proximal Femoral Nail 9.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 440mm 130°	1 Each
361E	Long Proximal Femoral Nail 10.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 440mm 130°	1 Each
361F	Long Proximal Femoral Nail 11.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 440mm 130°	1 Each
391D	Long Proximal Femoral Nail 9.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 440mm 135°	1 Each
391E	Long Proximal Femoral Nail 10.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 440mm 135°	1 Each
391F	Long Proximal Femoral Nail 11.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 440mm 135°	1 Each



- **358C** | XL Interlocking Nail Tibia 8.0mm - LENGTH - 285mm to 390mm (15mm Increment) - MATERIAL - Titanium / Stainless Steel
- **358D** | XL Interlocking Nail Tibia 9.0mm - LENGTH - 285mm to 390mm (15mm Increment) - MATERIAL - Titanium / Stainless Steel
- **358E** | XL Interlocking Nail Tibia 10.0mm - LENGTH -285mm to 390mm (15mm Increment) - MATERIAL - Titanium / Stainless Steel
- **358F** | XL Interlocking Nail Tibia 11.0mm - LENGTH -285mm to 390mm (15mm Increment) - MATERIAL - Titanium / Stainless Steel



- **369C** | XL Interlocking Nail (Distal Tip) 8.0mm - LENGTH - 275mm to 380mm (15mm Increment) MATERIAL - Titanium / Stainless Steel
- **369D** | XL Interlocking Nail (Distal Tip) 9.0mm - LENGTH - 275mm to 380mm (15mm Increment) MATERIAL - Titanium / Stainless Steel
- **369E** | XL Interlocking Nail (Distal Tip) 10.0mm - LENGTH - 275mm to 380mm (15mm Increment) MATERIAL - Titanium / Stainless Steel



- **359D** | XL Interlocking Nail Femoral 9.0mm - LENGTH - 345mm to 435mm (15mm Increment) MATERIAL - Titanium / Stainless Steel
- **359E** | XL Interlocking Nail Femoral 10.0mm - LENGTH - 345mm to 435mm (15mm Increment) MATERIAL - Titanium / Stainless Steel
- **359F** | XL Interlocking Nail Femoral 11.0mm - LENGTH - 345mm to 435mm (15mm Increment) MATERIAL - Titanium / Stainless Steel
- **359G** | XL Interlocking Nail Femoral 12.0mm - LENGTH - 345mm to 435mm (15mm Increment) MATERIAL - Titanium / Stainless Steel



- **355.02** | End Cap for XL Interlocking Nail Tibia & Femoral MATERIAL - Titanium / Stainless Steel



- **352** | 3.9mm Interlocking Bolt 18mm to 50mm (2mm Increment) MATERIAL - Titanium / Stainless Steel




- **353** | 4.9mm Interlocking Bolt 22mm to 60mm (2mm Increment) 65mm to 80mm (5mm Increment) MATERIAL - Titanium / Stainless Steel



• **805C.001** | XL Gig for Tibia And Femoral Nail



• **805D.001** | Proximal Bar for XL Tibia Nail



• **805E.001** | Proximal Bar for XL Femoral Nail



• **805A.002** | Conical Bolt for Tibia / Femoral XL Nail



• **805.003** | Reamer for Nail Neck 12mm



• **805A.003** | Protection Sleeve 10.0/8.0mm for Tibia / Femoral XL Nail



• **805A.004** | Drill Sleeve 8.0/3.2mm for Tibia / Femoral XL Nail



• **805A.005** | Drill Sleeve 8.0/4.0mm for Tibia / Femoral XL Nail



• **805A.006** | Pointer 8.0mm for Tibia / Femoral XL Nail



• **811.032** | Hammer Guide



• **811.033** | Connector For Hammer Guide



• **811.036** | Slotted Hammer with Handle



• **804A.008** | Depth Gauge for Tibia / Femoral XL Nail



- **804.009** | Hexagonal Screw Driver 2.5mm Tip Long for 3.9mm Bolt
- **804.010** | Hexagonal Screw Driver 3.5mm Tip Long for 4.9mm Bolt



• **804.011** | Medullary Tube



• **804.012** | Stopper Guide Wire 2.5mm Length 950mm



• 804.013 | Guide Wire 2.5mm Length 950mm



• 811.035 | Canulated Bone Awl



• 802.002 | Tissue Protector



• 602.011 | Drill Bits 3.2mm Length 250mm for Jacobs Chuck
• 602.030 | Drill Bits 4.0mm Length 300mm for Jacobs Chuck



• 805C.004 | Distal Bar for XL Tibia Nail (Left & Right)
• 805C.005 | Distal Bar for XL Femoral Nail (Left & Right)



• 805.009 | Holding Bolt for Distal locking Device



• 805.007 | 4mm Tommy Bar



• 724.001 | Flexible Reaming Shaft



- **724** | Bit for Flexible Reaming Shaft 8.0mm, 8.5mm, 9.0mm, 9.5mm, 10.0mm, 10.5mm, 11.0mm, 11.5mm, 12.0mm



- **811.040** | T Handle Hexagonal Wrench 8.0mm for universal



- **805.008** | Tibia / Femoral XL Nail Extractor Attachment



- **811.022** | Key For PFNA Blade With Combination Wrench 11mm

AGC.04 XL Interlocking Tibia Nail Set



Implants

- | | | |
|-------------|---|---------------|
| 358C | XL Interlocking Nail Tibia 8.0mm
Length : 285mm, 300mm, 315mm, 330mm, 345mm,
360mm, 375mm, 390mm | 1 Each |
| 358D | XL Interlocking Nail Tibia 9.0mm
Length : 285mm, 300mm, 315mm, 330mm, 345mm,
360mm, 375mm, 390mm | 1 Each |
| 358E | XL Interlocking Nail Tibia 10.0mm
Length : 285mm, 300mm, 315mm, 330mm, 345mm,
360mm, 375mm, 390mm | 1 Each |

AGC.05 XL Interlocking Femoral Nail Set



Implants

- | | | |
|-------------|--|---------------|
| 359D | XL Interlocking Nail Femoral 9.0mm
Length : 345mm, 360mm, 375mm, 390mm, 405mm,
420mm, 435mm | 1 Each |
| 359E | XL Interlocking Nail Femoral 10.0mm
Length : 345mm, 360mm, 375mm, 390mm, 405mm,
420mm, 435mm | 1 Each |
| 359F | XL Interlocking Nail Femoral 11.0mm
Length : 345mm, 360mm, 375mm, 390mm, 405mm,
420mm, 435mm | 1 Each |

AGC.56 Instrument for XL Interlocking Tibia & Femoral Nail with locking Bolt Set



Instruments

805C.001	XL Gig for Tibia and Femoral Nail	1 Each
805D.001	Proximal Bar for XL Tibia Nail	1 Each
805E.001	Proximal Bar for XL Femoral Nail	1 Each
805A.002	Conical Bolt for Tibia/Femoral XL Nail	1 Each
805.003	Reamer for Nail Neck 12mm	1 Each
805A.003	Protection Sleeve 10.0/8.0mm for Tibia/Femoral XL Nail	1 Each
805A.004	Drill Sleeve 8.0/3.2mm for Tibia/Femoral XL Nail	1 Each
805A.005	Drill Sleeve 8.0/4.0mm for Tibia/Femoral XL Nail	1 Each
805A.006	Pointer 8.0mm for Tibia/Femoral XL Nail	1 Each
811.032	Hammer Guide	1 Each
811.033	Connector For Hummer Guide	1 Each
811.036	Slotted Hammer with Handle	1 Each
804A.008	Depth Gauge for Tibia/Femoral XL Nail	1 Each
804.009	Hexagonal Screw Driver 2.5mm Tip Long for 3.9mm Bolt	1 Each
804.010	Hexagonal Screw Driver 3.5mm Tip Long for 4.9mm Bolt	1 Each
804.011	Medullary Tube	1 Each
804.012	Stopper Guide Wire 2.5mm Length 950mm	1 Each
804.013	Guide Wire 2.5mm Length 950mm	1 Each
811.035	Canulated Bone Awl	1 Each
802.002	Tissue Protector	1 Each

Instruments

602.011	Drill Bits 3.2mm Length 250mm for Jacobs Chuck	1 Each
602.030	Drill Bits 4.0mm Length 300mm for Jacobs Chuck	1 Each
805C.004	Distal Bar for XL Tibia Nail (Left & Right)	1 Each
805C.005	Distal Bar for XL Femoral Nail (Left & Right)	1 Each
805.009	Holding Bolt for Distal Locking Device	1 Each
805.007	4mm Tommy Bar	1 Each
724.001	Flexible Reaming Shaft	1 Each
724	Bit for Flexible Reaming Shaft 8.0mm, 8.5mm, 9.0mm, 9.5mm, 10.0mm, 10.5mm, 11.0mm, 11.5mm, 12.0mm	1 Each
811.040	T Handle Hexagonal Wrench 8.0mm for universal	1 Each
805.008	Tibia / Femoral XL Nail Extractor Attachment	1 Each
811.022	Key For Pfna Blade With Combination Wrench 11mm	1 Each

Implants

352	3.9mm Interlocking Bolt 18mm, 20mm, 22mm, 24mm, 26mm, 28mm, 30mm, 32mm, 34mm, 36mm, 38mm, 40mm, 42mm, 44mm, 46mm, 48mm, 50mm	2 Each
353	4.9mm Interlocking Bolt 22mm, 24mm, 26mm, 28mm, 30mm, 32mm, 34mm, 36mm, 38mm, 40mm, 42mm, 44mm, 46mm, 48mm, 50mm, 52mm, 54mm, 56mm, 58mm, 60mm, 65mm, 70mm, 75mm, 80mm	2 Each



• **363A** | Humeral Nail 6.0mm

MATERIAL - Titanium / Stainless Steel

LENGTH - 180mm to 300mm (20mm Increment)



• **363B** | Humeral Nail 7.0mm

MATERIAL - Titanium / Stainless Steel

LENGTH - 180mm to 300mm (20mm Increment)



• **363C** | Humeral Nail 8.0mm

MATERIAL - Titanium / Stainless Steel

LENGTH - 180mm to 300mm (20mm Increment)



• **350** | 2.9mm Interlocking Bolt

MATERIAL - Titanium / Stainless Steel

LENGTH - 20mm to 40mm (2mm Increment)



• **351** | 3.4mm Interlocking Bolt

MATERIAL - Titanium / Stainless Steel

LENGTH - 16mm to 50mm (2mm Increment)



• **355.05** | End Cap for Humeral Nail

MATERIAL - Titanium / Stainless Steel



• 809.001 | Proximal Gig for Humeral



• 809.002 | Threaded Conical bolt for Humeral Nail



• 809.003 | Sleeve For Compression Device



• 809.004 | Compression Device



• 809.005 | Starting Reamer for Humeral Nail with Q/C 8.0mm



• 809C.003 | Protection Sleeve 10.0/8.0mm



• 809C.006 | Drill Sleeve 8.0/2.7mm for Humeral
• 809C.007 | Drill Sleeve 8.0/2.5mm for Humeral




• 809C.004 | Pointer 8.0mm for Humeral



804.007

|

Extractor-Impactor Rod with Hammer (Set of Two)



804E.008

|

Depth Gauge for Humeral Nail



804.009

|

Hexagonal Screw Driver Long 2.5mm Tip for 2.9 & 3.4mm Bolt



809.008

|

Guide Wire 1.8mm Length 750mm



803.006

|

Bone Awl With Handle Small



602.006

|

Drill Bits 2.5mm x 180mm for Jacobs Chuck

602.024

|

Drill Bits 2.7mm x 230mm for Jacobs Chuck



895.019

|

Spanner 19mm



809.010

|

Head For Threaded Conical bolt



809.009

|

Starting Reamer for Humeral Nail Fix T Handal 8.0mm



Instruments

809.001	Proximal Gig for Humeral	1 Each
809.002	Threaded Conical bolt for Humeral Nail	1 Each
809.003	Sleeve For Compression Device	1 Each
809.004	Compression Device	1 Each
809.005	Starting Reamer for Humeral Nail with Q/C 8.0mm	1 Each
809C.003	Protection Sleeve 10.0/8.0mm for Humeral	1 Each
809C.006	Drill Sleeve 8.0/2.7mm for Humeral	1 Each
809C.007	Drill Sleeve 8.0/2.5mm for Humeral	1 Each
809C.004	Pointer 8.0mm for Humeral	1 Each
804.007	Extractor-Impactor Rod with Hammer (Set of Two)	1 Each
804E.008	Depth Gauge for Humeral Nail	1 Each
804.009	Hexagonal Screw Driver 2.5mm Tip Long for 3.9mm Bolt	1 Each
809.008	Guide Wire 1.8mm Length 750mm	1 Each
803.006	Small Bone Awl Diamond Pointed	1 Each
602.006	Drill Bits 2.5mm x 180mm for Jacobs Chuck	1 Each
602.024	Drill Bits 2.7mm x 230mm for Jacobs Chuck	1 Each
809.009	Starting Reamer for Humeral Nail Fix T Handal 8.0mm	1 Each
809.010	Head For Threaded Conical bolt	1 Each
895.019	Spanner 19mm	1 Each

Implats

363A	Humeral Nail 6.0mm x 180mm, 200mm, 220mm, 240mm, 260mm, 280mm, 300mm	1 Each
363B	Humeral Nail 7.0mm x 180mm, 200mm, 220mm, 240mm, 260mm, 280mm, 300mm	1 Each
363C	Humeral Nail 8.0mm x 180mm, 200mm, 220mm, 240mm, 260mm, 280mm, 300mm	1 Each
350	2.9mm Interlocking Bolt 20mm, 22mm, 24mm, 26mm, 28mm, 30mm, 32mm, 34mm, 36mm, 38mm, 40mm, 45mm, 50mm	2 Each
351	3.4mm Interlocking Bolt 16mm, 18mm, 20mm, 22mm, 24mm, 26mm, 28mm, 30mm, 32mm, 34mm, 36mm, 38mm, 40mm, 45mm, 50mm	2 Each




- **374** | Elastic Nail 2.0mm, 2.5mm, 3.0mm, 3.5mm, 4.0mm x 440mm

MATERIAL - Titanium



- **T355** | End Cap For Elastic Nail
SIZE - 2.0mm, 2.5mm, 3.0mm, 3.5mm, 4.0mm
MATERIAL - Titanium



• **810.001** | Elastic Nail Inserter/Extractor with SS Chuck



• **810.002** | Hammer Guide For Elastic Nail



• **810.003** | Slide Hammer



• **810.004** | F - Tool



• **810.005** | Pliers For Elastic Nail



• **810.006** | Pin Wrench



• **810.007** | Rod Cutter



• **810.008** | Oblique Temp

Instruments



• **810** | Standred Temp 2.5mm, 3.0mm, 3.5mm, 4.0mm



• **803.006** | Small Bone Awl Diamond Pointed



• **602.007** | Drill Bits 2.7mm x 95mm for Jacobs Chuck
• **602.009** | Drill Bits 3.2mm x 130mm for Jacobs Chuck
• **602.016** | Drill Bits 4.5mm x 130mm for Jacobs Chuck



• **605.045** | Drill And Tap Sleeve 4.5/3.2mm



Instruments

810.001	Elastic Nail Inserter/Extractor with SS Chuck	1 Each
810.002	Hammer Guide For Elastic Nail	1 Each
810.003	Slide Hammer	1 Each
810.004	F - Tool	1 Each
810.005	Pliers For Elastic Nail	1 Each
810.006	Pin Wrench	1 Each
810.007	Rod Cutter	1 Each
810.008	Oblique Temp	1 Each
810	Standred Temp 2.5mm, 3.0mm, 3.5mm, 4.0mm	1 Each
803.006	Small Bone Awl Diamond Pointed	1 Each
602.007	Drill Bits 2.7mm x 95mm for Jacobs Chuck	1 Each
602.009	Drill Bits 3.2mm x 130mm for Jacobs Chuck	1 Each
602.016	Drill Bits 4.5mm x 130mm for Jacobs Chuck	1 Each
605.045	Drill And Tap Sleeve 4.5/3.2mm	1 Each

Implant

374	Elastic Nail 2.0mm, 2.5mm, 3.0mm, 3.5mm, 4.0mm x 440mm	1 Each
355	End Cap For Elastic Nail 2.0mm, 2.5mm, 3.0mm, 3.5mm, 4.0mm	1Each



◦ 4002 | G-PFN Nail 130°

MATERIAL - Titanium / Stainless Steel
DIAMETER - 9mm, 10mm, 11mm, 12mm
LENGTH - 170mm, 200mm, 240mm



◦ 4003 | G-PFN Nail 135°

MATERIAL - Titanium / Stainless Steel
DIAMETER - 9mm, 10mm, 11mm, 12mm
LENGTH - 170mm, 200mm, 240mm



◦ 4005 | Long G-PFN Nail 130° (Left & Right)

MATERIAL - Titanium / Stainless Steel
DIAMETER - 9mm, 10mm, 11mm
LENGTH - 340mm, 360mm, 380mm, 400mm, 420mm



◦ 4006 | Long G-PFN Nail 135° (Left & Right)

MATERIAL - Titanium / Stainless Steel
DIAMETER - 9mm, 10mm, 11mm
LENGTH - 340mm, 360mm, 380mm, 400mm, 420mm



◦ 353 | 4.9mm Interlocking Bolt

MATERIAL - Titanium / Stainless Steel
LENGTH - 22mm to 60mm (2mm Increment)
60mm to 80mm (5mm Increment)



◦ 4004 | G-PFN Screw

MATERIAL - Titanium / Stainless Steel
LENGTH - 60mm to 120mm (5mm Increment)



◦ 355.14 | End Cap For G-PFN Nail

MATERIAL - Titanium



◦ 355A.14 | Lock Screw For G-PFN Nail

MATERIAL - Titanium



◦ 355B.14 | End Cap For G-PFN Nail Screw

MATERIAL - Titanium



• **813.001** | Insertion Handle For G-PFN



• **813.003** | Aiming Arm 130° For G-PFN

• **813.004** | Aiming Arm 135° For G-PFN



• **813.005** | Aiming Arm For Dynamic Locking For G-PFN



• **813.006** | Threaded Conical Bolt For G-PFN Nail



• **813.007** | Coupling Screw For G-PFN Screw



• **813.008** | Wrench For Attach G-PFN Screw



• **813.009** | Reamer 10mm For G-PFN Screw



• **813.010** | Protection Sleeve 16 x 10mm For G-PFN

• **813.011** | 10mm Pointer For G-PFN

• **813.012** | Guide Wire Sleeve 10 x 2.3mm For G-PFN

• **811.007** | Buttress Compression Nut For Pfna Blade



◦ 802.002 | Tissue Protector



◦ 804.010 | Hexagonal Screw Driver 3.5mm Tip
Long for 4.9mm Bolt



◦ 804.011 | Medullary Tube



◦ 807.004 | Protection Sleeve 10.0/8.0mm Long
◦ 807.007 | Drill Sleeve 8.0/4.0mm Long



◦ 807.014 | Protection Sleeve with Handle



◦ 807.022 | Radiographic Ruler For PFN Nail




◦ 811.015 | Sleeve 17 x 3.2mm For Guide Wire
◦ 811.039 | Pointer for Sleeve 17 x 3.2mm



◦ 811.019 | Pin Wrench L & Key Type



• **811.022** | Key For Pfna Blade With Combination Wrench 11mm



• **804C.008** | Depth Gauge for Proximal Femoral Nail



• **811.032** | Hammer Guide



• **811.033** | Connector For Hammer Guide



• **811.034** | 8mm Hexagonal Shaft With Spherical Head



• **811.035** | Cannulated Bone Awl



• **811.036** | Slotted Hammer with Handle



• **811.037** | T Handle Cannulated



• **811.038** | Connecting Bolt For PFN Aiming Arm



• **811.040** | T Handle Hexagonal Wrench 8.0mm for universal



• **811.025** | Starting Reamer 16.5mm For PFNA Nail



• **807.026** | Guide wire pointer 6.2mm



• **811.030** | Guide Wire 3mm For PFNA Blade



• **811.041** | PFN-2 Nail Extractor Attachment




• **602.030** | Drill Bits 4.0 mm Length 300mm for Jacobs Chuck



• **724.001** | Flexible Reaming Shaft




• **724** | Bit for Flexible Reaming Shaft 8.0mm, 8.5mm, 9.0mm, 9.5mm, 10.0mm, 10.5mm, 11.0mm, 11.5mm, 12.0mm



• **804.012** | Stopper Guide Wire 2.5mm Length 950mm



• **804.013** | Guide Wire 2.5mm Length 950mm



• **807.015** | Guide Wire 2.0mm X 500mm With Threaded Tip

AGC.77 Instrument for G-PFN Nail 130° & 135° with short G-PFN and locking Bolt Set



Instruments

813.001	Insertion Handle For G-PFN	1 Each
813.003	Aiming Arm 130° For G-PFN	1 Each
813.004	Aiming Arm 135° For G-PFN	1 Each
813.005	Aiming Arm For Dynamic Locking For G-PFN	1 Each
813.006	Threaded Conical Bolt For G-PFN Nail	2 Each
813.007	Cannulated Coupling Screw For G-PFN Screw	1 Each
813.008	Wrench For Attach G-PFN Screw	1 Each
813.009	Reamer 10mm For G-PFN Screw	1 Each
813.010	Protection Sleeve 16 x 11mm For G-PFN	1 Each
813.011	Trocar 3.2mm For G-Pfn	1 Each
813.012	Trocar Sleeve 11 x 3.2mm For G-PFN	1 Each
813.013	Fixation Sleeve G-PFN for Reamer 10mm	1 Each
813.014	Guide Wire Sleeve 11.0/2.0mm for 2.0mm guide wire	1 Each
813.015	T-Wrench For G-Pfn	1 Each
802.002	Tissue Protector	1 Each
804.010	Hexagonal Screw Driver 3.5mm Tip Long for 4.9mm Bolt	1 Each
804.011	Medullary Tube	1 Each
807.004	Protection Sleeve 10.0/8.0mm Long	1 Each
807.007	Drill Sleeve 8.0/4.0mm Long	1 Each
811.014	Protection Sleeve With Handle 20 x 17mm	1 Each
807.022	Radiographic Ruler For PFN Nail	1 Each
807.026	Guide wire pointer 6.2mm	1 Each
811.015	Sleeve 17 x 3.2mm For Guide Wire	1 Each
811.007	Buttress Compression Nut For Pfna Blade	1 Each
811.019	Pin Wrench L & Key Type	1 Each
811.022	Key For Pfna Blade With Combination Wrench 11mm	1 Each
804C.008	Depth Gauge for Proximal Femoral Nail	1 Each
811.032	Hammer Guide	1 Each
811.033	Connector For Hummer Guide	1 Each
811.034	8mm Hexagonal Shaft With Spherical Head	1 Each
811.035	Cannulated Bone Awl	1 Each
811.036	Slotted Hammer with Handle	1 Each
811.037	T Handle Cannulated	1 Each
811.038	Connecting Bolt For PFN Aiming Arm	1 Each
811.040	T Handle Hexagonal Wrench 8.0mm for universal	1 Each
811.025	Starting Reamer 16.5mm For Pfna Nail	1 Each
811.039	Pointer for Sleeve 17 x 3.2mm	1 Each
807.015	Guide Wire 2.0mm X 500mm With Threaded Tip	2 Each
811.041	PFN-2 Nail Extractor Attachment	1 Each
602.030	Drill Bits 4.0mm x 300mm for Jacobs Chuck	1 Each
804.012	Stopper Guide Wire 2.5mm Length 950mm	1 Each
804.013	Guide Wire 2.5mm Length 950mm	1 Each
724.001	Flexible Reaming Shaft	1 Each
724	Bit for Flexible Reaming Shaft 8.0mm, 8.5mm, 9.0mm, 9.5mm, 10.0mm, 10.5mm, 11.0mm, 11.5mm, 12mm	1 Each

Implants

4002	G-PFN Nail 9.0mm, 10mm, 11mm, 12mm x 170mm, 200mm, 240mm 130°	1 Each
4003	G-PFN Nail 9.0mm, 10mm, 11mm, 12mm x 170mm, 200mm, 240mm 135°	1 Each
353	4.9mm Interlocking Bolt 22mm, 24mm, 26mm, 28mm, 30mm, 32mm, 34mm, 36mm, 38mm, 40mm, 42mm, 44mm, 46mm, 48mm, 50mm, 52mm, 54mm, 56mm, 58mm, 60mm, 65mm, 70mm, 75mm, 80mm	2 Each
4004	G-PFN Screw 60mm, 65mm, 70mm, 75mm, 80mm, 85mm, 90mm, 95mm, 100mm, 105mm, 110mm, 115mm, 120mm	1 Each
355.14	End Cap For G-PFN Nail	5 Each
355A.14	Lock Screw For G-PFN Nail	5 Each

PHOTO NOT AVAILABLE

Implants

4005D	Long G-PFN Nail 9.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 130°	1 Each
4005E	Long G-PFN Nail 10.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 130°	1 Each
4005F	Long G-PFN Nail 11.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 130°	1 Each
4006D	Long G-PFN Nail 9.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 135°	1 Each
4006E	Long G-PFN Nail 10.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 135°	1 Each
4006F	Long G-PFN Nail 11.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 135°	1 Each

PHOTO NOT AVAILABLE

Implants

4005D	Long G-PFN Nail 9.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 130°	1 Each
4005E	Long G-PFN Nail 10.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 130°	1 Each
4005F	Long G-PFN Nail 11.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 130°	1 Each
4006D	Long G-PFN Nail 9.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 135°	1 Each
4006E	Long G-PFN Nail 10.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 135°	1 Each
4006F	Long G-PFN Nail 11.0mm x 340mm, 360mm, 380mm, 400mm, 420mm, 135°	1 Each



HARDIK[®]
INTERNATIONAL

House of Ortho Tools

Hardik International Pvt. Ltd.

Plot No. G-2041/42/A, Kishan Gate, Lodhika G.I.D.C.,
Kalawad Road, Metoda, RAJKOT-360 021, (INDIA).

info@hardikinternational.com

www.hardikinternational.com

☎ : +91 97370 47470/1/2/4 | 📞 : hardikinternational

Follow Us :     