

01

COMPANY PROFILE
公司简介

03

ADVANCED PRODUCTION EQUIPMENT
先进生产设备

05

QUALIFICATIONS AND HONOR
公司资质和荣誉

09

HIGH TORQUE HANDPIECE
力王手机

11

E-GENERATOR LED HANDPIECE
风光手机

13

FIBER OPTIC HANDPIECE
1:5 光纤手机

15

ORTHODONTICS CONTRA ANGLE
正畸弯机

17

ELECTRIC MICRO MOTOR
电动打磨机

19

DIAMOND BURS
金刚砂车针

21

SPINDLE
轴

23

HIGH SPEED HANDPIECE SERIES
高速手机系列

27

LOW SPEED HANDPIECE SERIES
低速手机系列

30

MICRO MOTOR SERIES
电动打磨机系列

31

CARTRIDGE OF HIGH SPEED HANDPIECE
高速手机机芯

34


ACCESSORIES SERIES
配件系列

39

BUR SERIES
车针系列


MARKING NOTES


标识说明

 标准头 /Standard

 大头 /Super torque


 迷你头 /Mini

 取针式 /Wrench chuck


 按压式 /Push button


 单点喷雾 /Single water spray


 多点喷雾 /Multiple water spray


 2孔连接 /2 hole Connector


 4孔连接 /4 hole Connector

 6孔连接 /6 hole Connector

 LED 照明 /LED light


 光纤照明 /Fiber optic


 外置水道 /External type


 内置水道 /Internal channel type


FUNCTIONAL OVERVIEWS


功能概述


 抗过热按钮
It reduces the risk of patient was burnt by the head of high speed handpiece to make sure the operating comfort and patient's safety 从而降低了高速手机头部灼伤患者的风险，确保您的操作舒适性和患者的安全。


 车针锁紧
The most accurate chucking and rotaing technology for burs to create the most stable operation. 最精确的车针锁紧和旋转驱动技术，创造出最佳平稳无震动的仪器。


 内置水道
Inner channel sprinkler system, 4 hole water, fully cooled to avoid gas blockage. 内水道喷水系统，4孔出水，全面冷却，避免气体堵塞。


 镀钛机身
Titanium surface treatment,prevent scratches and achieve durable use. 钛机身表面处理，坚固耐用，防止刮花。


 陶瓷轴承
Longer service life,higher accuracy,lower noise,less vibration. 更长的使用寿命，更高的精度，更低的噪音，更少震动。


 微型电机
The internal coil of the motor adopts a unique technique. When a set of coil fails, the other set of coils can continue to work, which is equivalent to the standby system to ensure the service life. 电机内部线圈采用独特技术，当一组线圈出现故障时另外一组线圈可以继续工作，相当于备用系统，确保使用寿命

 光纤照明
Sunlight fiber lighting, comfortable and not dazzling, helps to achieve accurate diagnosis, to achieve the desired therapeutic effect. 日光级光纤照明，舒适不刺眼，有助于实现准确的诊断达到理想的治疗效果。

 LED 照明
7000lux luminance,perfect and natural white light,it can be autoclaved repeatedly. 7000lux 亮度，完美的自然白光，可经受反复消毒。

 水气雾化
Water vapor separation forms 360 degrees atomization, efficient and comprehensive cooling cutting surface. 水气分离形成 360 度雾化，高效且全面冷却切割面。

 电动驱动
The speed is fast, the maximum speed can reach 200 thousand rpm, the motor turns around, and the bending machine turns 5 times. 转速快，最高转速可达 20 万转，马达转一圈，弯机转 5 圈。

 135°C 高温消毒
135°C Autoclavable 可在 135 摄氏度下高温高压消毒。



OMPANY PROFILE
公司概况

Wenjian Medical Instrument Co.,Ltd is a professional dental handpiece manufacturer that focusing on R&D produce, sales and after services. We are manufacturer of Dental High-speed & Low-speed handpiece, E-generator LED handpiece, cartridges, Micro Motor, 3 way syringe and other dental instruments. We are one of the best dental handpiece manufacturers in China. We are committing to supply competitive products to every customers so our products are selling all over the world. In the future, we will continue to do as our principle that is "Quality Safely, Technology Innovation". We will supply better products and services to satisfy different requirements of customers.

佛山市稳健医疗器械有限公司是一家集研发、制造、销售和售后服务于一体的专业牙科手机制造企业。产品包含高速涡轮手机，低速涡轮手机，三用枪，打磨机等及相关配件。近年来因生产专业技术不断提高，产品质量得到了广大客户朋友的认可，我们的产品已畅销世界 100 多个国家。我们将一直秉承“稳健制造，安全可靠，技术创新，按压不飞针”的宗旨，为客户提供更好的产品和服务。



ADVANCED PRODUCTION EQUIPMENT



EXCLUSIVE INTRODUCTION OF "STAR" "CITIZEN"
ADVANCED PRODUCTION EQUIPMENT TO CREATE SUPERIOR PRODUCTS

独家引进 STAR 野村 西铁城 先进生产设备，缔造卓越产品





QUALIFICATIONS AND HONOR 公司资质和荣誉

Wenjian Medical Instruments Co., Ltd, is the National Science and Technology Enterprises, obtained the Medical Device Manufacturing Enterprise License that included dental highspeed handpiece, dental lowspeed handpiece, diamond bur and so on. Until now, we obtained CE and ISO Certification, FDA Certification and Korean GMP Certification, National Practical New-type Patent, National Invention Patent, Guangdong Science and Technology Products Certification and so on.

稳健医疗器械有限公司，是国家高新技术企业单位，取得医疗器械生产企业许可证，车针/牙科高/低速手机等国家医疗器械注册证。稳健医疗产品进入市场至今为止，公司已荣获CE、美国FDA、德国TUV、韩国GMP等认证、国家实用新型专利、国家发明专利及广东省高新技术产品证书十余项。





IGH SPEED HANDPIECE SERIES
高速手机系列


HIGH TORQUE HANDPIECE 力王手机

Wenjian has been committed to research and develop more closed with the user requirements. High Torque Handpiece has the following advantages like comfortable, lightweight, imported titanium material, the unique coating process and the torque up to 40N.

稳健一直致力于研发更符合使用者要求的产品。力王手机，舒适、轻盈，进口镀钛材质，独特的涂层工艺，其扭力高达 40N。

STRENGTH PARAMETERS

力度参数



OPTIONAL

可选配



Water vapor separation forms 360 degrees atomization, efficient and comprehensive cooling cutting surface. 水气分离形成 360 度雾化，高效且全面冷却切割面。



135°C Autoclavable 可在 135 摄氏度下高温高压消毒。



Titanium surface treatment, prevent scratches and achieve durable use. 钛机身表面处理，坚固耐用，防止刮花。






E-GENERATOR LED HANDPIECE 风光手机

Wenjian independently research and develop new products, unique noise reduction technology, feel what is the peaceful world together. Accurate design, professional skills, create the unparalleled cutting performance, this design has been exclusive patent, only belong to you.

稳健自主研发的全新产品，独一无二的降噪工艺，带你见证安静的秘密。纤细设计，精湛做工，铸就无与伦比的切削性能，此产品已获得独家专利，属于你的专属定制手机。

OPTIONAL 可选配



微型电机

The internal coil of the motor adopts a unique technique. When a set of coil fails, the other set of coils can continue to work, which is equivalent to the standby system to ensure the service life.

电机内部线圈采用独特技术，当一组线圈出现故障时另外一组线圈可以继续工作，相当于备用系统，确保使用寿命



陶瓷轴承

Longer service life, higher accuracy, lower noise, less vibration.

更长的使用寿命，更高的精度，更低的噪音，更少震动。



LED照明

7000lux luminance, perfect and natural white light, it can be autoclaved repeatedly.

7000lux 亮度，完美的自然白光，可经受反复消毒。



高温消毒

135°C Autoclavable
可在 135 摄氏度下高温高压消毒。



LOW SPEED HANDPIECE SERIES
低速手机系列

FIBER OPTIC HANDPIECE 1: 5 光纤手机

Connected with the electric motor, it can replace the high-speed handpiece totally and focus on breaking teeth crown, teeth bone, teeth cavity and adjacent teeth repair.

与电动马达连接，可以取代高速手机，扭力大，转速稳定，用于破冠，去骨，窝洞及邻牙修复。



电动驱动

The speed is fast, the maximum speed can reach 200 thousand rpm, the motor turns around, and the bending machine turns 5 times. 转速快，最高转速可达 20 万转，马达转一圈，弯机转 5 圈。



内置水道

Inner channel sprinkler system, 4 hole water, fully cooled to avoid gas blockage. 内水道喷水系统，4 孔出水，全面冷却，避免气体堵塞。



光纤照明

Sunlight fiber lighting, comfortable and not dazzling, helps to achieve accurate diagnosis, to achieve the desired therapeutic effect. 日光级光纤照明，舒适不刺眼，有助实现准确的诊断达到理想的治疗效果。



OPTIONAL
可选配



Fiber





RTHODONTICS SERIES
正畸系列

ORTHODONTICS CONTRA ANGLE 正畸弯机

In orthodontic care process, the need to remove the demand of the adjacent surface enamel and tooth clearance gap, using traditional high-speed car needle, have certain limitations and difficulties of treatment. Orthodontic bending machine, which has different coarse line sand, can be fast and accurate with the principle of linear motion, complete the function such as the removal of the adjacent enamel and the opening of the gap between the teeth. Different type of sand, reduce the cost of purchasing and solve the problem of orthodontics.

在正畸医疗过程中，需要去除邻面釉质和牙间隙开缝的需求，使用传统的高速车针，对治疗有一定的局限性和难度。正畸弯机，配置不同粗线的砂片，以线形运动的原理，可以快速，准确，完成邻面釉质去除和牙间隙开缝等功能。配置不同型号的砂片，降低采购成本，解决正畸难题。



OPTIONAL

可选配

- 15 μ m- 预抛光
15micron - for pre-polishing
- 25 μ m- 精修
25 micron - for finishing
- 40 μ m- 仿型
40 micron - for contouring

- 60 μ m- 减速去釉
60 micron - for reduction
- 90 μ m- 颌面开齿
90 micron - for use as a saw for opening contacts





ICRO MOTOR SERIES
电动打磨机系列

ELECTRIC MICRO MOTOR

电动打磨机



Electric micro motor is more powerful, can replace air motor to solve the lack of air motor power, the difficult of rotation adjustment and so on.

电动打磨机，动力强劲，可代替气动马达，解决因气动马达力量不足，转速不易调节等问题。

E Type compatible connector, can be matched with contra angle 1:1, 16:1, 1:5.
E型通配接口，可搭配 1:1, 16:1, 1:5 弯机使用。

High rotation, high torque, low noise, no vibration, more comfortable operation
转速高，力度大，噪音低，无振动，操作更舒适。

Inner type LED lighting electric motor system.
内置式 LED 发光电动马达系统。



DB

UR SERIES
车针系列

DIAMOND BURS 金刚砂车针

Imported stainless steel.
进口不锈钢材质。



Using high quality diamond grains, through to complex, magnetic, spheronization, purified by adhesive, electroplating and metal materials and methods.
采用优质金刚砂颗粒，经过去杂，去磁，滚圆，提纯处理后，经电镀方法与金属材料粘接而成。

International standard length:
14~22mm
国际标准长度 :14~22mm

High hardness, strong cutting force, sharpness and wear resistance, excellent performance and better service life.
硬度高，切削力强，锋利及耐磨度上表现优异，性能和使用寿命更佳。



WJ diamond bur adopt imported stainless steel material.

稳健金刚砂车针，采用进口不锈钢材质，坚固，经久耐用，不易变形。



ACCESSORIES SERIES
配件系列

SPINDLE 轴

Spindle the soul of the handpiece; After thousands of testing, keep the ultra-high quality. 轴，手机的灵魂；经过数千次检测，铸造超高品质。

Adopt the imported stainless steel, the high precision technology of locked the bur and drive rotation. The concentricity is less than 0.01 mm. 采用进口不锈钢材质，更精确的车针锁紧和旋转驱动技术，精密度极高，同心度在 0.01 毫米以下。

Locked the diamond bur more tightly and make sure safely.

具有强大的夹持力，安全不飞针。



OPTIONAL

可选配



push shaft
按压轴



Key shaft
取针轴

WJ MEDICAL



Close type handpiece series 高速筒夹手机

Rotation speed(转速)	280.000-400.000 RPM
Air pressure(工作气压)	0.25Mpa-0.30Mpa
Torque(扭矩) W	大头 (T): 21/ 标准头 (S): 19/ 迷你头 (M): 17
Bur applicable (适用车针)	Φ1.595-Φ1.600mm
Noise(噪音)	≤ 65db

WJ MEDICAL



high speed color handpiece 高速彩虹手机

Rotation speed(转速)	280.000-400.000 RPM
Air pressure(工作气压)	0.25Mpa-0.30Mpa
Torque(扭矩) W	大头 (T): 21/ 标准头 (S): 19
Bur applicable (适用车针)	Φ1.595-Φ1.600mm
Noise(噪音)	≤ 65db
shell design (花式)	

WJ MEDICAL



Anti-retraction handpiece series 高速环保手机

Rotation speed(转速)	280.000-400.000 RPM
Air pressure(工作气压)	0.25Mpa-0.30Mpa
Torque(扭矩) W	大头 (T): 21/ 标准头 (S): 19
Bur applicable (适用车针)	Φ1.595-Φ1.600mm
Noise(噪音)	≤ 65db
shell design (花式)	

WJ MEDICAL



High-Torque type handpiece 力王手机

Rotation speed(转速)	280.000-400.000 RPM
Air pressure(工作气压)	0.25Mpa-0.30Mpa
Torque(扭矩) W	标准头 (S): 23
Bur applicable (适用车针)	Φ1.595-Φ1.600mm
Noise(噪音)	≤ 65db

WJ MEDICAL



E-generator LED handpiece 高速风光手机

Rotation speed(转速)	280.000-400.000 RPM
Air pressure(工作气压)	0.25Mpa-0.30Mpa
Torque(扭矩) W	大头 (T): 21/ 标准头 (S): 19
LED life(寿命)	10000hours
Noise(噪音)	≤ 65db
shell design (花式)	

WJ MEDICAL



F type E-generator LED handpiece 鱼纹风光手机

Rotation speed(转速)	280.000-400.000 RPM
Air pressure(工作气压)	0.25Mpa-0.30Mpa
Torque(扭矩) W	大头 (T): 21/ 标准头 (S)
Bur applicable (适用车针)	Φ1.595-Φ1.600mm
Noise(噪音)	≤ 65db
LED life(寿命)	10000hours

WJ MEDICAL | WENJIAN



U type E-generator LED handpiece 弧形风光手机

Rotation speed(转速)	280.000-400.000 RPM
Air pressure(工作气压)	0.25Mpa-0.30Mpa
Torque(扭矩) W	大头 (T): 21/ 标准头 (S): 19
Bur applicable (适用车针)	Φ 1.595-Φ 1.600mm
Noise(噪音)	≤ 65db
LED life(寿命)	10000hours

WJ MEDICAL | WENJIAN



LED handpiece with MSDG-L 高速开关风光手机

Rotation speed(转速)	280.000-400.000 RPM
Air pressure(工作气压)	0.25Mpa-0.30Mpa
Torque(扭矩) W	大头 (T): 21/ 标准头 (S): 19
Bur applicable (适用车针)	Φ 1.595-Φ 1.600mm
Noise(噪音)	≤ 65db
LED life(寿命)	10000hours

WJ MEDICAL | WENJIAN



Imitate Fiber handpiece 六孔假光纤手机

Rotation speed(转速)	280.000-400.000 RPM
Air pressure(工作气压)	0.25Mpa-0.30Mpa
Torque(扭矩) W	大头 (T): 21/ 标准头 (S): 19/ 迷你头 (M) : 17
LED life(寿命)	10000hours
Noise(噪音)	≤ 65db
shell design (花式)	

WJ MEDICAL | WENJIAN



Fiber optice handpiece 六孔光纤手机

Rotation speed(转速)	280.000-400.000 RPM
Air pressure(工作气压)	0.25Mpa-0.30Mpa
Torque(扭矩) W	大头 (T): 21/ 标准头 (S): 19/ 迷你头 (M) : 17
Bur applicable (适用车针)	Φ 1.595-Φ 1.600mm
Noise(噪音)	≤ 65db
Light intensity(光强度)	≥ 7000LUX
Notes(备注)	与 KAVO 快接通用, 机芯与 kavo 通用 (645、642、660、655、105) compatible with KAVO coupling ,compatible with KAVO cartridge

WJ MEDICAL | WENJIAN



45 degree LED handpiece 45 度阻生齿风光手机


Rotation speed(转速)	280.000-400.000 RPM
Air pressure(工作气压)	0.25Mpa-0.30Mpa
Torque(扭矩) W	标准头 (S): 23
Bur applicable (适用车针)	Φ 1.595-Φ 1.600mm
Noise(噪音)	≤ 65db
LED life(寿命)	10000hours
Notes(备注)	(采用后排气, 机头前端没有气, 不会造成手术交叉感染) No coolant Air and Rear Exhaust, avoid cross infection

WJ MEDICAL | WENJIAN




Nsk type quick coupling E-generator LED handpiece 高速风光 N 快插手机


Rotation speed(转速)	280.000-400.000 RPM
Air pressure(工作气压)	0.25Mpa-0.30Mpa
Torque(扭矩) W	大头 (T): 21/ 标准头 (S): 19
Bur applicable (适用车针)	Φ 1.595-Φ 1.600mm
Noise(噪音)	≤ 65db
shell design (花式)	
Notes	备注: (普通、风光可选, 与 NSK 快接通用) Notes: Optional : Nomal type , E-generator type, compatible with NSK




KAVO type quick coupling E-generator LED handpiece 高速风光 K 快插手机	
Rotation speed(转速)	280.000-400.000 RPM
Air pressure(工作气压)	0.25Mpa-0.30Mpa
Torque(扭矩) W	大头 (T): 21/ 标准头 (S): 19
LED life(寿命)	10000hours
Noise(噪音)	≤ 65db
shell design (花式)	<p>备注: (普通风光可选, 与 KAVO 快接通用) Notes: Optional: Nomal type, E-generator type, compatible with KAVO</p>




High speed handpiece set 3+1 套装	
Rotation speed(转速)	280.000-400.000 RPM
Air pressure(工作气压)	0.25Mpa-0.30Mpa
Torque(扭矩) W	大头 (T): 21/ 标准头 (S): 19 /45度 :23
LED life(寿命)	10000hours
Noise(噪音)	≤ 65db
shell design (花式)	<p>备注: 用于所有口腔手术, 配四转六转接头, 不用光纤牙椅 Notes: For all oral surgery, with 4 to 6 coupling, any chair can use</p>




External Low speed set series (straight/ contra angle / Airmotor) 外水道套套装 (直机 / 弯机 / 马达)	
Ratio(速比)	等速 1:1
Air pressure(工作气压)	0.3Mpa-0.35Mpa
Noise(噪音)	≤ 68db
Rotation speed(转速)	≥ 15.000 RPM
shell design (花式)	<p>旧款 新款</p>



Inner channel Low speed set series(straight/ contra angle / Airmotor) 内水道套装 (直机 / 弯机 / 马达)	
Ratio(速比)	等速 1:1
Air pressure(工作气压)	0.3Mpa-0.35Mpa
Noise(噪音)	≤ 68db
Rotation speed(转速)	≥ 15.000 RPM
Light intensity(光强度)	≥ 7000LUX
shell design (花式)	<p>备注: (普通、风光可选) Notes: Optional: Nomal type, E-generator type</p>



Reduction contra angle series 减速弯机 (10:1, 16:1, 64:1)	
Ratio(速比)	10:1、16: 1、64:1
Air pressure(工作气压)	0.3Mpa-0.35Mpa
Noise(噪音)	≤ 60db
Rotation speed(转速)	4000; 2500; 625rpm
Notes(备注)	内置冷却系统, 适合于 E 型接口的气动或电动马达 internal cooling system, compatible with Airmotor or Electric Micromoto



Implant contra angle 20:1 种植弯机	
Ratio(速比)	20:1
Air pressure(工作气压)	0.3Mpa-0.35Mpa
Noise(噪音)	≤ 60db
Rotation speed(转速)	2000r/min
Torque(扭矩)	≥ 55N.cm
LED life(寿命)	10000hours
Notes(备注)	内置冷却系统, 适合于 E 型接口的气动或电动马达, 普通、风光、光纤可选 internal cooling system, compatible with Airmotor or Electric Micromotor ; Optional: Nomal type, E-generator type, Fiber optic type

WJ MEDICAL | WENJIAN



Increasing contra angle 1:5 增速弯机

Ratio(速比)	1:5
Air pressure(工作气压)	0.3Mpa-0.35Mpa
Noise(噪音)	≤ 60db
Rotation speed(转速)	20.0000 RPM
LED life(寿命)	10000hours
Notes(备注)	内置冷却系统, 适合于E型接口的气动或电动马达, 普通、风光、光纤可选 internal cooling system ,compatible with Airmotor or Electric Micromotor ; Optional : Nomal type , E-generator type, Fiber optic type

WJ MEDICAL | WENJIAN



implant torsion handpiece 扭力弯机

Torque(扭矩)	20N/35N
Notes(备注)	又称种植扳手, 种牙时拧紧螺丝 Used to tighten the screws during planting

WJ MEDICAL | WENJIAN



orthodontics contra angle 正畸弯机

Noise(噪音)	≤ 60db
Air pressure(工作气压)	0.3Mpa
Torque(扭矩)	30N
Rotation speed(转速)	≥ 15.000 RPM

WJ MEDICAL | WENJIAN



Brushless Electric Micromotor 电动打磨机

Input power (输入电压)	DC24V
Torque(扭矩)	3.4Ncm
Noise(噪音)	≤ 40db
Rotation speed(转速)	2.000-40.000 RPM
Light(光源)	LED
Notes(备注)	E型ISO国际标准接口, 正反转可调, 与1:1/16:1/1:5等弯机通配 E-type ISOstandard ; Forward / Reverse Switchable ,compatible with 1:1/16:1/1:5 contra angle

WJ MEDICAL | WENJIAN



35,000RPM Brush Micromotor 3.5万转有碳刷打磨机

Input power (输入电压)	220V/110V
Rotation speed(转速)	3.5000 RPM
Noise(噪音)	≤ 68db
Notes(备注)	可选配 18/102/204 手柄 Select E-type motor / Lab motor

WJ MEDICAL | WENJIAN



50,000RPM Brushless Micromotor 5万转无碳刷打磨机

Input power (输入电压)	220V
Rotation speed(转速)	50.000 RPM
Noise(噪音)	≤ 60db

Push Button close type cartridge series 敞开式按压机芯



Bur extracting force (取针力度)	大头 (T) : $\leq 45N$, 标准 (S) : $\leq 40N$
Hardness (硬度)	60HRC
Clamping force (夹持力)	Push button(按压式):30N,Wrench chuck(取针式):50N
Accuracy of bur in rotation (旋转精度)	Max(最大):0.01mm
Notes(备注)	可定制机芯, 与 NSK,KAVO,W&H 等通用。 Customizable cartridge,general with NSK, KAVO, W&H.



Wrench chuck type open type cartridge series 敞开式取针机芯



Bur extracting force (取针力度)	大头 (T) : $\leq 45N$, 标准 (S) : $\leq 40N$
Hardness (硬度)	60HRC
Clamping force (夹持力)	Push button(按压式):30N,Wrench chuck(取针式):50N
Accuracy of bur in rotation (旋转精度)	Max(最大):0.01mm
Notes(备注)	可定制机芯, 与 NSK,KAVO,W&H 等通用。 Customizable cartridge,general with NSK, KAVO, W&H.



Anti-retraction open type cartridge series 敞开式环保机芯



Bur extracting force (取针力度)	大头 (T) : $\leq 45N$, 标准 (S) : $\leq 40N$
Hardness (硬度)	60HRC
Clamping force (夹持力)	Push button(按压式):30N,Wrench chuck(取针式):50N
Accuracy of bur in rotation (旋转精度)	Max(最大):0.01mm
Notes(备注)	可定制机芯, 与 NSK,KAVO,W&H 等通用。 Customizable cartridge,general with NSK, KAVO, W&H.



Push Button close type cartridge series 按压式筒夹机芯



Bur extracting force (取针力度)	大头 (T) : $\leq 45N$, 标准 (S) : $\leq 40N$
Hardness (硬度)	60HRC
Clamping force (夹持力)	Push button(按压式):30N,Wrench chuck(取针式):50N
Accuracy of bur in rotation (旋转精度)	Max(最大):0.01mm
Notes(备注)	可定制机芯, 与 NSK,KAVO,W&H 等通用。 Customizable cartridge,general with NSK, KAVO, W&H.



WJ MEDICAL | WENJIAN

Wrench chuck type close type cartridge series
取针式筒夹机芯



Bur extracting force (取针力度)	大头 (T) : $\leq 45N$, 标准 (S) : $\leq 40N$
Hardness (硬度)	60HRC
Clamping force (夹持力)	Push button(按压式):30N,Wrench chuck(取针式):50N
Accuracy of bur in rotation (旋转精度)	Max(最大):0.01mm
Notes(备注)	可定制机芯, 与 NSK,KAVO,W&H 等通用。 Customizable cartridge,general with NSK, KAVO, W&H.



WJ MEDICAL | WENJIAN



New type NSK coupling

新款 NSK 快接头

WJ MEDICAL | WENJIAN



Tube

一体管线接头

WJ MEDICAL | WENJIAN



3 way syringe

三用枪

WJ MEDICAL | WENJIAN



Disposable handpiece

一次性手机

WJ MEDICAL | WENJIAN



Lab Motor

204 手柄

WJ MEDICAL | WENJIAN



Lab Motor

102 手柄

WJ MEDICAL | WENJIAN



KAVO coupling

Kavo 快接头

WJ MEDICAL | WENJIAN



E-type motor

18 手柄



Head of Low speed contra angle

内 / 外水道 弯机头

WJ MEDICAL | WENJIAN



3 way syringe tips

三用枪嘴

WJ MEDICAL | WENJIAN



colorful 3 way syringe

彩色三用枪

WJ MEDICAL | WENJIAN



NSK coupling

NSK 快接头

WJ MEDICAL | WENJIAN



4 to 2 Adapter

4 转 2 转接头


WJ MEDICAL | WENJIAN



2 to 4 Adapter

2 转 4 转接头

WJ MEDICAL | WENJIAN



Wave gasket

波形垫片

WJ MEDICAL | WENJIAN



Gear

齿轮

WJ MEDICAL | WENJIAN



O-ring

O 型圈

WJ MEDICAL | WENJIAN



2 or 4 hole gasket

2/4 孔胶垫

WJ MEDICAL | WENJIAN



Wave gasket

平垫

WJ MEDICAL | WENJIAN



Cap

手机后盖

WJ MEDICAL | WENJIAN



push button/ Wrench
chuck spindle

按压 / 取针轴

WJ MEDICAL | WENJIAN



Shaft of E-generator

电机
传动轴

WJ MEDICAL | WENJIAN



Bearing

轴承

WJ MEDICAL | WENJIAN



Impeller

风轮

WJ MEDICAL | WENJIAN



Key

取针器

WJ MEDICAL | WENJIAN



Wrench

扳手

ISO 173/014 TF-11	ISO 173/014 TF-12	ISO 173/018 TF-13	ISO 172/023 TF-14	ISO 171/014 TF-20	ISO 171/016 TF-21	ISO 170/021 TF-22	ISO 170/018 TF-23	ISO 170/016 TF-31	ISO 170/014 TF-43	ISO 170/012 TF-42
TF Series										
ISO 169/011 TF-41	ISO 199/016 TR-11	ISO 199/016 TR-12	ISO 198/018 TR-13	ISO 198/022 TR-14	ISO 200/022 TR-15	ISO 200/024 TR-19	ISO 197/014 TR-20	ISO 197/016 TR-21	ISO 197/018 TR-24	ISO 199/016 TR-25
TR Series										
ISO 199/018 TR-26	ISO 160/016 TC-11	ISO 160/014 TC-21	ISO 160/012 TC-26	ISO 299/012 FO-11	ISO 298/014 FO-21	ISO 298/016 FO-22	ISO 257/018 FO-32	ISO 257/028 FO-25	ISO 257/032 FO-27	ISO 111/012 SF-11
TC Series			FO Series			SF Series				
ISO 111/014 SF-12	ISO 111/016 SF-13	ISO 110/014 SF-21	ISO 109/010 SF-41	ISO 288/012 SO-20	ISO 289/014 SO-21	ISO 141/012 SR-11	ISO 141/014 SR-12	ISO 141/016 SR-13	ISO 237/018 EX-20	ISO 237/012 EX-21
SF Series				SO Series		SR Series			EX Series	
ISO 234/014 EX-31	ISO 237/010 EX-41	ISO 220/017 EX-24	ISO 237/032 EX-26	ISO 039/032 EX-11	ISO 039/035 EX-12	ISO 001/018 BR-31	ISO 001/014 BR-41	ISO 001/012 BR-46	ISO 001/010 BR-45	ISO 001/008 BR-49
BR Series										
ISO 002/018 BC-31	ISO 002/014 BC-42	ISO 002/012 BC-43	ISO 019/010 DI-41	ISO 019/014 DI-42	ISO 010/016 SI-48	ISO 010/014 SI-47	ISO 010/012 SI-46	ISO 068/041 WR-13	ISO 171/014 TF-S20	ISO 171/016 TF-S21
BC Series			DI Series		SI Series		WR Series		S Series	
ISO 170/012 TF-S22	ISO 170/018 TF-S23	ISO 173/016 TF-S31	ISO 169/011 TF-S41	ISO 198/018 TR-S13	ISO 197/016 TR-S21	ISO 160/014 TC-S21	ISO 298/014 FO-S21	ISO 288/012 SO-S20	ISO 001/012 BR-S46	ISO 001/010 BR-S45
TF Series										
ISO 002/012 BC-S43	ISO 002/014 BC-S42	ISO 019/010 DI-S41	ISO 010/012 SI-S46	ISO 010/016 SI-S48	ISO 198/016 TF-SS31	ISO 197/016 TR-SS21	ISO 019/010 DI-SS41	ISO 173/018 TF-13C	ISO 198/018 TR-13C	ISO 197/018 TR-62C
SS Series								C Series		

ISO 160/016 TC-11C	ISO 190/018 FO-54C	ISO 068/042 WR-13C	ISO 173/016 TF-12F	ISO 171/016 TF-21F	ISO 199/016 TR-11F	ISO 198/018 TR-13F	ISO 197/016 TR-21F	ISO 199/016 TR-25F	ISO 199/018 TR-26F	ISO 160/016 TC-11F
F Series										
ISO 160/014 TC-21F	ISO 298/014 FO-21F	ISO 298/016 FO-22F	ISO 254/017 FO-30F	ISO 237/021 EX-21F	ISO 150/012 EX-18F	ISO 001/006 BR-48F	ISO 304/016 WR-31F	ISO 304/016 WR-S31F	ISO 173/016 TF-12EF	ISO 171/016 TF-21EF
EF Series										
ISO 199/016 TR-11EF	ISO 198/018 TR-13EF	ISO 197/016 TR-21EF	ISO 199/016 TR-25EF	ISO 199/018 TR-26EF	ISO 160/016 TC-11EF	ISO 160/011 TC-21F	ISO 248/014 FO-20EF	ISO 298/014 FO-21EF	ISO 298/016 FO-22EF	ISO 243/010 FO-40EF
FO Series										
ISO 248/012 FO-41EF	ISO 248/012 FO-42EF	ISO 237/021 EX-21EF	ISO 001/010 BR-40EF	ISO 196/020 CR-11F	ISO 196/019 CR-11EF	ISO 197/012 CR-12F	ISO 196/014 CR-12EF	ISO 197/012 CR-21F	ISO 289/013 CR-22F	ISO 001/006 CD-50F
FO Series						CR Series				CD Series
ISO 178/007 CD-51F	ISO 17/007 CD-52F	ISO 161/007 CD-53F	ISO 246/006 CD-54F	ISO 254/009 CD-55F	ISO 246/007 CD-56F	ISO 247/007 CD-57F	ISO 108/008 CD-58F	ISO 171/013 CD-59F		
CD Series										

Color Coding

Color code indication is provided at the end of working part according to the size of diamond grit conforming to the I.S.O. Specification.

Classification

BEFULL coding system reflects shapes of burs. The contour and sizes could be imagined by BEFULL nomenclature since I.S.O. numbers are also indicated. TF/Taper Flat End TR/Taper Round End TC/Taper Conical End FO/FlameOgivalEnd SF/Straight Flat End SO/Straight Ogival End SR/Straight Round End RS/Rounded Shoulder EX/Special (Extra)Shape BR/Ball round type BC/Ball Collar tapy DI/Double Inverted ConeSI/Single Inverted Cone WR/Wheel Round Edge CR/C Rinlay Preparation Dia-Burs. CD/Childrens Dentistry Dia-BursS/S hort Shank SS/Super Short Shank C/Coarse F/Fine EF/Extra Fine.

High Quality Control

100 Percent inspection with laser measurement machine cutting test and strength test for each lot are performed. All products comply with our own standards which are far more strict than I.S.O. standards.

Before use BE FULL DIA-BURS. please see the instruction as indicated below

For dental use only Non-sterile Insert the shank into handpiece according to the instruction of the handpiece manufacturer. Insertion should be the full length Before use, run the instruments to make sure there is no descentering. If the instruments become bent, be sure to remove them from any further use. Please always use spray. A plentiful water supply is necessary to maintain long-lived performance and also to prevent loading and pulpitis. Handle with care not to put excessive load. Pressure should be 20 to 50P for all instruments with diameter up to 0.23. Use below max allowable r.p.m. Of each instrument, see chart above, to avoid a possibility of being injured due to breakage. Avoid levering, tilting, twisting and jamming motion. Don't work with a pulling motion. Wear a safety glasses and/or face mask to protect from being injured. Practicable for using auto-clave, chemical vaporized-clave, dry heat and E.O.G. Sterilization. After using, follow appropriate treatment as the medical waste.



C E R T I F I C A T E

Full Quality Assurance System Medical Devices Directive 93/42/EEC Annex II (Excluding Section 4)

M.2020.106.13286-1 Design Examination Certificate Was Prepared for Class III Products Defined in This Certificate.

Company Name : Medico (Huaian) Co., Ltd

Company Address : No:9 South Guangzhou Road, Huaian, Jiangsu, China

Related Directives and Annex : 93/42/EEC Medical Devices Directive - Annex II (Excluding Section 4)

Product : - Sterile, Absorbable, Synthetic Braided Polyglycolic Acid (PGA)
Sutures With and Without Needle - Class III
- Sterile, Absorbable, Synthetic Braided Poly (Glycolide-co-L-lactide)
(PGLA) Sutures, With and Without Needle - Class III
- Sterile, Absorbable, Synthetic Polydioxanone Monofilament (PDO)
Sutures With and Without Needle - Class III

GMDN : 13908, 17471, 16584

Product Types are attached.

Certificate Number : M.2020.106.13286

Report Number : MD.3894.IB

Initial Assessment Date : 29.05.2019

Registration Date : 31.01.2020

Revision Date /No : -

Expiry Date : 27.05.2024


UDEM International Certification
Auditing Training Centre Industry
and Trade Inc. Co.

UDEM hereby declares that the requirements of Annex II, excluding section 4 of the 93/42/EEC Directive have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance audits, defined by Annex II, section 5 of the forementioned directive. According to Annex II, section 4 an EC design- examination certificate is required for placing the Class III devices on the market. UDEM's responsibility for class I devices covered by the EC certificate is limited to manufacturing issues related to safeguarding and maintaining sterile conditions, if the dev.ce is sterile; and manufacturing issues related to product's conformity with metrological requirements, if it has measurement function. This certificate remains as the property of UDEM International Certification Auditing Training Centre Industry and Trade Inc. Co. to whom it must be returned upon request. The above named company and UDEM must keep a copy of this certificate for 5 years from the registration of the certificate. Usage of the CE mark is under the responsibility of the manufacturer with the completion of EC Declaration of Conformity. The above mentioned company must notify all changes related with the approved product to UDEM. If UDEM will not renew the expiry date of this certificate in question, the mentioned company should stop placing the product on the market. The validity of the certificate can be checked through www.udem.com.tr.




Address: Mutlukent Mahallesi 2073 Sokak (Eski 93 Sokak) No:10 Çankaya – Ankara – TURKEY

Phone: +90 0312 443 03 90 **Fax:** +90 0312 443 03 76

E-mail: info@udemltd.com.tr www.udem.com.tr

模型图谱 Moulds



采玉牙科
CAIYU DENTAL



- 纳米塑钢牙引进世界专业高端制造设备和高科技成型工艺
Nano tech tooth introduced the world professional high-end manufacturing equipment and high-tech molding process
- 使用自主研发的纳米级互穿网络聚合物和超微填料
Use of independent research and development of nanoscale interpenetrating polymer networks and superfine packing
- 独特设计的前牙切端和后牙的颌面加厚釉质层，具有高冲击性和高耐磨性
Unique design of the front teeth cut terminal and the maxillofacial thickening of the back teeth enamel layer, high impact and high wear resistance
- 基层选用优质丙烯酸树脂材料组成，能够与基托树脂牢固结合
Acrylic resin base layer selection of quality materials, with the base resin solid combination
- 模仿天然牙的多层色结构、使层次达到完美逼真的效果
Imitation of natural teeth multilayer structure, make the levels reach perfect vivid effect

河南采玉牙科材料有限公司
CAI YU DENTAL MATERIALS CO., LTD.

地址：滑县产业集聚区标准化厂房35号
Add: no.35 standardized plant, industrial agglomeration area, hua county
电话/Tel: 0372-5501012
传真/Fax: 0372-5502012
网址/http: www.cyykcl.com www.zgcyyc.com.cn
邮箱/E-mail: zgcyyk@163.com caiYu.dental@gmail.com

纳米塑钢牙
NANO-TECH TEETH



前牙图谱 Anteriors



HSS3
L1 39.6
H1 9.3
L2 7.6
L3 31.5
H2 8.0
L4 4.9



HC3
L1 40.0
H1 9.5
L2 7.6
L3 33.0
H2 8.6
L4 5.1



HSS4
L1 41.9
H1 9.4
L2 8.0
L3 33.2
H2 8.5
L4 5.2



HC4
L1 43.6
H1 10.5
L2 8.0
L3 34.3
H2 8.9
L4 5.3



HSS5
L1 44.3
H1 9.9
L2 8.5
L3 35.5
H2 9.0
L4 5.6



HC5
L1 45.6
H1 10.9
L2 8.5
L3 36.3
H2 9.6
L4 5.6



HSS6
L1 47.5
H1 11.0
L2 9.0
L3 43.5
H2 9.5
L4 6.0



HC6
L1 48.0
H1 11.5
L2 9.0
L3 43.5
H2 10.0
L4 6.9

后牙图谱 Posteriors



M28
L5 25.0
L6 8.5
L7 30.5
L8 8.5



M30
L5 30.0
L6 9.0
L7 32.5
L8 9.0



M32
L5 32.0
L6 10.0
L7 34.5
L8 10.0



黑鹰/HEI YING (二层色/Two Layers)
颜色/Shade: Vita色系 A1、A2、A3、A3.5
包装/Packing: 全口/full set
28x1x3 84颗/盒/Pcs/box



红鹰/HONG YING (三层色/Three layers)
颜色/Shade: Vita色系 A1、A2、A3、A3.5
包装/Packing: 全口/full set
28x1x3 84颗/盒/Pcs/box



宝鹰/BAO YING (三层色/Three layers)
颜色/Shade: Vita色系 A1、A2、A3、A3.5
包装/Packing: 全口/full set
28x1x3 84颗/盒/Pcs/box



全球鹰/GUAN QIU YING (三层色/Three layers)
颜色/Shade: Vita色系 A1、A2、A3、A3.5
包装/Packing: 全口/full set
28x1 28颗/盒/Pcs/box



猎鹰/LIEYING 纳米瓷钢牙, 选用快速口原
材料和自主研发的超微填料注塑新工艺, 使产品
达到了高硬度、高韧性、耐磨耗等奇特效果。
包装/Packing: 全口/full set
28x1 28颗/盒/Pcs/box



前皇后钢牙/QIAN SU HOU GANG YA
瓷厂设计的新品种之一, 前牙采用有机复合填
料, 使牙齿内部组织致密, 韧性功能特征, 实
现了前牙与后牙之间生理性需求的完美结合。
包装/Packing: 全口/full set
28x1 28颗/盒/Pcs/box




前牙图谱 Anteriors


	390 L1 39.9 H1 9.4 L2 7.5		U25 L3 33.8 H2 8.9 L4 5.3
	419 L1 43.9 H1 9.4 L2 8.1		U33 L3 34.6 H2 10.1 L4 5.2
	425 L1 43.9 H1 10.1 L2 8.3		U34 L3 36.3 H2 8.6 L4 5.4
	426 L1 44.1 H1 9.0 L2 8.2		U37 L3 34.8 H2 9.1 L4 5.3
	438 L1 46.0 H1 10.6 L2 8.5		U38 L3 35.7 H2 9.2 L4 5.5
	447 L1 47.3 H1 11.1 L2 8.6		U67 L3 36.3 H2 10.0 L4 5.7
	448 L1 47.4 H1 10.0 L2 8.8		U88 L3 40.0 H2 11.4 L4 6.0
	468 L1 46.0 H1 10.9 L2 8.9		U90 L3 40.7 H2 8.8 L4 6.1
	501 L1 50.1 H1 13.2 L2 9.6		503 L1 52.8 H1 11.0 L2 9.6

后牙图谱 Posteriors


	28 L5 28.3 L6 8.9 L7 28.9 L8 7.5
	30 L5 30.0 L6 8.0 L7 31.7 L8 7.7
	32 L5 30.7 L6 8.8 L7 33.5 L8 8.2
	34 L5 32.5 L6 9.8 L7 35.3 L8 8.7
	36 L5 34.2 L6 10.0 L7 37.5 L8 9.3



玉唐/YU YING (二层色/Two Layers)
颜色/Shade: Vita色系 A1、A2、A3、A3.5
包装/Packing: 全口/Full set
28×1×3 84颗/盒/Pcs/box



银唐/YIN YING (三层色/Three layers)
颜色/Shade: Vita色系 A1、A2、A3、A3.5
包装/Packing: 全口/Full set
28×1×3 84颗/盒/Pcs/box



金唐/JIN YING (三层色/Three layers)
颜色/Shade: Vita色系 A1、A2、A3、A3.5
包装/Packing: 全口/Full set
28×1×3 84颗/盒/Pcs/box



玉鹰/YU YING (二层色/Two Layers)
型号/Model: 390, 419, 425, 426, 436, 447, 448, 466, 501, 503
颜色/Shade: Vita色系 A1, A2, A3, A3.5
包装/Packing: 前牙/Anteriors 6×1×16 96颗/盒/Pcs/box



玉鹰/YU YING (二层色/Two Layers)
型号/Model: 390, 419, 425, 426, 436, 447, 448, 466, 501, 503
颜色/Shade: Vita色系 A1, A2, A3, A3.5
包装/Packing: 前牙/Anteriors 2×1×16 64颗/盒/Pcs/box



玉鹰/YU YING (二层色/Two Layers)
型号/Model: 28, 30, 32, 34, 36
颜色/Shade: Vita色系 A1, A2, A3, A3.5
包装/Packing: 后牙/Posteriors 单六 6×1×12 72颗/盒/Pcs/box




玉鹰/YU YING (二层色/Two Layers)
型号/Model: 28, 30, 32, 34, 36
颜色/Shade: Vita色系 A1, A2, A3, A3.5
包装/Packing: 后牙/Posteriors 8×1×12 96颗/盒/Pcs/box



合成树脂牙/SYNTHETIC RESIN TEETH
二层色/Two Layers
型号/Model: 22, 23, 24, 25
颜色/Shade: Vita色系 A1, A2, A3, A3.5
包装/Packing: 全口/full set
28×1×10 280颗/盒/Pcs/box



予成牙列/Manufactured Artificial Dentition
二层色/Two Layers
型号/Model: 22, 23, 24
颜色/Shade: Vita色系 A1, A2, A3, A3.5
包装/Packing: 全口/full set
28×1×10 280颗/盒/Pcs/box



合成树脂牙/SYNTHETIC RESIN TEETH
单层色/Single layer color
型号/Model: 22, 23, 24
颜色/Shade: Vita色系 A1, A2, A3, A3.5
包装/Packing: 全口/full set
28×1×10 280颗/盒/Pcs/box

EC Certificate
Directive 93/42/EEC Annex V
Production Quality Assurance
Medical Devices

Registration No.: DD 60146058 0001

Report No.: 50224950 002

Manufacturer: CAI YU DENTAL MATERIALS CO.,LTD
No.35,standardization workshop of hua
county industrial gathering area,Hua county,
456400 Henan
P.R. China

Products: Synthetic Resin Teeth

Expiry Date: 2024-05-26

The Notified Body hereby declares that the requirements of Annex V of the directive 93/42/EEC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex V, section 4 of the aforementioned directive. For placing on the market of class IIb and class III devices covered by this certificate an EC type-examination certificate according to Annex III is required.

Effective Date: 2020-01-20

Date: 2020-01-20



TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg
TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 93/42/EEC concerning medical devices with the identification number 0197.



TÜVRheinland®
LGA

Precisely Right.

Business Stream Products
Certification Department

TÜV Rheinland LGA Products GmbH · 90431 Nürnberg

CAI YU DENTAL MATERIALS CO.,LTD
county industrial gathering area,Hu
456400 HENAN
P.R. CHINA

Contact

Tel. +49 911 655-5225
Mail service@de.tuv.com

Date January 20, 2020

Application for : QMS Produktion, Anhang V MDD
Certificate No. : DD 60146058 Sheet 0001
Device : Only for QM-System audit
Test requirement : Richtlinie 93/42/EWG

Dear Madame or Sir,

Your Quality Management System has been tested and found to be in accordance with the above mentioned requirements.

Enclosed please find the certificate
No. DD 60146058 0001.

Kind regards

Certification body

Jing Zhang

Test sample: no, documentation available

TÜV Rheinland
LGA Products GmbH

Tillystraße 2
90431 Nürnberg

Tel. +49 911 655-5225
Fax +49 911 655-5226
Mail service@de.tuv.com
Web www.tuv.com/safely

Board of Management

Dipl.-Ing.
Jörg Mähler, Spokesman

Dipl.-Kfm.
Dr. Jörg Schlösser

Chairman of the
Supervisory Board

Dipl.-Ing.
Ralf Scheller

Nuremberg HRB 26013
VAT No.: DE 611835490

Certificate

The Certification Body of
TÜV Rheinland LGA Products GmbH

hereby certifies that the organization

CAI YU DENTAL MATERIALS CO.,LTD
No.35,standardization workshop of hua
county industrial gathering area,Hua county,
456400 Henan
P.R. China

has established and applies a quality management system for medical devices
for the following scope:

Manufacture and Distribution of Synthetic Resin Teeth

Proof has been furnished that the requirements specified in

EN ISO 13485:2016

are fulfilled. The quality management system is subject to yearly surveillance.

Effective Date: 2020-01-20
Certificate Registration No.: SX 60146060 0001
An audit was performed. Report No.: 50224950 002
This Certificate is valid until: 2023-01-17

Certification Body



Date 2020-01-20



TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg
Tel.: +49 221 806-1371 Fax: +49 221 806-3935 e-mail:cert-validity@de.tuv.com <http://www.tuv.com/safety>



TÜVRheinland®
LGA

Precisely Right.

Business Stream Products
Certification Department

TÜV Rheinland LGA Products GmbH · 90431 Nürnberg

CAI YU DENTAL MATERIALS CO.,LTD
county industrial gathering area,Hu
456400 HENAN
P.R. CHINA

Contact

Tel. +49 911 655-5225
Mail service@de.tuv.com

Date January 20, 2020

Application for : QMS
Certificate No. : SX 60146060 Sheet 0001
Device : Only for QM-System audit
Test requirement : EN ISO 13485:2016

Dear Madame or Sir,

Your Quality Management System has been tested and found to be in accordance with the above mentioned requirements.

Enclosed please find the certificate
No. SX 60146060 0001.

Kind regards

Certification body


Jing Zhang

Test sample: no, documentation available

TÜV Rheinland
LGA Products GmbH

Tillystraße 2
90431 Nürnberg

Tel. +49 911 655-5225
Fax +49 911 655-5226
Mail service@de.tuv.com
Web www.tuv.com/safety

Board of Management

Dipl.-Ing.
Jörg Mähler, Spokesman

Dipl.-Kfm.
Dr. Jörg Schüssler

Chairman of the
Supervisory Board

Dipl.-Ing.
Ralf Scheller

Nuremberg HRB 26013
VAT No.: DE 811835490

EC Certificate
Directive 93/42/EEC Annex V
Production Quality Assurance
Medical Devices

Registration No.: DD 60149915 0001

Report No.: 16805222 008

Manufacturer: Shenzhen Dian Fong Abrasives
Company Limited
West 4th Floor, Block 2
No. 5, TianHua Rd.
PingHu, LongGang District
Shenzhen
518111 Guangdong
P.R. China

Products: Dental Burs

(see attachment for additional site included)

Repalces Approval, Registration No.: DD 60109667 0001

Expiry Date: 2024-05-26

The Notified Body hereby declares that the requirements of Annex V of the directive 93/42/EEC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex V, section 4 of the aforementioned directive. For placing on the market of class IIb and class III devices covered by this certificate an EC type-examination certificate according to Annex III is required.

Effective Date: 2020-06-30

Date: 2020-06-30



TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg
TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 93/42/EEC concerning medical devices with the identification number 0197.

Certificate

The Certification Body of
TÜV Rheinland LGA Products GmbH

hereby certifies that the organization
**Shenzhen Dian Fong Abrasives
Company Limited**
West 4th Floor, Block 2
No. 5, TianHua Rd.
PingHu, LongGang District
Shenzhen
518111 Guangdong
China

has established and applies a quality management system for medical devices
for the following scope:

Manufacture and Distribution of Dental Burs
(see attachment for additional site included)

Proof has been furnished that the requirements specified in

EN ISO 13485:2016

are fulfilled. The quality management system is subject to yearly surveillance.

Effective Date: 2019-11-01
Certificate Registration No.: SX 60135453 0001
An audit was performed. Report No.: 16805222 004
This Certificate is valid until: 2022-05-22

Certification Body



Date 2019-11-01



TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg
Tel.: +49 221 806-1371 Fax: +49 221 806-3935 e-mail: cert-validity@de.tuv.com http://www.tuv.com/safety



EC Certificate
Directive 93/42/EEC Annex V
Production Quality Assurance
Medical Devices

Registration No.: DD 60146058 0001

Report No.: 50224950 002

Manufacturer: CAI YU DENTAL MATERIALS CO.,LTD
No.35,standardization workshop of hua
county industrial gathering area,Hua county,
456400 Henan
P.R. China

Products: Synthetic Resin Teeth

Expiry Date: 2024-05-26

The Notified Body hereby declares that the requirements of Annex V of the directive 93/42/EEC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex V, section 4 of the aforementioned directive. For placing on the market of class IIb and class III devices covered by this certificate an EC type-examination certificate according to Annex III is required.

Effective Date: 2020-01-20

Date: 2020-01-20



TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg
TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 93/42/EEC concerning medical devices with the identification number 0197.



TÜVRheinland®
LGA

Precisely Right.

Business Stream Products
Certification Department

TÜV Rheinland LGA Products GmbH · 90431 Nürnberg

CAI YU DENTAL MATERIALS CO.,LTD
county industrial gathering area,Hu
456400 HENAN
P.R. CHINA

Contact

Tel. +49 911 655-5225
Mail service@de.tuv.com

Date January 20, 2020

Application for : QMS Produktion, Anhang V MDD
Certificate No. : DD 60146058 Sheet 0001
Device : Only for QM-System audit
Test requirement : Richtlinie 93/42/EWG

Dear Madame or Sir,

Your Quality Management System has been tested and found to be in accordance with the above mentioned requirements.

Enclosed please find the certificate
No. DD 60146058 0001.

Kind regards

Certification body

Jing Zhang

Test sample: no, documentation available

TÜV Rheinland
LGA Products GmbH

Tillystraße 2
90431 Nürnberg

Tel. +49 911 655-5225
Fax +49 911 655-5226
Mail service@de.tuv.com
Web www.tuv.com/safely

Board of Management


Dipl.-Ing.
Jörg Mähler, Spokesman

Dipl.-Kfm.
Dr. Jörg Schlösser

Chairman of the
Supervisory Board

Dipl.-Ing.
Ralf Scheller

Nuremberg HRB 26013
VAT No.: DE 611835490

 Kazan Medical Instruments Plant	Declaration of Conformity	DoC 2017 vs. 01 TF 4
		Page: 1of 6

DECLARATION OF CONFORMITY

1) Manufacturer: KAZAN MEDICAL INSTRUMENTS PLANT JSC (KMIP JSC)
Address: No.12, Salikh Saidashev str., Kazan, 420021, Republic of Tatarstan, RUSSIA

2) European authorized representative: SKIRGESA JSC
Address: ENERGETIKU STR.8, LT-52461, KAUNAS, LITHUANIA
On product labels printed as:
 SKIRGESA JSC, ENERGETIKU STR.8, LT-52461, KAUNAS, LITHUANIA
 FAX. +370 37458161, E-mail: info@skirgesa.lt

3) Product(s) (name, type or model/batch number, etc.):

Hard - alloy dental laboratory cutters TS 9433-103-05519988-2005E <i>see appendix</i>
--

4) The product(s) described above is in conformity with:

DIRECTIVES

General Applicable Directive:

Medical Device Directive: COUNCIL DIRECTIVE 93/42/EEC of 14 june 1993 Concerning medical device (MDD93/42/EEC)

Document No.	Edition / Date of issue
BS EN 62366	2008-04-30
EN ISO 14971:2012	2012-08-30
EN 1639:2009	2010-07-07
EN 1041:2008	2009-02 19
EN ISO 17664:2004	2005-09-30
EN ISO 17665-1:2006	2006-11-15

5) Additional information (conformity procedure, Registration nr., etc.):

Conformity assessment procedure for CE marking: Medical Device Directive, Annex VII

KAZAN, RUSSIA, 2017-05-19
 (Place & date of issue (yyyy-mm-dd))




Shokirov Nur-Khanzinovich, Director-General of KMIP JSC
 (name; function and signature of manufacturer)

KAUNAS, LITHUANIA, 2017-05-19
 (Place & date of issue (yyyy-mm-dd))

SKIRMANTAS AKELIS, DIRECTOR of SKIRGESA JSC
 (name; function and signature of authorized representative)



 Kazan Medical Instruments Plant	Declaration of Conformity	DoC 2017 vs. 01 TF 4
		Page: 2 of 6

Appendix

Date: 2017-05-19

List of devices.

Device name	Type/ model/ref number	Risk class / rule ¹	First date of CE-marking
Hard - alloy dental laboratory cutters	see page 3-6	Class I / Rule1	2006-07-17

¹ See risk classification in Medical Device Directive, annex IX

СЕРТИФИКАТ №.
CERTIFICATE N. 9124.KMIP

НАСТОЯЩИМ ПОДТВЕРЖДАЕМ, ЧТО СИСТЕМА МЕНЕДЖМЕНТА КАЧЕСТВА, ДЕЙСТВУЮЩАЯ В
WE HEREBY CERTIFY THAT THE QUALITY MANAGEMENT SYSTEM OPERATED BY

АКЦИОНЕРНОЕ ОБЩЕСТВО
«КАЗАНСКИЙ МЕДИКО-ИНСТРУМЕНТАЛЬНЫЙ ЗАВОД»

РОССИЙСКАЯ ФЕДЕРАЦИЯ, РЕСПУБЛИКА ТАТАРСТАН, 420021,
Г. КАЗАНЬ, УЛ. САЛИХА САЙДАШЕВА, Д. 12

KAZAN MEDICAL INSTRUMENTS PLANT
JOINT STOCK COMPANY

12, SALIKH SAIDASHEV STR., KAZAN, 420021,
THE REPUBLIC OF TATARSTAN, RUSSIAN FEDERATION

соответствует требованиям стандарта
IS IN COMPLIANCE WITH THE STANDARD

EN ISO 13485:2012

для следующей области деятельности
FOR THE FOLLOWING ACTIVITIES

Разработка и производство вращающихся стоматологических инструментов;
эндоканальных стоматологических инструментов; микрохирургических
инструментов; стоматологических наконечников и микромоторов
*Design and manufacture of dental rotating tools, dental endodontic instruments,
microsurgical instruments, dental handpieces and micromotors*

Дальнейшие разъяснения относительно применимости требований EN ISO 13485:2012 могут быть получены от организации
Further clarifications regarding the applicability of EN ISO 13485:2012 requirements may be obtained by consulting the organization

ИСПОЛЬЗОВАНИЕ И СРОК ДЕЙСТВИЯ СЕРТИФИКАТА ДОЛЖНЫ УДОВЛЕТВОРЯТЬ
ТРЕБОВАНИЯМ ПРАВИЛ ПО СЕРТИФИКАЦИИ СИСТЕМ МЕНЕДЖМЕНТА
THE USE AND THE VALIDITY OF THE CERTIFICATE SHALL SATISFY THE REQUIREMENTS
OF THE RULES FOR CERTIFICATION OF MANAGEMENT SYSTEMS

ДАТА: ПЕРВИЧНАЯ СЕРТИФИКАЦИЯ DATE: FIRST CERTIFICATION	ДАТА ВЫПУСКА CURRENT ISSUE	СРОК ДЕЙСТВИЯ EXPIRY DATE
2017-01-19	2017-01-19	2020-01-18

Организация обязана пройти процесс сертификации в соответствии с (требованиями)
ISO 13485:2016 до 28.02.2019; в противном случае срок действия данного сертификата истекает
The Organization shall obtain the certification according to ISO 13485:2016
within 2019/02/28; otherwise the validity of this certificate will expire



IMQ S.p.A. - VIA QUINTILIANO, 43 - 20138 MILANO
Management Systems Division - Flavio Ornago

Only the English text is authentic

CISQ è la Federazione Italiana di
Organismi di Certificazione dei
sistemi di gestione aziendale.

CISQ is the Italian Federation
of management system
Certification Bodies.

Certificate

The Certification Body of
TÜV Rheinland LGA Products GmbH

hereby certifies that the organization

CAI YU DENTAL MATERIALS CO.,LTD
No.35,standardization workshop of hua
county industrial gathering area,Hua county,
456400 Henan
P.R. China

has established and applies a quality management system for medical devices
for the following scope:

Manufacture and Distribution of Synthetic Resin Teeth

Proof has been furnished that the requirements specified in

EN ISO 13485:2016

are fulfilled. The quality management system is subject to yearly surveillance.

Effective Date: 2020-01-20
Certificate Registration No.: SX 60146060 0001
An audit was performed. Report No.: 50224950 002
This Certificate is valid until: 2023-01-17

Certification Body



Date 2020-01-20



TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg
Tel.: +49 221 806-1371 Fax: +49 221 806-3935 e-mail:cert-validity@de.tuv.com <http://www.tuv.com/safety>



TÜVRheinland®
LGA

Precisely Right.

Business Stream Products
Certification Department

TÜV Rheinland LGA Products GmbH · 90431 Nürnberg

CAI YU DENTAL MATERIALS CO.,LTD
county industrial gathering area,Hu
456400 HENAN
P.R. CHINA

Contact

Tel. +49 911 655-5225
Mail service@de.tuv.com

Date January 20, 2020

Application for : QMS
Certificate No. : SX 60146060 Sheet 0001
Device : Only for QM-System audit
Test requirement : EN ISO 13485:2016

Dear Madame or Sir,

Your Quality Management System has been tested and found to be in accordance with the above mentioned requirements.

Enclosed please find the certificate
No. SX 60146060 0001.

Kind regards

Certification body


Jing Zhang

Test sample: no, documentation available

TÜV Rheinland
LGA Products GmbH

Tillystraße 2
90431 Nürnberg

Tel. +49 911 655-5225
Fax +49 911 655-5226
Mail service@de.tuv.com
Web www.tuv.com/safety

Board of Management

Dipl.-Ing.
Jörg Mähler, Spokesman

Dipl.-Kfm.
Dr. Jörg Schüssler

Chairman of the
Supervisory Board

Dipl.-Ing.
Ralf Scheller

Nuremberg HRB 26013
VAT No.: DE 811835490



Certificate

ISO 13485 : 2016

Medico (Huaian) Co., Ltd.

No.9, South Guangzhou Road, Huaian, 223000 Jiangsu, CHINA

This certificate shows that the medical devices - quality management system (EN ISO 13485:2016) of the above company was approved by PCA Certification for the following scope, the validity of the certificate depends on the company's pass the annual surveillance audits and company's maintenance the related management system conditions according to international accreditation criteria

SCOPE

*Design & development, production and sales of absorbable sutures and non-absorbable sutures with or without needle,
Design & development, production and sales of lubricating jelly*

GROUP CODE

M01 - MDS7006

Certificate No : TC-75132
Registration Date : 15.05.2019
Reissue Date : 28.04.2020
Expiry Date : 14.05.2021
Certificate Period : 3 Years (Form the date of registration)
Exclusion : 7.5.3 / 7.5.4 / 7.5.10

IAS
ACCREDITED
Management
Systems
Certification Body
MSCB-103

PCA Certification Approval

PCA Sertifikasyon Hizmetleri Limited Şirketi
Orta Mah. Ordu Sk. İzpark C Blok No:26/23 Kartal / İSTANBUL
Tel: +90 216 510 63 48-49 Pbx Faks: +90 216 517 63 49
www.pca-tr.com info@pca-tr.com



EC Certificate
Directive 93/42/EEC Annex V
Production Quality Assurance
Medical Devices

Registration No.: DD 60149915 0001

Report No.: 16805222 008

Manufacturer: Shenzhen Dian Fong Abrasives
Company Limited
West 4th Floor, Block 2
No. 5, TianHua Rd.
PingHu, LongGang District
Shenzhen
518111 Guangdong
P.R. China

Products: Dental Burs

(see attachment for additional site included)

Repalces Approval, Registration No.: DD 60109667 0001

Expiry Date: 2024-05-26

The Notified Body hereby declares that the requirements of Annex V of the directive 93/42/EEC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex V, section 4 of the aforementioned directive. For placing on the market of class IIb and class III devices covered by this certificate an EC type-examination certificate according to Annex III is required.

Effective Date: 2020-06-30

Date: 2020-06-30



TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg
TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 93/42/EEC concerning medical devices with the identification number 0197.

Certificate

The Certification Body of
TÜV Rheinland LGA Products GmbH

hereby certifies that the organization
**Shenzhen Dian Fong Abrasives
Company Limited**
West 4th Floor, Block 2
No. 5, TianHua Rd.
PingHu, LongGang District
Shenzhen
518111 Guangdong
China

has established and applies a quality management system for medical devices
for the following scope:

Manufacture and Distribution of Dental Burs
(see attachment for additional site included)

Proof has been furnished that the requirements specified in

EN ISO 13485:2016

are fulfilled. The quality management system is subject to yearly surveillance.

Effective Date: 2019-11-01
Certificate Registration No.: SX 60135453 0001
An audit was performed. Report No.: 16805222 004
This Certificate is valid until: 2022-05-22

Certification Body



Date 2019-11-01



TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg
Tel.: +49 221 806-1371 Fax: +49 221 806-3935 e-mail: cert-validity@de.tuv.com http://www.tuv.com/safety

Dianfang ABRASIVES



ISO 13485
CE 0197

点锋研磨

第6版/Vol.6

Rotary things



laboratory

About Us , 关于我们

Dian Fong specialize in Rotary things of dental and Jewelry, who provide innovative and professional products and service.

Our products include DIAMOND bur and disc, TUNGSTEN CARBIDE bur or cutter, FLEXIBLE RUBBER bur and wheel, Polishing Brush, Mounted Stone, ZIRCONIA grinder and polish etc.

We make sure of the best quality and have certificates of ISO 13485 and CE0197. Whoever you are dentist or technician, we always make you relaxing, safe and efficient. We are your best friend in dentistry!

Content 目录 6.2

Laboratory, Technician's Use . 加工厂类

Diamond bur HP . 金刚石磨头	3
Fully sintered Diamond bur HP . 全砂烧结磨头	4
Diamond Disc , Plaster Trimming wheel. 金刚石切片,石膏模修整盘	6
Tungsten Carbide Bur, HP . 钨钢车针, 钨钢磨头	8
CAD/CAM Milling bur . 铣刀	10
Diamond stones and polishers for Zirconia . 锆牙磨头,抛光轮	12
Silicon, Rubber Polishing bur, wheel & point 硅胶磨头,胶轮胶粒	14
Mounted stone, Sandpaper bur, Mandrels . 车石, 砂纸棒,夹针	15
Polishing brush, Felt Bob & bur . 抛光刷,毛毡磨头	16
Technician or DIY Kit 加工厂, DIY套装	18
FlexSoft Partial Denture 隐形义齿胶	19
Samples Collection 样品册	5,8,15

一首小诗

我磨一根针，
足以穿风尘。
海内结知己，
笑淡了风云。

A Little Poem

A small bur which I made,
is sharp to teeth-grind ;
Anywhere, we make friends,
and smile as sky fade .



Diamond Bur HP
(Electroplated)
金刚石磨头
(电镀)



Sharp and Durable,
Grind or polish Ceramics
锋利,实用,磨削平稳
适合车瓷,不留黑印



Grit

- Black / Super Coarse
- Green / Coarse
- Blue / Medium
- Red / Fine

Fine, Red
细砂,红圈

	001 3.0 3.0	010 3.3 1.2	111 2.3 8.0	111 3.3 8.0	143 1.8 15.0	141 2.5 8.0	161 2.7 15.0	198 3.3 12.0	225 2.3 6.0	239 3.3 8.0	292 2.3 15.0	304 4.0 1.2
编号 No.	FD30	FK33	FA23	FA33	FC18	FC25	FM27	FL33	FN23	FP33	FQ23	FU40

For grind Zirconia 可用来车氧化锆

Medium Blue
标准砂,蓝圈

	001 1.0	1.4	1.8	2.3	2.7	3.5	4.0	010 1.6	2.1	2.7	3.7	4.7
编号 No.	MD10	MD14	MD18	MD23	MD27	MD35	MD40	MK16	MK21	MK27	MK37	MK47
	Round / 圆球						Inverted Cone / 倒锥					

6.0	8.0	10.0	041 6.0 1.0	068 6.0 2.0	111 1.6 8.0	112 2.3 10.0	3.3 10.0	5.0 11.0	143 1.8 15.0	141 2.5 10.0	3.0 10.0	4.0 10.0
MK60	MK80	MK100	MT60	MJ60	MA16	MA23	MA33	MA50	MC18	MC25	MC30	MC40
			Wheel / 平轮			Cylinder / 平柱			Round Cylinder / 圆头柱			

160 1.2 10.0	1.8 10.0	161 2.7 15.0	3.7 10.0	5.0 12.0	171 1.2 6.0	1.6 6.0	173 2.1 10.0	3.1 10.0	4.0 11.0	198 1.6 10.0	2.3 12.0	3.3 11.0
MM12	MM18	MM27	MM37	MM50	MB12	MB16	MB21	MB31	MB40	ML16	ML23	ML33
Needle / 尖针				Taper / 平锥								Round Taper / 圆头锥

4.0 12.0	5.0 13.0	225 1.8 5.0	2.3 6.0	239 3.3 8.0	250 1.6 10.0	2.3 10.0	5.0 12.0	243 2.3 5.0	266 6.0 18.0	292 2.3 14.0	304 5.0 1.5	10.0 2.5
ML40	ML50	MN18	MN23	MP33	MG16	MG23	MG50	MR23	MV60	MQ23	MU50	MU100
		Inverted Cone / 倒锥		Pearl	Flame / 火焰				R.Taper	Torpedo	Wheel / 火焰	

Super Coarse Black
特粗砂,黑圈

001 3.5	112 5.0 11.0	141 4.0 10.0	238 7.0 12.5	251 5.0 12.0	263 7.0 12.0	266 6.0 18.0	274 6.0 14.0					
编号 No.	SD35	SA50	SC40	SP70	SG50	SV70	SV60	SF60				

Coarse Green
粗砂,绿圈

001 2.3	3.5	068 6.0 2.0	112 5.0 11.0	143 1.8 15.0	141 4.0 10.0	173 3.1 10.0	198 4.0 12.0	292 2.3 14.0	266 6.0 18.0	274 6.0 14.0	
编号 No.	CD23	CD35	CJ60	CA50	CC18	CC40	CB31	CL40	CQ23	CV60	CF60

Fully-Sintered
Diamond. HP
全砂烧结
磨头



Extremely Durable,
Grind or polish smoothly
especial for metal
从里到外都是砂
极其耐用,磨削平稳
实惠,尤其适合车金属



Grit

- Black / Super Coarse
- Green / Coarse
- Blue / Medium
- Red / Fine
- Yellow / Extra Fine

Super Coarse
Black
特粗砂
黑圈
ISO 807.104.544



编号 No.	D060S	D070S	K080S	K100S	J070S	J100S	A050S	A060S	C050S	C060S	C070S	M050S
--------	-------	-------	-------	-------	-------	-------	-------	-------	-------	-------	-------	-------



编号 No.	M070S	B050S	B070S	L050S	L060S	L070S	G060S	F050S	F070S	R050S	R070S	
--------	-------	-------	-------	-------	-------	-------	-------	-------	-------	-------	-------	--

Coarse
Green
粗砂
绿圈
ISO 807.104.534



编号 No.	D050C	D060C	J060C	J070C	A037C	A050C	A060C	C037C	C050C	C060C	B037C	B050C
--------	-------	-------	-------	-------	-------	-------	-------	-------	-------	-------	-------	-------



编号 No.	B060C	L037C	L050C	L060C	N050C	G050C	F037C	F050C	R037C	R050C		
--------	-------	-------	-------	-------	-------	-------	-------	-------	-------	-------	--	--

Medium
Blue
标准砂
蓝圈
ISO 807.104.524



编号 No.	D031	D037	D050	D060	K031	K037	K050	K060	J060	J070	A037	A050
--------	------	------	------	------	------	------	------	------	------	------	------	------





编号 No.	C031	C037	C050	M031	M037	M050	B037	B050	L037	L050	L060	N037	G037
--------	------	------	------	------	------	------	------	------	------	------	------	------	------



编号 No.	G050	F037	F050	R037	R050								
--------	------	------	------	------	------	--	--	--	--	--	--	--	--

Fine Red
细砂红圈
ISO 807.104.514



Extra Fine Yellow
特细砂黄圈
ISO 807.104.504



Samples Collection

样品册

- Flip-Cart Plastic Pages ;
- Smart, Easily fill all burs ;
- Economic and Recyclable ;
- Nice-looking for showing to Customers ;
- Convenient to take to anywhere .



No. 171B5

Empty Collection Wallet
样品本(不带车针)



4 Colors Available ,
100 holes per page,
5 pages per set,
Smart to HP Burs



No. 1704D

Sintered Diamond Burs Collection
全砂烧结磨头样品本



115 specifications ,
Fully Sintered Diamond,
HP shanks




Diamond Disc

金刚石切片



How to order these descriptions:



Extra Thin, 0.15mm
Double sided
Extra Fine
超薄, 双面砂
ISO 806.900.504



编号 No.	A19D15	A22D15	Am13D15	Am16D15	Am19D15	Am22D15
--------	--------	--------	---------	---------	---------	---------

Thin, 0.20 x 19mm
Double sided
Fine Grit
较薄, 双面砂
ISO 806.900.504



编号 No.	A19D20	B19D20	C19D20	N19D20	Am19D20	Bm19D20	Em19D20	Hm19D20
--------	--------	--------	--------	--------	---------	---------	---------	---------

Thin, 0.20 x 22mm
Double sided
Fine Grit
较薄, 双面砂
ISO 806.900.514



编号 No.	Am22D20	Bm22D20	Cm22D20	Em22D20	Im22D20	Jm22D20	Hm22D20
--------	---------	---------	---------	---------	---------	---------	---------




A22D20	B22D20	C22D20	E22D20	F22D20	H22D20	J22D20	N22D20
--------	--------	--------	--------	--------	--------	--------	--------


Hard, 0.3 - 0.4mm
Double sided
Medium Grit
厚硬, 双面砂
ISO 806.900.524



编号 No.	Am19D30	Am22D30	Am22D40	A19D40	A22D40		
ISO 806.900.	19.0 x 0.3 mm	22.0 x 0.3 mm	22.0 x 0.4 mm	19.0 x 0.4 mm	22.0 x 0.4 mm		




Single-sided
仅单面镀砂




编号 No.	A22S20	Am22S20	Am22S30	Am19S50		
ISO 806.900.	22.0 x 0.2 mm, Fine	22.0 x 0.2 mm	22.0 x 0.3 mm, Medium	19.0 x 0.5 mm		


Fast Cutting
Plaster Disc
0.3mm, Double sided
石膏模切片
切口平直,
无偏摆



编号 No.	C30D30	A35D30	C35D30	J35D30
--------	--------	--------	--------	--------



A40D30	C40D30	J40D30
--------	--------	--------




Cb45D30	Ab45D30	Ch85D30
---------	---------	---------

Plaster Trimming wheel
石膏修整盘

Cuts Fast and last longer
Ensure immediate removal
Used dry or wet on all materials.

快速磨削,耐用;即装即用
干磨水磨都可以.



R10	R12	S10	C10
10' x 1/4' x 1' @ 250 x 6.3 x 25.4mm	12' x 1/4' x 1' @ 300 x 6.3 x 25.4mm	10' x 1/8' x 1' @ 250 x 3.0 x 25.4mm	10' x 1/4' x 1' @ 250 x 6.3 x 25.4mm
Resin Plastic base 塑料板, 轻便, 节能		Steel base 钢板	Silicon Carbide 碳化硅

Please notice Inner hole diameter ; 1.25' = 32 mm is also alternative .
重要提示: 请注意您设备转轴 的直径, 并选择相应的内孔
我们可提供 25mm 和 32mm 的内孔供您选择.

Tungsten Carbide Burs . HP

Cut sharply, precise, Durable

钨钢车针

很锋利,精准,耐用



Shape

- S1, ISO 001, Round
- S2, ISO 002, Inverted cone
- S21, S31, ISO 107, Cylinder
- S33, ISO 168, Taper
- S31R, ISO 137, Round, Cylinder

S1
球形/Round

001 0.5	001 0.6	001 0.7	001 0.8	001 0.9	001 1.0	001 1.2	001 1.4	001 1.6	001 1.8	001 2.1	001 2.3
1005	1006	1007	1008	1009	1010	1012	1014	1016	1018	1021	1023
HP1/4	HP1/2	HP1s	HP 1	HP 1b	HP 2	HP 3	HP 4	HP 5	HP 6	HP 7	HP 8

S2
倒锥/Inv. Cone

001 2.5	001 2.7	001 3.1	010 0.8	010 1.0	010 1.2	010 1.4	010 1.6	010 1.8	010 2.1	010 2.3	S21
1025	1027	1031	2008	2010	2012	2014	2016	2018	2021	2023	S21
HP 9	HP 10	HP 12	HP 34	HP 35	HP 36	HP 37	HP 38	HP 39	HP 40	HP 41	圆柱/Cylinder

S31
圆柱/Cylinder

107 0.8	107 0.9	107 1.0	107 1.2	107 1.4	107 1.6	107 1.8	107 2.1	107 2.3	107 0.8	107 0.9	107 1.0
21008	21009	21010	21012	21014	21016	21018	21021	21023	31008	31009	31010
HP55	HP 56	HP 57	HP 58	HP 59	HP 60	HP 61	HP 62	HP 63	HP 555	HP 556	HP 557

S33
平锥/Taper

107 1.2	107 1.4	107 1.6	107 1.8	107 2.1	107 2.3	168 0.8	168 0.9	168 1.0	168 1.2	168 1.6	168 2.1
31012	31014	31016	31018	31021	31023	33008	33009	33010	33012	33016	33021
HP 558	HP 559	HP 560	HP 561	HP 562	HP 563	HP 698	HP 699	HP 700	HP 701	HP 702	HP 703

S33L
平锥加长/Taper L

168 2.3	171 0.9	171 1.0	171 1.2	171 1.6	171 2.1	S31R					
33023	33L009	33L010	33L012	33L016	33L021	31R010	31R012	31R014			
HP 704	HP 699 L	HP 700 L	HP 701 L	HP 702 L	HP 703 L	HP 1557	HP 1558	HP 1559			

S31R
圆头柱 Round Cyl.

No. 1707T

Tungsten carbide Burs Collection
高品质钨钢磨头样品本

More than 100 specifications, High Quality TC burs, HP Shanks

No. 1705D


HP Diamond Burs, Mounted Stones
金刚石磨头, 车石样品本

More than 100 specifications, Electropolishing Diamond burs HP Shanks

Tungsten carbide
Cutter . HP
高品质
钨钢磨头



Cut sharply,
Precise and Durable
很锋利,精确耐用



Super Fine
Cross Cut
(SFE)

特细交叉纹
黄圈
ISO 500.104.110.



Fine
Cross Cut
(FE)


较细交叉纹
红圈
ISO 500.104.140.



编号 No.	0230.023	1210.016	0130.023	5730.040	7730.023	1030.023	7230.040	0430.023		0220.023	1220.016	0120.023
参考号 Ref.	C023SFE	L016SFE	L023SFE	L040SFE	P023SFE	G023SFE	F040SFE	Q023SFE		C023FE	L016FE	L023FE

Medium
Cross Cut
(CE)

标准交叉纹
蓝圈
ISO 500.104.190.



	5720.040	7720.023	6120.040	5320.060	5420.060	0420.023		7110.023	7110.030	7110.040	7110.050	0310.023	0210.023
	L040FE	P023FE	G040FE	G060FE	F060FE	Q023FE		D023CE	D030CE	D040CE	D050CE	A023CE	C023CE

	5910.060	0110.023	5710.040	5510.060	7710.023	5110.060	1010.023	6110.040	5310.060	7210.040	5410.060	0410.023	
	M060CE	L023CE	L040CE	L060CE	P023CE	P060CE	G023CE	G040CE	G060CE	F040CE	F060CE	Q023CE	

Medium
Plain Cut
(CN)

标准单纹
蓝圈
ISO 500.104.175.



Coarse
Cross Cut
(SCE)

粗糙交叉纹
绿圈
ISO 500.104.220



	7170.030	7170.040	7170.050	5770.040	5570.060	5970.060	5370.060	5470.060		5850.060	5750.040	5550.060
	D030CN	D040CN	D050CN	L040CN	L060CN	M060CN	G060CN	F060CN		C060SCE	L040SCE	L060SCE

Especial
特殊纹



	5150.060	5350.060	5450.060		6925.040	6055.070	7275.060	467.023	118.060			
	P060SCE	G060SCE	F060SCE		L040QSF	L070XCE	F060QSCN	AL023S4	ES801			


CAD / CAM Milling Burs.

铣刀



Advantages

1. Precisely match Various System
2. Working More Efficiently
3. Durable and Valuable
4. R & D Capabilities to New.
5. More Individual Designed



D1 : Cutting Edge Diameter (.mm) L1 : Cutting Edge Length (.mm)
D2 : Shank Diameter (.mm) L2 : Free Grinding to (.mm)
L3 : Total Length (.mm)

Coating

CrN : for PMMA
AL : for Metal
NC : for Zirconia
DC : for Zirconia,
More Durable.



Article No	Coating	D1	L1	L2	L3	D2
for Roland 50						
RL-50-03-NC	NC	0.3	1.5	6.0	50.0	4.0
RL-50-06-CRN	CRN	0.6	1.5	6.0	50.0	4.0
RL-50-06-NC	NC					
RL-50-10A-CRN	CRN	1.0	3.5	16.0	50.0	4.0
RL-50-10A-NC	NC					
RL-50-20A-CRN	CRN	2.0	6.0	16.0	50.0	4.0
RL-50-20A-NC	NC					
RL-50-20A-DC	DC					
RL-50-10B-CRN	CRN	1.0	3.5	20.0	50.0	4.0
RL-50-10B-NC	NC					
RL-50-20B-CRN	CRN	2.0	6.0	25.0	52.0	4.0
RL-50-20B-NC	NC					
RL-50-20B-DC	DC					

Article No	Coating	D1	L1	L2	L3	D2
for 3 M ESPE LAVA						
LV-00-10-CRN	CRN	1.0	3.5	16.0	38.0	3.0
LV-00-10-NC	NC					
LV-00-20-CRN	CRN	2.0	6.0	16.0	38.0	3.0
LV-00-20-NC	NC					
LV-00-20-DC	DC					
LV-01-10-CRN	CRN	1.0	2.0	25.0	50.0	3.0
LV-01-10-NC	NC					
LV-01-20-CRN	CRN	2.0	2.0	25.0	50.0	3.0
LV-01-20-NC	NC					
LV-01-20-DC	DC					

Article No	Coating	D1	L1	L2	L3	D2
for RK						
RK-00-10-DC	DC	1.0	3.5	16.0	47.0	3.0
RK-00-20-DC	DC	2.0	3.5	16.0	47.0	3.0

Article No	Coating	D1	L1	L2	L3	D2
for Imes-Icore 240, 250						
II-24-06-CRN	CrN	0.6	1.5	10.0	48.0	3.0
II-24-06-NC	NC					
II-24-10-CRN	CrN	1.0	5.0	14.0	48.0	3.0
II-24-10-NC	NC					
II-24-25-CRN	CrN	2.5	8.0	20.0	48.0	3.0
II-24-25-NC	NC					
II-24-25-DC	DC					

Article No	Coating	D1	L1	L2	L3	D2
for Imes-Icore 350, 550, 750						
II-75-06-CRN	CrN	0.6	2.5	12.0	53.0	6.0
II-75-06-NC	NC					
II-75-10-CRN	CrN	1.0	6.0	14.0	53.0	6.0
II-75-10-NC	NC					
II-75-25-CRN	CrN	2.5	10.0	20.0	53.0	6.0
II-75-25-NC	NC					
II-75-25-DC	DC					

Article No	Coating	D1	L1	L2	L3	D2
for Amann Girrbach						
AG-00-06-CRN	CrN	0.6	1.5	13.0	47.0	3.0
AG-00-06-NC	NC					
AG-00-10-CRN	CrN	1.0	3.5	16.0	47.0	3.0
AG-00-10-NC	NC					
AG-00-25-CRN	CrN	2.5	5.0	16.0	47.0	3.0
AG-00-25-NC	NC					
AG-00-25-DC	DC					

Article No	Coating	D1	L1	L2	L3	D2
for Amann Girrbach - Metal mill						
AG-01-10-AL	Alnova	1.0	5.0	7.0	43.0	3.0
AG-01-20-AL	Alnova	2.0	5.0	7.0	43.0	3.0

Article No	Coating	D1	L1	L2	L3	D2
for VHF K1, K4						
VH-K4-06-CRN	CRN	0.6	2.0	7.0	35.0	3.0
VH-K4-06-NC	NC					
VH-K4-10-CRN	CRN	1.0	3.5	16.0	35.0	3.0
VH-K4-10-NC	NC					
VH-K4-20-CRN	CRN	2.0	5.0	16.0	35.0	3.0
VH-K4-20-NC	NC					

Article No	Coating	D1	L1	L2	L3	D2
for VHF K5						
VH-K5-03-NC	CRN	0.3	1.2	13.0	40.0	3.0
VH-K5-06-CRN	CRN	0.6	1.5	7.0	40.0	3.0
VH-K5-06-NC	NC					
VH-K5-10-CRN	CRN	1.0	3.5	14.0	40.0	3.0
VH-K5-10-NC	NC					
VH-K5-20-CRN	CRN	2.0	4.5	16.0	40.0	3.0
VH-K5-20-NC	NC					
VH-K5-20-DC	DC					

Article No	Coating	D1	L1	L2	L3	D2
for Dentium						
DT-00-10-NC	NC	1.0	3.5	16.0	50.0	6.0
DT-00-20-NC	NC	2.0	3.5	16.0	50.0	6.0
DT-00-20-DC	DC					
DT-01-05-NC	NC	0.5	2.0	8.0	43.0	3.0
DT-01-10-CrN	CrN	1.0	3.5	16.0	50.0	3.0
DT-01-10-NC	NC					
DT-01-20-CrN	CrN	2.0	3.5	18.0	50.0	3.0
DT-01-20-NC	NC					
DT-01-20-DC	DC					

Article No	Coating	D1	L1	L2	L3	D2
for Zirkon Zahn M3						
ZZ-03-05-CrN	CrN	0.5	1.5	6.0	57.0	3.0
ZZ-03-05-NC	NC					
ZZ-03-10-CrN	CrN	1.0	6.0	12.0	57.0	3.0
ZZ-03-10-NC	NC					
ZZ-03-20-CrN	CrN	2.0	10.0	18.0	57.0	3.0
ZZ-03-20-NC	NC					
ZZ-03-20-DC	DC					

Article No	Coating	D1	L1	L2	L3	D2
for Zirkon Zahn M5, M6, etc						
ZZ-06-05-NC	NC	0.5	1.5	3.0	50.0	6.0
ZZ-06-10-NC	NC	1.0	6.0	12.0	50.0	6.0
ZZ-06-20-NC	NC	2.0	10.0	18.0	50.0	6.0
ZZ-06-20-DC	DC					

Article No	Coating	D1	L1	L2	L3	D2
for Sirona MX 5						
SR-MX-05-NC	NC	0.5	1.5	5.0	42.0	3.0
SR-MX-10-NC	NC	1.0	4.5	16.0	43.0	3.0
SR-MX-25-NC	NC	2.5	7.0	25.0	43.0	3.0
SR-MX-25-DC	DC					

Article No	Coating	D1	L1	L2	L3	D2
for Cercon						
CC-00-05-NC	NC	0.5	2.0	10.0	55.0	3.5
CC-00-10-NC	NC	1.0	8.0	30.0	55.0	3.5
CC-00-20-NC	NC	2.0	10.0	30.0	55.0	3.5

Article No	Coating	D1	L1	L2	L3	D2
for Arum						
AR-00-10-CrN	CrN	1.0	2.0	16.0	55.0	6.0
AR-00-10-NC	NC					
AR-00-20-CrN	CrN	2.0	6.0	20.0	55.0	6.0
AR-00-20-NC	NC					
AR-00-20-DC	DC					

Article No	Coating	D1	L1	L2	L3	D2
for Wieland Mini						
WL-MI-07-CrN	CrN	0.7	1.5	7.0	35.0	3.0
WL-MI-07-NC	NC					
WL-MI-10-CrN	CrN	1.0	3.5	16.0	35.0	3.0
WL-MI-10-NC	NC					
WL-MI-25-CrN	CrN	2.5	3.5	20.0	35.0	3.0
WL-MI-25-NC	NC					

Article No	Coating	D1	L1	L2	L3	D2
for Wieland S1						
WL-S1-07-CrN	CrN	0.7	1.5	7.0	40.0	3.0
WL-S1-07-NC	NC					
WL-S1-10-CrN	CrN	1.0	3.5	16.0	40.0	3.0
WL-S1-10-NC	NC					
WL-S1-25-CrN	CrN	2.5	5.0	20.0	40.0	3.0
WL-S1-25-NC	NC					
WL-S1-25-DC	DC					

Article No	Coating	D1	L1	L2	L3	D2
for Kavo						
KA-00-10-DC	DC	1.0	2.0	16.0	53.0	6.0
KA-00-20-DC	DC	2.0	6.0	19.5	53.0	6.0

Article No	Coating	D1	L1	L2	L3	D2
for Yena						
YN-00-10-NC	NC	1.0	3.5	16.0	45.0	4.0
YN-00-20-NC	NC	2.0	3.5	16.0	45.0	4.0
YN-00-20-DC	DC					

Article No	Coating	D1	L1	L2	L3	D2
for Arum - Metal Mill						
AR-01-10-AL	Alnova	1.0	2.0	10.0	50.0	6.0
AR-01-15-AL	Alnova	1.5	3.0	10.0	50.0	6.0
AR-01-20-AL	Alnova	2.0	4.0	12.0	50.0	6.0
AR-01-30-AL	Alnova	3.0	6.0	12.0	50.0	6.0

The following milling burs are matched for any China's Brand CAD/CAM

Article No	Coating	D1	L1	L2	L3	D2
XT-22-05-CrN	CrN	0.5	1.6	7.0	50.0	3.0
XT-22-05-NC	NC					
XT-22-10-CrN	CrN	1.0	3.5	16.0	50.0	3.0
XT-22-10-NC	NC					
XT-22-20-CrN	CrN	2.0	10.0	22.0	50.0	3.0
XT-22-20-NC	NC					
XT-22-20-DC	DC					

Article No	Coating	D1	L1	L2	L3	D2
SD-00-06-CrN	CrN	0.6	2.0	6.0	38.0	3.175
SD-00-06-NC	NC					
SD-00-10-CrN	CrN	1.0	3.5	16.0	38.0	3.175
SD-00-10-NC	NC					
SD-00-20-CrN	CrN	2.0	3.5	16.0	38.0	3.175
SD-00-20-NC	NC					
SD-00-20-DC	DC					
SD-01-10-CrN	CrN	1.0	3.5	20.0	45.0	3.175
SD-01-10-NC	NC					
SD-01-20-CrN	CrN	2.0	3.5	20.0	45.0	3.175
SD-01-20-NC	NC					
SD-01-20-DC	DC					

Article No	Coating	D1	L1	L2	L3	D2
CR-00-10-CrN	CrN	1.0	3.5	16.0	38.0	3.0
CR-00-10-NC	NC					
CR-00-20-CrN	CrN	2.0	3.5	16.0	38.0	3.0
CR-00-20-NC	NC					
CR-00-20-DC	DC					

Article No	Coating	D1	L1	L2	L3	D2
LN-00-10-CrN	CrN	1.0	3.5	16.0	50.0	6.0
LN-00-10-NC	NC					
LN-00-20-CrN	CrN	2.0	3.5	16.0	50.0	6.0
LN-00-20-NC	NC					
LN-00-20-DC	DC					




Tungsten Carbide Burs for simple Milling Machines


钨钢车针(精密研磨仪用)




Precisely match Various System in conformity with ISO Standards ;
Extra - long working life and high cutting capacity ;
Used for trimming and polishing implant teeth shoulders ;
milling precious or Non-precious metals or Titanium .



标准双纹, Standard Cross Cut



细齿单纹, Plain Cut Fine



粗齿双纹, Coarse Cross Cut



蜡型车针, For Wax mill


DuraCool
Diamond Stones for Zirconia and All-Ceramics

“酷耐”锆牙磨头




MORE COMFORTABLE at the minimal pressure;
HEATLESS no need water cooling ;
MORE DURABLE under our suggested PRM

1. 不产热, 无需水冷; 2. 更耐磨, 实惠;
3. 锋利, 手感好, 更轻松!




DuraCool 
Diamond Polishers for Zirconia and All-ceramics
“酷耐”硅胶金刚石抛光

Especially designed for polishing all types of Ceramics and Zirconia, with unique combination of Pure Diamond particles, Aluminum oxide, Glass etc. Create this spectacular range of abrasive instruments. Using the suggested RPM, guarante to cut and abrade materials with great ease.



Universal diamond polisher for zirconia, lithium disilicate ceramic, porcelain and composite. Suitable for large and inaccessible surface, enable users to polish different areas with just one instrument.

万能金刚石抛光轮,可用于氧化锆、玻璃陶瓷、烤瓷及树脂的快速抛光。适合颌面、缝隙、瓷贴面、树脂贴面等大范围或难以抛光的区域,快速高效。



Silicon Rubber, Rubber Polishing bur, Wheel & Point
胶轮胶粒, 硅胶磨头



C Series, Reduce, Smooth and polish Ceramics, C系列: 用于烤瓷牙抛光
 A Series, Polish Gold, Amalgam, Metal, Alloy etc, A系列: 用于金属, 合金抛光
 R Series, Fast reduce, finish or polish Acrylic. R系列: 用于塑料胶托抛光

Silicon Rubber Wheel & Point Ceramic Polishing
胶轮, 胶粒 (烤瓷牙抛光)



编号 No.	C 104 C	C 107 C	C 101 C	C 104 M	C 107 M	C 101 M	C 104 F	C 107 F	C 101 F
颜色, 粒度 Grit	White, Coarse; 白色. 粗糙			Pink, Medium; 粉色. 中粗			Blue, Fine; 蓝色. 较细		


Silicon Rubber Bur Ceramic Polishing

硅胶磨头 (烤瓷牙抛光)



编号 No.	C203C	C201C	C205C	C207C	C203M	C201M	C205M	C207M	C203F	C201F	C205F	C207F
颜色, 粒度 Grit	White, Coarse; 白色. 粗糙				Pink, Medium; 粉色. 中粗				Blue, Fine; 蓝色. 较细			

Durable Rubber Wheel & Point Metal, Alloy Polishing
胶轮, 胶粒 (金属, 合金抛光)



编号 No.	A102 C	A103C	A104 C	A102 M	A103 M	A104 M	A102 F	A103 F	A104 F
颜色, 粒度 Grit	Black, Coarse; 黑色. 粗糙			Brown, Medium; 棕色. 中粗			Green, Fine; 绿色. 较细		

Durable Rubber Bur Metal, Alloy Polishing

橡胶磨头 (金属, 合金抛光)



编号 No.	A203C	A201C	A206C	A207C	A203M	A201M	A206M	A207M	A203F	A201F	A206F	A207F
颜色, 粒度 Grit	Black, Coarse; 黑色. 粗糙				Brown, Medium; 棕色. 中粗				Green, Fine; 绿色. 较细			

Silicon Rubber Bur Acrylic Polishing Reinforced Shank

橡胶磨头 加硬柄 (塑料基托抛光)



编号 No.	Rh 101 C	Rh 103 C	Rh 106 C	Rh 101 M	Rh 103 M	Rh 106 M	Rh 101 F	Rh 103 F	Rh 106 F
颜色, 粒度 Grit	Green, Coarse; 绿色. 粗糙			Grey, Medium; 灰黑色. 中粗			Yellow, Fine; 黄色. 较细		

Silicon Rubber Bur Acrylic Polishing

橡胶磨头 (基托抛光)



编号 No.	R103C	R106C	R209C	R103M	R106M	R209M	R103F	R106F	R209F	R202F
颜色, 粒度 Grit	Green, Coarse; 绿色. 粗糙			Grey, Medium; 灰黑色. 中粗			Yellow, Fine; 黄色. 较细			

Green Stone Medium
Contour Ceramic
绿车石
调烤瓷

编号 No.	G-01	G-02	G-03	G-05	G-06	G-08	G-09	G-17	G-13	G-16	G-04
	6.0 18.0	13.0 2.0	5.0 13.0	13.0 3.0	7.0 13.0	5.0 12.0	9.0 18.0	6.0 12.0	3.0 7.0	4.0 9.0	13.0 5.5

Pink Stone Medium
for Metal, Alloy
粉车石
调金属,合金

编号 No.	P-01	P-02	P-03	P-05	P-06	P-08	P-09	P-17	P-13	P-16
	6.0 18.0	13.0 2.0	5.0 13.0	13.0 3.0	7.0 13.0	5.0 12.0	9.0 18.0	6.0 12.0	3.0 7.0	4.0 9.0

White Stone for Metal, Alloy Composite,
白车石
调金属,合金
复合材料

编号 No.	W-01	W-02	W-03	W-05	W-06	W-08	W-09	W-17	W-13	W-16
	6.0 18.0	13.0 2.0	5.0 13.0	13.0 3.0	7.0 13.0	5.0 12.0	9.0 18.0	6.0 12.0	3.0 7.0	4.0 9.0

Abrasive Sandpaper with mandrel
砂纸棒

编号,目数 No./Grit	240#	400#	600#	800#	1000#	1500#	2000#	2500#	3000#	5000#
----------------	------	------	------	------	-------	-------	-------	-------	-------	-------

Mandrels Stainless Reinforced
不锈钢夹针

编号 No.	32HS	51HS	52HS	81HS	25HS	53HS	55HS	D20HS	D30HS	52D
Head size	3.0mm	5.0mm	5.0mm	8.0mm	2.5mm	5.0mm	5.0mm	5.0mm	6.0mm	5.0mm

Mandrels Steel
普通夹针

编号 No.	51H	52H	72H	54H
Head size	5.0mm	5.0mm	7.0mm	5.0mm

No. 1708P
Polishing Burs, Mandrels Collection
抛光磨头,夹针样品本

More than 100 specifications,
Silicon, Rubber, Diamond burs
HP Shanks

Polishing brush
Wheel Shape

有柄平刷



碳化硅磨料丝	氧化铝磨料丝	曲钢丝	曲铜丝	黑猪鬃	灰羊毛	灰羊毛	灰羊毛	白羊毛
22mm								
编号 No.	01SiC22	01Al22	01A22	01B22	01D22	01E22	01F22	01G22
材质 Content	Silicon Carbide	Al. Oxide	Crimped steel wire	Crimped brass wire	Chungking Bristle	Grey Pony hair	Grey goat hair	White goat hair



百洁布红,硬	百洁布棕,中粗	百洁布紫,软	百洁布蓝,普通	黄棉布	白棉布	橘红羊皮	鹅黄羊皮双面绒	灰羊毛夹羊皮	
25mm			22mm						
编号 No.	01Mp25	01My25	01Ms25	01M22	01P22	01Q22	01R22	01Fr22	
材质 Content	Red.coarse	Brown. Medium	Purple.fine	Blue	Yellow cloth	White cloth	Chamois leather	Microfiber leather	Goat hair / chamois



羊绒	白棉线	黑猪鬃	灰羊毛	灰羊毛	白羊毛	黑猪鬃	灰羊毛	白羊毛	
22mm		19mm							
编号 No.	01S22	01T22	01D19	01E19	01F19	01G19	02D22	02E22	02G22
材质 Content	Cashmere	Cotton yarn	Chungking Bristle	Grey Pony hair	Grey goat hair	White goat hair	Chungking Bristle	Pony hair	White goat hair

Wheel Shape
without shank
无柄平刷



编号 No.	03A22	03B22	03D22	03E22	03F22	03G22	03Q22	03T22
材质 Content	Crimped steel wire	Brass wire	Chungking bristle	Pony hair	Grey goat hair	White goat hair	Muslin cloth	Cotton yarn

Bowl shape
12mm
杯型刷



Pen shape
5 mm
笔型刷




编号 No.	04A12	04B12	04D12	04F12	05SiC05	05Al05	05A05	05B05	05D05	05E05	05F05
材质 Content	Steel	Brass	Bristle	Goat Hair	SiC	Al. Oxide	Steel	Brass	Bristle	Pony Hair	Goat Hair

Felt Bobs & Burs 毛毡磨头



03N25	03Ns25	3P	2P	1P	1R
25 x 5 mm, 硬质	25 x 5 mm, 软质	15 x 25 mm	20 x 40 mm	22 x 50 mm	22 x 50 mm

FBB-01
Curved Brass Wire






6.0 15.0	9.0 18.0	6.0 15.0	9.0 18.0	22.0 4.0	19.0 4.0	22.0 5.0	22.0 5.0
G6	G9	A6	A9	N22	P19	P22	T22

Polishing brush & Wheel
in Bristle, Cotton Cloth
抛光刷




These brushes are made of bristle hair or White, yellow pure cotton cloth, assembled on polishing latch to finish or polish Plaster Moulds, Acrylic
台式抛光机配用的抛光刷/布轮,用于抛光塑料,树脂,模型等。



Polishing Brushes 抛光刷

Item	2A	2B	2C	6A	6B	6C	J03	J02	J01	P80	03B48
D	79	73	68	64	60	54	70	58	51	80	48
d	48	48	48	35	35	35	35	30	25	40	25
e	6.5	6.5	6.5	6.5	6.5	6.5	8.5/5	8/5	7/5	6.5	6.35
B	16	16	16	16	16	16	24	16.5	15	36	9
Rows	2	2	2	2	2	2	3	2	1	4	1




Converging

Straight

ZD 101

Popular HP diamond burs KIT
常用金刚石磨头套装



MD23, MK47, MJ60, MA23, MC18, MM27, MB21, ML40, MG16, MR23

10 Hot Spec. of HP Diamond, universally grind, remove, polish and finish of Alloy, metal, Jewellery, stones etc. It widely applied to Laboratory, DIY or Hobbies.

用途: 选取10支最常用的金刚石磨头, 形状全, 大小适中, 可以用来处理陶瓷, 金属, 树脂, 石头, 珠宝等 几乎所有材质。

ZD 102

Small & Tiny Cut, HP diamond burs KIT
金刚石磨头, 小号套装



MD10, MD14, MD23, MK16, MA16, MA23, MC25, MM12, MM18, MM27, MB12, MB16, MB21, ML16, ML23, MN18, MG16, MR23, MQ23, MU60.

20 Spec. small or sharp burs. To create tiny hole, slim curves or tiny, deep grind smartly. Used to Alloy, metal, ceramics, Jewellery, stones, hard wood etc. It widely applied to Laboratory, DIY or Hobbies.

用途: 20支金刚石磨头, 用来开小洞, 孔, 细边, 细槽, 或者小的磨削。很轻巧, 灵活。适用于各种金属, 合金, 陶瓷, 珠宝, 宝石, 木头等。

ZD 103

Fast Trimming, HP diamond burs KIT
金刚石磨头套装, 快速磨削



CQ23, CL40, SD35, SA50, SC40, SP70, SG50, SV70, SV60, SF60

20 Spec. HP Diamond for fast trimming, Grinding of Plastic, Jewellery, stones, Hard wood, etc. Save time when removing and work efficiently

用途: 10支常用的金刚石磨头, 头部相对比较大, 可以快速磨削塑料, 木头, 陶瓷和石头等等, 适用于齿科, 珠宝加工和DIY爱好者。

ZD 104

Hot Fully-Sintered diamond KITS
全砂烧结磨头套装, 常用



C 050S, F050S, L05S, B050C, A037, G050, C037F, F050F, M050EF, G050EF

10 spec. Hot Fully-sintered Diamond HP burs, universal grind, remove, polish and finish of Alloy, metal, Jewellery, stones etc. It widely applied to Lab, DIY or Hobbies.

用途: 10个常用全砂烧结磨头, 5种不同粗细。广泛适用于研磨, 抛光金属, 合金, 珠宝, 宝石。经久耐用, 方便实惠。

ZD 105

Fast Trimming Fully-Sintered diamond KIT
全砂烧结磨头套装, 大磨削




K080S, K100S, J100S, C070S, M070S, B070S, L070S, G060S, F070S, R070S

If workpieces are extremely hard and waste too much time when grinding, this fast trimming kit may help to work efficiently. Applied to Alloy, metal, Jewellery, stones, hard wood, etc

用途: 10个大头的全砂烧结磨头, 用于粗磨, 去块, 效率高, 经济实惠。

ZT 101

Popular HP Tungsten carbide burs KIT
常用低速钨钢针




HP 1/4, HP 1/2, HP 1, HP 2, HP 3, HP 4, HP 5, HP 8, HP 34, HP 36, HP 38, HP 41, HP 556, HP 557, HP 560, HP 563, HP 699, HP 701, HP 702, HP 704

20 Spec. HP Tungsten Carbide burs are used to create tiny hole, slim curves or tiny, deep grind very smartly. It is applied to Alloy, metal, ceramics, Jewellery, stones, hard wood etc.

用途: 20个型号钨钢车针, 用于雕牙花, 钻小孔, 洞, 开细槽, 细线等, 金属, 合金, 宝石, 坚木都可以用。

ZT 102

Hot Tungsten carbide cutters KIT
钨钢磨头套装, 常用



D050SP, L040SC, L060SC, F060SC, B023SC, L023SC, L040FC, C023FC, L040EFC, F060CC

10 spec. Tungsten carbide Cutter, universal grind, remove, polish and finish of Acrylic, Resin, stones, soft metal or wood.

用途: 10个常用钨钢磨头, 形状全, 大小适中。主要用来研磨, 抛光塑料, 树脂, 软金属, 石头等, 方便实惠。

ZF 101

Acrylic Polishing Kit
树脂基托抛光套装



F060SC, L040SC, U12 x 2, U13 x 2, U16 x 2

Simple and straight way to adjust Acrylic without the problem associated with pumice.

用途: 2个型号的钨钢磨头快速磨削, 三个不同粒度的硅胶磨头依次进行抛光。

ZA 101

Non-precise metal, Alloy Adjust and polish
非贵金属, 铬钴调整与抛光



W-01, W-09, W-08, W-16, W-03, W-06
A201C, A203C, A206C, A201F, A 203F

To adjust and polish Cobalt-chromium, nickel chromium Alloy. White stones for adjusting or grinding, Rubber burs for polishing.

用途: 试用于铬钴、镍铬合金等非贵金属的研磨和抛光 6支白车石用来研磨和调整, 5支硅胶磨头用来抛光。

ZP 100

Sandpaper Polishing burs Kit
砂纸棒套装




240#, 400#, 600#, 800#, 1000#, 1500#, 2000#, 2500#, 3500#, 5000#

10 different Grit Sandpaper burs are used for polishing Alloy, metal, Jewellery stones, hard woods etc., Working step by step, mirror surface in coming soon, water cooling or dry are acceptable.

用途: 10个不同粒度的砂纸棒, 从240目到5000目不等, 请根据需要, 选取合适的粒度, 耐心细致地使用。

Z 2001

Popular Felt Polishing bur Kit
毛毡抛光磨头套装



12 pieces of High quality Felt polishing burs are widely used to mould, Alloy, metal, Jewellery, stones, Hard woods etc.

12 pieces of High quality Felt polishing burs are widely used to mould, Alloy, metal, Jewellery, stones, Hard woods etc.

用途: 12支常用的毛毡磨头, 广泛应用于各种模具, 金属, 陶瓷, 珠宝, 宝石和木头等。

ZS 100

Hot HP Mounted stones mixed
车石套装, 常用。




G-01, G-02, G-03, G-06, G-08, G-17, P-01, P-02, P-03, P-06, P-08, P17
W-01, W-02, W-03, W-06, W-08, W-17, C-06, C-09 x 5 pcs of each

100 Pcs, Hot Mounted stones in Green, Pink, White, applied to Alloy, metal, ceramics, Jewellery, stones, hard wood etc.

用途: 100支常用的车石, 广泛用于磨削金属, 合金, 陶瓷, 宝石, 坚木等, 经济实惠。


隐形义齿胶 FlexSoft Partial Dentures



The alternative to conventional Partial Dentures ;
Flexible . Soft. Easy-to-finish. Easy-to-polish
弹性好，易成型，易抛亮，美观大方。



编号 No.	规格 Spec.
D2	Crystal ; 18 gs
C1	Ligth Pink ; 9 gs
C2	Ligth Pink ; 18 gs
B1	Standard Pink ; 9 gs
B2	Standard Pink ; 18 gs
A1	Deep Pink ; 9 gs
A2	Deep Pink ; 18 gs
For bulk pak of 1 kilo.	
C3	Ligth Pink ; 1000 gs
B3	Standard Pink ; 1000 gs
A3	Deep Pink ; 1000 gs

 Crystal (D)
(透明)

 Light pink (C)
(浅粉红)

 Standard pink (B)
(粉红)

 Deep pink (A)
(深粉红)



Instructions for Use

Suitable for: Partial and complete dentures, occlusal splints and nightguards.

Not Suitable for: Small unilateral restorations.

Warnings: Use eye protection and heat-resistant gloves when injecting.

Use proper ventilation, vacuum system and masks when grinding prosthodontic resins.

Precautions: Cartridge must be injected within 1 hour of opening the foil pouch.

I. PREPARING THE MASTER MODEL FOR DUPLICATION

All undercuts that will affect the path of insertion should be blocked-out with wax prior to duplication. Do not relieve necklines of teeth or areas desirable for retention. For free-end saddles, it is not necessary to remove the distal undercut of the last natural tooth, unless the area is too deep. For free-end saddle cases where anterior teeth are also being replaced, undercuts in the anterior should also be blocked out.

II. SET-UP

Duplicate blocked-out master model. Articulate the duplicate model and the master model (you will need the master model articulated at finishing steps to check occlusion). Set up the case on the duplicate model using only wax (no base plate).

Allow 1 mm of thickness under the teeth to flow to the labial side and into the retention diatorics in the teeth. Use a small round carbide bur to create mechanical retention in the teeth, from mesial to distal and connected to the center diatoric hole.

III. WAXING

Important: Wax the case with anatomical features to reduce grinding during the finishing stage.

Recommended thickness of wax-up:

1. Palatal thickness 1.2 - 1.5 mm
2. Buccal & Labial Flange 1.2 - 2.0 mm
3. Clasps 1.0 - 1.5 mm
4. Lingual Flange 1.2 - 2.0 mm
5. Lingual Major Connector (Mandibular) 1.7 - 2.0 mm

IV. INVESTING AND SPRUING

Invest the case in the center of the flask, aiming the posterior portion toward the sprue opening. After the stone (or plaster) has set, use two sprues from the flask opening to the most posterior portion of the wax on both sides (see diagram). Sprues should be approximately 7 mm. width. Do not create sharp angles with the sprues, as the heated resin must flow in a straight or slightly curved path. Apply a thin layer of petroleum jelly on the entire stone surface. Then invest top half of flask.

V. REMOVING WAX

After the stone has set, place flask in boiling water for 8 minutes. Open the flask and wash out the wax, making sure all wax residues are eliminated. Allow the flask to cool. Apply thin layers of separator on the entire surface of the stone and allow it to dry completely. Do not apply separator on teeth. When dry, close flask.

VI. INJECTING

Wear safety goggles and heat-protective gloves for this procedure.

Spray inside of the cartridge sleeve with mold release spray and place it into furnace. Turn on the furnace and allow it to reach approximately 550°F (290°C), preheating cartridge sleeve for approximately 15 minutes.

Remove cartridge from foil-lined bag (do not remove cartridge from bag until this step). Spray the cartridge and bronze disc with mold release spray. Insert the cartridge into the cartridge sleeve with the cap end first. Insert the bronze disc into the cartridge sleeve against the end of the cartridge. Use the thick bronze disc for medium tubes and the thin bronze disc for the large tubes. Heat cartridge for 12 minutes.

Position the flask with the sprue end up, directly under the shaft of the press. After 12 minutes, remove the cartridge sleeve from the furnace. Keep the cartridge sleeve in a horizontal position until resting on top of the flask, and then turn the cartridge sleeve into a vertical position. Immediately use the press to inject the heated cartridge into the flask using a steady fast motion. Continue applying pressure until tube bursts and continue with maximum downward pressure for 3 minutes.

Remove the cartridge sleeve from the flask with a slight twisting motion. After removing the cartridge sleeve, let the flask cool for at least 30 minutes. Do not put flask into a cold water bath. Use a plaster nipper to de-flask. Do not use a hammer. Remove the compressed cartridge from the cartridge sleeve with a knock-out punch. Important: Cartridge sleeve must be pre-heated for 10 minutes BEFORE the next injection.

VII. FINISHING AND POLISHING

Remove appliance from the model. Cut off the sprues with a cut-off wheel. Trim the edges of appliance with a coarse grinding wheel. Avoid overheating. Use rubber wheels to smooth scratches. Place the appliance on the articulated master model to check the fit, retention, thickness and occlusion. Clean the tissue surface with a small bristle brush.

To polish, use a coarse dental abrasive, such as pumice, with a wet rag wheel. Apply final shine with a hi-shine compound on low speed with a final buff wheel. Be careful not to overheat appliance.

铝桶

Empty Aluminum Cartridge

编号 No.	规格 D x H, mm
25085	25 x 85 mm
25100	25 x 100 mm
25CAP	24.2 x 16 mm, Cap
22065	22 x 65 mm
22085	22 x 85 mm
22100	22 x 100 mm
22120	22 x 120 mm
22CAP	21.2 x 16 mm, Cap



Dianfong ABRASIVES



ISO 13485 / CE 0197



阿里官网
Website on Alibaba



微信公众
WeChat Public

深圳点锋研磨材料有限公司
Shenzhen Dian Fong Abrasives Co.,Ltd

地址: 广东省深圳市龙岗区平湖天华路5号2栋4楼
Address: 4th Floor, Blk 2, No. 5, TianHua Rd, PingHu, LongGang Distr.,
Shenzhen, Guangdong Province, China 邮编: Post code: 518 111
网站(Web): www.dianfong.com; 邮箱(Email): dianfong@foxmail.com
电话(Tel.): 0755-2885 6248; 2834 4545; 传真(Fax): 0755-2885 6246
QQ: 27212989; Skype:welson881; WhatsApp:0086-13528829441

Dianfang ABRASIVES



ISO 13485
CE 0197

点锋研磨

第6版/Vol.6

Rotary things



dentistry

About Us , 关于我们

Dian Fong specialize in Rotary things of dental and Jewelry, who provide innovative and professional products and service.

Our products include DIAMOND bur and disc, TUNGSTEN CARBIDE bur or cutter, FLEXIBLE RUBBER bur and wheel, Polishing Brush, Mounted Stone, ZIRCONIA grinder and polish etc.

We make sure of the best quality and have certificates of ISO 13485 and CE0197. Whoever you are dentist or technician, we always make you relaxing, safe and efficient. We are your best friend in dentistry!

Content 目录6.1

Dentistry , Dentist's Use . 临床类	
Diamond bur . 金刚石车针	3
Zekrya Bur , Endo Z , Scaling Burs , Oral Surgical Burs 拔牙车针, 手术车针	7
Tungsten Carbide Bur . 钨钢车针	8
Diamond Polisher , Arkansas Stone , Green Stone , Diamond Stones 橡胶金刚石磨头 , 车石 , 陶瓷金刚石磨头	10
Silicon Polisher for Ceramic , Amalgam 硅胶磨头(用于烤瓷, 银汞抛光)	11
Polishing brush and Cup 抛光刷, 抛光杯	11
Mandrels, Blocks for Sterilizing 夹针, 消毒盒	12
Samples Collection 样品册	6, 12
Clinic Kit 套装	13, 14

一首小诗

我磨一根针，
足以穿风尘。
海内结知己，
笑淡了风云。

A Little Poem

A small bur which I made,
is sharp to teeth-grind ;
Anywhere, we make friends,
and smile as sky fade .



ISO 13485



CE 0197



Established in 2007


Diamond Bur

金刚砂车针



Advantages

1. Strictly comply with ISO Standard ;
2. Sharper and More Efficient ;
3. Durable and Valuable ;
4. Full Range of shapes and grits.
5. Fastly develop new items.
6. Individual Design and OEM Service































































Grit

- Black / Super Coarse
- Green / Coarse
- Blue / Medium
- Red / Fine
- Yellow / Extra Fine































































































001-008M Item No.
























































































Diamond @ Round - Preparation & Polish




































Most Frequently used for crown preparation , to cut depth guides, Occlusal rest preparation, endodontic access and cavity preparation.

 <p>Round 圆球</p>	 001/008 0.8/19.0 BR-49	 001/010 1.0/19.0 BR-45	 001/012 1.2/19.0 BR-46	 001/014 1.4/19.0 BR-41	 001/016 1.6/19.0 BR-40	 001/018 1.8/19.0 BR-31	 001/021 2.0/19.0 BR-28	 001/023 2.3/19.0 BR-26	 001/025 2.5/19.0 BR-25	 001/029 2.9/19.0 BR-23	 002/012 2.3/19.0 BC-43	 002/014 2.3/19.0 BC-42
	001-008M	001-010M	001-012M	001-014M	001-016M	001-018M	001-021M	001-023M	001-025M	001-029M	002-012M	002-014M
Medium grit										Round with Collar /BC		
 002/016 3.5/19.0 BC-32	 002/018 2.3/19.0 BC-31	 005/018 7.5/21.0 EX-28	 001/014 1.4/19.0 BR-41S	 001/016 1.6/19.0 BR-40S	 001/018 1.8/19.0 BR-31S	 001/021 2.0/19.0 BR-28S	 001/023 2.3/19.0 BR-26S	 001/025 2.5/19.0 BR-25S	 001/029 2.9/19.0 BR-23S	 001/014 1.4/19.0 BR-41C	 001/016 1.6/19.0 BR-40C	 001/018 1.8/19.0 BR-31C
002-016M	002-018M	005-018M	001-014S	001-016S	001-018S	001-021S	001-023S	001-025S	001-029S	001-014C	001-016C	001-018C
Super Coarse										Coarse		
 001/021 2.0/19.0 BR-28C	 001/023 2.3/19.0 BR-26C	 001/025 2.5/19.0 BR-25C	 001/029 2.9/19.0 BR-23C	 001/012 1.2/19.0 BR-46F	 001/014 1.4/19.0 BR-41F	 001/016 1.6/19.0 BR-40F	 001/018 1.8/19.0 BR-31F	 001/021 2.0/19.0 BR-28F	 001/023 2.3/19.0 BR-26F			
001-021C	001-023C	001-025C	001-029C	001-012F	001-014F	001-016F	001-018F	001-021F	001-023F			
Fine grit												
 001/012 1.2/19.0 BR-46EF	 001/014 1.4/19.0 BR-41EF	 001/016 1.6/19.0 BR-40EF	 001/018 1.8/19.0 BR-31EF	 001/021 2.0/19.0 BR-28EF	 001/023 2.3/19.0 BR-26EF							
001-012EF	001-014EF	001-016EF	001-018EF	001-021EF	001-023EF							
Extra Fine												
 25mm / 1 inch FG XL Surgical 001/010 1.0/25.0 BR-X45	 001/012 1.2/25.0 BR-X46	 001/014 1.4/25.0 BR-X41	 001/016 1.6/25.0 BR-X40	 001/018 1.8/25.0 BR-X31	 001/016 1.6/25.0 BR-X40S	 001/018 1.8/25.0 BR-X31S	 001/014 1.4/25.0 BR-X41C	 001/016 1.6/25.0 BR-X40C	 001/018 1.8/25.0 BR-X31C			
加长柄	001-X010M	001-X012M	001-X014M	001-X016M	001-X018M	001-X016S	001-X018S	001-X014C	001-X016C	001-X018C		
	Medium					Super Coarse			Coarse			
 22mm / 0.83 inch RA/CA Shank 001/012 1.2/22.0 BR-R46	 001/014 1.4/22.0 BR-R41	 001/016 1.6/22.0 BR-R40	 001/018 1.8/22.0 BR-R31	 001/023 2.3/22.0 BR-R26	 001/016 1.6/22.0 BR-R40F	 001/018 1.8/22.0 BR-R31F	 001/023 2.3/22.0 BR-R26F	 001/016 1.6/22.0 BR-R40EF	 001/018 1.8/22.0 BR-R31EF	 001/023 2.3/22.0 BR-R26EF		
弯机柄	001-R012M	001-R014M	001-R016M	001-R018M	001-R023M	001-R016F	001-R018F	001-R023F	001-R016EF	001-R018EF	001-R023EF	
	Medium					Fine			Extra Fine			



 010/010 1.0/19.0 SI-45	 010/012 1.2/19.0 SI-46	 010/014 1.6/19.0 SI-47	 010/016 1.8/19.0 SI-48	 019/010 2.3/19.0 DI-41	 019/014 3.2/19.0 DI-42	 038/032 5.0/19.0 EX-11	 039/035 7.0/20.0 EX-12	 041/025 0.5/19.0 WF-30	 055/045 2.0/19.0 WR-44	 068/042 1.8/19.0 WR-13	 109/010 4.0/19.0 SF-41	 109/012 4.0/19.0 SF-31
010-010M	010-012M	010-014M	010-016M	019-010M	019-014M	038-032M	039-035M	041-025M	055-045M	068-042M	109-010M	109-012M
Inverted Cone /SI		Inverted Cone with Collar /DI				Barrel / EX		Wheel / WF		Flat End Cylinder /SF		
 109/014 4.0/19.0 SF-42	 110/014 6.0/21.0 SF-21	 110/016 6.0/21.0 SF-15	 111/012 8.0/22.0 SF-11	 111/014 8.0/22.0 SF-12	 111/016 8.0/22.0 SF-13	 138/010 4.0/19.0 SR-41	 138/012 4.0/19.0 SR-31	 138/014 4.0/19.0 SR-42	 140/014 6.0/21.0 SR-21	 141/012 8.0/22.0 SR-11	 141/014 8.0/22.0 SR-12	 141/016 8.0/22.0 SR-13
109-014M	110-014M	110-016M	111-012M	111-014M	111-016M	138-010M	138-012M	138-014M	140-014M	141-012M	141-014M	141-016M
Round End Cylinder /SR												
 150/012 10.0/21.0 EX-18	 164/012 6.0/21.0 TC-26	 165/014 8.0/22.0 TC-21	 166/010 10.0/23.0 TC-27	 166/012 10.0/23.0 TC-28	 166/014 10.0/22.0 TC-10	 166/016 10.0/22.0 TC-11	 166/018 10.0/23.0 TC-16	 169/010 3.5/19.0 TF-41	 170/012 4.0/19.0 TF-42	 170/014 4.0/19.0 TF-43	 170/016 4.0/19.0 TF-31	 170/018 5.0/20.0 TF-23
150-012M	164-012M	165-014M	166-010M	166-012M	166-014M	166-016M	166-018M	169-010M	170-012M	170-014M	170-016M	170-018M
End Cut /EX		Needle /TC				Flat Cone /TF						
 170/021 5.0/20.0 TF-22	 171/014 6.5/21.0 TF-20	 171/016 6.5/21.0 TF-21	 172/010 8.0/22.0 TF-07	 172/012 8.0/22.0 TF-08	 172/014 8.0/22.0 TF-09	 172/016 8.0/22.0 TF-10	 172/018 8.0/22.0 TF-56	 172/023 9.0/22.0 TF-14	 173/014 10.0/24.0 TF-11	 173/016 10.0/22.5 TF-12	 173/018 10.0/22.5 TF-13	 196/010 4.0/19.0 CR-41
170-021M	171-014M	171-016M	172-010M	172-012M	172-014M	172-016M	172-018M	172-023M	173-014M	173-016M	173-018M	196-010M
 196/012 4.0/19.0 CR-31	 196/014 4.0/19.0 CR-42	 196/016 4.0/19.0 CR-12	 196/020 5.0/19.0 CR-11	 197/014 6.5/21.0 TR-20	 197/016 6.5/21.0 TR-21	 197/018 6.5/21.0 TR-62	 197/020 6.5/21.0 TR-63	 198/014 8.0/22.0 TR-25	 198/016 8.0/22.0 TR-24	 198/018 8.0/22.0 TR-23	 199/014 10.0/22.0 TR-11	 199/016 10.0/22.0 TR-12
196-012M	196-014M	196-016M	196-020M	197-014M	197-016M	197-018M	197-020M	198-014M	198-016M	198-018M	199-014M	199-016M
Round End Cone /CR, TR												
 199/018 9.0/22.0 TR-13	 199/022 9.0/22.0 TR-14	 200/022 11.0/24.0 TR-15	 200/024 11.0/24.0 TR-19	 233/010 2.8/19.0 EX-41	 234/014 4.0/19.0 EX-31	 237/018 5.0/19.0 EX-20	 237/021 5.0/19.0 EX-21	 237/032 5.0/19.0 EX-26	 248/012 6.0/21.0 FO-41	 248/014 6.0/21.0 FO-20	 249/014 8.0/22.0 FO-21	 249/016 8.0/22.0 FO-22
199-018M	199-022M	200-022M	200-024M	233-010M	234-014M	237-018M	237-021M	237-032M	248-012M	248-014M	249-014M	249-016M
				Pearl /EX				Flame /FO				
 250/012 10.0/24.0 FO-11	 250/014 10.0/24.0 FO-12	 250/016 10.0/24.0 FO-13	 257/018 3.6/19.0 FO-32	 257/023 5.0/19.0 FO-34	 257/028 5.0/19.0 FO-25	 257/020 5.0/19.0 FO-22	 277/030 5.0/19.0 TR-28	 288/012 6.0/21.0 SO-20	 289/014 8.0/22.0 SO-21	 298/016 8.0/22.0 878K16	 298/018 8.0/22.0 878K18	 299/016 10.0/24.0 879K16
250-012M	250-014M	250-016M	257-018M	257-023M	257-028M	257-020M	277-030M	288-012M	289-014M	298-016M	298-018M	299-016M
				Football /TR				Torpedo /SO				
 515/014 8.0/22.0 BC-11	 519/016 8.0/22.0 BC-12	 552/016 6.0/21.0 VP-20	 552/018 6.0/21.0 VP-21	 552/021 6.0/21.0 VP-22	 546/016 8.0/22.0 847KR16	 553/018 10.0/24.0 848KR18						
515-014M	519-016M	552-016M	552-018M	552-021M	546-016M	553-018M						
Gross Reduction /BC		Depth Marking		Taper Round edge								

Super Coarse 特粗砂	 010/016 1.8/19.0 SI-48S	 038/032 5.0/19.0 EX-11S	 039/035 7.0/20.0 EX-12S	 109/012 4.0/19.0 SF-31S	 109/014 4.0/19.0 SF-42S	 111/012 8.0/22.0 SF-11S	 111/014 8.0/22.0 SF-12S	 111/016 8.0/22.0 SF-13S	 166/016 10.0/22.0 TC-11S	 166/018 10.0/23.0 TC-16S	 171/016 6.5/21.0 TF-21S	 173/014 10.0/24.0 TF-11S
	Spec. FG-	010-016S	038-032S	039-035S	109-012S	109-014S	111-012S	111-014S	111-016S	166-016S	166-018S	171-016S
 173/016 10.0/22.0 TF-12S	 173/018 10.0/22.0 TF-13S	 197/016 6.5/21.0 TR-21S	 199/014 10.0/22.0 TR-11S	 199/016 10.0/22.0 TR-12S	 199/018 9.0/22.0 TR-13S	 237/018 5.0/19.0 EX-20S	 277/020 3.9/19.0 TR-22S					
173-016S	173-018S	197-016S	199-014S	199-016S	199-018S	237-018S	277-020S					
Coarse 粗砂	 010/014 1.6/19.0 SI-47C	 068/042 1.8/19.0 WR-13C	 111/014 8.0/22.0 SF-12C	 165/014 8.0/22.0 TC-21C	 173/018 10.0/22.0 TF-13C	 190/018 9.5/21.0 FO-54C	 197/018 6.5/21.0 TR-62C	 199/018 9.0/22.0 TR-13C	 257/018 3.6/19.0 FO-32C			
	Spec. FG-	010-014C	068-042C	111-014C	165-014C	173-018C	190-018C	197-018C	199-018C	257-018C		
Fine 细砂	 110/014 6.0/21.0 SF-21F	 150/012 10.0/21.0 EX-18F	 165/014 8.0/22.0 TC-21F	 166/016 10.0/22.0 TC-11F	 171/016 6.5/21.0 TF-21F	 173/014 10.0/24.0 TF-11F	 173/016 10.0/22.0 TF-12F	 173/018 10.0/22.0 TF-13F	 196/016 4.0/19.0 CR-12F	 196/020 5.0/19.0 CR-11F	 197/014 6.5/21.0 TR-20F	
	Spec. FG-	110-014F	150-012F	165-014F	166-016F	171-016F	173-014F	173-016F	173-018F	196-016F	196-020F	197-014F
 197/016 6.5/21.0 TR-21F	 197/018 6.5/21.0 TR-62F	 199/014 10.0/22.0 TR-11F	 199/016 10.0/22.0 TR-12F	 199/018 9.0/22.0 TR-13F	 234/014 4.0/19.0 EX-31F	 237/021 5.0/19.0 EX-21F	 249/014 8.0/22.0 FO-21F	 249/016 8.0/22.0 FO-22F	 254/016 3.2/19.0 FO-30F	 289/014 8.0/22.0 SO-21F		
197-016F	197-018F	199-014F	199-016F	199-018F	234-014F	237-021F	249-014F	249-016F	254-016F	289-014F		
Extra Fine 特细砂	 165/014 8.0/22.0 TC-21EF	 166/010 10.0/23.0 TC-27EF	 166/012 10.0/23.0 TC-28EF	 166/016 10.0/22.0 TC-11EF	 171/016 6.5/21.0 TF-21EF	 173/014 10.0/24.0 TF-11EF	 173/016 10.0/22.0 TF-12EF	 173/018 10.0/22.0 TF-13EF	 196/016 4.0/19.0 CR-12EF	 196/020 5.0/19.0 CR-11EF	 197/016 6.5/21.0 TR-21EF	 197/018 6.5/21.0 TR-62EF
	Spec. FG-	165-014EF	166-010EF	166-012EF	166-016EF	171-016EF	173-014EF	173-016EF	173-018EF	196-016EF	196-020EF	197-016EF
 198/016 8.0/22.0 TR-24EF	 199/014 10.0/22.0 TR-11EF	 199/016 10.0/22.0 TR-12EF	 199/018 9.0/22.0 TR-13EF	 237/021 5.0/19.0 EX-21EF	 243/010 4.3/23.5 FO-40EF	 248/012 6.0/21.0 FO-41EF	 248/014 6.0/21.0 FO-20EF	 249/014 8.0/22.0 FO-21EF	 249/016 8.0/22.0 FO-22EF	 254/016 3.2/19.0 FO-30EF	 257/023 5.0/19.0 FO-34EF	
198-016EF	199-014EF	199-016EF	199-018EF	237-021EF	243-010EF	248-012EF	248-014EF	249-014EF	249-016EF	254-016EF	257-023EF	
FG Short 短柄	 001/010 1.0/16.0 BR-S45	 001/014 1.4/16.0 BR-S41	 010/012 1.2/16.0 SI-S46	 010/016 1.8/16.0 SI-S48	 165/014 8.0/19.0 TC-S21	 169/010 3.5/16.0 TF-S41	 170/016 4.0/16.0 TF-S31	 170/018 5.0/18.0 TF-S23	 171/014 6.5/18.0 TF-S20	 171/016 6.5/18.0 TF-S21	 197/016 6.5/17.0 TR-S21	 249/014 8.0/19.0 FO-S21
	Spec. FG-	001-S010M	001-S014M	010-S012M	010-S016M	165-S014M	169-S010M	170-S016M	170-S018M	171-S014M	171-S016M	197-S016M
Medium Grit / Blue												

 FG XL Surgical	 010/012 1.2/25.0 SI-X46	 111/012 8.0/28.0 SF-X11	 111/014 8.0/28.0 SF-X12	 111/016 8.0/28.0 SF-X13	 166/016 10.0/29.0 TC-X11	 173/014 10.0/29.0 TF-X11	 173/016 10.0/29.0 TF-X12	 173/018 10.0/29.0 TF-X13	 199/016 10.0/29.0 TR-X12	 200/022 11.5/29.0 TR-X15	 249/014 8.0/28.0 FO-X21	 249/016 8.0/28.0 FO-X22
	加长柄	010-X012M	111-X012M	111-X014M	111-X016M	166-X016M	173-X014M	173-X016M	173-X018M	199-X016M	200-X022M	249-X014M
 RA/CA Shank	 010/010 1.0/22.0 SI-R45	 010/012 1.2/22.0 SI-R46	 010/014 1.6/22.0 SI-R47	 010/016 1.8/22.0 SI-R48	 109/010 4.0/22.0 SF-R41	 109/012 4.0/22.0 SF-R31	 109/014 4.0/22.0 SF-R42	 111/010 8.0/24.0 SF-R10	 111/012 8.0/24.0 SF-R11	 111/014 8.0/24.0 SF-R12	 141/012 8.0/24.0 SR-R11	 141/014 8.0/24.0 SR-R12
	弯机柄	010-R010M	010-R012M	010-R014M	010-R016M	109-R010M	109-R012M	109-R014M	111-R010M	111-R012M	111-R014M	141-R012M
 165/014 8.0/24.0 TC-R21	 166/016 10.0/26.0 TC-R11	 169/010 3.0/22.0 TF-R41	 170/012 4.0/22.0 TF-R42	 198/018 8.0/24.0 TR-R23	 199/016 10.0/27.0 TR-R12	 257/018 3.6/22.0 FO-R32	 165/014 8.0/24.0 TC-R21F	 257/018 3.6/24.0 FO-R32F				
165-R014M	166-R016M	169-R010M	170-R012M	198-R018M	199-R016M	257-R018M	165-R014F	257-R018F				

Samples Collection For Clinic Items

样品册



- Flip-Cart Plastic Pages ; Smart, Easily fill all burs ; Economic and Recyclable ; Nice-looking for showing to Customers ; Convenient to take to anywhere .



No. 170B5

Empty Collection Wallet
样品本(不带车针)



4 Colors Available ,
100 holes per page ,
5 pages per set,
Smart to FG or RA Burs



No. 1701D

Diamond Burs Collection Wallet
金刚石车针样品本



More than 300 specification ,
Diamond burs ,
FG, FG XL, FG 2XL,
FG Short, RA, RA Long





Four Colors Available

Brown

Black

Blue








Yellow



Sharp Burs for Cut teeth (Extraction) 拔牙车针

Applications: Sectioning an impacted wisdom tooth prior to its extraction / Separating roots
 Removing a broken root stump Cutting a tooth horizontally at neck level /
 Reworking temporary teeth of acrylic material / Preparing a shoulder on the vestibular side of the tooth










 <p>Zekrya</p>				 <p>EC</p>		
	11.0 23.0	11.0 25.0	11.0 28.0		7.0 25.0	7.0 28.0
	151	X151	XL151		EC25	EC28

Endo Z 安全端车针

Fast enter root Canal. the end is round and No-cut so it is safe and easily confirm the position.
 顶端圆洞无齿纹, 可以快速建立根管口通道, 避免髓壁穿孔, 安全可靠。














Scaling Bur 去结石车针











Remove Tartar, big scales when cleaning teeth
 洗牙的时候用来除去大块的牙结石, 牙垢

 <p>Endo Z</p>						
	9.0 21.0	9.0 23.0	9.0 25.0		6 sided, 6磨 3.0 21.0	6 sided, 6磨 5.5 21.0
	FG 152	FG X152	FG XL152		ET3	ET5

Oral Surgical Burs of Tungsten Carbide 手术用车针

Please select the following length of FG X Long shank of High speed or RA Long shank of Low speed.

 <p>25mm / 1 inch FG X Long</p>												
	0.8 25.0	1.0 25.0	1.2 25.0	1.4 25.0	1.6 25.0	1.8 25.0	2.3 25.0	4.0 25.0	4.0 25.0	4.0 25.0	4.0 25.0	4.0 25.0
	FG X1	FG X2	FG X3	FG X4	FG X5	FG X6	FG X8	FG X57	FG X58	FG X170	FG X171	FG X557


 <p>26mm / 1.02inch RA Long</p>									
	1.0 26.0	1.2 26.0	1.4 26.0	1.6 26.0	1.8 26.0	2.1 26.0	2.3 26.0	4.0 26.0	4.0 26.0
	RA X2	RA X3	RA X4	RA X5	RA X6	RA X7	RA X8	RA X557	RA X701

Tungsten Carbide Bur












钨钢车针

Advantages












1. Strictly comply with ISO Standard
2. Sharper and More Efficient
3. Durable and Valuable
4. FG, FG Short, FG Long, RA, RA Long etc




























Prepare Cavities 备洞												
		Ball	FG 1/4	FG 1/2	FG 1	FG 2	FG 3	FG 4	FG 5	FG 6		
				Plain Cut	001/005	001/006	001/008	001/010	001/012	001/014	001/016	001/018
FG 7	FG 8	FG 9			Invert Cone	FG 34	FG 35	FG 36	FG 37	FG 38	FG 39	FG 41
001/021	001/023	001/025			Plain Cut	010/008	010/010	010/012	010/014	010/016	010/018	010/023
		Cylinder	FG 56	FG 57	FG 58	FG 59	FG 60	Cylinder	FG 556	FG 557	FG 558	FG 559
		Plain Cut	107/009	107/010	107/012	107/014	107/016	Cross Cut	107/009	107/010	107/012	107/014
FG 560	FG 556L	FG 557L	FG 558L	FG 559L			Taper	FG 168	FG 169	FG 170	FG 171	FG 172
107/016	110/009	110/010	110/012	110/014			Plain Cut	168/008	168/009	168/010	168/012	168/016
FG 169L	FG 170L	FG 171L	Taper	FG 698	FG 699	FG 700	FG 701	FG 702	FG 699L	FG 700L	FG 701L	FG 702L
171/009	171/010	171/012	Cross Cut	168/008	168/009	168/010	168/012	168/016	171/009	171/010	171/012	171/016
Cutting of Crown 破冠												
		Round Cyl.	FG 1556	FG 1557	FG 1558	FG 1559	Crown Prep.	FG 1931	FG 1932			
				Cross Cut	137/009	137/010	137/012	137/014	Fine, Cross	237/010	237/012	
Crown Prep.	FG 1957	FG 1958	Crown Prep.	FG 2057	FG 2058							
Fine, Cross	137/010	137/012	Fine, Cross	107/010	107/012							
Preparation or Removal of Amalgam, metal 去金属												
		Pear	FG 245	FG 246	Egg	FG 329	FG 330	FG 331	FG 332	FG 331 L		
				Plain Cut	233/009	233/012	Plain Cut	232/006	232/008	232/010	232/012	232/010








Removal of Composite, Orthodontics 去树脂 / 正畸											
	Round Cyl. Plain Cut	FG 1155	FG 1156	FG 1157	FG 1158	FG 1159	FG 1160	Round Taper Plain Cut	FG 1169	FG 1170	FG 1171
		4.0 19.0	4.0 19.0	4.0 19.0	4.0 19.0	4.0 19.0	4.0 19.0		4.0 19.0	4.0 19.0	4.0 19.0
		137/008	137/009	137/010	137/012	137/014	137/016		194/009	194/010	194/012











				Shoulder Preparation 备肩台				Carving of tooth pattern (Morphology) 雕牙花		
FG 1172	Round Taper Cross Cut	FG 1701	FG 1702		End Cut	FG 957	FG 958		X-Mas Tree 3 Sided	FG 468
4.0 19.0		4.0 19.0	4.0 19.0			4.0 19.0	4.0 19.0		4.0 19.0	
194/016		194/012	194/016			150/010	150/012		3 Sided	468/014













16.5mm / 0.63 inch FG Short Shank 短柄											
	FG S1	FG S2	FG S3	FG S4	FG S5	FG S245	FG S330	FG S331	FG S556	FG S557	FG S1157
	0.8 16.5	1.0 16.5	1.2 16.5	1.4 16.5	1.6 16.5	2.5 16.5	1.8 16.5	2.0 16.5	4.0 16.5	4.0 16.5	4.0 16.5
	001/008	001/010	001/012	001/014	001/016	233/009	232/008	232/010	107/009	107/010	137/010

22mm / 0.83 inch RA / CA Shank 弯机柄												
	RA 1/2	RA 1	RA 2	RA 3	RA 4	RA 5	RA 6	RA 7	RA 8	RA 34	RA 35	RA 36
	0.6 22.0	0.8 22.0	1.0 22.0	1.2 22.0	1.4 22.0	1.6 22.0	1.8 22.0	2.1 22.0	2.3 22.0	0.8 22.0	1.0 22.0	1.2 22.0
	001/006	001/008	001/010	001/012	001/014	001/016	001/018	001/021	001/023	010/008	010/010	010/012


												
RA 37	RA 38	RA 39	RA 40	RA 41	RA 245	RA 556	RA 557	RA 558	RA 559	RA 560	RA 561	RA 562
1.4 22.0	1.6 22.0	1.8 22.0	2.1 22.0	2.3 22.0	2.5 22.0	4.0 22.0	4.0 22.0	4.0 22.0	4.0 22.0	4.0 22.0	4.0 22.0	4.0 22.0
010/014	010/016	010/018	010/021	010/023	233/009	107/009	107/010	107/012	107/014	107/016	107/018	107/021

												
RA 563	RA 699	RA 700	RA 701	RA 702	RA 703	RA 704						
4.0 22.0	4.0 22.0	4.0 22.0	4.0 22.0	4.0 22.0	4.0 22.0	4.0 22.0						
107/023	168/009	168/010	168/012	168/016	168/021	168/023						

Trimming & Finishing Gold TC bur 8-12 Blades 金色钨钢精修车针 (8-12刃)										
	Round	7004	7005	7006	Taper Cone	7204	7205	Taper Pointed	7611	7612
		1.4 19.0	1.6 19.0	1.8 19.0		9.0 23.0	9.0 23.0		4.2 19.0	6.0 21.0
	ISO 071	001/014	001/016	001/018	ISO 072	184/014	184/016	ISO 072	159/010	164/014


											
Egg	7404	7406	Taper T	7114	7214	7714	Sharp Flame	7901	7902	7903	7903L
	3.3 19.0	3.5 19.0		3.5 19.0	5.5 19.0	9.0 19.0		3.6 19.0	3.6 19.0	3.6 19.0	9.0 23.0
ISO 072	277/014	277/018	ISO 072	168/009	170/009	172/009	ISO 071	496/009	496/010	496/012	250/012

Diamond Polisher
旋风轮
 树脂抛光利器，也可以抛光瓷，氧化锆和玻璃陶瓷。



Coarse 粗
 Medium 中
 Fine 细













Sr204C **Sr204M** **Sr204F**













KIT CP 6001 套装











Arkansas Stone
白车石
 Finish, Polish Filling Composite 去树脂, 填充物及抛光。

											
2.8 6.8	2.8 7.2	2.5 6.2	2.5 6.2	3.0 3.0	3.4 2.0	2.8 6.8	2.8 7.2	2.5 6.2	2.5 6.2	3.0 3.0	3.4 2.0
FG小杆						RA弯头杆					
FL2-Fw	CN1-Fw	TC1-Fw	CY2-Fw	RD1-Fw	IC2-Fw	FL2-Rw	CN1-Rw	TC1-Rw	CY2-Rw	RD1-Rw	IC2-Rw

Green Stone
绿车石
 Grind Ceramics, Precious Metal 研磨烤瓷, 贵金属等。

									
2.8 6.8	2.8 7.2	2.5 6.2	2.5 6.2	3.0 3.0	2.8 6.8	2.8 7.2	2.5 6.2	2.5 6.2	3.0 3.0
FG小杆					RA弯头杆				
FL2-Fg	CN1-Fg	TC1-Fg	CY2-Fg	RD1-Fg	FL2-Rg	CN1-Rg	TC1-Rg	CY2-Rg	RD1-Rg


Diamond Stone
陶瓷金刚石磨头
 Grind, Finish All-Ceramic, Zirconia 全瓷, 氧化锆研磨

							
3.0 7.0	4.0 7.0	3.5 11.0	3.0 7.0	4.0 7.0	3.5 11.0	2.5 6.0	8.0 1.0
FG小杆			RA弯头杆				
DCF101	DCF102	DCF122	DCR101	DCR102	DCR122	DCR104	DCR105

Diamond Polisher . Polish All-Ceramic, Zirconia 橡胶金刚石磨头 (全瓷牙, 氧化锆牙抛光)


															
3.0 7.0	5.0 0.0	10.0 2.0	6.0 10.0	3.0 7.0	5.0 10.0	10.0 2.0	6.0 10.0	3.0 7.0	5.0 10.0	10.0 2.0	6.0 10.0	3.0 7.0	5.0 10.0	10.0 2.0	6.0 10.0
Blue, Medium Grit / 蓝色, 标准粒度				Light Grey, Fine Grit / 浅灰色, 较细				Yellow, Extra Fine / 黄色, 特细				Light Pink, Ultra Fine / 浅粉色, 超细, 高亮			
Sr101M	Sr103M	Sr105M	Sr106M	Sr101F	Sr103F	Sr105F	Sr106F	Sr101EF	Sr103EF	Sr105EF	Sr106EF	Sr101UF	Sr103UF	Sr105UF	Sr106UF

Silicon Polisher For Ceramic 硅胶磨头(烤瓷牙抛光)



3.0 7.0	5.0 10.0	10.0 1.0	10.0 2.0	7.0 11.0	11.0 9.0	3.0 7.0	5.0 10.0
Cr101C	Cr103C	Cr104C	Cr105C	Cr106C	Cr108C	Cf101C	Cf103C
White, Coarse / 白色, 粗糙					FG高速机用		
3.0 7.0	5.0 10.0	10.0 1.0	10.0 2.0	7.0 11.0	11.0 9.0	3.0 7.0	5.0 10.0
Cr101M	Cr103M	Cr104M	Cr105M	Cr106M	Cr108M	Cf101M	Cf103M
Pink, Medium / 粉色, 中粗					FG高速机用		
3.0 7.0	5.0 10.0	10.0 1.0	10.0 2.0	7.0 11.0	11.0 9.0	3.0 7.0	5.0 10.0
Cr101F	Cr103F	Cr104F	Cr105F	Cr106F	Cr108F	Cf101F	Cf103F
Blue, Fine / 蓝色, 较细					FG高速机用		

Silicon Polisher For Amalgam/Metal 硅胶磨头(银汞, 金属抛光)



3.0 7.0	5.0 10.0	10.0 1.0	10.0 2.0	7.0 11.0	11.0 9.0	3.0 7.0	5.0 10.0
Ar101C	Ar103C	Ar104C	Ar105C	Ar106C	Ar108C	Af101C	Af103C
Black, Coarse / 黑色, 粗糙					FG高速机用		
3.0 7.0	5.0 10.0	10.0 1.0	10.0 2.0	7.0 11.0	11.0 9.0	3.0 7.0	5.0 10.0
Ar101M	Ar103M	Ar104M	Ar105M	Ar106M	Ar108M	Af101M	Af103M
Brown, Medium / 棕色, 中粗					FG高速机用		
3.0 7.0	5.0 10.0	10.0 1.0	10.0 2.0	7.0 11.0	11.0 9.0	3.0 7.0	5.0 10.0
Ar101F	Ar103F	Ar104F	Ar105F	Ar106F	Ar108F	Af101F	Af103F
Light Green, Fine / 浅绿, 较细					FG高速机用		

Polishing Brush. RA 抛光刷

Shape of Cup 杯型 White, Hard 硬	White, 白	Clear, 透明	Blue, 蓝	Pink, 粉	Purple, 紫	Yellow, 黄	White, Soft 软	Shape of Pen 笔型 White, Hard 硬	White, 白	Clear, 透明		
R4Dw	R4Cw	R4Cc	R4Cb	R4Cr	R4Cp	R4Cy	R4G	R5Dw	R5Cw	R5Cc		
Bristle 猪鬃	Colorful Nylon 缤纷尼龙					Goat Hair 羊毛		Bristle 猪鬃	Colorful Nylon 缤纷尼龙			
Blue, 蓝	Pink, 粉	Purple, 紫	Yellow, 黄	White, Soft 软	Pointed Head 尖头 White, Hard 硬	Clear, 透明	White, Soft 软	Flat 平轮	Leather 毛皮	White, Soft 软		
R5Cb	R5Cr	R5Cp	R5cy	R5G	R7Dw	R7Cc	R7G	R1Ns	R1Ry	R1G	R1T	R9Cw
Goat Hair 羊毛			Bristle 猪鬃	Nylon 尼龙	Nylon 尼龙	Goat Hair 羊毛	Felt 毛毡	Leather 毛皮	Goat Hair 羊毛	Thread 棉线	Nylon 尼龙	
								Prophy Cup 抛光杯				
With Al. Oxide 氧化铝				With Silicon Carbide 碳化硅								
R4AI	R5AI	R7AI	R1AI	R4SiC	R5SiC	R7SiC	R1SiC					
Abrasive Fiber 磨料丝								4 Angles, Soft 软		6 Angles, Hard 硬		

Mandrels 夹针		
26RS		Adapter FG into RA
35RS		Snap
53RS		Pop-on
31RS		Head 3.0mm
31RSL		30 mm Long
52RS		Head 5 mm
51R		Head 5 mm Steel
55HS		Adapter FG into HP
32HS		HP, Head 3mm

Bur Blocks for Sterilizing 消毒盒		
No. 201 3 x 5 = 15 FG burs hole	No. 211 10 FG + 5 RA	No. 221 5 FG + 10 RA
No. 241 2 x 6 = 12 RA burs hole	No. 202 3 x 6 = 18 Root canal files	No. 212 18 HP Burs

RF Stainless Burs All Ball 碳钢车针,球形 ISO 330		
RA XL 30.0mm Implant use 种植用球钻 快速处理软组织	Ball	
	RA XL D10	1.0
	RA XL D20	2.0
	RA XL D30	3.0
	RA XL D40	4.0
	RA XL D50	5.0
HP XL 70.0mm Bone Surgical 骨科用球钻	Ball	
	HP XL D10	1.0
	HP XL D20	2.0
	HP XL D30	3.0
	HP XL D40	4.0
	HP XL D50	5.0

No. 1702T

Tungsten Carbide Collection Wallet
钨钢车针样品本




Nearly 200 specification, TC burs
FG, FG XL, FG Short, RA, RA Long
T & F Gold, Endo z etc.

No. 1703P

Polishing burs Collection Wallet
抛光器械样品本





















Nearly 150 specifications,
Silicon Rubber burs,
Diamond Abrasives & polishers,
Polishing brushes etc

Sharpening Stones 磨刀石



Check the instruments for contamination and clean, if required. Clean diamond instruments by Cleaning stones.

1W	1G	2W	2G
75 x 25 x 8 mm		100 x 15 x 15 mm	

	<p>Diamond burs of Ball shape Assorted in FG shank 圆球状金刚石车针混装(高速)</p>	<p>Diamond burs of Ball shape Assorted in FG LONG 圆球状金刚石车针混装(高速加长)</p>	<p>Diamond burs of Ball shape Assorted in RA shank 圆球状金刚石车针混装(低速)</p>
<p>KIT</p>	<p>BR1901 BR-49, BR-45, BR-46, BR-41, BR-40, BR-31, BR-28, BR-26, BR-25, BR-23.</p>	<p>BR1902 BR-X45, BR-X46, BR-X41, BR-X40, BR-X31, BR-X40S, BR-X31S, BR-X41C, BR-X40C, BR-X31C</p>	<p>BR1903 BR-R46, BR-R41, BR-R40, BR-R31, BR-R26, BR-R40F, BR-R31F, BR-R26F, BR-R40EF, BR-R26EF</p>
<p>Diamond burs of Inverted Cone Assorted in FG shank 常用倒锥金刚石车针混装(高速)</p>	<p>Regular Diamond burs of Taper Assorted in FG shank 常用锥形金刚石车针混装(高速)</p>	<p>Interproximal Enamel Reduction (IPR) for Ortho 邻面去釉套装</p>	<p>Safe End Diamond Preparation 前后牙预备(安全端金剛車針)</p>
			
<p>1804IV SI-45 x 2; SI-46 x 3, SI-47 x 3, SI-48 x 2.</p>	<p>1805TA TR-11 x 2, TR-12 x 2, TF-11 x 2; TF-12 x 1, TC-10 x 1, TC-27 x 1, TC-28 x 1</p>	<p>1806 TC-11EF, 820-016EF, 820-031F, 820-016F, TC-10F</p>	<p>4100 508T-014S, BR-41, 198-016S, 508T-016S, 508T-016F; 508T-014S, BR-540, 198-021S, 508T-021S, 508T-021F</p>
<p>Crown & Bridge Preparation for Anterior 前牙桥体/冠体预备</p>	<p>Crown & Bridge Preparation for Posterior 后牙桥体/冠体预备</p>	<p>Cavity Preparation, Restoratives 备洞</p>	<p>Inlay & Onlay Preparation 嵌体/高嵌体预备</p>
			
<p>4101 TF-21, TF-10, TF-12, TC-10, WR-13 SR-11, TR-11, EX-18, TR-12F, TC-11F</p>	<p>4102 EX-20S, SR-21, FO-34, TR-62C, TR-21 TR-24, SO-21, EX-18, TR-21F, EX-21F</p>	<p>4103 BR-45, BR-40, BC-32, EX-31 TF-20, SF-41, SF-21, SI-48</p>	<p>4104 SI-48, DI-42, TF-31, TF-22, SF-21 EX-18, TF-21F, SF-21F, FO-22EF, EX-21EF</p>
<p>Micro Invasive Surgery 微创备牙</p>	<p>Porcelain Veneers' Preparation 瓷贴面备牙套装</p>	<p>Children's Teeth Preparation 儿童牙预备</p>	<p>Ceramics / Composite Polishing 烤瓷牙/复合树脂抛光</p>
			
<p>4105 BR-49, SF-41, SR-41, TC-26, CR-41 EX-41, FO-40EF / FG1/4, FG1/2, FG 1</p>	<p>6000 BR-45, VP-20, VP-21, VP-22, SR-41, TR-25, SO-20, TR-13F, Cf101F</p>	<p>6002 CD-50F, CD-51F, CD-52F, CD-53F, CD-54F CD-55F, CD-56F, CD-57F, CD-58F, CD-59F</p>	<p>4106 BR-40EF, TC-11EF, TF-11EF, TF-12EF, TR-11EF, EX-21EF, FO-40EF, FO-21EF, FO-22EF, FO-30EF</p>
<p>Student or Beginner's Diamond FG 0.8mm 新手入门金刚石车针,高速</p>	<p>Student or Beginner's Diamond FG 0.9mm 新手入门金刚石车针,高速</p>	<p>Student or Beginner's Diamond RA 0.8mm 新手入门金刚石车针,慢速</p>	<p>Student or Beginner's Diamond RA 0.9mm 新手入门金刚石车针,慢速</p>
			
<p>FG008D BR-49, BR-49, SI-49, SF-49, TF-49</p>	<p>FG009D 001-009M, 001-009M, 010-009M, 109-009M, 170-009M</p>	<p>RA008D BR-R49, BR-R49, SI-R49, SF-R49, TF-R49</p>	<p>RA009D 001-R009M, 001-R009M, 010-R009M, 109-R009M, 170-R009M</p>

<p>TC Trimming & Finishing burs, FG 钨钢精修车针,FG</p>	<p>Tungsten carbide burs of Ball shape Assorted in FG shank 圆球状钨钢车针混装(高速)</p>	<p>Tungsten Carbide of FG shank 0.8mm for Students 初学者用的钨钢车针 (头径0.8mm,高速)</p>	<p>Tungsten Carbide of RA shank 0.8mm for Students 初学者用的钨钢车针 (头径0.8mm,低速)</p>
			
<p>TF7100 7004, 7404, 7214, 7714, 7903L</p>	<p>1900 T FG 1, FG 2, FG 3, FG 4, FG 5</p>	<p>FG008T FG 1, FG 1, FG 34, FG555, FG698</p>	<p>RA008T RA 1, RA 1, RA 34, RA 555, RA 698</p>
<p>FG White Stone Burs, to Finish Composite 白车石,高速</p>	<p>RA White Stone Burs, to Finish Composite 白车石,低速</p>	<p>FG Green Stone burs, Finish Ceramics 绿车石,高速</p>	<p>RA Green Stone burs, Finish Ceramics 绿车石,低速</p>
			
<p>7109FW FL2-Fw, CN1-Fw, RD1-Fw, RD2-Fw, TC1-Fw,</p>	<p>7109RW FL2-Rw, CN1-Rw, RD1-Rw, RD2-Rw, TC1-Rw,</p>	<p>7109FG FL2-Fg, CN1-Fg, RD1-Fg, RD2-Fg, TC1-Fg,</p>	<p>7109RG FL2-Rg, CN1-Rg, RD1-Rg, RD2-Rg, TC1-Rg,</p>
<p>Extraction Burs 拔阻生齿车针</p>	<p>Oral Surgery burs FG Long 手术车针(高速加长)</p>	<p>Porcelain / Natural Teeth Polishing RA 烤瓷牙,自然牙抛光,弯机柄</p>	<p>Amalgam / Metal Polishing RA 银汞合金,金属抛光,弯机柄</p>
			
<p>1904 X151, X1151, EC25, EC28, FG X6, EX-28, BR-X31, TR-X12, TR-X15, TF-X13</p>	<p>1905 FG X2, FG X4, FG X6, FG X558, FG X701, BR-X41, BR-X40, BR-X31, TF-X12, TR-X12</p>	<p>4112 Coarse : Cr101C, Cr103C, Cr104C, Cr106C Medium : Cr101M, Cr103M, Cr104M, Cr106M Fine : Cr101F, Cr103F, Cr104F, Cr106F</p>	<p>4114 Coarse : Ar101C, Ar103C, Ar104C, Ar106C Medium : Ar101M, Ar103M, Ar104M, Ar106M Fine : Ar101F, Ar103F, Ar104F, Ar106F</p>
<p>Composite, Ceramic, Zirconia Polishers . RA 树脂、烤瓷、氧化锆抛光</p>	<p>Composite, Ceramic, Zirconia Polishers . RA 树脂、烤瓷、氧化锆抛光</p>	<p>Diamond stones & polishers for grind & polish Zirconia / Lithium Disilicate 氧化锆,玻璃陶瓷研磨抛光</p>	<p>Diamond Polishers for Zirconia /Lithium Disilicate 氧化锆,玻璃陶瓷抛光</p>
			
<p>CP 6011 Medium : Sr 204 M x 3 Fine : Sr 204 F x 3</p>	<p>CP 6001 Coarse : Sr 204 C x 2 Medium : Sr 204 M x 2 Fine : Sr 204 F x 2</p>	<p>7115 Sr101UF, Sr103UF, Sr105UF, Sr106UF Sr101F, Sr103F, Sr105F, Sr106F DCR101, DCR102, DCR122, DCR105</p>	<p>4115 Medium : Sr103M, Sr105M, Sr106M Fine : Sr103F, Sr105F, Sr106F Extra Fine : Sr101EF, Sr103EF, Sr105EF, Sr106EF, Ultra Fine : Sr101UF, Sr103UF</p>
<p>Composite Finish and Polish RA 树脂成型和抛光,弯机柄</p>	<p>Composite Finish and Polish RA 树脂成型和抛光,弯机柄</p>	<p>Composite Finish and Polish RA 树脂成型和抛光,弯机柄</p>	
			
<p>4110V RA Stone Rw01, Rw02 x 2 RA Silicon Rubber Cr101M, Cr103M, Cr104M, Cr106M, Cr101F, Cr103F, Cr104F, Cr106F Cr10106M</p>	<p>1906 Cr103C, Cr104C, Cr106C, Cr103M, Cr106M, R4AL, R5AL, R1Ry, R1G, R1T</p>	<p>1907 Cr103C x 2, Cr103M, Cr104C, Cr106C R4AL, R5AL, R4Cb, R5Cb, R4Cr, R5Cr.</p>	

Dianfong ABRASIVES



ISO 13485 / CE 0197 / Certified



阿里官网
Website on Alibaba



微信公众
WeChat Public

深圳点锋研磨材料有限公司
Shenzhen Dian Fong Abrasives Co.,Ltd

地址: 广东省深圳市龙岗区平湖天华路5号2栋4楼
Address : 4th Floor, Blk 2, No. 5, TianHua Rd, PingHu, LongGang Distr.,
Shenzhen, Guangdong Province, China 邮编 : Post code : 518 111
网站(Web): www.dianfong.com ; 邮箱(Email): dianfong@foxmail.com
电话(Tel.): 0755-2885 6248 ; 2834 4545 ; 传真(Fax): 0755- 2885 6246
QQ: 27212989 ; Skype:welson881 ; WhatsApp/Wechat:0086-13528829441