



# **Disposable Hemostasis Clips & Delivery Systems**

## Intended use

Clamping the mucosa tissue of alimentary cavity under the guide of endoscope.

#### A) Hemostasis for

- 1) Mucosal / sub-mucosal defeats <3cm
- 2) Bleeding ulcers
- 3) Arteries <2mm
- 4) Surgery sites
- 5) Diverticula in colon

#### B) Closure of GI luminal perforations

#### **Features**

- · Accurate sync-rotation offers exact clip deployment.
- · Reinforced design makes the clip performance more steady.
- · Convenient re-load method.

### **Specifications**







### Single-Use Hemoclips

Model	Jaw angle	Arm length	Openning	Total closed length	Picture
SX-C-090	90°	5.0mm	12mm (max)	13.5mm	
SX-C-135	135°	5.7mm	12mm (max)	14.2mm	
					671.0

#### Re-Usable Delivery System(applicators)

Model	Working length	Diameter	Working channel	Sterile	Picture
SX-H-1650	1650mm	≤2.75mm	≥2.8mm	Auto-clavable / cold-immerse	
SX-H-1950	1950mm	≤2.75mm	≥2.8mm	Auto-clavable / cold-immerse	
SX-H-2300	2300mm	≤2.75mm	≥2.8mm	Auto-clavable / cold-immerse	

