


BOJIN MEDICAL www.bojin-medical.com



System 4200 Multi Functional Surgical Power Tools

For Joint surgery | Trauma surgery | Neurosurgery | Thoracic surgery



Auction no. 21075923 / ocds-b3wdp1-MD-1678893907083 Technical specifications

Electric Multifunctional Saw Drill BYJZ-I 42 (Lot 7)

Auction specifications	BOJIN Specifications		
Hand Piece - BJ4200			
Description of Function: The Battery Operate of attachment and accessories for use in general K-Wires.	ed Drill system is a battery driven power system with a range ral traumatology, including drilling, sawing, inserting/removing		
Electric Hand piece: Yes	Electric Hand Piece (Battery operated)		
Power: ≥ 100 W	Power: ≥ 100 W		
Weight: ≤ 950 g	Weight: ≤ 950 g		
Sterilization: Autoclavable 135°	Sterilization method: Autoclavable 135°		
Speed: ≥ 900 rpm	Adjustable speed: 0 - 1100 rpm		
Noice: ≤ 75 db	Noice: ≤ 75 db		
Blocking mechanism: Yes	Battery & Attachment blocking mechanism		
Quick Coupling: Yes	Quick Coupling Attachments		
Increasing temperature while working: ≤20°C	Increasing temperature while working: ≤ 20°C		
Sterilizable cover: Yes	Autoclavaible handpiece, Battery cover & Aseptic battery transfer kit		
	Forward and Reverse switch with safe mode		
	Pistol grip handpiece, ergonormic design		
Battery - BJ4001			
Quantity: 2 pcs., easily removable	2 pcs.		
	Insert/detachment method: telescopic, easily		
Run time: ≥ 100 min	Run time with fully charged battery: ≥ 30 min		
	Ni-Mh rechargeable battery, 14.4V/1800mAh		
	Non-memory		
	Non-autoclavaible		
	Shield to protect battery from the housing		
Battery Charger - BJ1002-4			
Simultaneosly charged batteries: 2 batteries	4 batteries		
Battery charge status: Yes	Light indicator		
	Input Voltage: 110-250V, 50/60Hz.		
	The output voltage: 17.5*4 (this is for four charger).		
	Output Current: 700mAh		
Drill Attachment - BJ4203B	A RIVANI		
	Adjustable speed: 0 800 rpm / 0 - 1100 rpm		
	Sterilization method: Autoclavable 135°		
	Output power, 3 Nm		
	Forward & Reverse mode		
	Dual mode function: drilling and K-Wire inserting		

Sagittal Saw Attachment - BJ4201				
360° rotating: Yes	Rotating: Up, down, left and right (360°)			
Sawing Speed: ≥ 1000 cpm	Adjustable speed: 0 - 15000 cycles/min			
Sterilization: Autoclavable	Sterilization method: Autoclavable 135°			
	Tool less blade mounting			
	Blade arc of excursion: 5 degree			
Reciprocating Saw Attachment - BJ4209	9			
Sawing Speed: ≥ 1000 cpm	Adjustable speed: 0 - 15000 cycles/min			
Sterilization: Autoclavable	Sterilization method: Autoclavable 135°			
Changing and attaching blades without a key: Yes	Tool less blade mounting			
Acetabulum Reaming Drill Attachment	- BJ4207B			
Speed: Adjustable	Adjustable speed: 0 - 300 rpm			
Sterilization: Autoclavable	Sterilization method: Autoclavable 135°			
	Forward & Reverse mode			
	Output power (to 250 rpm): 8 Nm			
K-Wire Drill Attachment (Ø 0.8mm - 2.0m	mm) - BJ4213A			
Ø 0.8mm - 4.0mm	Ø 0.8mm - 2.0mm			
	Adjustable speed: 0 - 1100 rpm			
	Forward & Reverse mode			
	Sterilization method: Autoclavable 135°			
K-Wire Drill Attachment (Ø 2.0mm - 3.0	mm) - BJ4213B			
Ø 0.8mm - 4.0mm	Ø 2.0mm - 3.0mm			
	Adjustable speed: 0 - 1100 rpm			
	Forward & Reverse mode			
	Sterilization method: Autoclavable 135°			
Other Attachment: Optional				
	High speed drill attachment BJ4202			
	High speed drill attachment BJ4202SAO			
	Cannulated drill attachment BJ4203D			
	Acetabulum reamer attachment BJ4207B			
Accessories				
Power cable:	European power cable (battery charger)			
Sagital saw blade 10 pcs.	3 different size (10 pcs.): BJ90221.0A, BJ90141.0A, BJ90251.0A			
Reciprocating saw blade 10 pcs.	2 different size (10 pcs.): WF70131.0B, 2WF70131.0B			
Drill bits 10 pcs	9 different size (10 pcs,): BJZC-1.5-35-115 BJZC-2.0-35-115 BJZC-2.5-35-115 BJZC-2.7-50-115 BJZC-3.0-50-115 BJZC-3.5-50-115 BJZC-4.0-50-115 BJZC-4.5-50-115			

Key 1 pc.	Jacob Chuck Key BJ1005 (2 pcs.)		
Sterilization container 1 pc.	Code: BJXDH4200.		
Case for storage and transportation: Yes	Container for sterilization, transportation and storage, with fixing accessories during sterilization.		
Kirschner wire:			
Ø 0.8mm - 50 pcs.	Ø 0.8mm x 150mm - 50 pcs.		
Ø 1.0mm - 100 pcs.	Ø 1.0mm x 150/225/300mm - 100 pcs.		
Ø 1.2mm - 100 pcs.	Ø 1.2mm x 150/225/300mm - 100 pcs.		
Ø 1.5mm - 100 pcs.	Ø 1.5mm x 225/250/300mm - 100 pcs.		
Ø 1.8mm - 100 pcs.	Ø 1.8mm x 300/400mm - 100 pcs.		
Ø 2.0mm - 200 pcs.	Ø 2.0mm x 300/400mm - 200 pcs.		
Ø 2.5mm - 50 pcs.	Ø 2.5mm x 300/400mm - 50 pcs.		
Ø 3.0mm - 100 pcs.	Ø 3.0mm x 300mm - 100 pcs.		
Other requirements			
Warranty: 2 years.	Handpiece, Attachments & Battery Charger: 2 years warranty.		

Date: 12.04.2023

Name: Anne Lee
Position: Deputy General Manager

BJ1002-4 Battery Charger (four chargers)



Fit for BJ1001, BJ4001, BJ6001 battery.

Input Voltage: 110-250V, 50/60Hz.

The output voltage: 17.5*4 (this is for four charger).

Output Current: 700mAh.

Battery charge status: Light indicator.

Size: 26(L)*26(W)*11cm(H)

Net weight: 1.4 kgs

Package size: 40(L)*28(W)*18(H)

GW: 2kgs

The manufacturers (sted in this document are not affiliated with Bojin and have not endorsed the Bojin products listed in this catalogue.

This documentation is intended exclusively for physicians and is not intended for laypersons. Internation on the products and procedures contained in this document is of a general nature and does not represent and does not constitute medical advice or recommendations. Because this information does not purport to constitute any diagnostic or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this document does not replace the need for such examination and/or advice in whole or in part.

The information and illustrations included in this catalogue are correct and true to the best of our knowledge. However, we cannot be regarded as binding warranty for audity. The user is responsible for checking the products prior to use whether they are suited for the intended purpose. We can only be held liable for these instruments pursuant to our terms of delivery and payment. Subject to technical alterations.





Shanghai Bojin Medical Instrument Co...Ltd

Company address: Room 606, No.3388, Gonghexin Road, Shanghai, China Fectory address: 1-D, 2-A, No.125, Longpon Road, Jioding District, Shanghai, China

7e0 +86-21-66308078 66308077 66526503 Fax:+86-21-66527863


Website: www.bajin-medical.com fmai: manager@bojin-medical.com



ORTHOPEDIC CONSUMABLES

OVERVIEW





Product code

BJ90181.0 A

Product code

BJ90181.2 A

Product code

BJ90251.0 A

Product code

BJ90251.2 A

Product code

BJ90251.27 A

50

50

Equivalent code

Equivalent code

Equivalent code

Equivalent code

Equivalent code

90 80 70 80



Width

18mm

Width

18mm

Width

25mm

Width

25mm

Width

25mm

Length

90mm

Length

90mm

Length

90mm

Length

90mm

Length

90mm

Thickness 1.2mm

A

Thickness

1.0mm

Thickness

1.2mm

Thickness

1.27mm

www.bojin-medical.com

A

Width Thickness Length 90mm 16mm 1.27mm

8 7 8 8 8 20

Equivalent code

Equivalent code Width Thickness. Length 80mm 34mm 0.8mm

8 8 8 8 8

Width Equivalent code Thickness Length 1.0mm 80mm 34mm

Length Width Thickness 80mm 34mm 1.2mm

70 60 50 40 Product code

BJ8034.127 A

Bojin'

Equivalent code

Width Length 80mm 34mm

Thickness. 1.27mm

www.bojin-medical.com


Α

A

A


Bojin'





Product code





Length

28mm



www.bojin-medical.com

Product code

WFSW2860.8 B

Equivalent code

Thickness


0.8mm

9


Width

6mm

ORTHOPEDIC DRILL BITS




Length(mm)	Connect/Shank Length(mm)	Effective Cutting Length(mm)	Mead Diameter(mm)
115	80	35	Φ1.5
115	80	35	Φ2,0
115	60	55	Φ2.5
115	60	55	Ф2.7
115	60	55	Ф3.2
115	60	55	Φ3.5
115	60	55	Φ4.0
115	60	\$5	Ф4.2
115	60	55	Φ4.5
115	60	55	Φ5.0



AO Drill bits (Synthes quick connect coupling)

Length(mm)	Connect/Shank Length(mm)	Effective Cutting Length(mm)	Head Diameter(mm)
125	100	25	Φ1.1
125	100	25	Φ1.5
125	100	25	Φ1.8
125	95	30	Ф2.0
125	95	30	Φ2.5
125	95	30	Ф2.7
125	85	40	Ф3.2
125	85.	40	Ф3.5
125	85	40	Φ4.0
125	85	40	Ф4.5

ORTHOPEDIC DRILL BITS



Length (mm)	Connect/Shank Length (mm)	Effective Cutting Length (mm)	Head Diameter (mm)	Inner cannulation (mm)
150	120	30	Φ2.7	Φ1.4
BJZCZ-2.7-30-150			<u> </u>	
150	115	35	Ф3.2	Ф1.7
BJZCZ-3.2-35-150				
150	100	50	Ф3.8	Φ1.7
BJ2CZ-3.8-50-150				
150	130	20	Φ2.0	Φ1.1
BJZCZ-2.0-20-150				
150	130	20	Φ2.2	Ф0.9
BJZCZ-2.2-20-150				
250	210	40	Φ4.5	Φ2.1
BJZCZ-4.5-40-250				
300	230	70	Φ6.0	Ф2.7
BJZCZ-6.0-70-300				