

Technical Data Sheet

Model No.: N61G01

Product Name: Surgical drape

Classification: Class Is, according to Annex VIII in MDR (EU) 2017/745

Description: Surgical drape is made of SMS. There are different

design styles: with fenestration or not, with tape or not. It could be used for different operations according to its intended use. It is intended to be used to isolate the dust and bacteria, prevent

the cross infection.

Intended use: Used to isolate the dust and bacteria, prevent the cross infection.


Material: SMS 30gsm

Approvals: CE, ISO 13485

Manufacturer: Leboo Healthcare Products Limited

Country of origin: China

Specification: Detail see specification as below.



Material	Weight (gsm)	Size(cm)	Measurement Location(cm)				Color	Stylo
			Α	В	С	D	Coloi	Style
SMS	30	200x160	200	160	1	1	Blue	Individually packaged; EO Sterile

Check: Date: Design: Approve:




⊚ H.O.: A-1506, Lar Valley International, No.168 Guang'anmen Wai St., Xicheng District, Beijing 100055, China.

Packaging:

- (1) Carton Requirements: Five layers corrugated
- (2) Packing Quantity: 1pc/pouch, 100pcs/ctn
- (3) 5 Layers Carton(Paper)





Date:

Storage:

- (1) Store in dry, clean conditions in original packaging.
- (2) Store away from direct sunlight, source of high temperature and solvent vapors.
- (3) Store with the temperature range -5 $^{\circ}$ C to +45 $^{\circ}$ C and with relative humidity below 65%.

Shelf Life:

Design:

Shelf life is 36 months from the date of manufacture when stored as stated above.

Cautions and disposal:

(1) Do not contact with heavy oils, sparks or flames;

Check:

- (2) Do not use in the environment with conditions of excessive heat;
- (3) Replace if damaged, heavily contaminated or in accordance with local practice;
- (4) Hand and dispose of contaminated products with care and in accordance with national regulations.

(€ 0197	CE Symbol	MD	Medical Device			
*	Keep away from sunlight	2	Do not reuse			
\bigcirc	Keep dry		Don't use when packing damaged and consult instructions for use			
45°C	Temperature limit	STERILE EO	EO Sterile			
		LEBOO HEALTHCARE PRODUCTS LIMITED				

Approve: