

Single-Use Cleaning Brushes – Combination Brushes



- Combination brush (channel-cleaning brush and channel-opening brush on one sheath)
- Stiffer and thicker bristles for higher cleaning efficiency
- High kink resistance and pushability due to new flexible sheath
- Plastic tip on distal end of brush wire to protect scope channel

2

Article number	Article Name	Quantity	Working Length	Instrument Channel ø	Endoscope Compatibility	Area Compatibility
N3640750	BW-411B	50	950 mm	2.0 - 3.2 mm	Channel-cleaning brush part: For Olympus bronchoscopes, choledochoscopes, nasopharyngoscopes, cystoscopes, rhinolaryngoscopes, and hysteroscopes, MAF-TM, MAF-TM2 Channel-opening brush part: For all Olympus endoscopes with an instrument channel diameter of 1.2–6.0 mm	ENT/SURG/URO
N3640850	BW-412T	50	2200 mm	2.0 - 4.2 mm	Channel-cleaning brush part: For all Olympus gastrointestinal endoscopes with an instrument channel working length of 1800 mm or less	