


INTRODUCERS

Comprehensive product line



- Excellent transitions
- Hemostasis valve
- Color-coded hubs
- Variety of sizes

BALTON
since
1980

BALTON®

CE 0197

INTRODUCERS

INTRODUCER

Set with hemostasis valve to introduce and exchange catheters and endocavitary electrodes

Set contains:

- sheath with a valve
- angiographic needle
- guide wire
- dilatator

Catalog Number	Sheath size / length	Colour identification	Dilatator size	Angiographic needle size	Guide wire size / length / type
INT3F	3F / 60 mm	white	3F	21G	.018" / 35 cm / STRAIGHT
INT4F	4F / 70 mm	red	4F	19G	.022" / 40 cm / "J"
INT5F	5F / 110 mm	grey	5F	18G	.035" / 40 cm / "J"
INT5F23	5F / 230 mm	grey	5F	18G	.035" / 60 cm / "J"
INT6F	6F / 110 mm	green	6F	18G	.035" / 40 cm / "J"
INT6F23	6F / 230 mm	green	6F	18G	.035" / 60 cm / "J"
INT7F	7F / 110 mm	orange	7F	18G	.035" / 40 cm / "J"
INT7F23	7F / 230 mm	orange	7F	18G	.035" / 60 cm / "J"
INT8F	8F / 110 mm	blue	8F	18G	.035" / 40 cm / "J"
INT8F23	8F / 230 mm	blue	8F	18G	.035" / 60 cm / "J"
INT9F	9F / 110 mm	black	9F	18G	.035" / 40 cm / "J"
INT9F23	9F / 230 mm	black	9F	18G	.035" / 60 cm / "J"
INT10F	10F / 110 mm	white	10F	18G	.035" / 40 cm / "J"
INT10F23	10F / 230 mm	white	10F	18G	.035" / 60 cm / "J"
INT11F	11F / 110 mm	blue	11F	18G	.035" / 40 cm / "J"
INT11F23	11F / 230 mm	blue	11F	18G	.035" / 60 cm / "J"

When ordering use Ref. Number, e.g.: INT6F.

Special order can contain introducer with radiopaque marker (available for sizes greater than 4F).

Model reference number: **INT8FM** – introducer with a marker.

Introducers are intended for introducing catheters and endocavitary electrodes.

The introducer's sheath is delivered with hemostasis valve, which prevents from losing blood during the surgery and it minimizes the embolism manifestation possibility.

The connection of the drain with sheath's head enables the infusion of fluids during the electrode placement procedure.

Caution:

This publication could include technical or other inaccuracies or typographical errors. Changes are periodically added to the information on herein; these changes will be incorporated in the new editions of the publication. Balton Ltd. may make improvements and/or changes in the product(s) described in this publication at any time. This publication is only a commercial information and cannot be handled the commercial offer in the meaning of the art. 66 §1 of Polish Civil Code.



294 Modlinska Str., 03-152 Warsaw, POLAND
ph.: (+48) 22 597 44 00, fax: (+48) 22 597 44 44
e-mail: balton@balton.pl, www.balton.pl