

MOBILE

- Compact design
- Robust metal housing
- Easy to move thanks to built-in wheels
- Optional trolley for a more fluid transport, especially for frequent movements or long distances and to fix the pre-treatment

EASY TO USE

- Simple access to settings via a 2-button user interface
- Low noise operation so it can be used beside the patient without any disturbance
- Integrated pre-filter with pressure reducer

PRODUCT RANGE OVERVIEW

HemoRO 3000

Reference codes	Description
02HRO3001N-UG	HemoRO 3000 (50 l/h) - 230 V / 50-60 Hz
02HRO3002N-UG	HemoRO 3000 (70 l/h) - 230 V / 50-60 Hz
02HRO3110N-UG	HemoRO 3000 (50 l/h) - 115 V / 60 Hz
02HRO3111N-UG	HemoRO 3000 (70 l/h) - 115 V / 60 Hz

HemoRO 3000 Pre-treatment

Reference codes	Description
S08BHH000-UG	Carriage for HemoRO 3000
09WWG0001-UG	Water guard
S07B059G8-UG	Softener
S08B09895A-UG	3 pre-filters: 10"-AC, 10"-20µm, 5"-5µm
EJMGE41015-UG	Spare membrane M4.6-11
S08B08490-UG	Citic acid set (5 pcs.)
20KSET110-UG	Sampling set for permeate (6 pcs)