

RADIFOCUS® TORQUE DEVICE



Accessories

The Radifocus Torque Device is a supportive device mainly intended for use when guiding a disposable guidewire, which is used for percutaneous angiography and treatment under fluoroscopy, to the target site.

Product Characteristics

Torque device is compatible with guidewires between 0.010" (0.26 mm) and 0.038" (0.97 mm) in size.

Item specifications

| Guidewire Compatibility - Diameter | Code |
|------------------------------------|----------|
| 0.010 in 0.038 in | XX*RF02M |

Please quote above item reference codes when placing an order