

Sales only to commercial customers. No sales to consumers within the meaning of §13 BGB.



Zoom in

Bacterial/Viral Filter, with HME, sterile

Bacterial/Viral Filter with HME and CO₂-port. The purpose of the filter is to protect the patient and to humidify the inhaled air. The filter is placed between the endotracheal tube/tracheal cannula and the Y-piece of the breathing system.

Shape (Ventilation Filter):	Round
ISO style (patient end):	Male/Female
Inner diameter (patient end):	20/15 mm
Sterile:	Yes

Country availability: Please note: not all products are available in all countries. Feel free to contact our international sales team about any questions relating to our sales and services.
Jutta Dorenberg
Director International Sales
Tel.: +49 2236 3913 - 405
eMail: jutta.dorenberg@dahlhausen.de

shorten text

Article no.	Application area	Particle retention	ISO style (machine end)	Inner diameter (machine end)	Max. tidal volume (Breathing Filter)	minimal quantity	
4300905707	Adult	99,999 %	Male/Female	22/13 mm	1.500 ml	50	