

TIP.DOC, 1

LEI: Noua Mii Trei Sute lei 00 bani

L.S.

CERTIFICAT
privind lipsa sau existența restanțelor față de bugetul public național

Nr.
№ **A2221379**

din
от **26.10.2022**

1. Destinația / Назначение

Pentru participarea la proceduri de achiziții publice

2. Date despre contribuabil / Информация о налогоплательщике

Denumirea Наименование	Codul fiscal / Numărul de identificare Фискальный код / Идентификационный номер
BIOSISTEM MLD S.R.L.	1010600028048
Adresa sediului de bază (strada, numărul) Адрес основного месторасположения (улица, номер)	Codul - Denumirea localității Код - Наименование населенного пункта
Albisoara nr.16 bl.1 of.7	0150-SEC.RISCANI

3. Atestarea lipsei sau existenței restanțelor conform datelor Sistemului Informațional Automatizat /

Подтверждение отсутствия или наличия недоимки согласно данных Информационной автоматизированной системы

La data emiterii prezentului certificat restanța față de bugetul public național constituie/ На дату выдачи данной справки недоимка перед национальным публичным бюджетом составляет:
0,00 lei/лей.

4. Valabil până la / Действителен до 10.11.2022

5. Autentificarea Serviciului Fiscal de Stat / Подтверждение Государственной налоговой службы

**Șef DDF Rîșcani
a DGDF**

Funcția/Должность

L.Ș/ М.П.

Executor:

Claudia GOJAN

Numele și prenumele/Фамилия и имя

Tel.(022)823102



Ana STOICOV

Numele și prenumele/Фамилия и имя

Este extras din Sistemul Informațional al SFS SIA „Contul curent al contribuabilului”// 26.10.2022 ora 13:13:21
cu aplicarea prevederilor pct. 82-83 Ordin IFPS nr.400 din 14.03.2014 (Monitorul Oficial 72-77/399, 28.03.2014)

NOTA (0,01)



BC "MOLDINDCONBANK" S.A.

Filiala "Invest"

Republica Moldova, MD-2068
mun. Chișinău, bd. Moscovei, 14/1
Tel. : (373-22) 43-44-81, 43-46-24
Fax : (373-22) 43-44-22
cod: MOLDMD2X329

Data 14. IAN. 2016
Nr. 03/2 - 19/23

Республика Молдова, MD-2068
мун. Кишинэу, бул. Московской, 14/1
Тел. : (373-22) 43-44-81, 43-46-24
Факс : (373-22) 43-44-22
код: MOLDMD2X329

Filiala „Invest” BC „Moldindconbank” SA confirmă existența contului curent în moneda națională al **“BIOSISTEM MLD” S.R.L. (c/f 1010600028048)**, cu **IBAN MD95ML000000002251429243**.

Codul băncii MOLDMD2X329.

Director

Nina Turcan

Director financiar



Nina Balmuş

Ex. Diana Brinza
Tel. 43-45-96

REPUBLICA



MOLDOVA

CERTIFICAT DE ÎNREGISTRARE

Societatea cu Răspundere Limitată "BIOSISTEM MLD"
— ESTE ÎNREGISTRATĂ LA CAMERA ÎNREGISTRĂRII DE STAT —

Numărul de identificare de stat - codul fiscal
1010600028048

Data înregistrării

12.08.2010

Data eliberării

12.08.2010

Svirepova Ludmila, registrator

*Funcția, numele, prenumele persoanei
care a eliberat certificatul*

L. Svirepova
semnătura

MD 0101250





I.P. "AGENȚIA SERVICII PUBLICE"

Departamentul înregistrare și licențiere a unităților de drept

EXTRAS

din Registrul de stat al persoanelor juridice

nr. 8506 din 28.04.2021

Denumirea completă: **Societatea cu Răspundere Limitată «BIOSISTEM MLD».**

Denumirea prescurtată: **«BIOSISTEM MLD» S.R.L.**

Forma juridică de organizare: **Societate cu Răspundere Limitată.**

Numărul de identificare de stat și codul fiscal: **1010600028048.**

Data înregistrării de stat: **12.08.2010.**

Sediul: **MD-2001, str. Albișoara, 16/1, ap.(of.) 7, mun. Chișinău, Republica Moldova.**

Obiectul principal de activitate:

- 1 Activitatea farmaceutică;**
- 2 Importul, fabricarea, comercializarea, asistența tehnică și (sau) reparația dispozitivelor medicale și (sau) a opticii;**
- 3 Acordarea asistenței medicale de către instituțiile medico-sanitare private;**
- 4 Comerțul cu ridicata al calculatoarelor, echipamentelor periferice și software-ului;**
- 5 Întreținerea și repararea mașinilor de birou și a tehnicii de calcul;**
- 6 Consultații în domeniul sistemelor de calcul.**

Capitalul social: **5400 lei.**

Administrator: POIATA VITALIE,

Asociați:

- 1. POIATA VITALIE 33,40 %**
- 2. NASEDCHIN ALEXANDR 33,30 %**
- 3. KOJEVNIKOV DMITRII 33,30 %.**

Prezentul extras este eliberat în temeiul art. 34 al Legii nr. 220-XVI din 19 octombrie 2007 privind înregistrarea de stat a persoanelor juridice și a întreprinzătorilor individuali și confirmă datele din Registrul de stat la data de: 28.04.2021.

Specialist coordonator
tel. 022-207-840



Lazari Aliona



EB 0358735

Lista fondatorilor Biosistem-mld SRL

Nr.	Nume, Prenume	IDNP
1.	Vitalie Poiata	0983103892591
2.	Alexandr Nasedchin	2002001070747
3.	Dmitrii Kojevnikov	0972305012362

Liga²®
Titanium Ligating Clips



TITANIUM

Reliable Ligation

Ligating Clips

LigaV[®]

LigaV[®] Ligating clips are implantable medical devices. They are widely used in surgery for marking and ligating purposes. Medical practitioners value the product for its hemostatic features.


LigaV[®] manufactured by Grena Ltd., comes with inner and outer serration which enhances its function. The outer serration helps to grip the clip in the applier's jaws; the inner serration enables its proper nourishment.

The product is offered in 4 standard sizes ensuring every user will be able to find appropriate clip for the given tissue. The LigaV[®] clips are delivered in colour coded cartridges to ensure easy and intuitive identification. A cartridge includes 4 or 6 clips that ensure economical and comfortable use. The product is disposable, sterilised with EO.

Clip appliers with proper sized jaws depending on the clip should be used.

Clip appliers are made of highest quality surgical grade steel. Handle rings are color coded matching cartridge hue.

Low density
High strength
Full biocompatibility
Safe for MRI



Ligating Clips

Liga²®

The LigaV® ligating clips are manufactured of a medical grade titanium which ensures proper level of biocompatibility. The applied technology guarantees the product meets the requirements of even the most demanding clinical institutes.



TITANIUM

Ligating Clips

Ligating Clips

**Cost-effective alternative
for simple procedures**



4 clips per cartridge


Key features

- 4 sizes to ligate tissue structures,
- made of apyrogenic titanium,
- easy loading system,
- clips delivered in color-coded cartridges for easy identification,
- 4 or 6 clips per cartridge to ensure economical and comfortable use of the product,
- V-shaped,
- color-coded applicator rings for easy identification



REF	Clip size	Compatible with appliers (jaws grip) (mm)	Product specification	Colour	Pcs per cartridge	Pcs cartridges per box
0301-01S	Small	0,59-0,75	V-shaped titanium ligating clips		6	50
0301-01M	Medium	0,84-1,00			6	50
0301-01ML	Medium-Large	1,16-1,32			6	20
0301-01ML04	Medium-Large	1,16-1,32			4	20
0301-01L	Large	1,26-1,42			6	20

Open Surgery Clip Appliers



Open Surgery Clip Appliers

**Ergonomic and effective approach
in endoscopic procedures**

Open Surgery Clip Appliers

Open Surgery Clip Appliers

Key features

- available for Small, Medium, Medium-Large and Large non-absorbable LigaV[®] clips to be used on large spectrum of tissue sizes,
- 15cm, 20cm and 28cm clip,
- non-detachable design for easier maintenance and durability, appliers with proper sized jaws made of high quality stainless steel,
- colour-coded rings matching cartridge hue,
- compatible with corresponding sizes of LigaV[®] titanium non-absorbable clips.



REF	Clip size	Applier length	Pcs per box
0301-02S15*	Small	15	1
0301-02S20	Small	20	1
0301-02S28*	Small	28	1
0301-02M15*	Medium	15	1
0301-02M20	Medium	20	1
0301-02M28*	Medium	28	1
0301-02ML20	Medium-Large	20	1
0301-02ML28	Medium-Large	28	1
0301-02L20	Large	20	1
0301-02L28	Large	28	1



* Manufactured on request.

Endoscopic Clip Appliers

Endoscopic Clip Appliers

Key features

- available for Medium, Medium-Large and Large non-absorbable LigaV[®] clips to be used on large spectrum of tissue sizes,
- jaws activated by front handle for better stability,
- 360° shaft rotation by index finger allows one-hand operation,
- flushing port facilitates easy cleaning, thus eliminating potential contamination,
- made of high quality stainless steel,
- colour coded size matching cartridge hue,
- compatible with corresponding sizes of LigaV[®] titanium non-absorbable clips.

REF 0301-02ME

Size	M
Colour	○
Shaft diameter (mm)	5
Trocar diameter (mm)	5
Working length (cm)	33
Loading	front - loaded
Box quantity (pcs)	1



REF 0301-02MLE

Size	ML
Colour	●
Shaft diameter (mm)	10
Trocar diameter (mm)	10
Working length (cm)	33
Loading	front - loaded
Box quantity (pcs)	1




REF 0301-02LE

Size	L
Colour	●
Shaft diameter (mm)	10
Trocar diameter (mm)	10
Working length (cm)	33
Loading	front - loaded
Box quantity (pcs)	1




Endoscopic Clip Appliers



**Ergonomic and effective approach
in endoscopic procedures**



Endoscopic Clip Appliers



M



ML



L



Bariatric Endoscopic Clip Appliers - 45 cm

Extra length for bariatric surgery

Endoscopic Clip Appliers

Bariatric Endoscopic Clip Appliers

Key features

- 45 cm long shaft,
- available for Medium, Medium-Large and Large non-absorbable LigaV[®] clips to be used on large spectrum of tissue sizes,
- jaws activated by front handle for better stability,
- 360° shaft rotation by index finger allows one-hand operation,
- flushing port facilitates easy cleaning, thus eliminating potential contamination,
- made of high quality stainless steel,
- colour coded size matching cartridge hue,
- compatible with corresponding sizes of LigaV[®] titanium non-absorbable clips.

REF 0301-02MEB

Size	M
Colour	○
Shaft diameter (mm)	5
Trocar diameter (mm)	5
Working length (cm)	45
Box quantity (pcs)	1



REF 0301-02MLEB

Size	ML
Colour	●
Shaft diameter (mm)	10
Trocar diameter (mm)	10
Working length (cm)	45
Box quantity (pcs)	1



REF 0301-02LEB


Size	L
Colour	●
Shaft diameter (mm)	10
Trocar diameter (mm)	10
Working length (cm)	45
Box quantity (pcs)	1



Accessories

Accessories

REF	Product	Clip size	Pcs per box
0301-05CB	Cartridge base	S, M, ML, L	1



New International Initiative



Grena forges Talents!



Who's next?

Visit us

<http://talents.grena.co.uk>

Brentford 23.04.2012

TECHNICAL PARAMETERS OF GRENA TITANIUM AND POLYMER LIGATING CLIPS

LigaV

REF	S	M	ML	L
Clip shape in the cartridge	V shaped	V shaped	V shaped	V shaped
Opened clip length (outer) after loading into applier	3,20 mm	4,70 mm	8,00 mm	11,90 mm
Opened clip width (internal) after loading into applier	2,10 mm	2,90 mm	5,20 mm	5,40 mm
Closed clip length	3,70 mm	5,40 mm	9,00 mm	12,30 mm
Recommended tissue thickness	0,30-1,50 mm	1,00-2,50 mm	2,50-4,00 mm	3,50-7,50 mm

Vclip

REF	XS	S	M	ML	L
Clip shape in the cartridge	Chevron shaped	Chevron shaped	Chevron shaped	Chevron shaped	Chevron shaped
Opened clip length (outer) after loading into applier	2,25 mm	3,20 mm	5,30 mm	8,10 mm	11,00 mm
Opened clip width (internal) after loading into applier	1,80 mm	2,00 mm	3,30 mm	5,30 mm	7,70 mm
Closed clip length	2,65 mm	3,70 mm	6,00 mm	9,10 mm	12,50 mm
Recommended tissue thickness	0,15-0,30 mm	0,30-1,50 mm	1,00-2,50 mm	2,50-4,00 mm	3,50-7,50 mm

Click'aV

REF	M	ML	L	XL
Clip shape in the cartridge	Divergent curved	Divergent curved	Divergent curved	Divergent curved
Outer leg length	7,00 mm	9,50 mm	13,00 mm	17,80 mm
Clip leg thickness	0,83 mm	0,83 mm	1,29 mm	1,29 mm
Retention bosses diameter	1,00 mm	1,10 mm	1,20 mm	1,30 mm
Recommended tissue thickness	2 - 7 mm	3 - 10 mm	5 - 13 mm	7 - 16 mm
Hinge type	Living hinge	Living hinge	Living hinge	Living hinge
Retention bosses distance in the applier	5,50 mm	7,50 mm	10,20 mm	14,00 mm

MERIL - BONEWAX™



MERIL
BONEWAX

MERIL - BONEWAX™

MERIL BONEWAX™ achieves local hemostasis of the bone by acting as a mechanical (tamponade) barrier. It does not act biochemically and is non-absorbable.

MERIL BONEWAX™ is intended to aid in the control of bleeding in bone injuries.

Product Overview

- Structure : Slab of 2 gm
- Colour : Colourless
- Chemical composition : Natural bees wax
- Sterilization : Gamma radiation

Description	FG Code	Size	Needle Profile Shape, Type & Length	Suture Length (cm)
OPEN REDUCTION KIT 1				
BONEWAX X 2 GM	BW810			76