RADIFOCUS® TORQUE DEVICE


Accessories



The Radifocus Torque Device is a supportive device that is intended mainly for use when guiding a disposable guidewire, which is used for percutaneous angiography and treatment under fluoroscopy, to the objective site.

Product Characteristics

• Torque device compatible with guidewires 0.010" (0.25 mm) to 0.038" (0.97 mm)



Please quote above item reference codes when placing an order



Accessories



A hemostasis valve prevents blood reflux when a catheter is inserted in a vessel for angiographic or other purposes.

Product Characteristics

- Cross-cut valve
- Compatible up to 7 Fr catheter
- With side tube and 3-Way Stopcock

Item Reference

RF*VA713M

Please quote above item reference codes when placing an order

