


# eclipse & copernic

Compliant and hypercompliant balloon catheters

Designed to temporarily stop or control blood flow, to treat vasospasms, and for embolization of aneurysms with balloons.



balloon size from 4 to 6 & 3 to 5mm



With courtesy of Dr Sukalyan Purkayastha, India

aneurysm treatment

## Navigation

Hydrophilic coating on both the catheter & the balloon for a smooth progression


Compatibility with hybrid1214D guidewire which ensures support & trackability

## Compliant & hypercompliant balloons

Made from a very supple elastomer\* to easily conform to the anatomy

Rapid inflation deflation with 1/2 contrast\*

Good visibility thanks to 2RX markers providing the balloon length




## Remodeling

Both Eclipse and Copernic can be used in parallel with a coiling microcatheter in a single 6F guiding catheter (FARGOMAX6F)


# occlusion catheters with balloon

## copernic




Reference	Total length (cm)	Balloon length (mm)	Balloon diameter (mm)	Maximum volume (mL)	Guidewire recommended
SCOPERNIC10	160	10	3 to 5	0.25	Compatible with .012" guidewire (packaging includes a HYBRID1214D)
SCOPERNIC15		15		0.35	
SCOPERNIC20		20		0.4	
SCOPERNIC30		30		0.5	
COPERNIC10	160	10	3 to 5	0.25	Compatible with .012" guidewire (not included in the package)
COPERNIC15		15		0.35	
COPERNIC20		20		0.4	
COPERNIC30		30		0.5	

## eclipse



Reference	Total length (cm)	Balloon length (mm)	Balloon diameter (mm)	Maximum volume (mL)	Guidewire recommended
SECLIPSE7	160	7	4 to 6	0.2	Compatible with .012" guidewire (packaging includes a HYBRID1214D)
SECLIPSE9		9		0.25	
SECLIPSE12		12		0.3	
SECLIPSE15		15		0.4	
SECLIPSE20		20		0.5	
ECLIPSE7	160	7	4 to 6	0.2	Compatible with .012" guidewire (not included in the package)
ECLIPSE9		9		0.25	
ECLIPSE12		12		0.3	
ECLIPSE15		15		0.4	
ECLIPSE20		20		0.5	

## baltacci



Reference	Total length (cm)	Balloon length (mm)	Balloon diameter (mm)	Maximum volume (mL)	Guidewire recommended
BALTACCIB1	163	6 to 9	4 to 6	0.2	/

## magic

Reference	Total length (cm)	Balloon length (mm)	Balloon diameter (mm)	Maximum volume (mL)	Guidewire recommended
MAGICB1	155	9	6	0.2	/
MAGICB2	155	18	6	0.5	/