

PAGINA 0

**CREATA PENTRU A
PUTEA SEMNA
ELECTRONIC
DOCUMENTUL**



HEALTHCARE PROFESSIONALS
Admiral™ Xtreme PTA Catheter
 Percutaneous Transluminal Angioplasty

Indications, Safety, and Warnings

RELATED PRODUCTS

- PTA Catheters
- Cardiovascular Products

PERIPHERAL CATALOGUE

Choose from a full portfolio of peripheral products for endovascular interventions – products that help you reach your PAD patient treatment goals.

DOWNLOAD CATALOGUE

REQUEST MORE INFORMATION

+1 763-514-4000

rs.ctechsupport@medtronic.com

OVERVIEW

Treat longer PAD lesions above the knee with Admiral™ Xtreme PTA Catheter. Offered also in 250 and 300 mm lengths, with 0.035" guidewire compatibility, featuring large inflation lumen for rapid inflation and deflation.

INDICATIONS

and supra-aortic vessels. The Admiral™ Xtreme PTA Catheter is contraindicated for cases where the guidewire cannot cross the target lesion.

PRODUCT DETAILS

LONG BALLOON LENGTHS

- 250 and 300 mm balloon lengths are available for single dilatation of long lesions.
- Proprietary technology prevents banana-shaped balloon.

COAXIAL DUAL-LUMEN SHAFT DESIGN

- Strong pushability and kink resistance are combined with flexibility.
- Offers support for contralateral maneuvers.
- Large inflation lumen enable fast inflation and deflation.

LOW BALLOON PROFILE

- Shows good refolding behavior.
- Cross lesions with confidence.

HYDROPHILIC COATING

- Facilitates crossing.

FLEXITEC™ XTREME BALLOON MATERIAL

- Balloon wall thickness is reduced.
- Flexitec™ offers resistance against scratching.
- Provides high rated burst pressure.

MANUALS AND TECHNICAL GUIDES

Instructions for Use

Find this technical manual in the product labeling supplied with each device.

MODEL SPECIFICATIONS¹

TECHNICAL SPECIFICATIONS

Type
Catheter Design
Balloon Material
Balloon Coating
Balloon Marker
Shaft Diameter
Usable Shaft Lengths
Introducer Sheath Compatibility
Maximum Guidewire Compatibility

ORDERING INFORMATION

Reference Number Usable Length 80 cm	Reference Number Usable Length 130 cm	Reference Number Usable Length 150 cm	Balloon Diameter (mm)	Balloon Length (mm)
SBI 030 020 080	SBI 030 020 130		3	20
SBI 030 040 080	SBI 030 040 130		3	40
SBI 030 080 080	SBI 030 080 130		3	80
SBI 030 100 080	SBI 030 100 130		3	100
SBI 030 120 080	SBI 030 120 130		3	120
SBI 040 020 080	SBI 040 020 130	SBI 040 020 150	4	20
SBI 040 040 080	SBI 040 040 130	SBI 040 040 150	4	40
SBI 040 060 080	SBI 040 060 130	SBI 040 060 150	4	60
SBI 040 080 080	SBI 040 080 130	SBI 040 080 150	4	80
SBI 040 100 080	SBI 040 100 130	SBI 040 100 150	4	100
SBI 040 120 080	SBI 040 120 130	SBI 040 120 150	4	120
SBI 040 150 080	SBI 040 150 130	SBI 040 150 150	4	150
SBI 040 200 080	SBI 040 200 130	SBI 040 200 150	4	200
SBI 040 250 080	SBI 040 250 130		4	250
SBI 040 300 080	SBI 040 300 130		4	300
SBI 050 020 080	SBI 050 020 130	SBI 050 020 150	5	20
SBI 050 040 080	SBI 050 040 130	SBI 050 040 150	5	40
SBI 050 060 080	SBI 050 060 130	SBI 050 060 150	5	60
SBI 050 080 080	SBI 050 080 130	SBI 050 080 150	5	80
SBI 050 100 080	SBI 050 100 130	SBI 050 100 150	5	100
SBI 050 120 080	SBI 050 120 130	SBI 050 120 150	5	120
SBI 050 150 080	SBI 050 150 130	SBI 050 150 150	5	150
SBI 050 200 080	SBI 050 200 130	SBI 050 200 150	5	200
SBI 050 250 080	SBI 050 250 130		5	250
SBI 050 300 080L	SBI 050 300 130L		5	300
SBI 060 020 080	SBI 060 020 130	SBI 060 020 150	6	20
SBI 060 040 080	SBI 060 040 130	SBI 060 040 150	6	40
SBI 060 060 080	SBI 060 060 130	SBI 060 060 150	6	60
SBI 060 080 080	SBI 060 080 130	SBI 060 080 150	6	80
SBI 060 100 080	SBI 060 100 130	SBI 060 100 150	6	100
SBI 060 120 080	SBI 060 120 130	SBI 060 120 150	6	120
SBI 060 150 080	SBI 060 150 130	SBI 060 150 150	6	150
SBI 060 200 080L	SBI 060 200 130L	SBI 060 200 150	6	200
SBI 060 250 080L	SBI 060 250 130L		6	250
SBI 060 300 080L	SBI 060 300 130L		6	300
SBI 070 020 080	SBI 070 020 130	SBI 070 020 150	7	20
SBI 070 040 080	SBI 070 040 130	SBI 070 040 150	7	40
SBI 070 060 080	SBI 070 060 130	SBI 070 060 150	7	60
SBI 070 080 080	SBI 070 080 130	SBI 070 080 150	7	80
SBI 070 100 080	SBI 070 100 130	SBI 070 100 150	7	100
SBI 070 120 080L	SBI 070 120 130L	SBI 070 120 150	7	120
SBI 070 150 080L	SBI 070 150 130L	SBI 070 150 150	7	150
SBI 070 200 080L	SBI 070 200 130L	SBI 070 200 150	7	200
SBI 070 250 080	SBI 070 250 130		7	250
SBI 080 020 080	SBI 080 020 130		8	20
SBI 080 040 080	SBI 080 040 130		8	40
SBI 080 060 080	SBI 080 060 130		8	60
SBI 080 080 080	SBI 080 080 130		8	80
SBI 090 020 080	SBI 090 020 130		9	20
SBI 090 040 080	SBI 090 040 130		9	40
SBI 090 060 080	SBI 090 060 130		9	60
SBI 090 080 080	SBI 090 080 130		9	80
SBI 100 020 080	SBI 100 020 130		10	20
SBI 100 040 080	SBI 100 040 130		10	40
SBI 120 020 080	SBI 120 020 130		12	20
SBI 120 040 080	SBI 120 040 130		12	40

See product catalogue for complete, detailed product information.

This page may include information about products that may not be available in your region or country. Please consult the approved indications for use. Content on specific Medtronic products is not intended for users in markets that do not have authorisation for use.

¹ Bench Test Data on File with reference PPSR-10642317DOC at Medtronic, Inc. Test data not indicative of clinical performance.

HEALTHCARE PROFESSIONALS

- Overview
- Products
- Therapies & Procedures
- Integrated Health Solutions

PATIENTS

Overview

OUR COMPANY

- Who We Are
- Key Facts
- Leadership
- History
- Mission
- Governance
- Locations
- Careers

OUR IMPACT

- Our Impact
- Healthcare access
- Inclusion, Diversity & Equity
- Protecting our planet
- ESG Report



HEALTHCARE PROFESSIONALS
Admiral™ Xtreme PTA Catheter
 Percutaneous Transluminal Angioplasty

▲ Indications, Safety, and Warnings

RELATED PRODUCTS

- PTA Catheters
- Cardiovascular Products

PERIPHERAL CATALOGUE

Choose from a full portfolio of peripheral products for endovascular interventions – products that help you reach your PAD patient treatment goals.

DOWNLOAD CATALOGUE

REQUEST MORE INFORMATION

+1 763-514-4000

rs.ctechsupport@medtronic.com



OVERVIEW

Treat longer PAD lesions above the knee with Admiral™ Xtreme PTA Catheter. Offered also in 250 and 300 mm lengths, with 0.035" guidewire compatibility, featuring large inflation lumen for rapid inflation and deflation.

INDICATIONS

and supra-aortic vessels. The Admiral™ Xtreme PTA Catheter is contraindicated for cases where the guidewire cannot cross the target lesion.

PRODUCT DETAILS

LONG BALLOON LENGTHS

- 250 and 300 mm balloon lengths are available for single dilatation of long lesions.
- Proprietary technology prevents banana-shaped balloon.

COAXIAL DUAL-LUMEN SHAFT DESIGN

- Strong pushability and kink resistance are combined with flexibility.
- Offers support for contralateral maneuvers.
- Large inflation lumen enable fast inflation and deflation.

LOW BALLOON PROFILE

- Shows good refolding behavior.
- Cross lesions with confidence.

HYDROPHILIC COATING

- Facilitates crossing.

FLEXITEC™ XTREME BALLOON MATERIAL

- Balloon wall thickness is reduced.
- Flexitec™ offers resistance against scratching.
- Provides high rated burst pressure.

MANUALS AND TECHNICAL GUIDES

Instructions for Use

Find this technical manual in the product labeling supplied with each device.

MODEL SPECIFICATIONS¹

TECHNICAL SPECIFICATIONS

Specification
Over the wire (OTW)
Flexitec™ Xtreme
LFC Hydrophilic
2 swaged (zero profile) Platinum Iridium
5F-5.3F
80, 130, 150 cm
5 F - 6 F - 7 F
0.035"

ORDERING INFORMATION

Device Length (cm)	Reference Number Usable Length 150 cm	Balloon Diameter (mm)	Balloon Length (mm)	Recommended Introducer Sheath (F)	RBP ² (atm)
0 020 130		3	20	5	18
0 040 130		3	40	5	18
0 080 130		3	80	5	18
0 100 130		3	100	5	18
0 120 130		3	120	5	18
0 020 130	SBI 040 020 150	4	20	5	18
0 040 130	SBI 040 040 150	4	40	5	18
0 060 130	SBI 040 060 150	4	60	5	18
0 080 130	SBI 040 080 150	4	80	5	18
0 100 130	SBI 040 100 150	4	100	5	18
0 120 130	SBI 040 120 150	4	120	5	18
0 150 130	SBI 040 150 150	4	150	5	14
0 200 130	SBI 040 200 150	4	200	5	14
0 250 130		4	250	5	14
0 300 130		4	300	5	14
0 020 130	SBI 050 020 150	5	20	5	17
0 040 130	SBI 050 040 150	5	40	5	17
0 060 130	SBI 050 060 150	5	60	5	17
0 080 130	SBI 050 080 150	5	80	5	15
0 100 130	SBI 050 100 150	5	100	5	15
0 120 130	SBI 050 120 150	5	120	5	15
0 150 130	SBI 050 150 150	5	150	5	14
0 200 130	SBI 050 200 150	5	200	5	14
0 250 130		5	250	5	14
0 300 130L		5	300	6	14
0 020 130	SBI 060 020 150	6	20	5	17
0 040 130	SBI 060 040 150	6	40	5	17
0 060 130	SBI 060 060 150	6	60	5	17
0 080 130	SBI 060 080 150	6	80	5	15
0 100 130	SBI 060 100 150	6	100	5	15
0 120 130	SBI 060 120 150	6	120	5	15
0 150 130	SBI 060 150 150	6	150	5	12
0 200 130L	SBI 060 200 150	6	200	6	12
0 250 130L		6	250	6	12
0 300 130L		6	300	6	12
0 020 130	SBI 070 020 150	7	20	5	16
0 040 130	SBI 070 040 150	7	40	5	16
0 060 130	SBI 070 060 150	7	60	5	14
0 080 130	SBI 070 080 150	7	80	5	14
0 100 130	SBI 070 100 150	7	100	5	12
0 120 130L	SBI 070 120 150	7	120	6	12
0 150 130L	SBI 070 150 150	7	150	6	12
0 200 130L	SBI 070 200 150	7	200	6	12
0 250 130		7	250	6	12
0 020 130		8	20	6	14
0 040 130		8	40	6	14
0 060 130		8	60	6	11
0 080 130		8	80	6	11
0 020 130		9	20	6	14
0 040 130		9	40	6	14
0 060 130		9	60	6	11
0 080 130		9	80	6	11
0 020 130		10	20	6	11
0 040 130		10	40	6	11
0 020 130		12	20	7	11
0 040 130		12	40	7	11

See product catalogue for complete, detailed product information.

This page may include information about products that may not be available in your region or country. Please consult the approved indications for use. Content on specific Medtronic products is not intended for users in markets that do not have authorisation for use.

¹ Bench Test Data on File with reference PPSR-10642317DOC at Medtronic, Inc. Test data not indicative of clinical performance.

HEALTHCARE PROFESSIONALS

- Overview
- Products
- Therapies & Procedures
- Integrated Health Solutions

PATIENTS

Overview

OUR COMPANY

- Who We Are
- Key Facts
- Leadership
- History
- Mission
- Locations
- Governance
- Careers

OUR IMPACT

- Our Impact
- Healthcare access
- Inclusion, Diversity & Equity
- Protecting our planet
- ESG Report