

## Abbreviation

AEC	: Automatic Exposure Control
APR	: Anatomical Program
HU	: Heat Unit
IBS	: Image Brightness Stabilizer
I.I.	: Image Intensifier
kV	: kilovolt
LED	: Light Emitting Diode
mA	: milliamperes
mAs	: milliamperes second
min	: minute
ms	: millisecond
sec	: second
RIS	: Radiology Information System
CR	: Computed Radiology

## Options

The "Options" described in this manual may be provided as standard, or may be unavailable, depending on the sales format.

## Related Documents

M516-E009 High Speed Rotation Starter SA-60 Operation Manual

## System Activation

The equipment needs activation during installation.

The activation work needs authorization by Shimadzu.

\*This is available for X-ray High Voltage Generator shipped after January, 2020.