



CATALOGUE 2021

CARDIAC CATHETERS  
AND CANNULAE



# Retrograde cardioplegia cannulae



Size	Balloon diameter	Reference
14 Fr	18 mm	RCC01


## PRODUCT INFORMATION

Andocor retrograde cardioplegia cannulae are intended to be connected to the cardioplegia line in order to infuse cardioplegic solution and blood into the patient's coronary sinus during open heart procedures.

**Q** Length: 31 cm  
Qty/box: 20 Pcs

**+** Self-inflating  
Silicone balloon  
Malleable stylet

## RETROGRADE CARDIOPLEGIA CANNULA PRESSURE DROP VS FLOW



# Ostial perfusion cannulae

CORONARY ARTERY OSTIAL CANNULA, BASKET TIP WITH FLANGE



Tip size	45° angled	90° angled
10 Fr	OC1045	OC1090
12 Fr	OC1245	OC1290
14 Fr	OC1445	OC1490


## PRODUCT INFORMATION

Andocor ostial perfusion cannulae are intended for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.

**Q** Connection: Female luer lock  
 Length: 21 cm  
 Qty/box: 20 Pcs

**+** Soft basket tip  
 Different configurations  
 Stainless steel shaft

## CORONARY ARTERY OSTIAL CANNULA PRESSURE DROP VS FLOW



# Pericardial sumps



Type A



Type B



Type C

Description	Size	Tubing	Tip length	Tubing length	Reference
Type A Metal tip	18 Fr	Flexible	2 cm	30 cm	PS40M
	18 Fr	Hard	2 cm	30 cm	PS40MH
	18 Fr	Flexible	2 cm	79 cm	PS40ML
	18 Fr	Hard	2 cm	79 cm	PS40MLH
Type B Perforated pool tip	23 Fr	Flexible	3,5 cm	30 cm	PS40M1
Type C Perforated pool tip	15 Fr	Flexible	2 cm	30 cm	PS40M2
	15 Fr	Hard	2 cm	79 cm	PS40MLH2

## PRODUCT INFORMATION

Andocor suckers are used for the removal of excess fluid from the surgical field.



Connection: 1/4"  
Qty/box: 20 Pcs



Smooth design  
Weighted tip  
Stainless steel tip

## WIREWOUND TIP

Description	Size	Tubing	Tip length	Tubing length	Reference
Wirewound tip	16 Fr	Flexible	4 cm	30 cm	PS40

