



ZERONE CO., LTD.

(Shinil IT UTO Bldg., Dangjeong-dong) 8F, LS-ro 13, Gunpo-si, Gyeonggi-do, 15843, Republic of Korea
Tel: 82-31-689-5610, Fax: 82-31-453-5444, E-mail: overseas@01zeus.com, Website: www.01zeus.com

Specification Statement

Subject: Clarification Regarding Self-Test Functionality on ZEUS VISION

Ref.No. ZOEC260514-A01

To Whom It May Concern

This letter serves to provide clarification regarding the self-test functionality of our ZEUS VISION. We understand there may be questions concerning the absence of a visual "Self Test" display upon device boot-up.

Please be advised that the loading bar during power-on the ZEUS VISION indicates that a self-test is in progress.

We wish to assure you that the self-test function is fully operational and executes internally each time the device is powered on. This internal process rigorously assesses the device's critical components and ensures its proper functionality.

We are committed to providing reliable and high-quality products. The absence of a direct "Self Test" display in no way compromises the integrity or effectiveness of this essential diagnostic feature.

ZERONE Co., Ltd.

ZERONE CO., LTD

YANG BAE, KIM / PRESIDENT

Yang Bae, Kim / President