

Pacific™ Plus

PTA Catheter OTW 0.018"

Versatility for everyday and beyond


- Short and long balloons, small and large diameters
- Provides versatility to treat both ATK and BTK lesions
- Proprietary balloon technology allows for great crossability*
- Improved shaft design for faster deflation*

Technical specifications

Catheter design	Over-the-wire (OTW), coaxial shaft
Balloon coating	Hydrophilic
Balloon marker	2 or 4 swaged, platinum-iridium
Usable shaft lengths	90, 130, 150, and 180 cm
Introducer sheath compatibility	4.0 and 5.0 F
Guidewire compatibility	0.018"



Tapered tip



New dual middle marker bands provide enhanced visibility for precise placement[‡]



*Six samples of each brand tested, 5 mm x 40 mm. Medtronic data on file. CP230TP01R Bench test results may not be indicative of clinical performance.
[‡]Accurate with 200-300 mm balloon lengths

Pacific™ Plus

PTA Catheter OTW 0.018"

Order information

Product catalogue number				Balloon diameter (mm)	Balloon length (mm)	Recommended introducer sheath (F)	Nominal pressure (atm)	RBP (atm)
Usable length 90 cm	Usable length 130 cm	Usable length 150 cm	Usable length 180 cm					
PCE 020 020 090	PCE 020 020 130	-	-	2.00	20	4	8	22
PCE 020 040 090	PCE 020 040 130	-	PCE 020 040 180	2.00	40	4	8	22
PCE 020 060 090	PCE 020 060 130	-	-	2.00	60	4	8	22
PCE 020 080 090	PCE 020 080 130	-	PCE 020 080 180	2.00	80	4	8	22
PCE 020 120 090	PCE 020 120 130	-	PCE 020 120 180	2.00	120	4	8	22
PCE 020 150 090	PCE 020 150 130	-	PCE 020 150 180	2.00	150	4	8	22
PCE 025 020 090	PCE 025 020 130	-	-	2.50	20	4	8	16
PCE 025 040 090	PCE 025 040 130	-	PCE 025 040 180	2.50	40	4	8	16
PCE 025 060 090	PCE 025 060 130	-	-	2.50	60	4	8	16
PCE 025 080 090	PCE 025 080 130	-	PCE 025 080 180	2.50	80	4	8	16
PCE 025 120 090	PCE 025 120 130	-	PCE 025 120 180	2.50	120	4	8	16
PCE 025 150 090	PCE 025 150 130	-	PCE 025 150 180	2.50	150	4	8	16
PCE 030 020 090	PCE 030 020 130	-	-	3.00	20	4	8	16
PCE 030 040 090	PCE 030 040 130	-	PCE 030 040 180	3.00	40	4	8	16
PCE 030 060 090	PCE 030 060 130	-	-	3.00	60	4	8	16
PCE 030 080 090	PCE 030 080 130	-	PCE 030 080 180	3.00	80	4	8	16
PCE 030 120 090	PCE 030 120 130	-	PCE 030 120 180	3.00	120	4	8	16
PCE 030 150 090	PCE 030 150 130	-	PCE 030 150 180	3.00	150	4	8	16
PCE 035 020 090	PCE 035 020 130	-	-	3.50	20	4	8	16
PCE 035 040 090	PCE 035 040 130	-	-	3.50	40	4	8	16
PCE 035 060 090	PCE 035 060 130	-	-	3.50	60	4	8	16
PCE 035 080 090	PCE 035 080 130	-	-	3.50	80	4	8	16
PCE 035 120 090	PCE 035 120 130	-	-	3.50	120	4	8	16
PCE 035 150 090	PCE 035 150 130	-	-	3.50	150	4	8	16

Order information

Product catalogue number				Balloon diameter (mm)	Balloon length (mm)	Recommended introducer sheath (F)	Nominal pressure (atm)	RBP (atm)
Usable length 90 cm	Usable length 130 cm	Usable length 150 cm	Usable length 180 cm					
PCE040020090	-	PCE040020150	-	4.00	20	4	8	14
PCE040040090	-	PCE040040150	-	4.00	40	4	8	14
PCE040060090	-	PCE040060150	-	4.00	60	4	8	14
PCE040080090	-	PCE040080150	-	4.00	80	4	8	14
PCE040100090	-	PCE040100150	-	4.00	100	4	8	14
PCE040120090	-	PCE040120150	-	4.00	120	4	8	14
PCE040150090	-	PCE040150150	-	4.00	150	4	8	14
PCE040200090	-	PCE040200150	-	4.00	200	4	8	14
PCE040250090	-	PCE040250150	-	4.00	250	4	8	14
PCE040300090	-	PCE040300150	-	4.00	300	4	8	14
PCE050020090	-	PCE050020150	-	5.00	20	4	8	14
PCE050040090	-	PCE050040150	-	5.00	40	4	8	14
PCE050060090	-	PCE050060150	-	5.00	60	4	8	14
PCE050080090	-	PCE050080150	-	5.00	80	4	8	14
PCE050100090	-	PCE050100150	-	5.00	100	4	8	14
PCE050120090	-	PCE050120150	-	5.00	120	4	8	14
PCE050150090	-	PCE050150150	-	5.00	150	4	8	14
PCE050200090	-	PCE050200150	-	5.00	200	4	8	14
PCE050250090	-	PCE050250150	-	5.00	250	5	8	14
PCE050300090	-	PCE050300150	-	5.00	300	5	8	14
PCE060020090	-	PCE060020150	-	6.00	20	4	8	14
PCE060040090	-	PCE060040150	-	6.00	40	4	8	14
PCE060060090	-	PCE060060150	-	6.00	60	4	8	14
PCE060080090	-	PCE060080150	-	6.00	80	4	8	14

Pacific™ Plus

PTA Catheter OTW 0.018"

Order information

Product catalogue number				Balloon diameter (mm)	Balloon length (mm)	Recommended introducer sheath (F)	Nominal pressure (atm)	RBP (atm)
Usable length 90 cm	Usable length 130 cm	Usable length 150 cm	Usable length 180 cm					
PCE060100090	-	PCE060100150	-	6.00	100	4	8	14
PCE060120090	-	PCE060120150	-	6.00	120	4	8	14
PCE060150090	-	PCE060150150	-	6.00	150	5	8	12
PCE060200090	-	PCE060200150	-	6.00	200	5	8	12
PCE060250090	-	PCE060250150	-	6.00	250	5	8	12
PCE060300090	-	PCE060300150	-	6.00	300	5	8	12
PCE070020090	-	PCE070020150	-	7.00	20	4	8	12
PCE070040090	-	PCE070040150	-	7.00	40	4	8	12
PCE070060090	-	PCE070060150	-	7.00	60	4	8	12
PCE070080090	-	PCE070080150	-	7.00	80	5	8	12
PCE070100090	-	PCE070100150	-	7.00	100	5	8	12
PCE070120090	-	PCE070120150	-	7.00	120	5	8	12
PCE070150090	-	PCE070150150	-	7.00	150	5	8	12
PCE070200090	-	PCE070200150	-	7.00	200	5	8	12
PCE070250090	-	PCE070250150	-	7.00	250	5	8	12