


# NUCLEOSHELL® RP 18

The separation of 13  $\beta$ -lactam antibiotics illustrates how time of analysis can be shortened to a fractional part by using core-shell particles all without loss of resolution at moderate back pressure.



## Ordering information

NUCLEOSHELL® RP 18				
Analytical EC columns NUCLEOSHELL® RP 18 (pack of 1)				
Length (mm)	ID (mm)	Particle size ( $\mu$ m)	REF	Guard columns*
150	4.6	2.7	763136.46	763138.30
150	4	2.7	763136.40	763138.30
150	3	2.7	763136.30	763138.30
150	2	2.7	763136.20	763138.20
100	4.6	2.7	763134.46	763138.30
100	4	2.7	763134.40	763138.30
100	3	2.7	763134.30	763138.30
100	2	2.7	763134.20	763138.20
50	3	2.7	763132.30	763138.30
50	2	2.7	763132.20	763138.20
250	4.6	5	763157.46	763158.30
250	4	5	763157.40	763158.30
250	3	5	763157.30	763158.30
150	4.6	5	763156.46	763158.30
150	4	5	763156.40	763158.30
150	3	5	763156.30	763158.30
100	4.6	5	763154.46	763158.30
100	3	5	763154.30	763158.30
100	2	5	763154.20	763158.20
50	4.6	5	763152.46	763158.30

\* Pack of 3, EC guard columns require column protection system REF 718966. For more information, see page 90.

For more products  
and information

Or visit [www.mn-net.com](http://www.mn-net.com)

