



citek
Bone Wax

Description

- Mixture of beeswax, paraffin, isopropyl palmitate.
- Single use
- Sterile: Gamma Radiation (R)
- Shelf life: 5 years

Advantages:

- Soft and easy to shape
- Practical
- Controls bleeding
- Achieve local hemostasis of bone by acting as a mechanical (tamponade) barrier
- Compatible with body tissues
- Does not act biochemically

Indications:

D-tek Mesh is indicated for use in traumatology and orthopedics, sternum and ribs surgery, ion maxillofacial surgery, general, plastic surgery and neurosurgery (trepanation), dental, oral and jaw surgery.

Properties:

- Individually wrapped, 2,5grams
- Double peel-packed