CARDIOVASCULAR





FEATURES


- Rapid & easy to operate
- · Visually read results
- · Instrument-free
- · High sensitivity & specificity
- Accelerate clinical diagnosis

SPECIFICATIONS

- Shelf life: 24 months
- · Storage: 2 30 °C
- · Internal control included

Heart and blood vessel disorders are collectively referred to as cardiovascular disease (CVD). The development of fatty deposits inside the arteries (atherosclerosis) and a higher risk of blood clots are typically connected with it. Damage to the arteries in several organs, including the kidneys, eyes, heart, and brain, can also be linked to CVD. A healthy lifestyle can frequently greatly reduce the risk of developing it. The four main CVD diseases are coronary heart disease, strokes and transient ischemic attacks (TIAs), peripheral arterial disease, and aortic disease. CVD is one of the leading causes of death and disability worldwide.

SIMPLE INTERPRETATION





CARDIOVASCULAR

NAME	SPECIMEN	ITEM CODE	FORMAT	CUT-OFF	иом	REGISTRATION
CK-MB Test	S/P	GDCKM-302a	Cassette	5 ng/mL	25 T/Kit	CE marked
	WB/S/P	GDCKM-402a	Cassette	5 ng/mL	25 T/Kit	CE marked
CRP C-Reactive Protein Semi-Quantitative Test	WB/S/P	GDCRP-402a	Cassette	1~3~10 ng/mL	25 T/Kit	CE marked
		GDCRP-T402a	Cassette	10~40~80 ng/mL	25 T/Kit	CE marked
D-Dimer Test	WB/S/P	GDDDI-402b	Cassette	500 ng/mL	25 T/Kit	CE marked
Myoglobin Test	WB/S/P	GDMYO-402a	Cassette	50 ng/mL	25 T/Kit	CE marked
Procalcitonin Test	WB/S/P	GDPCT-T402a	Cassette	0.5~2~10 ng/mL	25 T/Kit	CE marked
Troponin I Test	S/P	GDTRO-302a	Cassette	0.5 ng/mL	25 T/Kit	CE marked
	WB/S/P	GDTRO-402a	Cassette	0.5 ng/mL	25 T/Kit	CE marked
		GDTRO-402b	Cassette	0.5 ng/mL	25 T/Kit	CE marked
Cardiac Myoglobin/ CK-MB/cTnl Combo Test	S/P	GDCAR-335a	Cassette	50/5/0.5 ng/mL	25 T/Kit	CE marked
	WB/S/P	GDCAR-435a	Cassette	50/5/0.5 ng/mL	25 T/Kit	CE marked
		GDCAR-W435a	Cassette	50/5/0.5 ng/mL	20 T/Kit	CE marked