



Sterile First Aid  
**BURN**  **GEL**  
and **DRESSING**



Fire



Thermal &  
Scalding



Chemical



Radiation



Electrical



Sunburn

**EMERGENCY BURN CARE**

**CE**  
0425

Manufactured in a facility with ISO 13485:2016 QMS.



**GELS**

**DRESSINGS**

**RETAIL**

**KITS**

Sterile First Aid  
**BURN GEL**  
and **DRESSING**





## ABOUT US

Since 2000, we have been reaching to more than 128 countries of the world from Turkey by manufacturing and exporting **Medical Products, Disinfectant and Antiseptic Products, Non-prescription Products (OTC Products), Sexual Health Products, Dermocosmetic and Personal Care Products**. We own globally known brands; **Konix, KLY, Konicare, Mediseptica, Gynosed, Gynotal and Joydrops**. At the same time we have a huge experience in Private Label and **OEM production**. Our experience of Private Label and OEM production is more than 20 years. We are a leading manufacturer in the medical sector, representing our country successfully in the global arena.

We have achieved number 2 at 'Turkey's Exporters' with our export performance covering 70% of our production capacity in the medical products. We have **GMP (Good Manufacturing Practice) Certificate, indicating that we have international quality standards in all processes of the product, ISO 13485: 2016 certificate certifying that we produce Class I, CE certified Class IIA, Class IIB and Class III products and Halal Certificate** registering that we manufacture according to Islamic rules. In the medical sector, we are one of the leading manufacturers in Europe with production standards at Class III level.

We make world-class production and Research and development (R&D) activities in our 18,000 square meter closed area in Istanbul, which is operated in accordance with ISO 9001 standards. **Our R&D center, approved by Turkish Ministry of Science and Technology**, has been one of the pioneers of the sector by developing many patented products. Changing world trends are examined thoroughly in our R&D center and our laboratories equipped with modern devices. We continue to lead the sector by introducing our innovative products to the market with the implementation of the international production standards.



More than 200  
Employees



20,000 m2  
Production Area



Approved R&D  
Center



More Than 200  
Products



Capacity of 130 Million  
Units per Year



Export to More  
Than 128 Countries

**Konix® Sterile First Aid Gel** is an antiseptic first aid product group designed for first aid in all types of burn injuries, which eliminates inflammation on the skin, prevents the burn from spreading further on the skin and accelerates the healing process of the burn.

**Konix® Sterile First Aid Gel** cools the burn rapidly in order to prevent further tissue damage and assist in the healing process, and provides the patient with a rapid and continuous reduction in sensitivity.

**Konix® Sterile First Aid Gel** cools the burned area four times faster than water and instantly prevents the burn from deepening, effectively relieves pain, moisturizes and soothes the skin, helps prevent infection and accelerates the healing process.

### TYPES OF BURNS USED

- Open Fire Burns
- Sunburns
- Hot Water Burns  
Steam Burns
- Chemical Burns  
Gas Burns
- Electrical Burns
- Radioactive Burns

- It helps reduce sensitivity and pain in the burn area in a short time after application.
- It cools and moisturizes the burn area.
- It prevents the deepening of the burn by stopping the heat from passing into the lower layers of the skin.
- It completely covers the burn and prevents infection caused by the external environment.
- It minimizes the risk of infection, is non-toxic and non-irritating.
- It protects the burn area against abrasion, friction and drying.
- Konix® Sterile First Aid Gel and Burn Dressing can be easily dissolved with normal water and leaves no residue after the wounded area is washed.
- It is safe for both children and adults.

## Konix® Sterile First Aid Gel



### MELALEUCA (TEA TREE) OIL

Tea tree oil is an oil with strong antibacterial and antifungal properties. It helps reduce the risk of infection in the burn area. It has a soothing effect on the skin and reduces inflammation.



### VISCOUS HYDROGEL TECHNOLOGY

Combines a viscous hydrogel that is perfect for cooling the wound and protecting it from infection.




### ALOE VERA

Aloe Vera is one of the most commonly used natural ingredients in the treatment of burns. It moisturizes the skin, has anti-inflammatory properties and promotes rapid healing. Aloe vera also relieves skin pain and reduces the discomfort of the burn with its cooling effect.

**%7,5**  
HYDROGEL

**%1**  
TEA TREE  
OIL

**%0,5**  
ALOE  
VERA



**% 91**  
WATER

## Konix® Sterile First Aid Burn Dressing

Konix® Sterile First Aid Gel and Burn Dressing consists of sterile gel-impregnated burn dressings and gel bottles and bags used in emergency first aid and burn treatment.

- Sterile
- Minimizes trauma and skin damage
- Moisturizes, cools and soothes the burn
- Non-stick
- Non-toxic
- Non-irritating
- Safe for use on children
- Perfect for travel or home use

# Konix® Sterile First Aid Gel

## Areas of Use

- Sunburn
- Minor burns and scalds
- Hard-to-reach areas: e.g. behind the ears, wrinkles in the face and groin area
- Remoistening the Burnshield dressing

## HOW TO APPLY?

1

Apply Konix® Sterile First Aid Gel liberally to the burn area.

2

Remember to extend the application beyond the burn area.

3

Cover with tape or a non-stick bandage.

### Product Code

### Product Name

### Packaging Type

Product Code	Product Name	Packaging Type
Y0000.133.0006	KONIX 50 ML STERILE FIRST AID BURN GEL	Plastic Bottle
Y0000.133.0007	KONIX 100 ML STERILE FIRST AID BURN GEL	Plastic Bottle
Y0000.133.0008	KONIX 118 ML STERILE FIRST AID BURN GEL	Plastic Bottle
Y0000.133.0023	KONIX 59 ML STERILE FIRST AID BURN GEL	Plastic Bottle
Y0000.133.0028	KONIX 125 ML STERILE FIRST AID BURN SPRAY	Spray Bottle
Y0000.133.0009	KONIX 3,5 GR STERILE FIRST AID BURN GEL	Sachet
Y0000.133.0010	KONIX 5 GR STERILE FIRST AID BURN GEL	Sachet
Y0000.133.0030	KONIX 3,5 GR STERILE FIRST AID BURN GEL - 6 PACK	Sachet
Y0000.133.0031	KONIX 3,5 GR STERILE FIRST AID BURN GEL - 25 PACK	Sachet
Y0000.133.0011	KONIX 10 GR STERILE FIRST AID BURN GEL	Sachet
Y0000.133.0012	KONIX 20 GR STERILE FIRST AID BURN GEL	Sachet



# Konix®Sterile First Aid Burn Dressing

## Areas of Use and Benefits

Burn dressing is a form of burn gel applied to the entire surface of a special non-sticky, fiber-free sponge and allowed to absorb. The product covers the burn surface naturally.

- Can be used in all types of burns
- Optimal heat distribution
- Natural wound dressing
- Protects and moistens the burn area
- Adapts to the wound area
- Gel retention
- Limits exposure of bacteria to the wound

## HOW TO APPLY?



Clean the burn area and leave it exposed.



Place Konix® Sterile Burn Dressing directly on the burn surface.



Secure the dressing if necessary and keep it in your first aid kit.

\*See a healthcare provider for further treatment.

### Product Code

### Product Name

### Barcode

Product Code	Product Name	Barcode
Y0000.133.0015	KONIX 10X10 STERILE FIRST AID BURN DRESSING	8681715349169
Y0000.133.0016	KONIX 20X20 STERILE FIRST AID BURN DRESSING	8681715349176
Y0000.133.0017	KONIX 30X30 STERILE FIRST AID BURN DRESSING	8681715349183
Y0000.133.0018	KONIX 10X40 STERILE FIRST AID BURN DRESSING	8681715349190
Y0000.133.0019	KONIX 20X45 STERILE FIRST AID BURN DRESSING (FACE MASK)	8681715349206
Y0000.133.0020	KONIX 30X40 STERILE FIRST AID BURN DRESSING	8681715349213
Y0000.133.0021	KONIX 40X40 STERILE FIRST AID BURN DRESSING	8681715349220
Y0000.133.0022	KONIX 40X60 STERILE FIRST AID BURN DRESSING	8681715349237
Y0000.133.0025	KONIX 40X60 STERILE FIRST AID BURN DRESSING (FACE MASK)	8683229602068



## Konix® Sterile Burn Dressing Retail Sales

### Product Code

### Product Name

### Barcode

Product Code	Product Name	Barcode
Y0000.133.0013	KONIX 5X5 STERILE FIRST AID BURN DRESSING	8681715349145
Y0000.133.0014	KONIX 5X15 STERILE FIRST AID BURN DRESSING	8681715349152
Y0000.133.0015	KONIX 10X10 STERILE FIRST AID BURN DRESSING	8681715349169
Y0000.133.0030	KONIX 3,5 GR STERILE FIRST AID BURN GEL - 6 PACK	8681715349107
Y0000.133.0031	KONIX 3,5 GR STERILE FIRST AID BURN GEL - 25 PACK	868171534910



## Konix® Sterile First Aid Gel

### İçerik:

- x2 Konix Burn Gel Dressing 100 mm x 100 mm
- x1 Konix Burn Gel 118 gr
- x6 Konix Burn Gel Sachet 3,5 ml
- x1 Elastic Bandage 8 cm\*120 cm
- x1 Cohesive Bandage 7,6 cm\*4,5 m
- x1 Medical Cloth Flaster 2,5 cm\*500 cm
- x1 Pair of Gloves
- x1 First Aid Scissors

**\*14 Item**

### İçerik:

- x1 Konix Burngel Dressing 30 mm x 30 mm
- x2 Konix Burngel Dressing 200 mm x 200 mm
- x4 Konix Burngel Dressing 100 mm x 100 mm
- x2 Konix Burngel Hydrogel 118 gr
- x6 Konix Burn Gel Hydrogel Sachet 3,5 ml
- x2 Elastic Bandage 8 cm\*120 cm
- x1 Medical Cloth Flaster 2,5 cm\*500 cm
- x1 Cohesive Bandage 7,6 cm\*4,5 m
- x4 Pair of Gloves
- x1 First Aid Scissors

**\*24 Item**



Product Code: Y0000.133.0032



Product Code: Y0000.133.0033

### Contraindications

Konix® sterile First Aid Gel and Burn Dressing is contraindicated in the situations described below; It should not be used by people who are sensitive and allergic to the ingredients, especially Melaleuca (Tea Tree) oil.

- Gel retention
- Limits exposure of bacteria to the wound

### Side Effect

- In very rare cases, local sensitivity reactions such as dryness or burning may occur.
- Propylene glycol can cause skin inflammation.
- Melaleuca (Tea Tree) oil can cause allergic reactions.

### Precautions

- Konix® Sterile First Aid Gel and Burn Dressing is suitable for children, adults and infants over 1 month old.
- It should be used under physician control in pregnant and breastfeeding mothers.
- The product is for human use only. It should not be used for animals.
- The products are sterile and single-use.
- Keep out of the sight and reach of children.
- Do not use after the expiry date indicated on the packaging and bottle.
- If an unexpected effect is seen or irritation occurs, stop using it and consult your doctor.
- İçinde bulunan maddelerden herhangi birine karşı aşırı hassasiyet olduğu zaman kullanmayınız.
- Do not use when there is extreme sensitivity to any of the ingredients. Is sterilized with radiation after the packaging.
- In case the package damaged, do not use since it will lose the sterility.