




HONJIN



Chimb	≤3%	
Open Force	≤35%	
Sterilization Force	Sterilization about 30min under 121±0.2°C, not change color, and will not crack	
Fastness of aluminium coat	Brush the aluminum Surface by 80% Ethanol about 30 seconds, then brush the surface by 70% isopropanol about 30 seconds, the aluminum coat not appear abrasion	
Material	Aluminium-plastics combination caps	
Size	Inner Diameter (B)	14.2±0.1 mm
	Thickness (F)	0.20±0.01 mm
	Inner Height (D)	6.2±0.1 mm
Produce basis: YBB00082005-2015		

Shanxi Hongjin Medical Glass Co.,Ltd			Title: Aluminum-pastic combination caps (flip-off caps) ZB13		
		UNIT:	GEN. TOLERANCE:		
DRAWING BY:	DEPARTMANTION	NAME - SURNAME	SIGN	DATE	SIZE: A4 REV:0
Hou Mingzhen	FACTORY MANAGER	Yang Wen			DWG NO: HJ-ZB13
Date:	QA / RA	Pei Zhuo			FILE:
2022/08/25					SCALE:1/1 PAGE:1