




SPECIFICATION

DISPOSABLE GRASPING FORCEPS

Item No	O. D. (mm)	Working Length (mm)	Coating	Prong/Jaw Shape
PTG-2316-(Z3/Z4/Z5)	2.3	1600/1800/2300	/	3/4/5 prongs
PTG-18(16/18/23)C	1.8	1600/1800/2300	No	Rat Tooth With Alligator Jaws
PTG-23(16/18/23)S	2.3	1600/1800/2300	No	Rat Tooth Type
PTG-23(16/18/23)SB	2.3	1600/1800/2300	Yes	Rat Tooth Type
PTG-23(16/18/23)E	2.3	1600/1800/2300	No	Alligator Type
PTG-23(16/18/23)EB	2.3	1600/1800/2300	Yes	Alligator Type
PTG-23(16/18/23)C	2.3	1600/1800/2300	No	Rat Tooth With Alligator Jaws
PTG-23(16/18/23)CE	2.3	1600/1800/2300	Yes	Rat Tooth With Alligator Jaws
PTG-23(16/18/23)T	2.3	1600/1800/2300	No	Pelican Type
PTG-23(16/18/23)TB	2,3	1600/1800/2300	Yes	Pelican Type

REUSABLE GRASPING FORCEPS

Item No.	O. D. (mm)	Working Length (mm)	Prong/Jaw Shape
PTRG-18(16/18/23)C	1.8	1600/1800/2300	Rat Tooth With Alligator Jaws
PTRG-23(16/18/23)S	2.3	1600/1800/2300	Rat Tooth Type
PTRG-23(16/18/23)E	2.3	1600/1800/2300	Alligator Type
PTRG-23(16/18/23)C	2.3	1600/1800/2300	Rat Tooth With Alligator Jaws
PTRG-23(16/18/23)T	2.3	1600/1800/2300	Pelican Type

This device is used for grasping the cut polyps,lesion tissue or mass foreign body.

RETRIEVAL NET

Item No.	O. D. (mm)	Working Length (mm)	Opening Diameter (mm)
PTG-C-(20/25/30/35/40)-26(160/180/230)	2.6	1600/1800/2300	20/25/30/35/40

This device can accurately grasp the foreign body firmly with a net design.

DISPOSABLE STONE EXTRACTION BASKETS

Item No	O. D (mm)	Working Length (mm)	Basket Width (mm)	Basket Shape	Endoscope Channel (mm)
PTQ-18(07/10/18)-(20/25/30)-(A/B/C)	1.8	700/1000/1800	20/25/30	Diamond(A) Oval橢(B) Spiral(C)	2.0
PTQ-23(16/18/23)-(20/25/30)-(A/B/C)	2.3	1600/1800/2300	20/25/30	Diamond(A) Oval(B) Spiral(C)	2.8

This device is used to grasp the stones,polyps after removal of polyps or foreign bodies under endoscope.

Rotatable handle can make the basket a 360 rotation,which effectively captures small stones and crushed stones.

Medical-use stainless steel tightens the basket well,which makes it uneasy to fall off after taking.