

fixLOCK Plate, 3.5 mm

Item Code: 845

Indications

Temporary internal fixation and stabilization of radius, ulna, fibula, small sized humerus shaft fractures

- Periprosthetic fractures
- Comminuted fractures
- Fractures in osteopenic bone
- Non-unions
- Mal-unions



SMALL FRAGMENT



Features

- Used over the tensile surface of bones to work on tension band principle
- Uniform stiffness of all segments, prolonging fatigue life of the implant
- Smaller area of post-fixation avascularity
- The combination of conventional and locking screws offers optimum fixation regardless of bone density
- Dynamic compression unit offers 2 mm of axial compression
- Available in both Titanium & Stainless Steel

Sizes Available

No. of Holes	Code	
	S.S.	Titanium
5	845.05	845.05T
6	845.06	845.06T
7	845.07	845.07T
8	845.08	845.08T
9	845.09	845.09T
10	845.10	845.10T
11	845.11	845.11T
14	845.14	845.14T