

Imported by  
**Beybi®**

**STERIL BISTURİ UCU  
STERILE SURGICAL BLADES**

**20**



**KARBON ÇELİK 100 ADET  
CARBON STEEL 100 PCS**

**CE**  
0197



Tek Kullanımlık





Üretici / Manufactured by : Huaian Angel Medical Instruments Co., Ltd  
Huaian Angel Medical Instruments Co., Ltd  
19 East Zhuhai Road, Huaian, Jiangsu 223001, P.R.C.

EC REP

Avrupa Temsilcisi / European Representative :  
Shanghai International Holding Corp. GmbH (Europe)  
Add: Eiffestrasse 80, 20537 Hamburg , Germany  
Tel: +49-40-2513175 Fax: +49-40-255726

İthalatçı / Imported by: **BEYBİ PLASTİK FABRİKASI SAN. A.Ş.**

Yukarı Dudullu Organize Sanayi Bölgesi 1.Cad. No:27 Ümraniye, İstanbul TURKEY

Tel : +90 216 420 11 01 (6 Lines) Made in P.R.C.



# serra's

Classic

## YATAK KORUYUCU ÖRTÜ UNDERPAD


SERRAS Yatak Koruyucu Örtü  
60 x 90 cm 30'lu

- \* Yumuşak üst yüzey,  
Soft top surface
- \* Sıvının sızmasını önleyen kenarlar,  
Liquid-leaking edges
- \* Sıvının yatağa geçmesini önleyen, koruyucu dış yüzey  
Special outer surface for protecting your mattress
- \* Yırtılmalara karşı dayanıklı, esneyen dış yapı  
Tear resistant, stretching outer structure

**30** ADET  
PIECES  
عدد  
60x90 cm











# THREE-WAY STOPCOCK

- E.O. Sterilized    ISO:13485 CERTIFIED
- Pyrogen free
- Do not use if package is damaged or opened
- Discard after single use    Made In China

LOT:2022

MFG:2022

EXP:2027



CCCSA

ORAL  
NASAL

ID 8.0

OD 10.7

18 cm

20

22

24

26

28

30

32

34

36

38

40

42

44

46

48

50

52


54

56

58

60

62





# CE

# ENDOTRACHEAL TUBE



### Directions:

- Prior to intubation, deflate the cuff completely.
- After intubation, inflate the cuff using the minimum volume of air required to provide an effective seal.
- Immediately after cuff inflation, auscultate both lung fields. If breath sounds diminished over one lung field or absent over one or both fields, adjust the tube as required.
- Endotracheal tube placement should be confirmed by viewing the position of the tube tip with a chest radiograph.

### Warning:

- Single use only.
- Sterile if package is unopened or undamaged.
- Do not resterilize.
- Do not expose to temperatures above 49 °C.

Type: cuffed  
I.D.: 8.0 mm  
O.D.: 10.7 mm  
Lot No.: 20210510  
Mfg. Date: 10/05/2021  
Exp. Date: 09/05/2026





# 100% SILICONE FOLEY CATHETER



STERILE EO

AVOID FREEZING AND EXCESSIVE HEAT

**DIRECTIONS:** After aseptic catheterization, inflate the retention balloon by inserting the syringe tip into the valve and introduce the necessary amount of sterile water.

**WARNING:** Do not use ointments or lubricants having a petroleum base. They will damage the catheter and cause the balloon to burst.

Size: 3-way 20Fr (6.7mm) 30ml

Lot No.: 20230210

Mfg.Date: 20230210

Exp.Date: 20280209

CE 0123

Sterile if package is unopened and undamaged.