


## COM 30

### Combined dressing fabric

A three-layer combined covering consisting of a polyester mesh, PUR foam and polyamide knitted fabric.

#### INDICATIONS:

COM 30 is used as a temporary covering on wound surfaces after necrectomies in case of third-degree burns, infected third-degree burns, covering of donor sites after dermal grafts have been taken, for treatment of decubitus ulcer, but also on other skin and subcutaneous defects, and calf-ulcers. COM 30 shall be applied with the mesh (textured layer) to the defect.



#### The fabric

- is easy to remove
- is sterilized by irradiation
- is supplied in unit packaging
- expiration period is 5 years

## Sizes and catalogue codes:

### COM 30 – Combined dressing fabric

| Size [mm] | Catalogue code |
|-----------|----------------|
| 100 x 75  | 140-1075       |
| 150 x 100 | 140-1510       |
| 200 x 150 | 140-2015       |
| 300 x 200 | 140-3020       |
| 300 x 400 | 140-3040       |

*Other dimensions are produced on the basis of customer specification.*

### VUP Medical, a.s.

Sujanovo nam. 3 | 602 00 Brno

Czech Republic | [www.vup.cz](http://www.vup.cz)

Phone: +420 543 160 444

e-mail: [sales@vup.cz](mailto:sales@vup.cz)