

Part No	OEM No	Description		Connector	Patient group
EKG-10012A	21714 401	Suction Electrode use with DIN 3.0 connector; ball diameter: 26mm; metal electrode diameter: 21 mm;			adult
EKG-10012B		Suction Electrode use with banana 4.0 connector; ball diameter: 26mm; metal electrode diameter: 21 mm;			adult size
EKG-10012C	compatible with Philips 40490 E	suction ECG electrode use with both DIN 3.0 and banana 4.0 connector; ball diameter: 26mm; metal electrode diameter: 21 mm;			adult size
EKG-10012D	20298 92-001; Philips 40418 A	Multi-Use Suction Electrode use for snap/banana/DIN/clip connectors; ball diameter: 26mm; metal electrode diameter: 21 mm; AgCl surfacing;			adult size
EKG-10012E		Suction Electrode use with DIN 3.0 connector; ball diameter: 23mm; metal electrode diameter: 17mm;			pediatric size
EKG-10012F		Suction Electrode use with banana 4.0 connector; ball diameter: 23mm; metal electrode diameter: 17mm;			pediatric size
EKG-10012G		suction ECG electrode use with both DIN 3.0 and banana 4.0 connector; ball diameter: 23mm; metal electrode diameter: 17mm;			pediatric size
EKG-10012H		Multi-Use Suction Electrode use for snap/banana/DIN/clip connectors; ball diameter: 23mm; metal electrode diameter: 17mm; AgCl surfacing;			pediatric size

Limb ECG Electrode

Part No	OEM No	Description		Connector	Compatibility
---------	--------	-------------	--	-----------	---------------