

# Carotid WALLSTENT™ Monorail™

Carotid Stenting Device

Order Number	Size of Carotid WALLSTENT™ implanted in Vessel					Compatibility: Minimum ID			
	Fully Open Diameter (mm)	Vessel Diam. (mm)	Stent Length (mm)	Vessel Diam. (mm)	Stent Length (mm)	System OD (F / mm)	Working Length (cm)	Guiding Catheter (F / in / mm)	Guiding Sheath (F / in / mm)
H965SCH647010	6	5	30	4	36	5.0 / 1.75	135	7 / 0.073 / 1.85	5 / 0.073 / 1.85
H965SCH647070	8	7	30	6	36	5.0 / 1.75	135	7 / 0.073 / 1.85	5 / 0.073 / 1.85
H965SCH647080	8	7	40	6	48	5.0 / 1.75	135	7 / 0.073 / 1.85	5 / 0.073 / 1.85
H965SCH647090	8	7	50	6	62	5.0 / 1.75	135	7 / 0.073 / 1.85	5 / 0.073 / 1.85
H965SCH647120	10	9	30	8	36	5.9 / 2.1	135	8 / 0.086 / 2.18	6 / 0.086 / 2.18
H965SCH647130	10	9	40	8	49	5.9 / 2.1	135	8 / 0.086 / 2.18	6 / 0.086 / 2.18
H965SCH647140	10	9	50	8	59	5.9 / 2.1	135	8 / 0.086 / 2.18	6 / 0.086 / 2.18

Recommended guidewire: 0.014 in (0.36 mm)



- 2  
Balloon Dilatation Catheters
- 3  
Specialty Balloons
- 4  
Catheters
- 5  
Accessories
- 6  
Stents
- 7  
Drug-Eluting Technologies