



AESCULAP® POWER SYSTEMS

CATALOG 2022

ELAN 4 and Acculan 4

AESCULAP® ELAN 4

PRECISION. PERFORMANCE. PLEASURE.

With its technical concept ELAN 4 addresses important requirements of the surgeon and all parties involved in handling power systems in the hospital.

- Convenient and reliable grip
- Powerful, precise and reliable work
- Easy to set up and simple to use



ELAN 4 electro BASIC COMPONENTS


GA800

ELAN 4 electro CONTROL UNIT

with irrigation pump; touch screen;
power cord not included
pump flow 0 - 65 ml/min

100 - 240 V AC; 50/60 Hz;
380 x 330 x 200 mm; 9,500 g

consisting parts:
GD412804 bottle holder



GA808

ELAN 4 electro FOOT CONTROL

protection type IPX8; anaesthetic
proof class AP

5.0 m;
245 x 144 x 76 mm; 2,240 g



GA810

ELAN 4 electro WIRELESS FOOT CONTROL

protection type IPX8; anaesthetic
proof class AP

271 x 266 x 85 mm; 3,800 g

Requires control unit with software
version V3.00.



GA806

GB073R

ELAN 4 electro MOTOR CABLE

for use with ELAN 4 electro applied
parts

4.0 m; 270 g

consisting parts:
GB600860 Sterilit® spray adapter
ELAN 4 electro



AESCULAP® ELAN 4

ELAN 4 electro

BASIC COMPONENTS, PERFORATOR DRIVER & LOWSPEED MOTOR

GA804

GB073R

ELAN 4 electro MOTOR CABLE WITH HAND SWITCH

for use with ELAN 4 electro applied parts

4.0 m; 374 g

consisting parts:
GB600860 Sterilit® spray adapter
ELAN 4 electro



GA805

GB073R

ELAN 4 electro MOTOR CABLE WITH HAND CONTROL

for use with ELAN 4 electro applied parts

4.0 m; 369 g

consisting parts:
GB600860 Sterilit® spray adapter
ELAN 4 electro



GA822

GB076R

ELAN 4 electro PERFORATOR DRIVER

for use with cranial perforators with Hudson shank

180 W; 0 - 1,200 rpm; 250 Ncm;
166 x 35 mm; 410 g



GA824


GB077R

ELAN 4 electro LOWSPEED MOTOR

with intra coupling

180 W; 0 - 40,000* rpm;
10 Ncm; 123 x 23 mm; 150 g

* Requires control unit with software version V3.00. With lower versions the motor runs max. 20,000 rpm.



ELAN 4 electro CRANIOTOME & ATTACHMENTS

GA849

GB085R

ELAN 4 electro CRANIOTOME AND MULTIFUNCTION HANDPIECE

2-ring coding; for use with
dura guards and holding sleeves
GB941R - GB947R; spray nozzle
GA259SU (not included)

140 W; ↻ ↺; 0 - 80,000 rpm; 2.2 Ncm;
156 x 17 mm; 📖 121 g



GB941R

GB719R

ELAN 4 DURA GUARD PEDIATRIC

fixed; 1 dot coding;
for craniotome and multifunction
handpieces GA749, GA849



GB942R

GB719R

ELAN 4 DURA GUARD STANDARD

fixed; 2 dot coding;
for craniotome and multifunction
handpieces GA749, GA849



GB943R

GB719R

ELAN 4 DURA GUARD LONG

fixed; 3 dot coding;
for craniotome and multifunction
handpieces GA749, GA849



GB945R

GB719R

ELAN 4 HOLDING SLEEVE

for craniotome and multifunction
handpieces GA749, GA849



GB947R

GB719R

ELAN 4 DURA GUARD STANDARD

turnable; 2 dot coding;
for craniotome and multifunction
handpieces GA749, GA849



AESCULAP® ELAN 4

ELAN 4 electro

1-RING HIGHSPEED HANDPIECES

GA861

GB083R

ELAN 4 electro 1-RING HANDPIECE
L4

1-ring coding;
spray nozzle GA261SU (not included)
140 W; ↗ ↘; 0 - 80,000 rpm; 2.2 Ncm;
152 x 16 mm; 80 g



GA862

GB083R

ELAN 4 electro 1-RING HANDPIECE
L7

1-ring coding;
spray nozzle GA262SU (not included)
140 W; ↗ ↘; 0 - 80,000 rpm; 2.2 Ncm;
177 x 16 mm; 83 g



GA863

GB083R

ELAN 4 electro 1-RING HANDPIECE
L10

1-ring coding;
spray nozzle GA263SU (not included)
140 W; ↗ ↘; 0 - 80,000 rpm; 2.2 Ncm;
206 x 16 mm; 87 g



GA864

GB083R

ELAN 4 electro 1-RING HANDPIECE
L13

1-ring coding;
spray nozzle GA264SU (not included)
140 W; ↗ ↘; 0 - 80,000 rpm; 2.2 Ncm;
237 x 16 mm; 91 g



AESCLAP® ELAN 4

ELAN 4 electro


SAWS



GA836

GB086R

ELAN 4 electro MICRO SAGITTAL SAW
spray nozzle GA258SU (not included)
180 W; 0 - 20,000 rpm; 161 x 23 mm;
 \sphericalangle 7°; \updownarrow 190 g



GA837

GB087R

ELAN 4 electro
MICRO RECIPROCATING SAW
spray nozzle GA258SU (not included)
180 W; 0 - 20,000 rpm; 184 x 23 mm;
 \updownarrow 1.3 mm; \updownarrow 198 g



GA831

GB088R

ELAN 4 electro SAGITTAL SAW
spray nozzle GA258SU (not included)
180 W; 0 - 20,000 rpm; 129 x 23 mm;
 \sphericalangle 4.5°; \updownarrow 220 g



GA832

GB087R


ELAN 4 electro RECIPROCATING SAW
spray nozzle GA258SU (not included)
180 W; 0 - 20,000 rpm; 194 x 23 mm;
 \updownarrow 2.5 mm; \updownarrow 238 g



GA833

GB088R

ELAN 4 electro TRANSVERSAL SAW
spray nozzle GA258SU (not included)
180 W; \sphericalangle 10°; 0 - 20,000 rpm; 10 Ncm;
130 x 23 mm; \sphericalangle 10°; \updownarrow 185 g
to be used with saw blade shafts
GP580R, GP582R, GP586R and GP590R



AESCLAP[®] ELAN 4

ELAN 4 ACCESSORIES

GB861P

GB727R

ELAN 4 LENGTH REDUCTION SLEEVE
for use with 1-ring handpieces



GB863P

GB727R

ELAN 4 LENGTH REDUCTION SLEEVE
for use with 2-ring handpieces



GB862P

GB727R

ELAN 4 MIS SLEEVE
for use with 1-ring handpieces



GB866R

ELAN 4 LATERAL BURR GUARD
for use with 1-ring handpieces and
burrs up to 4.5 mm diameter



TE780

POWER CORD
Europe,  1.5 m

TE730

POWER CORD
Europe,  5.0 m

TE734

POWER CORD
UK,  5.0 m

TE735


POWER CORD
USA, Canada, Japan,  3.5 m



GA395SU

sterile  

ELAN 4 electro DISPOSABLE TUBE SET



ELAN 4 SPRAY NOZZELS

sterile 

GA259SU

for craniotome attachments

GA261SU

for 1-ring handpiece L4

GA262SU

for 1-ring handpiece L7

GA263SU

for 1-ring handpiece L10

GA264SU

for 1-ring handpiece L13

GA265SU

for 2-ring handpiece L7

GA266SU

for 2-ring handpiece L10

GA267SU

for 2-ring handpiece L13

GA268SU

for 2-ring handpiece L17

GA269SU

for 2-ring handpiece L22

GB796SU

for MIS handpiece shafts L10

GB797SU

for MIS handpiece shafts L13



GA258SU

for saws



AESULAP[®] MECHANICAL PROCESSING

ECCOS[®] SETS

ELAN 4

GB261R

ELAN 4 electro ECCOS[®] SET

 56 mm

consisting parts:

GB073R

GB076R


GB083R

GB085R

GB719R

GB072R

JF212R



GB262R

ELAN 4 electro ECCOS[®] SET FOR DRILL GA844

 106 mm

consisting parts:

GB072R

GB073R

GB482R


GB715R

GB716R




























2x GB690R

JG342

JG344









ECCOS®
ELAN 4 electro

<p>GB073R</p> <p>ECCOS® HOLDER for motor cable GA804, GA805, GA806 (set consisting of 2 pieces)</p>  56 mm	<p>GB076R</p> <p>ECCOS® HOLDER for perforator driver GA822</p>  56 mm	<p>GB077R</p> <p>ECCOS® HOLDER for intra motor GA824</p>  56 mm
 	 	 
<p>GB082R</p> <p>ECCOS® HOLDER for 1 handpiece</p>  56 mm	<p>GB083R</p> <p>ECCOS® HOLDER for 3 handpieces</p>  56 mm	<p>GB085R</p> <p>ECCOS® HOLDER for craniotome GA849</p>  56 mm
 	 	 
<p>GB086R</p> <p>ECCOS® HOLDER for micro sagittal saw GA836</p>  56 mm	<p>GB087R</p> <p>ECCOS® HOLDER for reciprocating saws GA832, GA837</p>  56 mm	<p>GB088R</p> <p>ECCOS® HOLDER for saws GA831, GA833</p>  56 mm
 	 	 

AESULAP[®] MECHANICAL PROCESSING

ECCOS[®]

ELAN 4 electro

GB482R	GB715R	GB716R
<p>ECCOS[®] HOLDER for ELAN 4 electro drill GA844</p>  106 mm	<p>ECCOS[®] HOLDER for 1 attachment GB870R, GB871R, GB872R, GB878R, GB881R, GB884R, GB892R, GB896R, GB897R</p>  106 mm	<p>ECCOS[®] HOLDER for 3 attachments GB870R, GB871R, GB872R, GB878R, GB881R, GB884R, GB892R, GB896R, GB897R</p>  106 mm
 	 	 

ECCOS® & RINSING DEVICES


ELAN 4 electro

GB084R

ECCOS® HOLDER

for MIS handpiece GA860

 56 mm



GB723R

ECCOS® HOLDER

for MIS handpiece shafts GB920R - GB927R

 56 mm




GB692R

RINSING DEVICE

for applied parts, incl. ceramic filter GB693

(not suitable for autoclaving)





GB679R

RINSING DEVICE

for 4 MIS shafts, incl. ceramic filter GB693

(not suitable for autoclaving)



GB698R

RINSING ADAPTER

for applied parts

(for manual cleaning)



GB699R

RINSING ADAPTER

for MIS shafts

























(for manual cleaning)



AESCULAP® MECHANICAL PROCESSING

ECCOS®

ELAN 4 ACCESSORIES

<p>GB718R</p>	<p>GB719R</p>	<p>GB720R</p>
<p>ECCOS® HOLDER for 12 1- and 2-ring burrs  56 mm</p>	<p>ECCOS® HOLDER for 3 craniotome attachments  56 mm</p>	<p>ECCOS® UNIVERSAL HOLDER for saw blades  56 mm</p>
		
		
<p>GB729R</p>	<p>GB730R</p>	<p>GB072R</p>
<p>ECCOS® HOLDER for 2 reciprocating saw shafts  56 mm</p>	<p>ECCOS® HOLDER for 2 transversal saw shafts  56 mm</p>	<p>ECCOS® CABLE/HOSE POS. GUIDE for cable and hoses (set of 4 pieces)  56 mm</p>
		
		
<p>GB727R</p>	<p>GB572R</p>	
<p>ECCOS® HOLDER for handpiece holding-sleeves GB861P, GB862P, GB863P  56 mm</p>	<p>ECCOS® FIXATION for GB300R - GB304R  56 mm</p>	
		
		

AESULAP® BURRS & BLADES













ELAN 4 BURRS – 1-RING











ELAN 4 Fräser – 1-Ring







Fraises ELAN 4 – 1 anneau






Fresa ELAN 4 – 1 anillos

Fresa ELAN 4 – 1 anelli


Ø	0.6 mm	0.8 mm	1.0 mm	1.0 mm	1.4 mm	1.4 mm	1.8 mm	1.8 mm	2.3 mm	2.3 mm		
												sterile
I	GP111R	GP112R	GP113R	GP114R	GP115R	GP116R	GP117R	GP118R	GP120R	GP121R		Rosen












Ø	2.7 mm	2.7 mm	3.0 mm	3.5 mm	4.0 mm	4.5 mm	5.0 mm	6.0 mm				
												sterile
I	GP122R	GP123R	GP124R	GP125R	GP126R	GP127R	GP128R	GP129R				Rosen

Ø	3.0 mm	4.0 mm	5.0 mm	6.0 mm	
Soft cut					
					
I	GP133R	GP134R	GP135R	GP136R	Rosen

Ø	4.0 mm	5.0 mm	6.0 mm	7.0 mm	
					
I	GP184R	GP185R	GP186R	GP187R	










Twin-Cut burr
Twin-Cut Fräser
Fraise Twin-Cut
Fresa Twin-Cut
Fresa Twin-Cut



Ø	0.6 mm	0.8 mm	1.0 mm	1.0 mm	1.4 mm	1.4 mm	1.8 mm	1.8 mm	2.3 mm	2.3 mm		
												
I	GP141R	GP142R	GP143R	GP144R	GP145R	GP146R	GP147R	GP148R	GP149R	GP150R		

 sterile

Diamond burr
 Diamantfräser
 Fraise à diamant
 Fresa diamante
 Fresa a diamante

Ø	2.7 mm	2.7 mm	3.0 mm	3.5 mm	4.0 mm	4.5 mm	5.0 mm	6.0 mm				
												
I	GP151R	GP152R	GP153R	GP154R	GP155R	GP156R	GP157R	GP158R				

 sterile

Diamond burr
 Diamantfräser
 Fraise à diamant
 Fresa diamante
 Fresa a diamante

ELAN 4



AESULAP® BURRS & BLADES












ELAN 4 BURRS – 1-RING








ELAN 4 Fräser – 1-Ring


Fraises ELAN 4 – 1 anneau

Fresa ELAN 4 – 1 anillos

Fresa ELAN 4 – 1 anelli

Ø	2.0 mm	2.0 mm	2.3 mm	2.3 mm	3.0 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm			
												
I	GP161R	GP162R	GP163R	GP164R	GP165R	GP166R	GP168R	GP169R	GP170R			
												 sterile Diamond burr coarse, Diamantfräser grob, Fraise diamantée rugueuse, Fresa diamantada de grano grueso, Frese a diamante grosse

Ø	3.0 mm	4.0 mm	5.0 mm	6.0 mm	7.0 mm							
												
I	GP173R	GP174R	GP175R	GP176R	GP177R							
												 sterile Diamond burr extra coarse, Diamantfräser extra grob, Fraise diamantée très rugueuse, Fresa diamantada de grano extra grueso, Frese a diamante extra-grosse



Ø	1.5 mm	2.0 mm	2.5 mm	3.0 mm	3.0 mm
US Type					
I	GP201R	GP202R	GP203R	GP204R	GP205R

sterile

Neuro cutter
Olivenfräser
Fraise à olive
Fresa oliva
Fresa ad oliva

Ø	2.0 mm	2.5 mm	3.0 mm
AESCLAP® Type			
I	GP208R	GP209R	GP210R

sterile

Neuro cutter
Olivenfräser
Fraise à olive
Fresa oliva
Fresa ad oliva

Ø	1.5 mm	1.5 mm	2.0 mm	2.0 mm	2.5 mm	2.5 mm	3.0 mm	3.0 mm	4.0 mm	4.0 mm
I	GP211R	GP212R	GP213R	GP214R	GP215R	GP216R	GP217R	GP218R	GP219R	GP220R

sterile

Neuro cutter diamond
Olivenfräser diamantiert
Fraise à olive diamantée
Fresa oliva diamantada
Fresa ad oliva diamant

Ø	4.0 mm	5.0 mm	6.0 mm
Standard			
I	GP224R	GP225R	GP226R

sterile

Barrel burr
Walzenfräser
Fraise cylindrique
Fresa cilíndrica
Fresa tangenziale

Ø	4.0 mm	5.0 mm	6.0 mm
Soft cut			
I	GP228R	GP229R	GP230R

sterile

Barrel burr
Walzenfräser
Fraise cylindrique
Fresa cilíndrica
Fresa tangenziale

Ø	4.0 mm	5.0 mm	6.0 mm
I	GP234R	GP235R	GP236R

sterile

Cone burr
Kegelfräser
Fraise conique
Fresa cónica
Fresa conica



AESCULAP® BURRS & BLADES



ELAN 4 BURRS – 1-RING

ELAN 4 Fräser – 1-Ring


Fraises ELAN 4 – 1 anneau

Fresa ELAN 4 – 1 anillos


Fresa ELAN 4 – 1 anelli




Acorn burr
Flammenfräser
Fraise flamme
Fresa flama
Fresa di flamma




Oval burr
Ellipsenfräser
Fraise elliptique
Fresa elíptica
Fresa elíptica




Reverse taper burr coarse diamond
Osteophytenfräser
Fraise pour ostéophytes
Fresa per resezione osteofiti
Fresa de diamante para osteofitos




Rosen




Neuro cutter
Olivenfräser
Fraise à olive
Fresa oliva
Fresa ad oliva



Lindemann




Twist drill
Spiralbohrer
Fraise spirale
Broca espiral
Punta elicoidale



Pin cutter
Pinfräser
Fraise pointue
Fresa punta
Fresa puntuale

ELAN 4 BURRS – 2-RING

ELAN 4 Fräser – 2-Ring
 Fraises ELAN 4 – 2 anneaux
 Fresa ELAN 4 – 2 anillos
 Fresa ELAN 4 – 2 anelli



Ø	2.3 mm	3.0 mm	3.5 mm	4.0 mm	4.5 mm	5.0 mm	6.0 mm	7.0 mm	8.0 mm			
II	GP301R	GP302R	GP303R	GP304R	GP305R	GP306R	GP307R	GP308R	GP309R			1 sterile Rosen

Ø	2.3 mm	3.0 mm	3.5 mm	4.0 mm	4.5 mm	5.0 mm	6.0 mm					
II	GP311R	GP312R	GP313R	GP314R	GP315R	GP316R	GP317R					1 sterile Diamond burr Diamantfräser Fraise à diamant Fresa diamante Fresa a diamante

Ø	2.3 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm							
II	GP321R	GP322R	GP323R	GP324R	GP325R							1 sterile Diamond burr coarse, Diamantfräser grob, Fraise diamantée rugueuse, Fresa diamantada de grano grueso, Frese a diamante grosse

Ø	4.0 mm	5.0 mm	6.0 mm									
II	GP328R	GP329R	GP330R									1 sterile Diamond burr extra coarse, Diamantfräser extra grob, Fraise dia- mantée très rugueuse, Fresa diamantada de grano extra grueso, Frese a diamante extra-grosse

ELAN 4

AESCULAP® BURRS & BLADES


ELAN 4 BURRS – 2-RING

ELAN 4 Fräser – 2-Ring


Fraises ELAN 4 – 2 anneaux

Fresa ELAN 4 – 2 anillos


Fresa ELAN 4 – 2 anelli





Neuro cutter
Olivenfräser
Fraise à olive
Fresa oliva
Fresa ad oliva





Neuro cutter
Olivenfräser
Fraise à olive
Fresa oliva
Fresa ad oliva




Neuro cutter diamond
Olivenfräser diamantiert
Fraise à olive diamantée
Fresa oliva diamantada
Fresa ad oliva diamant




Barrel burr
Walzenfräser
Fraise cylindrique
Fresa cilíndrica
Fresa tangenziale





Barrel burr
Walzenfräser
Fraise cylindrique
Fresa cilíndrica
Fresa tangenziale




Cone burr
Kegelfräser
Fraise conique
Fresa cónica
Fresa conica





Acorn burr
Flammenfräser
Fraise flamme
Fresa flama
Fresa di flamma



Oval burr
Ellipsenfräser
Fraise elliptique
Fresa elíptica
Fresa elíptica




Reverse taper burr coarse diamond
Osteophytenfräser
Fraise pour ostéophytes
Fresa per resezione osteofiti
Fresa de diamante para osteofitos




1 sterile

Craniotome cutter – spiral type
 Kraniotomfräser – spiralgig
 Fraise craniotome – en spirale
 Fresa cranéotomo – helicoidal
 Fresa di craniotomo – a spirale




1 sterile

Craniotome cutter – straight type
 Kraniotomfräser – geradverzahnt
 Fraise craniotome – dents droites
 Fresa cranéotomo – dientes rectos
 Fresa di craniotomo – denti dritti




1 sterile

Lindemann




1 sterile

Twist drill
 Spiralbohrer
 Fraise spirale
 Broca espiral
 Punta elicoidale




1 sterile

Twist drill
 Spiralbohrer
 Fraise spirale
 Broca espiral
 Punta elicoidale



1 sterile

Pin cutter
 Pinfräser
 Fraise pointue
 Fresa punta
 Fresa puntuale



1 sterile

Twin-Cut burr
 Twin-Cut Fräser
 Fraise Twin-Cut
 Fresa Twin-Cut
 Fresa Twin-Cut

AESCULAP® BURRS & BLADES


ELAN 4 BURRS – 2-RING

ELAN 4 Fräser – 2-Ring

Fraises ELAN 4 – 2 anneaux


Fresa ELAN 4 – 2 anillos

Fresa ELAN 4 – 2 anelli



Tungsten Carbide Burrs

Ø	2.3 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm	7.0 mm	8.0 mm
II	GP301TC	GP302TC	GP304TC	GP306TC	GP307TC	GP308TC	GP309TC



Rosen

Ø	1.2 mm	1.6 mm
II	GP395TC	GP396TC



Taper burr
Konusfräser
Fraise conique
Fresa cónica
Fresa conica

Ø	3.0 mm
II	GP367TC



Neuro cutter
Olivenfräser
Fraise à olive
Fresa oliva
Fresa ad oliva

Burrs for metal dissection

Ø	2.3 mm
II (2)	GP399TC-SU



Side cutter
Zylinderfräser
Fraise cylindrique
Fresa cilíndrica
Fresa cilíndrica

Ø	25 mm
II (2)	GP398SU



Diamond disc
Diamant-Trennscheibe
Disque cou pant diamant
Disco cortante de diamante
Disco tagliante diamante

ELAN 4 BURRS – MIS

ELAN 4 Fräser – MIS
Fraises ELAN 4 – MIS
Fresa ELAN 4 – MIS
Fresa ELAN 4 – MIS

GA860



GB920R



GB921R




GB922R



L10

GB925R



GB926R



GB927R






L13

AESCU LAP[®] BURRS & BLADES




ELAN 4 BURRS – MIS

ELAN 4 Fräser – MIS
 Fraises ELAN 4 – MIS
 Fresa ELAN 4 – MIS
 Fresa ELAN 4 – MIS

Ø	3.0 mm	4.0 mm		
				
L10 (2)	GP422SU	GP423SU		
L13 (2)	GP442SU	GP443SU		







Rosen

Ø	2.0 mm	3.0 mm		
				
L10 (2)	GP420SU	GP421SU		
L13 (2)	GP440SU	GP441SU		







Neuro cutter
 Olivenfräser
 Fraise à olive
 Fresa oliva
 Fresa ad oliva

Ø	3.0 mm	4.0 mm	5.0 mm	
				
L10 (2)	GP432SU	GP433SU	GP434SU	
L13 (2)	GP452SU	GP453SU	GP454SU	








Diamond burr
 Diamantfräser
 Fraise diamantée
 Fresa diamantada
 Frese a diamante

Ø	3.0 mm	4.0 mm	5.0 mm	
				
L10 (2)	GP435SU	GP436SU	GP437SU	
L13 (2)	GP455SU	GP456SU	GP457SU	







Diamond burr extra coarse
 Diamantfräser extra grob
 Fraise diamantée très ruguese
 Fresa diamantada de grano extra grueso
 Frese a diamante extra-grosse

Ø	2.0 mm	3.0 mm	4.0 mm	5.0 mm
				
L10 (2)	GP424SU	GP425SU	GP426SU	GP427SU
L13 (2)	GP444SU	GP445SU	GP446SU	GP447SU



Neuro cutter diamond
 Olivenfräser diamantiert
 Fraise à olive diamantée
 Fresa oliva diamantada
 Fresa ad oliva diamante

Ø	3.0 mm	4.0 mm	5.0 mm	
				
L10 (2)	GP428SU	GP429SU	GP430SU	
L13 (2)	GP448SU	GP449SU	GP450SU	




Neuro cutter diamond extra coarse
 Olivenfräser diamantiert extra grob
 Fraise à olive diamantée très ruguese
 Fresa oliva diamantada de grano extra grueso
 Fresa ad oliva diamante extra-grosse

AESCULAP® BURRS & BLADES

SKULL PERFORATORS, CRANIOTOME BURRS


Schädeltrepane, Kraniotomfräser
 Trépans crâniens, fraises de craniotome
 Trépanos y fresas para craneótomo
 Perforatori cranici, frese per craniotomia



Ø/mm


GB300R	6/9	4.0 mm*
GB302R	9/12	5.0 mm*
GB304R	12/15	5.4 mm*

Hudson



Ø/mm

TE561	6/9	4.0 mm*
TE562	9/12	5.0 mm*
TE563	12/15	5.4 mm*



* Minimal cranial bone thickness
 Mindestdicke des Schädelknochens
 Épaisseur minimal de l'os crânien
 Grosor mínimo del hueso craneal
 Spessore minimo dell'osso cranico


* Minimal cranial bone thickness
 Mindestdicke des Schädelknochens
 Épaisseur minimal de l'os crânien
 Grosor mínimo del hueso craneal
 Spessore minimo dell'osso cranico




Ø

GC303R	20 mm	1.6 mm
GC307R	20 mm	2.0 mm

GC303R



GC307R



» GB265R, GB268R



GC304R


» GB292

GC304R



CARE AND MAINTENANCE

Pflege und Aufbereitung
 Entretien et maintenance
 Conservación y mantenimiento
 Cura e manutenzione



GB600860
 ELAN 4 electro



GB600870
 ELAN 4 MIS shafts,
 microspeed uni XS



GB600880
 Acculan 4 small drill



GB600840
 Acculan 4, Acculan 3Ti



GB600820
 HiLAN, HiLAN XS, Hi-Line
 (XS, XXS), Microspeed EC,
 microspeed uni, microLAN,
 miniLAN



GB600830
 mini-Line



GB600810
 micro-Line



GB600 
 Sterilit® Power Systems
 consisting parts:
 3x GB600890 spray tube

etc.




Project Overview

Description of Similarity:

equal – the same article

like – nearly equal but with some slight differences

maybe – only as an alternative

Pos No.	Item Code	Description / Comment	Qty.
360	GB600	<p>(GB600) Sterilit Power Systems</p>  <p>Customer Reference: 1.30</p>	1
370	JK442	<p>(JK442) Bottom, without base perforation, Standard 1/1, outer length: 592 mm, outer width: 274 mm, outer height: 135 mm, inner length: 544 mm, inner width: 258 mm, inner height: 120 mm</p>  <p>Customer Reference: 1.31</p>	1
380	JF221R	<p>(JF221R) Perforated basket, Standard 1/1, 540 x 253 x 36 mm, with feet</p>  <p>Customer Reference: 1.31</p>	1





Project Overview

Description of Similarity:

equal – the same article

like – nearly equal but with some slight differences

maybe – only as an alternative

Pos No.	Item Code	Description / Comment	Qty.
390	JK489	<p>(JK489) Inner Lid, basis, Standard 1/1, silver, outer length: 582 mm, outer width: 291 mm, outer heigth: 36 mm</p>  <p>Customer Reference: 1.31</p>	1
400	JG785B	<p>(JG785B) Accessories Container, identification label, standard, aluminium, red, outer length: 56 mm, outer width: 18 mm</p>  <p>Customer Reference: 1.31</p>	4
410	JK090	<p>(JK090) Reusable Filter, package of 10 pieces, diam. 190 mm</p>  <p>Customer Reference: 1.31</p>	1
420	JG744	<p>(JG744) Indicator Label, package of 1000 pieces, outer length: 80 mm</p>  <p>Customer Reference: 1.31</p>	1