


UNDERPAD
REF.NO:800.003C-1.02
PRODUCT DATA SHEET

medcare

**PRODUCT DATA SHEET
UNDERPAD-NON STERILE
800.003C-1.02**

1) PRODUCT COMPONENTS



ITEM:1
Underpad
90x180 cm
Absorbent Area
50x80 cm

2) DETAILS OF PRODUCT COMPONENTS

ITEM	QTY	PRODUCT	PRODUCTION TYPE	DIMENSION	COLOUR
1	1	Underpad	Total 60 gr (±3)	90x180 cm (±3)	White
			Backsheet – PE Film 18 gsm		
			Topsheet 8 gsm		
			Pulp-Absorbent Area-25.52 gr	50x80 cm (±3)	

3) TEST RESULT (ISO 11948-1)

ABSORPTION CAPACITY (ML)*	891
----------------------------------	-----

**PRODUCT DATA SHEET
UNDERPAD-NON STERILE
800.003C-1.02**

4) GENERAL INFORMATION

Product Description	UNDERPAD Non-Sterile
Certificates	EN ISO 13485:2016
Standards	EN ISO 14971:2019 EN ISO 15223-1:2021
Intended Use	Underpad product is used as a mattress protector in health institutions such as hospitals, polyclinics, intensive care units, and home care patients.
Instructions For Usage	Cut the products bag open with a scissor from the bottom side of the packaging. Doing so will present you a better place to hold onto the pad as you take it out of the packaging. Take the underpad out and place it in a folded-up state, onto the surface you'll be using it on. Unfold the pad. Take the edge of the product and unfurl it outward, away from you. Allow it to unfold completely. Place the pad down on the surface, with the non plastic side up. The soft side absorbs the urine, while the plastic-looking side will prevent any liquid from going through and onto the surface. The person will lie down on the softer side.

5) PACKAGING

Packaging

Bag in 30 pcs