



CATALOGUE 2021

CARDIAC CATHETERS
AND CANNULAE



Retrograde cardioplegia cannulae



Size	Balloon diameter	Reference
14 Fr	18 mm	RCC01


PRODUCT INFORMATION

Andocor retrograde cardioplegia cannulae are intended to be connected to the cardioplegia line in order to infuse cardioplegic solution and blood into the patient's coronary sinus during open heart procedures.

Q Length: 31 cm
Qty/box: 20 Pcs

+ Self-inflating
Silicone balloon
Malleable stylet

RETROGRADE CARDIOPLEGIA CANNULA PRESSURE DROP VS FLOW



Ostial perfusion cannulae

CORONARY ARTERY OSTIAL CANNULA, BASKET TIP WITH FLANGE



Tip size	45° angled	90° angled
10 Fr	OC1045	OC1090
12 Fr	OC1245	OC1290
14 Fr	OC1445	OC1490

PRODUCT INFORMATION

Andocor ostial perfusion cannulae are intended for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.

Q Connection: Female luer lock
 Length: 21 cm
 Qty/box: 20 Pcs

+ Soft basket tip
 Different configurations
 Stainless steel shaft

CORONARY ARTERY OSTIAL CANNULA PRESSURE DROP VS FLOW

