


Esophageal Balloon Dilator

Used to dilate strictures of the gastrointestinal tract including strictures of the esophagus through the scope

High Radial Vector Force
Atraumatic Tip
Two radiopaque markers inside the balloon



Art. No.	Shaft Ø (mm)	Guidewire	Balloon Ø (mm)	Balloon Length (mm)	Working Length (cm)
EBD-0630	2.5	0.035"	6	30	180
EBD-0830	2.5	0.035"	8	30	180
EBD-1040	2.5	0.035"	10	40	180
EBD-1280	2.5	0.035"	12	80	180
EBD-1480	2.5	0.035"	14	80	180
EBD-1680	2.5	0.035"	16	80	180
EBD-18230	2.5	0.035"	18	80	230
EBD-20230	2.5	0.035"	20	80	230


Other sizes on request

Class IIa device – for single use only

Achalasia-Cardia Balloon Dilator

Used for dilation of the lower esophageal sphincter in patient with achalasia (stomach entry doesn't open properly)

High Radial Vector Force
Atraumatic Tip
Two radiopaque markers inside the balloon




Art. No.	Shaft Ø (mm)	Guidewire	Balloon Ø (mm)	Balloon Length (mm)	Working Length (cm)
SABD-25	5	0.038"	25	80	110
SABD-30	5	0.038"	30	80	110
SABD-35	5	0.038"	35	80	110
SABD-40	5	0.038"	40	80	110

Class IIa device – for single use only

NEG Nitinol Guidewire

Used to gain access through difficult biliary and pancreatic strictures and to facilitate exchange of compatible wire guided accessories to the desired biliary or pancreatic ducts.



Art. No.	Guidewire Ø / Tip	Length (cm)	Art. No.	Guidewire Ø / Tip	Length (cm)
NEG-025-S	0.025" Straight	450	NEG-DBE-025-S	0.025" Straight*	525
NEG-025-A	0.025" Angled	450	NEG-DBE-025-A	0.025" Angled*	525
NEG-025-J	0.025" J Tip	450	NEG-DBE-025-J	0.025" J Tip*	525
NEG-032-S	0.032" Straight	450	NEG-DBE-032-S	0.032" Straight*	525
NEG-032-A	0.032" Angled	450	NEG-DBE-032-A	0.032" Angled*	525
NEG-032-J	0.032" J Tip	450	NEG-DBE-032-J	0.032" J Tip*	525
NEG-035-S	0.035" Straight	450	NEG-DBE-035-S	0.035" Straight*	525
NEG-035-A	0.035" Angled	450	NEG-DBE-035-A	0.035" Angled*	525
NEG-035-J	0.035" J Tip	450	NEG-DBE-035-J	0.035" J Tip*	525
NEG-038-S	0.038" Straight	450	NEG-DBE-038-S	0.038" Straight*	525
NEG-038-A	0.038" Angled	450	NEG-DBE-038-A	0.038" Angled*	525
NEG-038-J	0.038" J Tip	450	NEG-DBE-038-J	0.038" J Tip*	525

For other length please add to Art. No. -150 or -260


*for Double Balloon Enteroscope

Class IIa device – sterile – for single use only

Hydra NEG Nitinol Guidewire

Used to gain access through difficult biliary and pancreatic strictures and to facilitate exchange of compatible wire guided accessories to the desired biliary or pancreatic ducts

With Extra Filiform Wire



Art. No.	Guidewire Ø	Tip	Length (cm)
HYNEG-032-S	0.032"	Straight	450
HYNEG-032-A	0.032"	Angled	450
HYNEG-032-J	0.032"	J Tip	450
HYNEG-035-S	0.035"	Straight	450
HYNEG-035-A	0.035"	Angled	450
HYNEG-035-J	0.035"	J Tip	450
HYNEG-038-S	0.038"	Straight	450
HYNEG-038-A	0.038"	Angled	450
HYNEG-038-J	0.038"	J Tip	450

For other length please add to Art. No. -150 or -260

Class IIa device – sterile – for single use only