



JF SERIES

- Black PEEK plastic feet for stacking
- Fold-down handles on the side
- Mat surface
- 4 mm base perforation / 2 mm side perforation
- Can be combined with Aesculap InOS®
- Sloped edges


1/2

Top basket


Basket in combination with:

		Alone	JF111R	JF111R with lid*	JF112R	JF112R with lid*	JF113R	JF113R with lid*	JF114R	JF114R with lid*	
Bottom basket	JF111R	Container bottom	JK340	JK341	JK341	JK341	JK342	JK342	JK344	JK344	
	JF111R with lid*		JN340	JN341	JN341	JN342	JN342	JN344	JN344	JN344	JN344
	JF112R		JK340	JK341	JK342	JK342	JK342	JK344	JK344	JK344	JK346
	JF112R with lid*		JN341	JN342	JN342	JN344	JN344				
	JF113R		JK341	JK342	JK344	JK344	JK344	JK346	JK346	JK346	JK346
	JF113R with lid*		JK341	JN344	JN344	JN344	JN344	JK344	JK346	JK346	JK346
	JF114R		JK342	JK344	JK344	JK344	JK344	JK346	JK346	JK346	JK346
	JF114R with lid*		JN344	JN344	JN344	JN344	JN344	JN346	JN346	JN346	JN346
JF117R lid	JK344	JK344	JK344	JK344	JK344	JK346	JK346	JK346	JK346		

*JF117R lid



- 1 EXAMPLE**
- Basket alone: JF112R (without lid)
 - ➔ suitable container bottom: JK340 or JN341



- 2 EXAMPLE**
- Top basket: JF112R (without lid)
 - Bottom basket: JF113R with lid
 - ➔ suitable container bottom: JK344 or JN344

JK = Solid bottom
(fractionated vacuum process)











JN = Perforated bottom
(gravity process)



NOTE: ■ According to EN 868-8, the load must not exceed a total of 5 kg in order to avoid excessive condensate formation.
 ■ Caution! The inner height and therefore the appropriate basket changes depending on whether JK or JN bottoms are used.









BASKET WITH FEET



		Outside dim. Inside dim.		Outside dim. Inside dim.		Outside dim. Inside dim.		Outside dim. Inside dim.
 36	JF111R	243 x 253 x 36 241 x 251 x 22	JF251R	406 x 253 x 36 404 x 251 x 22	JF211R	485 x 253 x 36 483 x 251 x 22	JF221R	540 x 253 x 36 538 x 251 x 22
 56	JF112R	243 x 253 x 56 241 x 251 x 42	JF252R	406 x 253 x 56 404 x 251 x 42	JF212R	485 x 253 x 56 483 x 251 x 42	JF222R	540 x 253 x 56 538 x 251 x 42
 76	JF113R	243 x 253 x 76 241 x 251 x 62	JF253R	406 x 253 x 76 404 x 251 x 62	JF213R	485 x 253 x 76 483 x 251 x 62	JF223R	540 x 253 x 76 538 x 251 x 62
 106	JF114R	243 x 253 x 106 241 x 251 x 92	JF254R	406 x 253 x 106 404 x 251 x 92	JF214R	485 x 253 x 106 483 x 251 x 92	JF224R	540 x 253 x 106 538 x 251 x 92

BASKET WITHOUT FEET



























		Outside dim. Inside dim.		Outside dim. Inside dim.		Outside dim. Inside dim.		Outside dim. Inside dim.
 24	JG111R	243 x 253 x 24 241 x 251 x 22	JG251R	406 x 253 x 24 404 x 251 x 22	JG211R	485 x 253 x 24 483 x 251 x 22	JG221R	540 x 253 x 24 538 x 251 x 22
 44	JG112R	243 x 253 x 44 241 x 251 x 42	JG252R	406 x 253 x 44 404 x 251 x 42	JG212R	485 x 253 x 44 483 x 251 x 42	JG222R	540 x 253 x 44 538 x 251 x 42
 64	JG113R	243 x 253 x 64 241 x 251 x 62	JG253R	406 x 253 x 64 404 x 251 x 62	JG213R	485 x 253 x 64 483 x 251 x 62	JG223R	540 x 253 x 64 538 x 251 x 62
 94	JG114R	243 x 253 x 94 241 x 251 x 92	JG254R	406 x 253 x 94 404 x 251 x 92	JG214R	485 x 253 x 94 483 x 251 x 92	JG224R	540 x 253 x 94 538 x 251 x 92

AESCULAP® MECHANICAL PROCESSING

ECCOS®

ELAN 4 ACCESSORIES

<p>GB718R</p>	<p>GB719R</p>	<p>GB720R</p>
<p>ECCOS® HOLDER for 12 1- and 2-ring burrs  56 mm</p>	<p>ECCOS® HOLDER for 3 craniotome attachments  56 mm</p>	<p>ECCOS® UNIVERSAL HOLDER for saw blades  56 mm</p>
 	 	 
<p>GB729R</p>	<p>GB730R</p>	<p>GB072R</p>
<p>ECCOS® HOLDER for 2 reciprocating saw shafts  56 mm</p>	<p>ECCOS® HOLDER for 2 transversal saw shafts  56 mm</p>	<p>ECCOS® CABLE/HOSE POS. GUIDE for cable and hoses (set of 4 pieces)  56 mm</p>
 	 	 
<p>GB727R</p>	<p>GB572R</p>	
<p>ECCOS® HOLDER for handpiece holding-sleeves GB861P, GB862P, GB863P  56 mm</p>	<p>ECCOS® FIXATION for GB300R - GB304R  56 mm</p>	
 	 	

ECCOS® SETS

Acculan 4

GB257R

ECCOS® SET

for Acculan 4 small drill GA344

 106 mm

consisting parts:

GB482R

GB487R

GB488R


GB715R

GB716R

2x GB690R

JG342

JG344



GB243R

ECCOS® SET

for 2 Acculan 4 machines

 106 mm

consisting parts:

GB481R

GB487R

GB488R

GB495R


GB496R

GB497R

2x GB690R

JG342

JG344



GB256R

ECCOS® SET

for Acculan 4 dermatome

 106 mm


consisting parts:

GB487R

GB488R

GB498R



















JJ464R



AESCALAP® MECHANICAL PROCESSING

ECCOS®


Acculan 4

<p>GB482R</p> <p>ECCOS® HOLDER for Acculan 4 small drill GA344</p>  106 mm	<p>GB715R</p> <p>ECCOS® HOLDER for 1 attachment GB870R, GB871R, GB872R, GB878R, GB881R, GB884R, GB892R, GB896R, GB897R</p>  106 mm	<p>GB716R</p> <p>ECCOS® HOLDER for 3 attachments GB870R, GB871R, GB872R, GB878R, GB881R, GB884R, GB892R, GB896R, GB897R</p>  106 mm
 	 	 
<p>GB495R</p> <p>ECCOS® HOLDER for GA330, GA331, GA334</p>  106 mm	<p>GB496R</p> <p>ECCOS® HOLDER for GB620R, GB621R, GB623R, GB639R, GB641R, GB645R GB660R, GB663R, GB664R, GB665R, GB657R, GB667R, GB668R, GB669R, GB670R</p>  106 mm	<p>GB497R</p> <p>ECCOS® HOLDER for GB620R, GB621R, GB623R, GB639R, GB641R, GB645R GB660R, GB663R, GB664R, GB665R, GB657R, GB667R, GB668R, GB669R, GB670R</p>  106 mm
 	 	 

GB487R


ECCOS® HOLDER
for GA675

 56 mm

**GB488R**

ECCOS® HOLDER
for GA678

 76 mm

**GB489R**


ECCOS® HOLDER
for GA679

 106 mm

**GB498R**


ECCOS® HOLDER
for GA336, GA337, GA340, GA341

 106 mm

**GB686R**

ECCOS® HOLDER
for GB436R















 56 mm



AESCALAP® MECHANICAL PROCESSING





ECCOS® HOLDER

Acculan 4

<p>GB485R</p>	<p>GB486R</p>	<p>GB688R</p>
<p>ECCOS® HOLDER for GA875  56 mm</p>	<p>ECCOS® HOLDER for GA878  56 mm</p>	<p>ECCOS® HOLDER for BA720R  106 mm</p>
		
		
<p>GB481R</p>	<p>JF148R</p>	<p>JJ589R</p>
<p>ECCOS® HOLDER for saw blades and other accessories  56 mm</p>	<p>MINI TRAY WITH LID outer size: 145 x 44 x 25 mm inner size: 140 x 39 x 18 mm for storage of Acculan 4 Dermatome rod and flaps during automatic cleaning process</p>	<p>SOCKET WRENCH AF7</p>
 		

ECCOS® HOLDER

Universal Tools

<p>GB459R</p> <p>ECCOS® MOUNTING KIT consisting of: socket wrench AF10 and hex-wrench</p>	<p>GB690R</p> <p>ECCOS® HOLDER 5 mm for drilling wire guide GA672207, GA344244</p>	<p>JG303R </p> <p>REPLACEMENT SCREWS-SET screw nut size AF7</p>
		
<p>JG342</p> <p>SILICONE CLAMP 8.5 mm screw JG303R to be ordered separately</p>	<p>JG344</p> <p>SILICONE CLAMP 12 mm screw JG303R to be ordered separately</p>	
